

ENG

Product Data Sheet:

**HUMAN DHN-MA LIPID**

**PEROXIDATION ENZYME ELISA**

Catalogue number:  
**RA19023R**

**For research use only!**

Example Version

1.	DHN-MA LIPID PEROXIDATION ELISA	3
2.	PRECAUTION FOR USE	4
3.	BACKGROUND	5
4.	PRINCIPLE OF THE ASSAY	7
5.	MATERIAL REQUIRED BUT NOT PROVIDED	8
6.	SAMPLE COLLECTION AND PREPARATION	8
7.	REAGENT PREPARATION	9
8.	ASSAY PROCEDURE	11
9.	DATA ANALYSIS	14
10.	ACCEPTABLE RANGE	14
11.	TYPICAL RESULTS	15
12.	ASSAY VALIDATION AND CHARACTERISTICS	16
13.	ASSAY TROUBLE SHOOTING	22
14.	BIBLIOGRAPHY	23
15.	EXPLANATION OF SYMBOLS	24

## 1. DHN-MA LIPID PEROXIDATION ELISA

96 wells

Storage: -20°C

Expiry date: stated on the package

This kit contains:

Reagents	Colour code	Quantity	Form
Antibody Coated Microtiter Strips	Blister with zip	1	-
Conjugate Solution (DHN-MA Tracer)	Green	1	lyophilized
DHN-MA Antiserum	Red	1	lyophilized
DHN-MA Standard	Blue with red septum	2	lyophilized
Quality Control	Green with red septum	2	lyophilized
Dilution Buffer (EIA buffer)	Blue	1	lyophilized
Wash Solution	Silver	1	liquid
Substrate Solution (Ellman's Reagent)	Black with red septum	2	lyophilized
Tween 20	Transparent	1	liquid
Cover Sheet	-		-

Each kit contains sufficient reagents for 96 wells. This allows for the construction of one standard curve in duplicate and the assay of 34 samples in duplicate.

If you want to use the kit in two times, we provide one additional vial of Standard, one of Quality Control and one of Substrate Solution.

## 2. PRECAUTION FOR USE

Users are recommended to carefully read all instructions for use before starting work.

Each time a new pipette tip is used, aspirate a sample or reagent and expel it back into the same vessel. Repeat this operation two or three times before distribution in order to equilibrate the pipette tip.

- For research laboratory use only.
- Not for human diagnostic use.
- Do not pipet liquids by mouth.
- Do not use kit components beyond the expiration date.
- Do not eat, drink or smoke in area in which kit reagents are handled.
- Avoid splashing.

The total amount of reagents contains less than 100 µg of sodium azide. Flush the drains thoroughly to prevent the production of explosive metal azides.

Wearing gloves, laboratory coat and glasses is recommended when assaying kit materials and samples.

### **Temperature:**

Unless otherwise specified, all the experiments are done at room temperature (RT), that is around 20°C. Working at 25°C or more affects the assay and decreases its efficiency.

### 3. BACKGROUND

#### Acetylcholinesterase AChE<sup>®</sup> Technology

Acetylcholinesterase (AChE<sup>®</sup>), the enzymatic label for EIA, has the fastest turnover rate of any enzymatic label. This specific AChE is extracted from the electric organ of the electric eel, *Electrophorus electricus*, and it's capable of massive catalytic turnover during the generation of the electrochemical discharges. The use of AChE as enzymatic label for EIA has been patented by the French academic research Institute CEA [1, 2, 3].

AChE<sup>®</sup> assays are revealed with Substrate Solution (Ellman's reagent), which contains acetylthiocholine as a substrate. The final product of the enzymatic reaction (5-thio-2-nitrobenzoic acid) is bright yellow and can be read at 405-414 nm. AChE<sup>®</sup> offers several advantages compared to enzymes conventionally used in EIAs:

##### Kinetic superiority and high sensitivity

AChE<sup>®</sup> shows true first-order kinetics with a turnover of 64,000 sec<sup>-1</sup>. That is nearly 3 times faster than Horse Radish Peroxidase (HRP) or alkaline phosphatase. AChE<sup>®</sup> allows a greater sensitivity than other labeling enzymes.

##### Low background

non-enzymatic hydrolysis of acetylthiocholine in buffer is essentially absent. So, AChE<sup>®</sup> allows a very low background and an increased signal/noise ratio compared to other substrate of enzymes which is inherently unstable.

##### Wide dynamic range

AChE<sup>®</sup> is a stable enzyme and its activity remains constant for many hours as, unlike other enzymes, its substrate is not suicidal. This permits simultaneous assays of high diluted and very concentrated samples.

##### Versatility

AChE<sup>®</sup> is a completely stable enzyme, unlike peroxidase which is suicidal. Thus, if a plate is accidentally dropped after dispatch of the AChE<sup>®</sup> Substrate Solution (Ellman's Reagent) or if it needs to be revealed again, one only needs to wash the plate, add fresh Substrate Solution and proceed with a new development. Otherwise, the plate can be stored at 4°C with Wash Buffer in wals while waiting for technical advice from the Bioreagent Department.

#### DHN-MA and lipid peroxidation

Reactive Oxygene Species (ROS) play an important role in pathogenicity of several diseases (cardiovascular diseases such as atherosclerosis, cerebral or heart ischemia-reperfusion injury, neurodegenerative diseases, diabetes, infl ammation and cancer) but are also involved in cell signaling. Consequently there is an increasing need in assays to monitor those biomarkers.


For decade, 8-Isoprostane (8-isoPGF<sub>2</sub>α, 8-isoprostaglandin F<sub>2</sub> α, 8-epiPGF<sub>2</sub>α) has been used as the biomarker of lipid peroxidation. This biomarker is diffi cult to measure and needs an extraction whatever the method used. Also very popular are beta-cleavage products of polyunsaturated fatty acids (PUFA), such as alkanes, ketones or aldehydes.

Two well-known aldehydes formed during the lipid peroxidation process, namely malondialdehyde (MDA) and 4-hydroxy-2-nonenal (HNE), have already been used as lipid peroxidation biomarkers for decades [4]. MDA assays using thiobarbituric reaction are now named TBARS assay (ThioBarbituric Acid Reactive Substances assay) due to the lack of specificity of the reaction. Measurement of HNE and protein/HNE adducts in tissues makes HNE more attractive in the field of clinical and experimental studies [8, 9, 10].

Moreover, **1,4-dihydroxynonane mercapturic acid (DHNMA)**, the major urinary metabolite of HNE, is present at a physiological level in rat and human urine. It was demonstrated that DHN-MA measured without extraction in rat urines treated with BrCCl<sub>3</sub>, which induces lipid peroxidation, mainly correlates with MDA and 8-Isoprostane, which were measured concomitantly [7].

HNE and thus DHN-MA are generated from PUFA (alimentation), conclusions on lipid peroxidation should be considered accordingly [5].

DHN-MA Structure



## 4. PRINCIPLE OF THE ASSAY

This Enzyme Immunoassay (EIA) is based on the competition between unlabelled DHN-MA and acetylcholinesterase (AChE) labelled DHN-MA (Tracer) for limited specific rabbit anti-DHN-MA antiserum sites.

The complex rabbit antiserum - DHN-MA (free DHN-MA or Tracer) binds to the mouse monoclonal anti-rabbit antibody coated in the well.

The plate is washed to remove any unbound reagent, and Substrate Solution (Ellman's Reagent - enzymatic substrate for AChE and chromogen) is added to the wells.


AChE tracer acts on the Substrate Solution (Ellman's Reagent) to form a yellow compound that strongly absorbs at 414 nm.

The intensity of the colour, determined by spectrophotometry, is proportional to the amount of tracer bound to the well and is inversely proportional to the amount of free DHN-MA present in the well during the immunological incubation.

The kit has been developed and validated for urine samples.

For any other sample it is the responsibility of the user to check for potential interferences (see our web site or contact our technical support).

The principle of the assay is summarised below:



## 5. MATERIAL REQUIRED BUT NOT PROVIDED

In addition to standard laboratory equipment, the following material is required:

- Precision micropipettes (20 to 1000  $\mu\text{L}$ )
- Multichannel pipette 100  $\mu\text{L}$  or 200  $\mu\text{L}$  and disposable tips
- Spectrophotometer plate reader (405 or 414 nm filter)
- Microplate washer (or washbottles)
- Microplate shaker
- Magnetic stirring bar
- UltraPure water
- Polypropylene tubes

**Water used to prepare all EIA reagents and buffers must be UltraPure** (deionized & free from organic contaminants traces).

Otherwise, organic contamination can significantly affect the enzymatic activity of the tracer Acetylcholinesterase (AChE). Do not use distilled water, HPLC-grade water or sterile water.

UltraPure water may be purchased from BioVendor #S0001.

## 6. SAMPLE COLLECTION AND PREPARATION

Urine samples collected in tubes are to be used immediately or stored at  $-20^{\circ}\text{C}$  for long term storage.

Before use, urine samples must be centrifuged at 1500 g for 5 minutes at  $4^{\circ}\text{C}$ .

Samples are diluted in Dilution Buffer (recommended between 1/20 and 1/200). We advise to assay each sample in duplicate, in at least 3 dilutions. The urine sample minimal dilution is 1/5.

Avoid thawing samples more than three times.



## 7. REAGENT PREPARATION

Each kit contains sufficient reagents for 96 wells. This allows for the construction of one standard curve in duplicate and the assay of 34 samples in duplicate.

If you want to use the kit in two times, we provide one additional vial of Standard, one of Quality Control and one of Substrate Solution.

All reagents need to be brought to room temperature, around 20°C, prior to the assay.

### Dilution Buffer

Reconstitute the vial of Dilution Buffer with 50 mL of UltraPure water.

Allow it to stand 5 minutes until completely dissolved and then mix thoroughly by gentle inversion.

Stability at 4°C: 1 month

### DHN-MA Standard

Reconstitute one vial of DHN-MA Standard with 1 mL of UltraPure water.

Allow it to stand 5 minutes until completely dissolved and then mix thoroughly by gentle inversion.

The concentration of the first standard **S1** is 1000 pg/mL. Prepare seven propylene tubes for the other standards and add 500 µL of Dilution Buffer into each tube. Then prepare the standards by serial dilutions as follows:

Standard	Volume of Standard	Volume of Dilution Buffer	Standard concentration
S1	-	-	1000 pg/mL
S2	500 µL of S1	500 µL	500 pg/mL
S3	500 µL of S2	500 µL	250 pg/mL
S4	500 µL of S3	500 µL	125 pg/mL
S5	500 µL of S4	500 µL	62.5 pg/mL
S6	500 µL of S5	500 µL	31.3 pg/mL
S7	500 µL of S6	500 µL	15.7 pg/mL
S8	500 µL of S7	500 µL	7.8 pg/mL

Stability at 4°C: 48 hours

### DHN-MA Quality Control

Reconstitute one Quality Control vial with 1 mL UltraPure water. Allow it to stand 5 minutes until completely dissolved and then mix thoroughly by gentle inversion.

Stability at +4°C: 48 hours

### DHN-MA Conjugate Solution

Reconstitute one vial with 5 mL of Dilution Buffer. Allow it to stand 5 minutes until completely dissolved and then mix thoroughly by gentle inversion.

Stability at 4°C: 1 week

### DHN-MA Antiserum

Reconstitute the vial with 5 mL of Dilution Buffer. Allow it to stand 5 minutes until completely dissolved and then mix thoroughly by gentle inversion.

Stability at 4°C: 1 week

### Wash Buffer

Dilute 2 mL of concentrated Wash Buffer with 800 mL of UltraPure water. Add 400 µL of Tween 20. Use a magnetic stirring bar to mix the content.

Stability at 4°C: 1 week

### Substrate Solution (Ellman's Reagent)

**5 minutes before use** (development of the plate), reconstitute one vial of Substrate Solution with 50 mL of UltraPure water. The tube content should be thoroughly mixed.

Stability at 4°C and in the dark: 24 hours

Example Version

## 8. ASSAY PROCEDURE

It is recommended to perform the assays in duplicate and to follow the instructions hereafter.

### Plate preparation

Prepare the Wash Buffer as indicated in the reagent preparation section.

Open the plate packet and select the sufficient strips for your assay and place the unused strips back in the packet.

Stability at 4°C: 1 month

Rinse each well 4 times with Wash Buffer (300 µL/well).

Just before distributing the reagents and samples, remove the buffer from the wells by inverting the plate and shaking out the last drops on a paper towel.

### Plate set-up

A plate set-up is suggested hereafter.

The contents of each well may be recorded on the template sheet provided at the end of this technical booklet.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Bk	S8	S4	*	*	*	*	*	*	*	*	*
B	Bk	S8	S4	*	*	*	*	*	*	*	*	*
C	Bk	S7	S3	*	*	*	*	*	*	*	*	*
D	NSB	S7	S3	*	*	*	*	*	*	*	*	*
E	NSB	S6	S2	*	*	*	*	*	*	*	*	*
F	NSB	S6	S2	*	*	*	*	*	*	*	*	*
G	B0	S5	S1	*	*	*	*	*	*	*	*	QC
H	B0	S5	S1	*	*	*	*	*	*	*	*	QC

Bk : Blank

NSB : Non Specific Binding

B0: Maximum Binding

S1-S8: Standards 1-8

\*: Samples

QC: Quality Controls

## Pipetting the reagents

All samples and reagents must reach room temperature prior to performing the assay.

Use different tips to pipette the buffer, standard, sample, conjugate, antiserum and other reagents.

Before pipetting, equilibrate the pipette tips in each reagent. Do not touch the liquid already in the well when expelling with the pipette tip.

### Dilution Buffer

Dispense 100  $\mu\text{L}$  to Non Specific Binding (NSB) wells and 50  $\mu\text{L}$  to Maximum Binding (B0) wells.

### DHN-MA Standards

Dispense 50  $\mu\text{L}$  of each of the eight standards S1 to S8 in duplicate to appropriate wells. Start with the lowest concentration standard S8 and equilibrate the tip in the next higher standard before pipetting.

### Quality Control and Samples

Dispense 50  $\mu\text{L}$  in duplicate to appropriate wells. Highly concentrated samples may be diluted in Dilution Buffer.

### Conjugate Solution

Dispense 50  $\mu\text{L}$  to each well, **except** Blank (Bk) wells.

### DHN-MA Antiserum

Dispense 50  $\mu\text{L}$  to each well **except** Blank (Bk) wells and Non Specific Binding (NSB) wells.

## Incubating the plate

Cover the plate with the cover sheet and incubate over night at 4°C.

## Developing and reading the plate

- Reconstitute Substrate Solution as mentioned in the Reagent preparation section.
- Empty the plate by turning over. Rinse each well 4 times with 300  $\mu$ L Wash Buffer. At the end of the last washing step, empty the plate and blot the plate on a paper towel to discard any trace of liquid.
- Add 200  $\mu$ L of Substrate Solution to each well. Cover the plate with an aluminium sheet and incubate in the dark at room temperature. Optimal development is obtained using an orbital shaker.
- Wipe the bottom of the plate with a paper towel, and make sure that no liquid has splashed outside the wells.
- Read the plate at a wavelength between 405 and 414nm (yellow colour).

**After addition of Substrate Solution the absorbance has to be checked periodically (every 30 minutes) until the maximum absorbance has reached a minimum of 200 mAU (blank subtracted).**

Enzyme Immunoassay Protocol (volumes are in $\mu$ L)						
Volume/ Wells	Blank	NSB	B0	Standard	QC	Sample
Dilution Buffer	-	100	50	-	-	-
Standard	-	-	-	50	-	-
QC	-	-	-	-	500	-
Sample	-	-	-	-	-	50
Conjugate Sol.	-	-	-	50		
Antiserum	-	-	-	50		
Cover plate, incubate overnight at 4°C						
Wash plate 4 times, & discard liquid from the wells						
Substrate Sol.	-	-	-	200		
Incubate with an orbital shaker in the dark at RT						
Read the plate between 405 and 414 nm						

## 9. DATA ANALYSIS

Make sure that your plate reader has subtracted the absorbance readings of the blank well (absorbance of Substrate Solution alone) from the absorbance readings of the rest of the plate. If not, do it now.

- Subtract the average absorbance of NSB for each B0, standards, quality control and samples.
- Calculate the average absorbance for each B0, standard, quality control and samples.
- Calculate the B/B0 (%) for each standard, QC and sample (average absorbance of standards, QC or sample divided by average absorbance of B0) & multiplied by 100.
- Using a semi-log graph paper for each standard point, plot the B/B0 (%) on **y** axis versus the concentration (pg/mL) on **x** axis. Draw a best-fit line through the points.
- To determine the concentration of your sample, the corresponding B/B0 (%) value has to be comprised between 20% and 80%. Find the B/B0 (%) value on the **y** axis. Read the corresponding value on the **x** axis which is the concentration of your unknown sample.
- Diluted samples which concentration determined on standard curve is greater than 1000 pg/mL should be re-assayed after appropriated dilution in Dilution buffer.
- Samples with a concentration greater than 1000 pg/mL should be re-assayed after dilution in Dilution Buffer.
- Most plate readers are supplied with curve-fitting software capable of graphing these data (4-parameter logistic fit 4PL). If you have this type of software, we recommend using it. Refer to it for further information.

**Two vials of Quality Control are provided with this kit.**

**Your standard curve is validated only if the calculated concentration of the Quality Control obtained with the assay is +/- 25% of the expected concentration (see the label of QC vial).**

## 10. ACCEPTABLE RANGE

- B0 absorbance > 200 mAU blank subtracted in the conditions indicated above
- NSB absorbance < 35 mAU
- IC50: 65 to 85 pg/mL (mean: 75 pg/mL)
- QC sample:  $\pm 25\%$  of the expected concentration (see the label of QC vial)


## 11. TYPICAL RESULTS

The following data are for demonstration purpose only. Your data may be different and still correct.

These data were obtained using all reagents as supplied in this kit under the following conditions: 120 minutes developing at 20°C, reading at 414 nm. A 5-parameter curve-fitting was used to determine the concentrations.

	DHN-MA pg/mL	mAU	B/B0 (%)
Standard S1	1000	41	8.9
Standard S2	5000	63	13.7
Standard S3	250	102	22.4
Standard S4	125	172	37.7
Standard S5	62.5	270	59.2
Standard S6	31.3	342	74.8
Standard S7	15.7	392	85.8
Standard S8	7.8	433	94.8
B0	0	456	100

Typical DHN-MA standard curve



## 12. ASSAY VALIDATION AND CHARACTERISTICS

The Enzyme Immunoassay of DHN-MA has been validated for its use in human urine.

For additional information regarding the validation of immunoassay for protein biomarkers in biological samples, please refer to bibliography [11, 12].

The **Limit of Detection (LOD)** of DHN-MA corresponding to the B0 average minus three standard deviations is around 10 pg/mL.

The **IC50** is the concentration in DHN-MA corresponding to 50 % of the maximum Binding is around 75 pg/mL.

### Intra-assay variation in Dilution Buffer

QC level	QC1	QC2	QC3	QC4
QC dilution (1/x)	200	100	50	20
Mean of calculated concentration (pg/mL)	30.7	67.6	137	340
Mean concentration x dilution (pg/mL)	6133	6762	6842	6800
CV %	12.1 %	6.6 %	2.7 %	2.7 %

### Inter-assay variation in Dilution Buffer

QC level	QC1	QC2	QC3	QC4
QC dilution (1/x)	200	100	50	20
Mean of calculated concentration (pg/mL)	31.9	60.8	131	332
Mean concentration x dilution (pg/mL)	6370	6083	6567	6633
CV %	17.9 %	16.2 %	17.7 %	13.2 %

Due to the endogenous presence of DHN-MA in urine, the intraassay and inter-assay variations were performed on a pool of human urines at four dilutions (1/20, 1/50, 1/100 and 1/200) corresponding for each at a level of QC.

For the **intra-assay** validation, the number of replicates (n) is equal to 6 for each level of QC. The four validation levels were analysed along with the calibration curve for a unique experiment.

For the **inter-assay** validation, the number of replicates (n) is equal to 6 for each level of QC. The four validation levels were analysed along with the calibration curve for a total of 6 independent runs.



### Matrix variability

Matrix	1	2	3
Mean endogenous conc. x dil factor* (pg/mL)	2033	2193	2733
CV% (n=3)	17.6 %	13.7 %	10.2 %
Mean spiked <sup>§</sup> matrix con. x dil factor* (pg/mL)	3120	1260	1737
CV% (n=3)	2.6 %	9.6 %	10.5 %
Back calculated spiked concentration <sup>§</sup> (pg/mL)	1087	933	997
Back calculated spiked concentration <sup>§</sup> (recovery)	8.67 %	-6.67 %	-0.33 %

\* Matrix dilution factor: 1/100

<sup>§</sup> Matrix spiked with DHN-MA 1000 pg/mL

Three individual urine samples were spiked or not with DHN-MA 1000 pg/mL. Each sample (spiked or not) was evaluated 100 fold diluted in triplicate and analysed against a calibration curve.

### Cross-reactivity [6]

DHN-MA (racemic)	100 %
MA(N-acetyl-cysteine)	< 0.1 %
3-methylindole-MA	< 0.1 %
4-HNE hemiacetal-MA	< 0.1 %
HNA-lactone-MA	< 0.1 %
HNA-MA	30 %
butan-1-ol-MA	< 0.1 %
1-hexanol-MA	< 0.1 %
1-nonanol-MA	< 0.1 %
DHN	< 0.1 %
HNA	7 %
1,4-dihydroxyoctane MA	< 0.1 %
1,4-hydroxyundecane MA	< 0.1 %
1,3,4-trihydroxynonane MA	< 0.1 %
1,4-dihydroxyhexane MA	0.6 %

## Linearity

Matrix	Endogenous calculated conc. (pg/mL)	Spiked DHN-MA (pg/mL)	Dilution (1/x)	Mean calculated conc. (pg/mL) n=3	Mean conc. X dil (pg/mL) n=3	CV %	Back calculated spiked conc. (pg/mL)	Back calculated spiked conc. % RE
1	7 096	300 000	50	ND	ND	ND	ND	ND
			100	ND	ND	ND	ND	ND
			200	ND	ND	ND	ND	ND
			500	504	252 167	4.89 %	245 071	-18.3 %
			1000	239	239 333	3.60 %	232 237	-22.6 %
			2000	123	246 000	5.86 %	238 904	-20.4 %
2	5 432	300 000	50	ND	ND	ND	ND	ND
			100	ND	ND	ND	ND	ND
			200	1 260	252 067	3.39 %	246 635	-17.8 %
			500	548	274 167	9.30 %	268 735	-10.4 %
			1000	261	260 667	15.70 %	255 235	-14.9 %
			2000	133	266 667	5.43 %	261 235	-12.9 %

ND: outside range of standard curve

Two individual human urine samples were spiked or not with DHN-MA 300 ng/mL.

Each sample (spiked or not) was evaluated at 6 serial dilutions in order to be detected for at least 3 dilutions in the range of the standard curve.

Each dilution was tested in triplicate and analysed against a calibration curve.

## Parallelism

Matrix	Dilution factor (1/x)	Back calculated DHN-MA (pg/mL)	DHN-MA corrected dilution factor (pg/mL)	CV %
1	5	488	2440	15.1%
	10	253	2530	
	20	108	2160	
	50	35.6	1780	
	100	ND	ND	
	200	ND	ND	
	5	605	3025	18.8%
	10	248	2480	
	20	125	2500	
	50	37.7	1885	
	100	ND	ND	
	200	ND	ND	
2	5	661	3305	19.8%
	10	263	2630	
	20	107	2140	
	50	46.2	2310	
	100	ND	ND	
	200	ND	ND	
	5	ND	ND	11.7%
	10	586	5860	
	20	318	6360	
	50	137	6950	
	100	69.1	6910	
	200	40	8000	
	5	1305	6525	4.8%
	10	700	7000	
	20	326	6520	
	50	147	7350	
	100	71.1	7110	
	200	34.9	6980	
5	ND	ND	5.8%	
10	659	6590		
20	326	6520		

Matrix	Dilution factor (1/x)	Back calculated DHN-MA (pg/mL)	DHN-MA corrected dilution factor (pg/mL)	CV %
	50	148	7400	
	100	72.9	7290	
	200	35.3	7060	
3	5	360	1800	1.5%
	10	177	1770	
	20	91.1	1822	
	50	ND	ND	
	100	ND	ND	
	200	ND	ND	
	5	306	1530	
	10	160	1600	
3	20	82.3	1646	25.1%
	50	50	2500	
	100	ND	ND	
	200	ND	ND	
	5	352	1760	
	10	182	1820	
3	20	80.2	1604	7.7%
	50	38.7	1935	
	100	ND	ND	
	200	ND	ND	
	200	ND	ND	

ND: outside range of standard curve

Three individual human urine samples were diluted between 1:5 and 1:200 by serial dilution in order to be detected for at least 3 dilutions in the range of the standard curve.

Each dilution was tested in triplicate and analysed against a calibration curve.

### Stability tests (freezing, thawing and 24 hours at 4°C or 20-25°C)

QC level	QC condition	Mean concentration (pg/mL) n = 3	CV %
LQ	Run QC	57.3	1.7 %
	Freeze stability (3	67.1	8.0 %
	24h at RT stability	58.9	0.6 %
	24h at 4 °C	53.2	3.7 %
HQ	Run QC	345.3	7.4 %
	Freeze stability (3	358.7	1.6 %
	24h at RT stability	345.7	3.4 %
	24h at 4 °C	330.3	1.3 %

The stability tests were performed on a pool of human urines at two dilutions (1/50 and 1/200) corresponding to Low QC and High QC respectively.

QCs were prepared, frozen at -20°C, and then stored 24h at room temperature or 4°C, or frozen/thawed 3 times. QCs were then analysed against the calibration curve.

### Long term stability tests

QC level	QC condition	Mean concentration (pg/mL) n = 3	CV %
LQ	Run QC	62.5	0.96 %
	1 month stability	60.2	4.07 %
HQ	Run QC	333	0.35 %
	1 month stability	335	2.37 %

QC level	QC condition	Mean concentration (pg/mL) n = 3	CV %
LQ	Run QC	68.0	7.39 %
	3 month stability	72.2	14.30 %
HQ	Run QC	353	15.80 %
	3 month stability	326	1.69 %

QC level	QC condition	Mean concentration (pg/mL) n = 3	CV %
LQ	Run QC	66.9	4.01 %
	6 month stability	73.8	3.42 %
HQ	Run QC	324	2.94 %
	6 month stability	359	10.70 %

The long term stability tests were performed on a pool of human urines at two dilutions (1/50 and 1/200) corresponding to Low QC and High QC respectively.

QCs were prepared then frozen at -20°C before to be assayed at different times. QCs were then analysed against the calibration curve.

### 13. ASSAY TROUBLE SHOOTING

#### Absorbance values too low:

- organic contamination of water,
- one reagent has not been dispensed,
- incorrect preparation/dilution,
- assay performed before reagents reached room temperature,
- reading time not long enough.

#### High signal and background in all wells:

- inefficient washing,
- overdeveloping (incubation time should be reduced),
- high ambient temperature.

#### High dispersion of duplicates:

- poor pipetting technique,
- irregular plate washing.

#### If a plate is accidentally dropped after dispatch of the AChE® substrate (Substrate Solution) or if it needs to be revealed again:










- one only needs to wash the plate, add fresh Substrate Solution and proceed with a new development.
- otherwise, the plate can be stored at 4°C with Wash Buffer in wells while waiting for technical advice from the Bioreagent Department.

These are a few examples of troubleshooting that may occur.


## 14. BIBLIOGRAPHY

1. Grassi J., Pradelles P.  
Compounds labelled by the acetylcholinesterase of *Electrophorus Electricus*. Its preparation process and its use as a tracer or marker in enzyme-immunological determinations. *United States patent, N° 1,047,330. September 10, 1991*
2. Grassi J., Pradelles P.  
The use of Acetylcholinesterase as a Universal marker in Enzyme-Immunoassays. *Proceedings of the Third International Meeting on Cholinesterases, American Chemical Society (1991)*
3. Pradelles P., Grassi J., Maclouf J.  
Enzyme Immunoassays of Eicosanoids Using Acetylcholinesterase. *Methods in Enzymology (1990), vol. 187, 24-34*
4. Spickett CM.  
The lipid peroxidation product 4-hydroxy-2-nonenal: Advances in chemistry and analysis. *Redox Biol. 2013 Jan 21;1(1):145-152. eCollection 2013.*
5. Pierre F., Peiro G., Taché S. et al.  
New marker of colon cancer risk associated with heme intake: 1,4-dihydroxynonane mercapturic acid. *Cancer Epidemiol Biomarkers Prev. 2006 Nov;15(11):2274-9.*
6. Guéraud F, Peiro G, Bernard H. et al.  
Enzyme immunoassay for a urinary metabolite of 4-hydroxynonanal as a marker of lipid peroxidation *Free Radic Biol Med. 2006 Jan 1;40(1):54-62. Epub 2005 Sep 2*
7. Peiro G., Alary J., Cravedi JP. et al.  
Dihydroxynonane mercapturic acid, a urinary metabolite of 4-hydroxynonanal, as a biomarker of lipid peroxidation. *Biofactors. 2005;24(1-4):89-96.*
8. Zarkovic N.  
4-hydroxynonanal as a bioactive marker of pathophysiological processes. *Mol Aspects Med. 2003 Aug-Oct;24(4-5):281-91. Review.*
9. Uchida K., Szweda LI., Chae HZ. et al.  
Immunochemical detection of 4-hydroxynonanal protein adducts in oxidized hepatocytes. *Proc Natl Acad Sci U S A. 1993 Sep 15;90(18):8742-6.*
10. Esterbauer H., Schaur RJ., Zollner H.  
Chemistry and biochemistry of 4-hydroxynonanal, malonaldehyde and related aldehydes. *Free Radic Biol Med. 1991;11(1):81-128. Review.*
11. Valentin MA., Ma S, Zhao A. et al.  
Validation of immunoassay for protein biomarkers: Bioanalytical study plan implementation to support pre-clinical and clinical studies. *J Pharm Biomed Anal. (2011) 55(5) : 869-877*
12. European Medicines Agency  
*Guideline on bioanalytical method validation, 21 July 2011*

## 15. EXPLANATION OF SYMBOLS

	Catalogue number
	Batch code
	Caution
	Use by date
	Temperature limit
	Manufacturer
 <a href="http://www.biovendor.com">www.biovendor.com</a>	Read electronic instructions for use - eIFU
	The content is sufficient for 96 tests
	Biological risks







BioVendor – Laboratorní medicína a.s.  
Karásek 1767/1, 621 00 Brno, Czech Republic  
+420 549 124 185  
[info@biovendor.com](mailto:info@biovendor.com)  
[sales@biovendor.com](mailto:sales@biovendor.com)  
[www.biovendor.com](http://www.biovendor.com)

Example Version

