

# NOVEL IMMUNOASSAY TOOL FOR QUANTIFICATION OF STROKE RELATED miRNA BIOMARKERS

Tereza Batkova<sup>1,2,3</sup>, Iveta Krepelkova<sup>2,5</sup>, Barbora Dvorakova<sup>2</sup>, Eva Bace<sup>2</sup>, Robert Mikulik<sup>1,3,4</sup>, Viktor Ruzicka<sup>2</sup>

<sup>1</sup>Stroke Program, International clinical research center (FNUSA-ICRC), Brno; <sup>2</sup>BioVendor-laboratorní medicína a.s., Brno; <sup>3</sup>Faculty of medicine, Masaryk University, Brno; <sup>4</sup>Department of Neurology, St Anne's University Hospital, Brno; <sup>5</sup>Department of Oncology, Faculty of medicine, Masaryk University, Brno

## INTRODUCTION

- miRNAs are promising stroke biomarkers, but methods for their detection have limitations especially reverse transcription step in the most widely used RT-qPCR principle [1].
- We present novel method for quantitative miRNA measurement based on immunoassay – miREIA.
- Analytical characteristics and correlation with RT-qPCR (ID3EAL miRNA qPCR Assays; MIRXES Life Science) were evaluated on 3 stroke related miRNAs hsa-miR-223-3p, hsa-miR-486-5p and hsa-miR-let-7b-5p [2-4].

## METHODS

**miREIA** - an enzyme immunoassay for miRNA quantification involves

- hybridization of miRNA isolated from a patient sample to complementary biotinylated DNA probe (Fig. 1a).
- the DNA/RNA hybrids are transferred into microplate wells pre-coated with monoclonal antibody specific to perfectly matched DNA/miRNA hybrids (Fig. 1b).
- after washing, the solid phase is incubated with streptavidin-HRP conjugate followed by another washing step. The resulting complexes are visualized by chromogenic substrate (Fig. 1b).
- the absorbance is proportional to the concentration of target miRNA [5].




Fig. 1: miREIA workflow (a – hybridization; b – immunoreaction and detection)

## RESULTS

**miREIA analytical performance** determined on the assay hsa-miR-223-3p was following:

- calibration range 12,5 - 0,39 amol/ul (Fig.2)
- limit of detection 0,13 amol/ul
- intra-assay (n=8) CV= 8,0 %
- spiking recovery 98,1 %
- dilution linearity 101,6 %
- repeatability CV = 12 % on average (Fig. 3)




Fig. 2: Typical miREIA standard curve – similar to ELISA principle the absorbance is proportional to the absolute concentration of target miRNA






Fig.3: Repeatability – measured in 20 consecutive experiments on 4 different concentrations showed low variation what supports miREIA as a good validation tool and potentially applicable method in clinical practice.

## CORRELATION


RNA from 38 samples (brain tissue, PBMC, whole blood and plasma) was isolated and analyzed by novel miREIA method and RT-qPCR. We found high correlation between both methods (Graphs 1-3).



Graph 1: correlation of hsa-miR-223-3p



Graph 2: correlation of hsa-miR-486-5p



Graph 3: correlation of hsa-miR-let-7b-5p

## CONCLUSION

- miREIA has favourable analytical characteristics.
- Great correlation with RT-qPCR technique measuring 3 stroke related miRNAs.
- Method is well compatible with classic clinical workflow and could help with validation of stroke biomarkers.

## REFERENCES

- [1] STÄHLBERG, Anders, et al. Properties of the reverse transcription reaction in mRNA quantification. *Clinical chemistry*, 2004, 50.3: 509-515.
- [2] CHEN, Yajing, et al. Increased circulating exosomal miRNA-223 is associated with acute ischemic stroke. *Frontiers in neurology*, 2017, 8: 57.
- [3] LONG, G., et al. Circulating miR-30a, miR-126 and let-7b as biomarker for ischemic stroke in humans. *BMC neurology*, 2013, 13.1: 178.
- [4] ZHANG, Yu, et al. Clinical predictor and circulating microRNA profile expression in patients with early onset post-stroke depression. *Journal of affective disorders*, 2016, 193: 51-58.
- [5] www.biovendor.com/mirna