

ENG

Product Data Sheet:

HUMAN APRIL ELISA

Catalogue number:

RAG023R

For research use only!

Example Version

1. INTENDED USE	3
2. HANDLING, STORAGE	3
3. INTRODUCTION	3
4. TEST PRINCIPLE	4
5. TECHNICAL HINTS	4
6. REAGENT SUPPLIED	5
7. MATERIAL REQUIRED BUT NOT SUPPLIED	5
8. PREPARATION OF REAGENTS	6
9. PREPARATION OF SAMPLES	8
10. ASSAY PROCEDURE	9
11. CALCULATIONS	10
12. PERFORMANCE CHARACTERISTICS	11
13. TROUBLESHOOTING	13
14. REFERENCES	14
15. EXPLANATION OF SYMBOLS	15

1. INTENDED USE

The Human APRIL ELISA Kit is to be used for the in vitro quantitative determination of human APRIL in cell culture supernatants, serum and plasma. This ELISA Kit is for research use only.

2. HANDLING, STORAGE

- Reagent must be stored at 2-8°C when not in use.
- Plate and reagents should be at room temperature before use.
- Do not expose reagents to temperatures greater than 25°C.

3. INTRODUCTION

The B cell-stimulating molecules, BAFF (B cell activating factor) and APRIL (a proliferation-inducing ligand), are critical factors in the maintenance of the B cell pool and humoral immunity (1). APRIL binds to transmembrane activator and CAML interactor (TACI), B cell maturation antigen (BCMA), and heparan sulfate proteoglycans (HSPG) within the extracellular matrix or on the surface of cells such as plasma cells (2). APRIL maintains B cell homeostasis by acting at a later stage, modulating the function and survival of antigen-experienced B cells. APRIL (as well as BAFF) stimulates class-switch recombination (CSR), hence contributes to shaping humoral effector mechanisms. With regards to humoral memory, APRIL is involved in the establishment and survival of the long-lived plasma cell (LLPC) pool in the bone marrow (BM) (3).

APRIL is expressed by a number of myeloid-derived cell types including BM granulocytes, megakaryocytes, eosinophils and osteoclasts and by dendritic cells following exposure to IFN α , IFN γ or CD40L. APRIL expression is induced during hematopoiesis in the bone marrow. APRIL expression is not limited to cells of myeloid origin, but can also be found in epithelial cells of the gut, tonsil, breast and skin. Finally, APRIL is expressed in tumor cell lines and human cancer cells of colon, thyroid and lymphoid origin (1).

The B cell-stimulating molecules, BAFF (B cell activating factor) and APRIL (a proliferation-inducing ligand), are implicated in several human autoimmune diseases with autoreactive B cell involvement, including systemic lupus erythematosus (SLE) (4), Sjögren's syndrome (SS) (5), IgA nephropathy (IgAN) (6) and rheumatoid arthritis (RA) (7). APRIL might also function in enhancing proliferation of some tumor cells, especially B-cell malignancies (8, 9).

4. TEST PRINCIPLE

This assay is a sandwich Enzyme Linked-Immunosorbent Assay (ELISA) for quantitative determination of human APRIL in cell culture supernatants, serum and plasma (EDTA, heparin or citrate). A monoclonal antibody specific for APRIL has been precoated onto the 96-well microtiter plate. Standards (STD) and samples are pipetted into the wells for binding to the coated antibody. After extensive washing to remove unbound compounds, APRIL is recognized by the addition of a biotinylated monoclonal antibody specific for APRIL (DET). After removal of excess biotinylated antibody, streptavidin-peroxidase (STREP-HRP) is added. Following a final washing, peroxidase activity is quantified using the substrate 3,3',5,5'-tetramethylbenzidine (TMB). The intensity of the color reaction is measured at 450 nm after acidification and is directly proportional to the concentration of APRIL in the samples.

5. TECHNICAL HINTS

- It is recommended that all standards and samples be run in duplicate.
- Do not combine leftover reagents with those reserved for additional wells.
- Reagents from the kit with a volume less than 100 μ l should be centrifuged.
- Residual wash liquid should be drained from the wells after last wash by tapping the plate on absorbent paper.
- Crystals could appear in the 10X solution due to high salt concentration in the stock solutions. Crystals are readily dissolved at room temperature or at 37°C before dilution of the buffer solutions.
- Once reagents have been added to the 16-well strips, DO NOT let the strips DRY at any time during the assay.
- Keep Substrate Solution protected from light.
- The Stop Solution consists of sulfuric acid. Although diluted, the Stop Solution should be handled with gloves, eye protection and protective clothing.

6. REAGENT SUPPLIED

Kit Components	Quantity
1 plate coated with APRIL Antibody	6x16-well strips
2 bottles Wash Buffer 10X	30 ml
2 bottle ELISA Buffer 10X	30 ml
1 vial Detection Antibody	20 μ l
1 vial HRP Labeled Streptavidin (lyophilized)	2 μ g
1 vial human APRIL Standard (lyophilized)	100 ng
1 bottle TMB Substrate Solution	12 ml
1 bottle Stop Solution	12 ml
2 plate sealers (plastic film)	
2 silica Gel Minibags	

7. MATERIAL REQUIRED BUT NOT SUPPLIED

- Microtiterplate reader at 450 nm
- Calibrated precision single and multi-channel pipettes. Disposable pipette tips
- Deionized water
- Microtubes or equivalent for preparing dilutions
- Disposable plastic containers for preparing working buffers
- Plate washer: automated or manual
- Glass or plastic tubes for diluting and aliquoting standard

8. PREPARATION OF REAGENTS

NOTE: Prepare just the appropriate amount of the buffers necessary for the assay.

Wash Buffer 10X

has to be diluted with deionized water 1:10 before use (e.g. 30 ml Wash Buffer 10X + 270 ml water) to obtain Wash Buffer 1X.

ELISA Buffer 10X

has to be diluted with deionized water 1:10 before use (e.g. 10 ml ELISA Buffer 10X + 90 ml water) to obtain ELISA Buffer 1X.

Detection Antibody (DET)

has to be diluted 1:1000 in ELISA Buffer 1X (10 μ l DET + 10 ml ELISA Buffer 1X).

NOTE: The diluted Detection Antibody is not stable and cannot be stored!

HRP Labeled Streptavidin (STREP-HRP)

has to be reconstituted with 100 μ l of ELISA Buffer 1X.

After reconstitution of STREP-HRP, prepare aliquots and store them at -20°C. Avoid freeze/thaw cycles

Dilute the reconstituted STREP-HRP to the working concentration by adding 50 μ l in 10 ml of ELISA Buffer 1X (1:200).

NOTE: The diluted STREP-HRP is not stable and cannot be stored!

Human APRIL Standard (STD)

has to be reconstituted with 100 µl of ELISA Buffer 1X.

This reconstitution produces a stock solution of 1 µg/ml. Mix the standard to ensure complete reconstitution and allow the standard to sit for a minimum of 15 minutes at room temperature. Mix well prior to making dilutions.

NOTE: The reconstituted standard is aliquoted and stored at -20°C.

Dilute the standard protein concentrate (STD) (1 µg/ml) in ELISA Buffer 1X. A seven-point standard curve using 2-fold serial dilutions in ELISA Buffer 1X is recommended.

Suggested standard points are:

250, 125, 62.5., 31.125, 15.5625, 7.78125, 3.90625 and 0 pg/ml.

Start with the dilution of the concentrate (STD):

To obtain	Add	Into
10 ng/ml	10µl of APRIL (STD) (1 µg/ml)	990 µl of ELISA Buffer 1X
1 ng/ml	100µl of APRIL (10 ng/ml)	900 µl of ELISA Buffer 1X

Dilute further for the standard curve:

To obtain	Add	Into
250 pg/ml	300 µl of APRIL (1 ng/ml)	900 µl of ELISA Buffer 1X
125 pg/ml	300 µl of APRIL (250 pg/ml)	300 µl of ELISA Buffer 1X
62.5 pg/ml	300 µl of APRIL (125 pg/ml)	300 µl of ELISA Buffer 1X
31.25 pg/ml	300 µl of APRIL (62.5 pg/ml)	300 µl of ELISA Buffer 1X
15.625 pg/ml	300 µl of APRIL (31.25 pg/ml)	300 µl of ELISA Buffer 1X
7.781 pg/ml	300 µl of APRIL (15.625 pg/ml)	300 µl of ELISA Buffer 1X
3.906 pg/ml	300 µl of APRIL (7.781 pg/ml)	300 µl of ELISA Buffer 1X
0 ng/ml	300 µl of ELISA Buffer 1X	Empty tube

9. PREPARATION OF SAMPLES

Serum

Use a serum separator tube. Let samples clot at room temperature for 30 minutes before centrifugation for 20 minutes at 1,000xg. Assay freshly prepared serum or store serum in aliquot at $\leq -20^{\circ}\text{C}$ for later use. Avoid repeated freeze/thaw cycles.

Plasma

Collect plasma using heparin, citrate or EDTA as an anticoagulant. Centrifuge for 15 minutes at 1000xg within 30 minutes of collection. Assay freshly prepared plasma or store plasma sample in aliquot at $\leq -80^{\circ}\text{C}$ for later use. Avoid repeated freeze/ thaw cycles.

Serum, Plasma and Cell Culture Supernatant

have to be diluted in ELISA Buffer 1X. Samples containing visible precipitates must be clarified before use.

NOTE: As a starting point, 1/10 dilution of serum or 1/4 dilution of plasma is recommended! If sample values fall outside the detection range of the assay, a lower or higher dilution may be required!

Example Version

10. ASSAY PROCEDURE

1

Determine the number of 16-well strips needed for the assay and insert them in the frame for current use. The extra strips are left in the bag with 2 silica gel minibags and stored at 4°C.

NOTE: Remaining 16-well strips coated with APRIL antibody when opened can be stored in the presence of 2 silica gel minibags at 4°C for up to 1 month.

2.

Add 100 µl of the different standards into the appropriate wells in duplicate! At the same time, add 100 µl of diluted serum, plasma or cell culture supernatant samples in duplicate to the wells (**see 8.1. Preparation and Storage of Reagents and 8.2 Preparation of Samples**).

3.

Cover the plate with plastic film and incubate for **2 hours at Room Temperature**.

4.

Aspirate the coated wells and add 300 µl of Wash Buffer 1X using a multichannel pipette or auto-washer. Repeat the process for a total of five washes. After the last wash, complete removal of liquid is essential for good performance.

5.

Add 100 µl to each well of the diluted Detection Antibody (**DET**) (**see 8.1 Preparation and Storage of Reagents**).

6.

Cover the plate with plastic film and incubate for **1 hour at room temperature**.

7.

Aspirate the coated wells and add 300 µl of Wash Buffer 1X using a multichannel pipette or auto-washer. Repeat the process for a total of five washes. After the last wash, complete removal of liquid is essential for good performance.

8.

Add 100 µl to each well of the diluted HRP Labeled Streptavidin (**STREP-HRP**) (**see 8.1. Preparation and Storage of Reagents**).

9.

Cover the plate with plastic film and incubate for **30 min at room temperature**.

10.

Aspirate the coated wells and add 300 µl of Wash Buffer 1X using a multichannel pipette or auto-washer. Repeat the process for a total of five washes. After the last wash, complete removal of liquid is essential for good performance.

11.

Add 100 µl to each well of TMB substrate solution (**TMB**).

12.

Allow the color reaction to develop **at room temperature in the dark for 10 minutes**. Do not cover the plate.

13.

Stop the reaction by adding 100 µl of Stop Solution (**STOP**). Tap the plate gently to ensure thorough mixing. The substrate reaction yields a blue solution that turns yellow when Stop Solution (**STOP**) is added.

! CAUTION: CORROSIVE SOLUTION !


14.

Measure the OD at 450 nm in an ELISA reader.

11. CALCULATIONS

- Average the duplicate readings for each standard and sample and subtract the average blank value (obtained with the 0 pg/ml point).
- Generate the standard curve by plotting the average absorbance obtained for each standard concentration on the horizontal (X) axis vs. the corresponding APRIL concentration (pg/ml) on the vertical axis (see TYPICAL DATA).
- Calculate the APRIL concentrations of samples by interpolation of the regression curve formula as shown above in a form of a quadratic equation
- If the test sample was diluted, multiply the interpolated value by the dilution factor to calculate the concentration of human APRIL in the sample

The following data are obtained using the different concentrations of standard as described in this protocol:



Standard APRIL (pg/ml)	Optical Density (mean)
250	2.694
125	1.580
62.50	0.951
31.25	0.549
15.63	0.297
7.81	0.188
3.91	0.118
0	0.062

Figure: Standard curve

12. PERFORMANCE CHARACTERISTICS

Sensitivity (Limit of detection)

The lowest level of human APRIL that can be detected by this assay is 1 pg/ml.

NOTE: The Limit of detection was measured by adding three standard deviations to the mean value of 50 zero standard.

Assay Range

3.90 pg/ml – 250 pg/ml

Specificity

This ELISA is specific for the measurement of natural and recombinant human APRIL. Detection of APRIL (human) in biological fluids by this ELISA kit is abolished in the presence of a APRIL receptor such as TACI:Fc.




Figure: Specific quantitation of APRIL in human serum.

Method: Serum from a healthy patient is left untreated or treated with 1µg/ml of the APRIL receptor, TACI (h):Fc (h) . APRIL levels were measured using the Human APRIL ELISA Kit.

Precision:

Intra-assay:

Four samples of known concentrations of human APRIL were assayed in replicates 8 times to test precision within an assay.

Samples	Means (ng/ml)	SD	CV (%)	n
A1	2.86	0.15	5.22	8
A2	2,60	0.15	5.70	8
A3	1.62	0.10	6.34	8
A4	0.96	0.05	5.39	8

Inter-assay

Four samples of known concentrations of human APRIL were assayed in 5 separate assays to test precision between assays.

Samples	Means (ng/ml)	SD	CV (%)	n
B1	2.792	0.193	6.94	5
B2	0.155	0.015	9.37	5
B3	0.100	0.008	8.17	5
B4	0.775	0.076	9.80	5

Recovery

When samples are spiked with known concentrations of human APRIL, the recovery averages range from 92% to 109%.

Linearity

Different samples containing human APRIL were diluted several fold (1/2 to 1/8 for plasmas and 1/20 to 1/80 for sera) and the measured recoveries ranged from 93% to 109%.

Expected values:

Human APRIL levels range in serum from **0.1 ng/ml to >20 ng/ml** and in plasma from **20 pg/ml to > 1 ng/ml**. Levels of APRIL detected in plasma are lower than in serum.










13. TROUBLESHOOTING


PROBLEM	POSSIBLE CAUSES	SOLUTIONS
No signal or weak signal	Omission of key reagent	Check that all reagents have been added in the correct order.
	Washes too stringent	Use an automated plate washer if possible.
	Incubation times inadequate	Incubation times should be followed as indicated in the manual.
	Plate reader settings not optimal	Verify the wavelength and filter setting in the plate reader.
	Incorrect assay temperature	Use recommended incubation temperature. Bring substrates to room temperature before use.
High background	Concentration of STREP-HRP too high	Use recommended dilution factor.
	Inadequate washing	Ensure all wells are filling wash buffer and are aspirated completely.
Poor standard curve	Wells not completely aspirated	Completely aspirate wells between steps.
	Reagents poorly mixed	Be sure that reagents are thoroughly mixed.
Unexpected results	Omission of reagents	Be sure that reagents were prepared correctly and added in the correct order.
	Dilution error	Check pipetting technique and double-check calculations.

14. REFERENCES

1. Targeting BAFF and APRIL in systemic lupus erythematosus and other antibody-associated diseases: E. Samy, et al.; *Int. Rev. Immunol.* 36, 3 (2017)
2. Identification of proteoglycans as the APRIL-specific binding partners: K. Ingold, et al.; *J. Exp. Med.* 201, 1375 (2005)
3. Factors of the bone marrow microniche that support human plasma cell survival and immunoglobulin secretion: D.C. Nguyen, et al.; *Nat. Commun.* 12, 3698 (2018)
4. Raised serum APRIL levels in patients with systemic lupus erythematosus: T. Koyama, et al.; *Ann. Rheum. Dis.* 64, 1065 (2005)
5. The expression of APRIL in Sjogren's syndrome: Aberrant expression of APRIL in the salivary gland: J.L. Vosters, et al.; *Rheumatology (Oxford)* 51, 1557 (2012)
6. Increased APRIL Expression Induces IgA1 Aberrant Glycosylation in IgA Nephropathy: Y.L. Zhai, et al.; *Medicine (Baltimore)* 95, e3099 (2016)
7. The BAFF/APRIL system: an important player in systemic rheumatic diseases: F. Mackay, et al.; *Curr. Dir. Autoimmun.* 8, 243 (2005)
8. APRIL, a New Ligand of the Tumor Necrosis Factor Family, Stimulates Tumor Cell Growth: M. Hahne, et al.; *J. Exp. Med.* 188, 1185 (1998)
9. In situ detection of APRIL-rich niches for plasma-cell survival and their contribution to B-cell lymphoma development: M. Burjanadze, et al.; *Histol. Histopathol.* 24, 1061 (2009)

15. EXPLANATION OF SYMBOLS

	Catalogue number
	Batch code
	Caution
	Use by date
	Temperature limit
	Manufacturer
 www.biovendor.com	Read electronic instructions for use - eIFU
	The content is sufficient for 96 tests
	Biological risks





BioVendor – Laboratorní medicína a.s.
Karásek 1767/1, 621 00 Brno, Czech Republic
+420 549 124 185
info@biovendor.com
sales@biovendor.com
www.biovendor.com

Example Version

