DCM001-13  
Ed. 09/2022

# CORTISOL ELISA

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### 18. BIBLIOGRAFIA

1. Ramamoorthy, S. & Cidlowski, J. A. Corticosteroids: Mechanisms of Action in Health and Disease. *Rheumatic diseases clinics of North America* 42, 15-31, doi:10.1016/j.rdc.2015.08.002 (2016).
2. Wild, D. *The immunoassay handbook*. (Elsevier, 2005).
3. Cain, D.W. & Cidlowski, J. A. Specificity and sensitivity of glucocorticoid signaling in health and disease. *Best practice & research. Clinical endocrinology & metabolism* 29, 545-556, 2015.04.007 (2015).
4. Schacke, H., Docke, W. D. & Asadullah, K. Mechanisms involved in the side effects of glucocorticoids. *Pharmacology & therapeutics* 96, 23-43 (2002).
5. Charmandari, E., Nicolaides, N. C. & Chrousos, G. P. Adrenal insufficiency *Lancet* 383, 2152-2167, (13)61684-0 (2014).
6. Bancos, I., Hahner, S., Tomlinson, J. & Arlt, W. Diagnosis and management of adrenal insufficiency. *The lancet. Diabetes & endocrinology* 3, 216-226, (14)70142-1 (2015).
7. Arlt, W. & Allolio, B. Adrenal insufficiency. *Lancet* 361, 1881-1893, (03)13492-7 (2003).
8. Lacroix, A., Feelders, R. A., Stratakis, C. A. & Nieman, L. K. Cushing's syndrome. *Lancet* 386, 913-927, (14)61375-1 (2015).
9. Arnaldi, G. et al. Diagnosis and complications of Cushing's syndrome: a consensus statement. *The Journal of clinical endocrinology and metabolism* 88, 5593-5602 (2003).
10. Yaneva, M. et al. Midnight salivary cortisol for the initial diagnosis of Cushing's syndrome of various causes. *The Journal of clinical endocrinology and metabolism* 89, 3345-3351, (2004).
11. Krieger, D. T., Allen, W., Rizzo, F. & Krieger, H. P. Characterization of the normal temporal pattern of plasma corticosteroid levels. *The Journal of clinical endocrinology and metabolism* 32, 266-284, (1971).
12. Ceccato, F. et al. Performance of salivary cortisol in the diagnosis of Cushing's syndrome, adrenal incidentaloma, and adrenal insufficiency. *European journal of endocrinology / European Federation of Endocrine Societies* 169, 31-36, (2013).
13. Basic QC Practices On-line Course; <http://www.Westgard.com>.
14. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

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# CORTISOL ELISA

for routine analysis

Direct immunoenzymatic determination of Cortisol in human serum or plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO001

## 1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use

For Laboratory Professional Use

Cortisol ELISA is a manual in vitro diagnostic device intended for the quantitative determination of cortisol in human serum or plasma.

## 2. CLINICAL SIGNIFICANCE

Cortisol is a glucocorticoid steroid hormone synthesised from cholesterol in the zona fasciculata of the adrenal cortex. Glucocorticoids are predominantly stress hormones that regulate a variety of physiological processes by binding to glucocorticoid receptors expressed throughout the body<sup>1</sup>. Endogenous glucocorticoids like cortisol play a crucial role in establishing the body's homeostasis after physical or psychological stress. Additionally, cortisol plays a role in the metabolism of carbohydrates, fat and protein, the maintenance of myocardial function<sup>2</sup> and is involved in the development and function of the central nervous system, the immune system and the gastrointestinal system<sup>3</sup>.

Approximately 80-90% of cortisol in plasma is protein-bound to an  $\alpha$ 2-globulin, cortisol-binding-globulin (CBG), also known as transcortin. The remaining 5 to 15% is bound to albumin, with only 5% in a free form. The non-protein bound free cortisol is the biologically active hormone, with higher levels of CBG protecting against high cortisol concentrations<sup>2,3</sup>. Increased levels of cortisol in the circulation feedback to the hypothalamus and pituitary gland, decreasing the levels of ACTH and corticotropin-releasing hormone in a classic negative feedback loop<sup>4</sup>. The hypothalamic-pituitary-adrenal axis (HPA) and consequently cortisol, is under the control of three major factors including stress levels, the negative feedback loop and a circadian rhythm<sup>2</sup>. Cortisol determinations are frequently used in the assessment of adrenocortical function and other disturbances of the hypothalamic-pituitary-adrenal axis. These disturbances can result in under- or over- production of glucocorticoids, most commonly referred to as Addison's disease and Cushing's syndrome respectively<sup>3</sup>. Addison's disease is a failure of the adrenal gland to produce cortisol due to a primary adrenal insufficiency, most frequently caused by autoimmune adrenalitis. Secondary adrenal insufficiency is due to an impaired regulation of the adrenal gland by the HPA axis. Whereas patients with primary

adrenal insufficiency can be diagnosed based on a combined measurement of early morning serum cortisol and plasma ACTH levels, secondary adrenal insufficiency most commonly requires dynamic testing such as the insulin tolerance test or the overnight metyrapone test<sup>5-7</sup>. Cushing's syndrome (CS) is defined as a large group of signs and symptoms to prolonged and inappropriately high exposure to glucocorticoids<sup>8</sup>. The most common cause is iatrogenic due to use of prescribed corticosteroids. Endogenous CS may be caused by excess ACTH production (80-85% of cases) or ACTH-independent due to excess secretion of cortisol by the adrenal gland<sup>9</sup>.

Under normal physiological conditions, the HPA axis has a characteristic circadian rhythm with highest cortisol levels found in the morning, reaching a nadir in the late evening<sup>10,11</sup>. As cortisol secretion is pulsatile, assessment of levels in a single sample may not provide adequate diagnostic accuracy. Absence of the normal diurnal rhythm is a good indication of the presence of CS<sup>12</sup>. Whereas in healthy individuals late night levels are undetectable, patients with increased function of the adrenal lose the normal circadian rhythm and have detectable levels of cortisol at midnight<sup>9</sup>. The time of collection is therefore an important consideration in the assessment of cortisol levels. Additionally, concentrations should be referred to assay-specific normal reference range profiles and clinical cut-off levels.

## 3. PRINCIPLE OF THE METHOD

The Cortisol ELISA is a competitive enzyme immunometric assay (ELISA) where Cortisol (antigen) in the sample competes with the antigenic Cortisol conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti-Cortisol coated on the microplate (solid phase).

After the incubation, the bound/free separation is performed by a simple solid phase washing. Then, the enzyme HRP in the bound fraction reacts with the Substrate ( $H_2O_2$ ) and the TMB Substrate and develops a blue colour that changes into yellow when the Stop Solution ( $H_2SO_4$ ) is added. The colour intensity is inversely proportional to the Cortisol concentration in the sample.

Cortisol concentration in the sample is calculated through a calibration curve.

## 4. REAGENTS, MATERIALS AND INSTRUMENTATION

### 4.1. Reagents and materials supplied in the kit

#### 1. Calibrators (5 vials, 1 mL each)

CAL0	<b>REF</b> DCE002/0106-0
CAL1	<b>REF</b> DCE002/0107-0
CAL2	<b>REF</b> DCE002/0108-0
CAL3	<b>REF</b> DCE002/0109-0
CAL4	<b>REF</b> DCE002/0110-0

#### 2. Control (1 vial, 1 mL)

See Control concentration on the Certificate of Analysis  
**REF** DCE045/0103-0

#### 3. Conjugate (1 vial, 21 mL)

Cortisol conjugated with Horseradish peroxidase (HRP)  
**REF** DCE002/0102-0

#### 4. Coated Microplate (1 breakable microplate)

Microplate coated with anti-Cortisol antibodies  
**REF** DCE002/0103-0

#### 5. TMB Substrate (1 vial, 15 mL)

H<sub>2</sub>O<sub>2</sub>-TMB 0.26 g/L (*avoid any skin contact*)  
**REF** DCE004-0

#### 6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (*avoid any skin contact*)  
**REF** DCE005-0

#### 7. 10X Conc. Wash Solution (1 vial, 50 mL)

Phosphate buffer 0.2M pH 7.4 **REF** DCE054-0

### 4.2. Materials required but not provided

Distilled water

### 4.3. Auxiliary materials and instrumentation

Automatic dispenser  
Precision Pipetting Devices  
Microplate reader (450 nm, 620-630 nm)

## 5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of ProClin™ 300 as preservative. Avoid contact with skin or mucosa.  
Classification under CLP  
Skin Sens. 1

#### Hazard statements

H317: May cause an allergic skin reaction

#### Precautionary statements

P261: Avoid breathing dust / fumes / gas / mist / vapours / spray.

P272: Contaminated work clothing should not be allowed out of the workplace.

P280: Wear protective gloves / protective clothing / eye protection / face protection.

P302/352 IF ON SKIN: Wash with plenty of water / soap and water.

P321: Specific treatment (see instructions on the label)  
P333/313: If skin irritation or rash occurs: get medical attention.

- The TMB Substrate contains an irritant, which harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to direct sunlight, metals or oxidants. Do not freeze the solution.

## 6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the kit is stable at 2 – 8°C for 5 months.
- The diluted wash solution is stable for 30 days at 2-8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

## 8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin, potassium EDTA) samples.

Sample Storage	Duration
2 – 8 °C	96 hours
Freeze/thaw cycles	3 cycles

## 9. PROCEDURE

### 9.1. Preparation of Calibrators and Controls

The Calibrators are ready for use and have the following concentration of Cortisol:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>
ng/mL	0	10	50	150	500

The calibrator is ready to use; the concentration of the Calibrator is printed on the label.

The Control is ready for use.

### 9.2. Preparation of the Conjugate

The conjugate is ready to use. Mix gently, for 5 minutes, on a roller mixer.

### 9.3. Preparation of the Wash Solution

Dilute the content of the vial "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

### 9.4. Preparation of Samples

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Avoid repetitive freezing and thawing of samples.

### 9.5. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2-8°C: avoiding long exposure to room temperature.

- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Calibrator C <sub>0</sub> -C <sub>4</sub>	20 µL		
Sample/Control		20 µL	
Conjugate	200 µL	200 µL	

Incubate 1 hour at 37°C.

Remove the contents from each well. Wash the wells 3 times with 300 µL of diluted wash solution.

**Important note:** during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

**Automatic washer:** in case you use an automatic washer, it is advised to do 6 washing steps.

TMB Substrate	100 µL	100 µL	100 µL
Stop Solution	100 µL	100 µL	100 µL

Incubate 15 minutes in the dark at room temperature (22±28°C).

Shake the microplate gently.

Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

## 10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit control provided in the kit should be tested as unknown and are intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis

relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories<sup>13</sup>.

### 11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, including Calibrator 0 is required.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

#### Conversion of units

To convert results to SI units:  
 $\text{nmol/L} = \text{ng/mL} \times 2.7586$

To convert results to mass units:  
 $\text{ng/mL} = \text{nmol/L} \times 0.3625$

### 12. MEASURING RANGE

The assay measuring range (AMR) is 7.16 – 500 ng/mL (19.75 – 1379.3 nmol/L).

Any value that reads below 7.16 ng/mL (19.75 nmol/L) should be reported as “< 7.16 ng/mL (19.75 nmol/L)”. Any value that reads above 500 ng/mL (1379.3 nmol/L) should be reported as “> 500 ng/mL (1379.3 nmol/L)”.

### 13. METROLOGY AND TRACEABILITY

The calibrators of this kit are traceable to a commercial cortisol assay which in turn is traceable to the Cortisol reference serum panel [ERM DA541\_IFCC].

### 14. EXPECTED VALUES

The following ranges were determined using the Cortisol ELISA and are provided for information only. The 95 % reference interval for apparently healthy adults were calculated by a non-parametric method following guidance from CLSI C28-A3 “Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory”.

	Morning 6 – 10 am	Evening 4 – 8 pm
No. of subjects	41	51
Mean (ng/mL)	108.41	66.81
SD (ng/mL)	42.16	25.06
Median (ng/mL)	112.64	61.95
Reference Interval (ng/mL)	47.1 - 187	31.4 – 114.2

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

### 15. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

#### 15.1. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, “Protocols for Determination of Limits of Detection and Limits of Quantitation” using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	3.47 ng/mL
Limit of Detection (LoD)	5.79 ng/mL
Limit of Quantitation (LoQ)	7.16 ng/mL

#### 15.2. Trueness

Trueness has been demonstrated through comparison of the Cortisol ELISA to a commercially available cortisol assay, which in turn is traceable to ID GC-MS, with a total of 34 individual donor native serum samples.

n	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
34	1.03 [0.98 – 1.10]	-0.13 [-0.62 – 0.36]	0.99

#### 15.3. Precision

Precision of the Cortisol ELISA was determined by performing a complex precision study.

**Repeatability:** A total of 5 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Data from one representative lot is shown below:

Sample	n	Mean Conc. (ng/mL)	Within run (Repeatability)	
			SD	CV%
1	75	38.38	3.43	8.9%
2	75	88.34	3.92	4.4%
3	75	156.08	13.08	8.4%
4	75	257.48	21.89	8.5%
5	75	379.41	31.18	8.4%
6	75	393.23	37.78	9.5%

**Reproducibility:** A total of 5 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Results for the combined data from two lots is shown below:

Sample	n	Mean Conc. (ng/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	150	37.93	4.63	12.2%
2	150	88.06	7.06	8.0%
3	150	156.98	18.26	11.6%
4	150	261.15	33.44	12.8%
5	150	385.80	52.37	13.6%
6	150	400.14	56.95	14.2%

#### 15.4. Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For cortisol concentration by Cortisol ELISA, the measurement procedure shows linearity for the interval from 7.16 to 500 ng/mL (19.75 to 1379.3 nmol/L) within the allowable deviation of linearity (ADL) of  $\pm 15\%$ .

#### 15.5. Method comparison

The Cortisol ELISA was compared against a commercially available quantitative automated assay, following CLSI EP-9A, "Method Comparison and Bias Estimation Using Patient Samples". A total of 100 samples, selected to represent a wide range of cortisol concentrations, was assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

n	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
100	1.06 [1.01 – 1.10]	-4.33 [-8.33 – 0.44]	0.98

#### 15.6. Analytical Specificity

The specificity was assessed with the following cross-reactants.

Cross-reactant	Concentration tested (unit)	Mean %Cross reactivity
Prednisolone	200 ng/ml	79%
11-Deoxycortisol	100 ng/ml	67%
21-Deoxycortisol	1 000 ng/ml	14%
Cortisone	1 000 ng/ml	14%
Cortiscosterone	1 000 ng/ml	11%
6-b-Hydroxycortisol	1 000 ng/ml	3%
17 $\alpha$ – OH Progesterone	10 000 ng/ml	2%
11 - Deoxycorticosterone	10 000 ng/ml	1%
Aldosterone	10 000 ng/ml	1%
Dexamethasone	10 000 ng/ml	1%
Prednisone	10 000 ng/ml	3%
Progesterone	10 000 ng/ml	1%

The following substances do not interfere with a bias of  $\geq \pm 15\%$  in the Cortisol ELISA assay when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	40 mg/dL
Bilirubin, unconjugated	40 mg/dL
Haemoglobin	200 mg/dL
Total Protein	15 g/dL
Triglyceride	1500 mg/dL

#### 15.7. Serum-plasma study

The Cortisol ELISA matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP9-A3 guidelines. A total of 20 samples (16 native, 4 spiked) to cover the assay range were evaluated. Passing-Bablok regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
SST	0.96 [0.91 – 1.04]	2.72 [-2.38 – 9.67]	1.0
Lithium Heparin	1.03 [0.91 – 1.18]	-2.87 [-13.10 – 5.81]	0.99
Sodium Heparin	0.99 [0.90 – 1.20]	0.14 [-13.69 – 10.67]	0.99
EDTA	1.00 [0.90 – 1.06]	-1.08 [-5.44 – 6.58]	0.99

### 15.8. Dilution recovery

For samples exceeding the analytical measurement interval (500 ng/mL), a dilution recovery study was performed. A total of 5 samples were prepared by spiking with a high concentration stock solution. The samples were then diluted and run in duplicate with the Cortisol ELISA.

Samples with analyte concentrations greater than 500 ng/mL can be diluted 1:4 with Calibrator 0. The assay result must then be corrected by applying the appropriate dilution factor.

### 16. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays<sup>14</sup>. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

### 17. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

### 18. BIBLIOGRAPHY

1. Ramamoorthy, S. & Cidlowski, J. A. Corticosteroids: Mechanisms of Action in Health and Disease. *Rheumatic diseases clinics of North America* 42, 15-31, doi:10.1016/j.rdc.2015.08.002 (2016).
2. Wild, D. *The immunoassay handbook*. (Elsevier, 2005).
3. Cain, D.W. & Cidlowski, J. A. Specificity and sensitivity of glucocorticoid signaling in health and disease. *Best practice & research. Clinical endocrinology & metabolism* 29, 545-556, 2015.04.007 (2015).
4. Schacke, H., Docke, W. D. & Asadullah, K. Mechanisms involved in the side effects of glucocorticoids. *Pharmacology & therapeutics* 96, 23-43 (2002).

5. Charmandari, E., Nicolaidis, N. C. & Chrousos, G. P. Adrenal insufficiency *Lancet* 383, 2152-2167, (13)61684-0 (2014).
6. Bancos, I., Hahner, S., Tomlinson, J. & Arlt, W. Diagnosis and management of adrenal insufficiency. *The lancet. Diabetes & endocrinology* 3, 216-226, (14)70142-1 (2015).
7. Arlt, W. & Allolio, B. Adrenal insufficiency. *Lancet* 361, 1881-1893, (03)13492-7 (2003).
8. Lacroix, A., Feelders, R. A., Stratakis, C. A. & Nieman, L. K. Cushing's syndrome. *Lancet* 386, 913-927, (14)61375-1 (2015).
9. Arnaldi, G. et al. Diagnosis and complications of Cushing's syndrome: a consensus statement. *The Journal of clinical endocrinology and metabolism* 88, 5593-5602 (2003).
10. Yaneva, M. et al. Midnight salivary cortisol for the initial diagnosis of Cushing's syndrome of various causes. *The Journal of clinical endocrinology and metabolism* 89, 3345-3351, (2004).
11. Krieger, D. T., Allen, W., Rizzo, F. & Krieger, H. P. Characterization of the normal temporal pattern of plasma corticosteroid levels. *The Journal of clinical endocrinology and metabolism* 32, 266-284, (1971).
12. Ceccato, F. et al. Performance of salivary cortisol in the diagnosis of Cushing's syndrome, adrenal incidentaloma, and adrenal insufficiency. *European journal of endocrinology / European Federation of Endocrine Societies* 169, 31-36, (2013).
13. Basic QC Practices On-line Course; <http://www.Westgard.com>.
14. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

### 19. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

### 20. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: [www.diametra.com](http://www.diametra.com).

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1. Ramamoorthy, S. & Cidlowski, J. A. Corticosteroids: Mechanisms of Action in Health and Disease. Rheumatic diseases clinics of North America 42, 15-31, doi:10.1016/j.rdc.2015.08.002 (2016).
2. Wild, D. The immunoassay handbook. (Elsevier, 2005).
3. Cain, D.W. & Cidlowski, J. A. Specificity and sensitivity of glucocorticoid signaling in health and disease. Best practice & research. Clinical

- endocrinology & metabolism 29, 545-556, 2015.04.007 (2015).
4. Schacke, H., Docke, W. D. & Asadullah, K. Mechanisms involved in the side effects of glucocorticoids. *Pharmacology & therapeutics* 96, 23-43 (2002).
  5. Charmandari, E., Nicolaidis, N. C. & Chrousos, G. P. Adrenal insufficiency *Lancet* 383, 2152-2167, (13)61684-0 (2014).
  6. Bancos, I., Hahner, S., Tomlinson, J. & Arlt, W. Diagnosis and management of adrenal insufficiency. *The lancet. Diabetes & endocrinology* 3, 216-226, (14)70142-1 (2015).
  7. Arlt, W. & Allolio, B. Adrenal insufficiency. *Lancet* 361, 1881-1893, (03)13492-7 (2003).
  8. Lacroix, A., Feelders, R. A., Stratakis, C. A. & Nieman, L. K. Cushing's syndrome. *Lancet* 386, 913-927, (14)61375-1 (2015).
  9. Arnaldi, G. et al. Diagnosis and complications of Cushing's syndrome: a consensus statement. *The Journal of clinical endocrinology and metabolism* 88, 5593-5602 (2003).
  10. Yaneva, M. et al. Midnight salivary cortisol for the initial diagnosis of Cushing's syndrome of various causes. *The Journal of clinical endocrinology and metabolism* 89, 3345-3351, (2004).
  11. Krieger, D. T., Allen, W., Rizzo, F. & Krieger, H. P. Characterization of the normal temporal pattern of plasma corticosteroid levels. *The Journal of clinical endocrinology and metabolism* 32, 266-284, (1971).
  12. Ceccato, F. et al. Performance of salivary cortisol in the diagnosis of Cushing's syndrome, adrenal incidentaloma, and adrenal insufficiency. *European journal of endocrinology / European Federation of Endocrine Societies* 169, 31-36, (2013).
  13. Basic QC Practices On-line Course; <http://www.Westgard.com>.

14. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

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**ERROR POSSIBLE CAUSES / SUGGESTIONS****No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

**Too low reaction (too low ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

**Too high reaction (too high ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

**Unexplainable outliers**

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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