DCM053-12  
Ed. 01/2021

# ALDOSTERONE ELISA

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# ALDOSTERONE ELISA

for routine analysis

Direct immunoenzymatic determination of Aldosterone in human serum, human plasma or urine.

IVD

LOT

See external label

2°C 8°C

$\Sigma = 96$  tests

REF DKO053

## INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of Aldosterone concentration in human serum, human plasma or urine.

Aldosterone ELISA kit is intended for laboratory use only.

## 1. CLINICAL SIGNIFICANCE

Aldosterone is a steroid hormone produced by the adrenal cortex in the adrenal gland, is the most potent mineralocorticoid in humans, it regulate sodium and potassium balance in the blood.

Aldosterone secretion appears to be stimulated primarily through the renin-angiotensin system.

Acting on mineralocorticoid receptors (MR) on principal cells in the collecting ducts of the kidneys, it increases the permeability of their apical (luminal) membrane to potassium and sodium and activates their basolateral Na<sup>+</sup>/K<sup>+</sup> pumps, stimulating ATP hydrolysis, reabsorbing sodium (Na<sup>+</sup>) ions and water into the blood, and excreting potassium (K<sup>+</sup>) ions into the urine. Aldosterone regulate plasma bicarbonate (HCO<sub>3</sub><sup>-</sup>) levels and its acid/base balance.

Aldosterone is responsible for the reabsorption of about 2% of filtered sodium in the kidneys.

Plasma aldosterone levels normally vary with body position (upright>supine) and salt intake. Overall plasma aldosterone levels show a circadian rhythm which is similar to but less marked than cortisol, with peak levels in the early morning; about 75% of the daily production is secreted between 04:00 am and 10:00 am each day. Age-related levels tend to decline from fetal through adult life.

Abnormally high plasma aldosterone concentrations can occur in adenomas, glucocorticoid-responsive hyperaldosteronism, idiopathic.

Abnormally low aldosterone secretion occurs in a number of conditions including salt-wasting forms of congenital adrenal hyperplasia, nephropathy, and renal tubular acidosis.

## 2. PRINCIPLE

The principle of this enzyme immunoassay test follow the typical competitive binding scenario. Competition occurs between an unlabeled antigen (present in calibrators, control and patient samples) and an

enzyme-labelled antigen (conjugate) for a limited number of antibody binding site on the microwell plate. The washing steps remove unbound materials. After the washing step, the enzyme substrate (TMB) is added. The enzymatic reaction is terminated by addition of the Stop Solution. The intensity of the colour is inversely proportional to the concentration of Aldosterone in the sample. The absorbance is measured on a microtiter plate reader. A set of calibrators is used to plot a calibration curve from which the amount of Aldosterone in patient samples and controls can be directly read.

## 3. REAGENTS, MATERIALS AND INSTRUMENTATION

### 3.1. Reagents and materials supplied in the kit

#### 1. Calibrators (6 vials, 1 mL each)

CAL0	REF DCE002/5306-0
CAL1	REF DCE002/5307-0
CAL2	REF DCE002/5308-0
CAL3	REF DCE002/5309-0
CAL4	REF DCE002/5310-0
CAL5	REF DCE002/5311-0

#### 2. Controls (2 vial, 1 mL each)

Control A	REF DCE045/5303A-0
Control B	REF DCE045/5303B-0

Control concentration is indicated on the Certificate of Analysis

#### 3. Conjugate (1 vial, 15 mL)

Aldosterone conjugated with horseradish peroxidase (HRP) REF DCE002/5302-0

#### 4. Coated Microplate (1 breakable microplate)

Microplate coated with anti aldosterone antibody REF DCE002/5303-0

#### 5. TMB Substrate (1 vial, 15 mL)

H<sub>2</sub>O<sub>2</sub>-TMB (0.26 g/L) (avoid any skin contact) REF DCE004-0

#### 6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (avoid any skin contact) REF DCE005-0

#### 7. 10X Conc. Wash Solution (1 vial, 50 mL)

0.2M Phosphate buffer pH 7,4 REF DCE054-0

### 3.2. Reagents necessary not supplied

Distilled water.

### 3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

#### Note

*Store all reagents between 2-8°C in the dark.*

*Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until expiry date of the kit.*

## 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of Proclin 300<sup>R</sup> as preservative. Avoid the contact with skin or mucosa.

#### **Hazard statement (Classification under CLP)**

H317 (Skin Sens. 1): may cause an allergic skin reaction

#### **Precautionary statements**

P261: avoid breathing dust / fumes / gas / mist / vapours / spray.

P272: contaminated work clothing should not be allowed out of the workplace.

P280: wear protective gloves / protective clothing / eye protection / face protection.

P302/352 IF ON SKIN: Wash with plenty of water, soap and water.

P321: specific treatment (see instructions on this label)

P333/313: if skin irritation or rash occurs, get medical attention.

- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of Aldosterone from 41,5 pg/mL (LOD) to 2000 pg/mL
- The treatment with natural or synthetic steroids can affect blood levels of Aldosterone.

## 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 6. STORAGE AND STABILITY

Store the kit at 2-8°C; the kit is stable until the expiry date claimed on the kit label and in the Certificate of Analysis.

Do not use the kit or its components after the expiry date.

## 7. PROCEDURE

### 7.1. Preparation of Calibrators and Controls

Leave on a rotating mixer for at least 5 minutes before using.

The Calibrators are ready to use and have the following concentration:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
pg/mL	0	20	80	300	800	2000

The Controls are ready to use.

Once opened, Calibrators and Controls are stable 6 months at 2-8°C.

### 7.2. Preparation of 10X Wash Solution

Dilute the content of the vial "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2÷8°C. In concentrated wash solution is possible to observe the presence of crystals, in this case mix at room temperature until complete dissolution of crystals; for greater accuracy dilute the whole bottle of concentrated wash solution to 500 mL on taking care also transfer the crystals, then mix until crystals are completely dissolved.

### 7.3. Preparation of Conjugate

The Conjugate is ready to use.

### 7.4. Preparation of the Sample

The determination of Aldosterone can be performed in human serum, human plasma or in urine.

#### – serum samples

SST tubes ("Serum Separation Tube") can be used to obtain serum without any interference with the assay.

#### – plasma samples

Plasma samples can be obtained with EDTA, Litium heparin and Sodium heparin without any interference with the assay.

#### – urine samples

To obtain urine samples use the following procedure:

1. pipet 250 µL of each urine sample into a polypropylene tube;
2. add 25 µL of 3.2 N HCl to every tube and mix thoroughly;  
Cap and incubate for 24h at room temperature (22-28°C) in the dark;
3. add 25 µL of 3.2N NaOH to every tube and mix thoroughly;
4. dilute 100 µL of the neutralized samples with 1400 µL of saline solution NaCl 0.9%;
5. use 50 µL of this solution to perform the assay as stated at par. 7.5.

Note: this pretreatment leads to a 1:18 dilution. The dilution factor 18 has to be taken into account for the calculation of the final concentration of the urine sample. Please refer to par. 9.

Store the sample at -20°C if the determination is not performed on the same day of the sample connection. Samples with concentration greater than 2000 pg/mL has not to be diluted; such samples has to be reported as "> 2000 pg/mL".

### 7.5. Procedure

- Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes. At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Samples/Control	Blank
Calibrator C <sub>0</sub> -C <sub>5</sub>	50 µL		
Samples/ Controls		50 µL	
Conjugate	100 µL	100 µL	
Incubate at +37°C for 1 hour. Remove the contents from each well; wash the wells 3 times with 300 µL of diluted wash solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22÷28°C) for 20 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 8. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Aldosterone for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of

kit reagents. Fresh reagents should be used to determine the reason for the variations.

## 9. RESULTS

### 9.1. Mean Absorbance

Calculate the mean of the absorbance (Em) for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>) and of each sample.

### 9.2. Calibration curve

Plot the values of absorbance (Em) of the Calibrators (C<sub>0</sub>-C<sub>5</sub>) against concentration. Draw the best-fit curve through the plotted points (es: Four Parameter Logistic).

### 9.3. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in pg/mL.

**Important note:** for urine samples the concentration in pg/ml read from the calibration curve has to be multiplied to the dilution factor 18.

To convert the Aldosterone concentration in “µg Aldosterone/24h”, calculate as above and correct for total volume of urine collected in 24 hours:

$$\text{pg/mL} \times \text{Vol (mL) urine 24 h} / 1.000.000 = \mu\text{g Aldosterone/24h}$$

## 10. REFERENCE VALUES

To determine the normal range for serum / plasma samples, EDTA plasma samples from 150 apparently healthy male and female adults have been tested. The posture of the patients, whether they were upright or supine before the collection of the specimen, is unknown.

Result:

Serum / plasma normal range
< 14,6 - 174 pg/mL

To determine the normal range for urine samples, 60 apparently healthy male and female adults have been tested; all the samples have been pretreated using the protocol stated at par. 7.4.

Result:

Urine normal range (24h)
2 - 23 µg/24h

Please pay attention to the fact that the determination of a range of expected values for a “normal” population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

## 11. PERFORMANCE AND CHARACTERISTICS

### 11.1. Analytical sensitivity

Analytical sensitivity has been investigated through the LOB (Limit of Blank), LOD (Limit of Detection), LOQ (limit of Quantitation) and Analytical Sensitivity (A.S.).

The table below shows the criteria of the study and the results obtained.

	Criteria	Results (pg/mL)
LOB	60 replicates of Cal 0, used as "Blank", have been investigated in 9 different sessions over 5 days	15,3
LOD	6 low serum samples have been analyzed in 10 different sessions over 5 days	41,5
LOQ	6 low serum samples have been analyzed in 10 different sessions over 5 days	79,1
A.S.	20 replicates of Cal 0 and 5 replicates of Cal 2 have been assayed. A.S. has been calculated by linear regression.	14,6

### 11.2. Precision and reproducibility (Complex Precision)

To determine precision and reproducibility 5 different serum samples have been used.

The table below shows the Within Run and Total CV%.

Sample	n°	Mean (pg/mL)	Within Run CV%	Total CV%
PS1	20	1410,102	5,1%	7,2%
PS2	20	824,239	5,2%	8,4%
PS3	20	483,272	5,7%	8,4%
PS4	20	257,171	6,8%	12,3%
PS5	20	127,102	10,4%	13,9%

### 11.3. Sample type study

Sample type study has been evaluated by performing the assay in parallel with serum and plasma belonging to the same subject.

20 different subjects have been used.

Serum obtained with SST tubes, and plasma obtained with EDTA, lithium-heparin and sodium-heparin have been investigated.

The interference has been evaluated as "significant" if it causes a concentration (bias) >10% between the sample and the reference (EDTA plasma).

The table below shows the results obtained:

Sample	Bias	Interference
EDTA plasma	reference	/
Serum	9,1%	No
Serum (SST pipes)	6,4%	No
Lithium heparin plasma	5,1%	No
Sodium heparin plasma	7,6%	No

Conclusion: following the study, serum and plasma samples can be used indifferently with this kit.

Urine samples can be used according to par. 7.4.

## 11.4. Analytical specificity

### 11.4.1. Interfering substances

Interference by Conjugated Bilirubin, Hemoglobin and Triglycerides has been investigated by adding the interfering substance to the serum sample and by comparing its concentration to the unspiked sample. Samples with low and high concentration of Aldosterone have been analyzed.

The interference has been evaluated as "significant" if it causes a concentration bias >10% between spiked and unspiked sample.

The table below shows the results obtained:

Substance	Conc. assayed	Interference
Triglycerides	600 mg/dL	No
Bilirubin, conjugated	33,1 mg/dL	No
Hemoglobin	6,6 mg/dL	No

Conclusion: following the study, there is no significant interference by Conjugated Bilirubin, Hemoglobin and Triglycerides at the concentration assayed.

### 11.4.2. Cross-reactivity

The following cross-reactants have been spiked in serum samples with low and high Aldosterone concentration. The respective low and high unspiked samples were used as reference.

The interference has been evaluated as "significant" if it causes a concentration bias >10% between spiked and unspiked sample.

The table below shows the results obtained:

Cross-reactant	Conc. assayed	Interference
Corticosterone	0,1 µg/mL	No
Androsterone	5 µg/mL	No
Cortisone	0,1 µg/mL	No
Diidrotosterone (DHT)	0,5 µg/mL	No
Estradiol	10 µg/mL	No
Estrone	10 µg/mL	No
Estriol	10 µg/mL	No
Testosterone	0,2 µg/mL	No
DHEA-S	10 µg/mL	No

## 11.5. Accuracy

### 11.5.1. Linearity

Three tests have been performed using 3 pair of serum (with low and high Aldosterone concentration); 11 dilutions spanning the range low to high have been analyzed.

The assay is linear over the analyzed range (approximately 85,8 - 1520,9 pg/mL).

### 11.5.2. Recovery

Three tests have been performed using 3 serum samples with low aldosterone concentration, each one spiked with the Cal 5 in different percentages.

The mean recovery in the three samples has been 109%, 102% and 111%.

## 11.6. Correlation

### – Plasma samples

126 EDTA plasma samples have been tested with the Dia.Metra Srl Aldosterone ELISA kit and with the IDS - iSYS Aldosterone assay (reference method).

The linear regression curve is:

$$Y = 1,106 \cdot X - 36,06$$

$$r^2 = 0,94$$

### – Urine samples

48 urine samples have been tested with the Dia.Metra Srl Aldosterone ELISA kit and with a commercially available Aldosterone assay (reference method).

The results obtained from each assay were multiplied to the respective dilution factor and converted to ng/mL and a Passing Bablok has been performed; the final result is:

$$Y = 1,133 \cdot X + 0,8645$$

$$r^2 = 0,87$$

## 12. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

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**ERROR POSSIBLE CAUSES / SUGGESTIONS****No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

**Too low reaction (too low ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

**Too high reaction (too high ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

**Unexplainable outliers**

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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