DCM048-8
Ed. 09/2018

THYROGLOBULIN ELISA

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BIBLIOGRAFIA

1. Beever K, Bradbury J, Phillips D, et al, "Highly sensitive assays of autoantibodies to Thyroglobulin and Thyroid Peroxidase", *Clin Chem*, **35**, 1949-1954 (1989).
2. Ladenson PW, "Optimal laboratory testing for diagnosis and monitoring of thyroid nodules, goiter, and thyroid cancer", *Clin Chem*, **42**, 183-187 (1996).
3. Mayo Medical Laboratories: test Catalog, Rochester, MN (1997).
4. Spencer CA, Takeuchi M, Kazarosyn M, "Current status and performance goals for serum thyroglobulin assays", *Clin Chem*, **42**, 164-173 (1996).
5. Tietz N. Ed: *Clinical Guide to Laboratory Tests*. 3rd Ed. Philadelphia. Saunders (1995).
6. Surks, MI, Chopra, IJ, Mariash, CN, "American Thyroid Association guidelines for use of laboratory tests in thyroid disorders", *JAMA*, **263**, 1529-1532 (1990).
7. Ng, M., Rajna, A, Khalid, B, "Enzyme immunoassay for simultaneous measurements of autoantibodies against thyroglobulin and thyroid

- microsomes in serum”, *Clin Chem*, **33**, 2286-2288 (1987).
8. Spencer, CA, Takeuchi M, Kazarosyn M, Wang CC, Guttler RB, Singer PA, et al, “Serum thyroglobulin autoantibodies; prevalence, influence on serum thyroglobulin measurements, and prognostic significance in patients with differentiated thyroid carcinoma”, *J Clin Endocrinol Metab*, **83**, 1121-27 (1998).
 9. Spencer CA, LoPresti JS, Fatemi S, Nicoloff JT, “Detection of residual and recurrent differentiated thyroid carcinoma by serum thyroglobulin measurements”, *Thyroid*, **9**, 435-41 (1999).
 10. Schlumberger M, Baudin E, “Serum thyroglobulin determinations in the follow up of patients with differentiated thyroid carcinoma”, *Eur J. Endocrinol*, **138**, 249-252 (1998).

Ed. 09/2018

DCM048-8

Dia.Metra S.r.l.

Via Pozzuolo 14, 06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: info@diametra.com

Example Version

DCM048-8
Ed. 09/2018

THYROGLOBULIN ELISA

for routine analysis

Determination of Thyroglobulin in human serum.

IVD

LOT

See external label

Σ = 96 tests

REF DKO048

INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of Thyroglobulin concentration in human serum.

Thyroglobulin ELISA kit is intended for laboratory use only.

1. CLINICAL APPLICATIONS

Thyroglobulin (TG), a glycoprotein with a molecular weight of about 660,000 Daltons, is the thyroid's main iodine protein and the most important compound of follicular colloid. Thyroglobulin is the form under which the active hormones T₃ and T₄ together with their immediate forerunners MIT and DIT are laid inside the thyroid gland. The clinical applications of the TG dosage seem to originate from its specificity for the thyroid and related cells.

The dosage of TG can be used as support to scintigraphies or other techniques for studying pathogenesis, making a diagnosis and analyzing the course of thyroid disorders.

The dosage of TG before and after replacement treatment with L-Thyroxin cannot be established in cases of hypothyroidism due to thyroid agenesis. In cases of secondary hypothyroidism with a dysglandular goiter or ectopic thyroid, the levels of TG are normal or high. The circulating levels of TG tend to increase in several thyroid disorders such as toxic and atoxic goiter, subacute thyroiditis, Basedow's disease and carcinoma. In Basedow's disease the TG dosage is a potentially interesting index of normalization of hyperthyroidism in patients treated with anti-thyroid drugs. In the oncology field and more specifically for differentiated thyroid carcinoma, there are very promising applications linked to the ability of thyroid tumors tissues to concentrate iodine and synthesize TG as a normal thyroid. Basically the dosage of TG can be used as follows:

a. Pre-operating diagnosis of thyroid tumors.

This application does not allow the differentiated diagnosis of the tumor as the values of TG seen in malignant and benign nodules are superimposable.

b. Post-operation monitoring

In patients treated surgically or with radiotherapy, long lasting TG levels suggest the presence of a residual carcinoma and/or carcinoma with

metastasis.

c. Monitoring of totally thyroidectomized patients


The use of circulating TG as an indicator of recurrent tumors (metastasis marker) has an established clinical value: the increase of Thyroglobulinaemia indicates the need to undergo further analysis for confirming the diagnosis. Interesting advantages can come from: a) a reduced use of scintigraphic diagnostic techniques as they imply regular suspension of replacement treatment and frequent exposure to radiation, b) and complete completion of the information obtained via scintigraphy

2. PRINCIPLE

The essential reagents required for an immunoenzymometric assay include high affinity and specificity antibodies (enzyme and immobilized), with different and distinct epitope recognition, in excess, and native antigen. In this procedure, the immobilization takes place during the assay at the surface of a microplate well through the interaction of the streptavidin coated on the well and exogenously added biotinylated monoclonal Thyroglobulin antibody.

When monoclonal biotinylated antibody is mixed with a serum containing the Thyroglobulin antigen and the antibody, an antibody-antigen complex is formed.

The interaction is illustrated in the following equation:

BtAb(m) = biotinylated monoclonal antibody (excess quantity)

Ag(Tg) = native antigen (variable quantity)

Ag(Tg)-BtAb(m) = antigen-antibody complex (variable quantity)

k_a = rate constant of association

k_{-a} = rate constant of dissociation

Simultaneously, the complex is deposited into the well through the high affinity reaction of streptavidin and biotinylated antibody.

This interaction is illustrated below:

Ag(Tg)-BtAb(m) + Streptavidin C.W. → Immobilized complex

Streptavidin C.W. = Streptavidin immobilized on well
Immobilized complex = sandwich complex bound to the well

After a suitable incubation period, the antibody-antigen bound fraction is separated from unbound antigen by a washing step. Another antibody (directed at a different epitope) labelled with an enzyme is added. Another interaction occurs to form and enzyme labelled antibody-antigen-biotinylated-antibody complex on the surface of the wells.

Excess enzyme is washed off via a wash step. A suitable substrate is added to produce color measurable with the use of a microplate spectrophotometer.

The enzyme activity in the well is directly proportional to the native antigen concentration. By utilizing several different serum references of known antigen concentration, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

3. REAGENTS, MATERIALS AND INSTRUMENTATIONS

3.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials, 1 mL each)

CAL0	REF DCE002/4806-0
CAL1	REF DCE002/4807-0
CAL2	REF DCE002/4808-0
CAL3	REF DCE002/4809-0
CAL4	REF DCE002/4810-0
CAL5	REF DCE002/4811-0

2. Biotin Reagent (1 vial, 13 mL)

Biotinylated monoclonal anti Tg antibody
REF DCE019/4819-0

3. Enzyme Reagent (1 vial, 13 mL)

Anti Tg IgG antibody labeled with horseradish peroxidase (HRP)
REF DCE002/4802-0

4. Coated Microplate (1 breakable microplate)

Microplate coated with streptavidin
REF DCE002/4803-0

5. TMB Substrate (1 vial, 12 mL)

H₂O₂-TMB 0,26 g/L (avoid any skin contact)
REF DCE004/4804-0

6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0,15 mol/L (avoid any skin contact)
REF DCE005-0

7. 50X Conc. Wash Solution (1 vial, 20 mL)

Buffered solution
REF DCE006/4806-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

Note

Store all reagents between 20-8°C in the dark.

Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable 2 months at 2-8°C.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the reagents should be handled in the same manner as potentially infectious material.
- Some reagents contain small amounts of Proclin 300^R as preservatives. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.

- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems, it is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibrators (C₀...C₅)

There is no known, internationally accepted thyroglobulin standard available. The thyroglobulin used in the serum based calibrators is a highly purified human Tg preparation that is calibrated gravimetrically against the reference material obtained from Community Bureau of Reference CRM 457.

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
ng/mL	0	2	10	40	100	250

Once opened, the Calibrators are stable 2 months at 2-8°C.

6.2. Preparation of the Sample

The assay can be performed in human serum. The usual precautions in the collection of venipuncture samples should be observed. For accurate comparison to established normal values, a fasting morning serum sample should be obtained. The blood should be collected in a plain red-top venipuncture tube without additives or gel

barrier. Allow the blood to clot, centrifuge the specimen to separate the serum from the cells. Samples may be refrigerated at 2-8°C for a maximum period of 5 days. If the specimens cannot be assayed within this time, they may be stored at temperatures of -20°C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, 0.100 mL of the specimen is required.

6.3. Preparation of Wash Solution

Dilute the content of each vial of the "50X Conc. Wash Solution" with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 60 days at room temperature (22-28°C).

6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each sample, one for Blank.

Reagent	Calibrator	Sample	Blank
Calibrator C ₀ -C ₅	50 µL		
Sample		50 µL	
Biotin Reagent	100 µL	100 µL	

Swirl the microplate gently for 20-30 seconds to mix. Cover with a plastic wrap or microplate cover. Incubate at 37°C for 90 minutes. Remove the contents from each well; wash the wells 3 times with 0.350 mL of diluted wash solution.

Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: if you use automated equipment, wash the wells at least 5 times.

Enzyme Reagent	100 µL	100 µL	
<p>Do not shake the plate after enzyme addition. Cover with a plastic wrap or microplate cover. Incubate at 37°C for 60 minutes. Remove the contents from each well; wash the wells 3 times with 0.350 mL of diluted wash solution.</p> <p>Washing: follow the same indications of the previous point.</p>			
TMB Substrate	100 µL	100 µL	100 µL
<p>Do not shake the plate after substrate addition. Cover with a plastic wrap or microplate cover. Incubate at room temperature (22-28°C) for 15 minutes in the dark.</p>			
Stop Solution	100 µL	100 µL	100 µL
<p>Shake the microplate gently. Read the absorbance (E) at 450 nm using a reference wavelength of 620-630 nm or against Blank within 30 minutes.</p>			

7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of hypothyroid, euthyroid and hyperthyroid for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. CALCULATION OF RESULTS

A dose response curve is used to ascertain the concentration of human thyroglobulin (Tg) in unknown specimens.

- Record the absorbance obtained from the printout of the microplate reader.
- Plot the absorbance for each duplicate serum reference versus the corresponding Tg concentration in ng/mL on linear graph paper (do not average the duplicates of the serum references before plotting)
- Draw the best fit curve through the plotted points.
- To determine the concentration of Tg for an unknown locate the average absorbance of the duplicates for each unknown on the vertical axis

of the graph, find the intersecting point on the curve, and read the concentration (in ng/mL) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance (0.424) intersects the dose response curve at 25.2 ng/mL Tg concentration.

Note: computer data reduction software designed for ELISA assays may also be used for the data reduction. If such software is utilized, the validation of the software should be ascertained.

8.1. Example of Calculation

The values shown below must be considered as an example and must not be used in place of experimental data.

Calibrator /Sample	OD 450 nm	Mean OD	Value Tg ng/mL
Cal 0	0.047	0.047	0
	0.047		
Cal 1	0.093	0.091	2
	0.090		
Cal 2	0.221	0.217	10
	0.214		
Cal 3	0.612	0.625	40
	0.634		
Cal 4	1.343	1.339	100
	1.335		
Cal 5	2.596	2.577	250
	2.557		
Sample	0.426	0.424	25.2
	0.422		

9. QUALITY CONTROL PARAMETERS

In order for the assay results to be considered valid the following criteria should be met.

- The absorbance (OD) of Calibrator 0 should be 0.10.
- The absorbance (OD) of Calibrator 5 should be 1.3.
- Respect the assigned ranges indicated on the Certificate of Analysis.

10. REFERENCE VALUES

Based on clinical data gathered in concordance with the published literature a normal range was established:

Thyroglobulin (ng/mL)	
Normal (adult)	3.5 - 56 ng/mL

Thyroglobulin is found to be elevated in patients with thyroid follicular and papillary carcinoma, thyroid adenoma, subacute thyroiditis, Hashimoto's thyroiditis and graves' disease. Low levels of Tg are an indication of thyrotoxicosis factitia.

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent

on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

11. PERFORMANCE AND CHARACTERISTICS

11.1. Precision

The intra and inter assay precisions were determined by analyses on three different levels of pool control sera:

Intra assay				
Serum	Mean ng/mL	SD	% CV	Replicates
1	6.2	0.41	6.6 %	22
2	64.4	2.23	3.6 %	22
3	194.1	8.17	4.2 %	22

Inter assay				
Serum	Mean ng/mL	SD	% CV	Replicates
1	5.8	0.52	9.0 %	10
2	62.2	3.82	6.1 %	10
3	192.3	10.90	5.7 %	10

11.2. Sensitivity

The analytical sensitivity (detection limit) was calculated by determining the variability of the 0 ng/mL serum calibrator and using 2 standard deviation (95% certainty) statistics to calculate the minimum dose. The assay sensitivity was found to be 0.44 ng/mL.

11.3. Accuracy

Dia.Metra Thyroglobulin ELISA assay was compared with a reference coated tube radioimmunoassay (IRMA). Biological specimens from population (symptomatic and asymptomatic) were used.

The linear regression curve is:

$$Y = 0.908 \cdot X + 2.55$$

$$(R^2=0.975)$$

11.4. Specificity

The cross-reactivity of Dia.Metra Thyroglobulin assay to selected substances was evaluated by adding the interfering substance to a serum matrix at the following concentrations. The cross-reactivity was calculated by deriving a ratio between dose of interfering substance to dose of Thyroglobulin needed to produce the same absorbance.

Substance	Conc.	Cross-reactivity
Thyroglobulin	100 ng/mL	100 %
Triiodothyronine	1000 ng/dL	ND
Thyroxine	1000 ng/mL	ND
TBG	100 ng/mL	ND

11.5. High Dose Effect

Since the assay is sequential in design, high concentrations of Tg do not show the hook effect. Samples with concentration over 50000 ng/mL demonstrated extremely high levels of absorbance.

12. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

1. Beever K, Bradbury J, Phillips D, et al, "Highly sensitive assays of autoantibodies to Thyroglobulin and Thyroid Peroxidase", *Clin Chem*, **35**, 1949-1954 (1989).
2. Ladenson PW, "Optimal laboratory testing for diagnosis and monitoring of thyroid nodules, goiter, and thyroid cancer", *Clin Chem*, **42**, 183-187 (1996).
3. Mayo Medical Laboratories: test Catalog, Rochester, MN (1997).
4. Spencer CA, Takeucho M, Kazarosyn M, "Current status and performance goals for serum thyroglobulin assays", *Clin Chem*, **42**, 164-173 (1996).
5. Tietz N. Ed: *Clinical Guide to Laboratory Tests*. 3rd Ed. Philadelphia. Saunders (1995).
6. Surks, MI, Chopra, IJ, Mariash, CN, "American Thyroid Association guidelines for use of laboratory tests in thyroid disorders", *JAMA*, **263**, 1529-1532 (1990).
7. Ng, M., Rajna, A, Khalid, B, "Enzyme immunoassay for simultaneous measurements of autoantibodies against thyroglobulin and thyroid microsomes in serum", *Clin Chem*, **33**, 2286-2288 (1987).
8. Spencer, CA, Takeucho M, Kazarosyn M, Wang CC, Guttler RB, Singer PA, et al, "Serum thyroglobulin autoantibodies; prevalence, influence on serum thyroglobulin measurements, and prognostic significance in patients with differentiated thyroid carcinoma", *J Clin Endocrinol Metab*, **83**, 1121-27 (1998).
9. Spencer CA, LoPresti JS, Fatemi S, Nicoloff JT, "Detection of residual and recurrent differentiated thyroid carcinoma by serum thyroglobulin measurements", *Thyroid*, **9**, 435-41 (1999).
10. Schlumberger M, Baudin E, "Serum thyroglobulin determinations in the follow up of patients with differentiated thyroid carcinoma", *Eur J. Endocrinol*, **138**, 249-252 (1998).

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ERROR POSSIBLE CAUSES / SUGGESTIONS**No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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