Total PSA

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14. BIBLIOGRAFIA

- Ornstein DK, Pruthi RS. Prostate-specific antigen. Expert Opin Pharmacother. 2000 Dec;1(7):1399-411.
- Kostova MB, Brennen WN, Lopez D, Anthony L, *et al.* PSA-alpha-2-macroglobulin complex is enzymatically active in the serum of patients with advanced prostate cancer and can degrade circulating peptide hormones. Prostate. 2018 Aug;78(11):819-829. doi: 10.1002/pros.23539. Epub 2018 Apr 16. PMID: 29659051.
- Stamey TA, Yang N, Hay AR, McNeal JE, *et al.* Prostate-specific antigen as a serum marker for adenocarcinoma of the prostate. N Engl J Med. 1987 Oct 8;317(15):909-16.
- Albertsen PC. Prostate cancer screening with prostate-specific antigen: Where are we going? Cancer. 2018 Feb 1;124(3):453-455. doi: 10.1002/cncr.31140. Epub 2017 Dec 12. PMID: 29231972.
- Thompson IM, Pauler DK, Goodman PJ, Tangen CM, *et al.* Prevalence of prostate cancer among men with a prostate-specific antigen level < or =4.0 ng per milliliter. N Engl J Med. 2004 May 27;350(22):2239-46. doi:

- 10.1056/NEJMoa031918. Erratum in: N Engl J Med. 2004 Sep 30;351(14):1470. PMID: 15163773.
6. Schröder FH, Hugosson J, Roobol MJ, Tammela TL, *et al*; ESRPC Investigators. Prostate-cancer mortality at 11 years of follow-up. N Engl J Med. 2012 Mar 15;366(11):981-90. doi: 10.1056/NEJMoa1113135. Erratum in: N Engl J Med. 2012 May 31;366(22):2137. PMID: 22417251; PMCID: PMC6027585.
7. Moyer VA; U.S. Preventive Services Task Force. Screening for prostate cancer: U.S. Preventive Services Task Force recommendation statement. Ann Intern Med. 2012 Jul 17;157(2):120-34. doi: 10.7326/0003-4819-157-2-201207170-00459. PMID: 22801674.
8. Li B. Prostate cancer targeted therapy. In: Schwab M, ed. *Encyclopedia of Cancer*. Berlin, Germany: Springer; 2012:3062–3064.

Ed. 09/2022

DCM137-6

Legal Manufacturer:

Dia.Metra Srl
Via Pozzuolo 14
06038 SPELLO (PG) Italy
Tel. +39-0742-24851
Fax +39-0742-316197
E-mail: info@diametra.com

Example Version

Total PSA

for routine analysis

Enzyme immunoassay for the quantitative determination of total PSA (Prostate Specific Antigen) in human serum and plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 test

REF DKO137

INTENDED USE

Diametra Total PSA ELISA is used for the quantitative determination of total Prostate Specific Antigen (t-PSA) in human serum or plasma (EDTA, lithium heparin or citrate plasma) samples. The determination of total PSA levels is used to help estimate the risk of prostate carcinoma in men in conjunction with digital rectal examination (DRE) or to help monitor the effectiveness of prostate carcinoma treatment in patients.

1. CLINICAL SIGNIFICANCE

PSA is a serine protease with chymotrypsin-like activity, whose primary function is seminal clot liquification. The production of PSA is regulated by androgens and occurs almost exclusively in the prostatic epithelial cells¹. PSA is found in circulation in a number of different forms which can be classified into 2 general groups – complexed PSA (bound to serum protease inhibitors) and free PSA (unbound and inactive)². The predominant form of PSA is that bound to the serum protease inhibitor alpha-1 anti-chymotrypsin and comprises 65 – 95% of the circulating level²; the remaining 5 – 35% of PSA in circulation is free PSA and that complexed with alpha-2 macroglobulin.

PSA is produced by normal, as well as malignant, cells of the prostate gland and a number of benign (not cancerous) conditions can cause a man's PSA level to rise. The most frequent benign prostate conditions that cause an elevation in PSA level are prostatitis (inflammation of the prostate) and benign prostatic hyperplasia (BPH) (enlargement of the prostate). However, prostate cancer is most commonly associated with elevated levels of PSA and since the correlation between PSA and prostate cancer volume was first published in 1987³, PSA has become the most widely used tumour marker for prostate cancer⁴.

Previously digital rectal examination (DRE) was frequently used as only diagnostic modality for the detection of early stages of prostate cancer. More recently assessment of serum PSA in conjunction with digital rectal examination (DRE) has been found to increase the chance of early detection of prostate cancer improving the diagnostic specificity of DRE. In male serum, the normal concentration of PSA is < 4 ng/mL with levels > 4 ng/mL indicative of carcinoma⁵.

The determination of PSA serum levels is not only important for the screening of patients for prostate cancer, but also for monitoring patients that have been treated for this disease. Regular PSA measurements are an important tool to examine the potential and actual effectiveness of surgery or other

therapies. An increase of PSA in patients after radical prostatectomy or radiotherapy may allow an earlier discovery of residual or recurrent carcinoma⁶⁻⁸.

2. PRINCIPLE OF THE METHOD

The Total PSA assay is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the sandwich principle. The microtiter wells are coated with a monoclonal antibody, directed towards a specific epitope on the PSA molecule. An aliquot of patient serum is incubated in the coated plate; an enzyme (horseradish peroxidase – HRP) labelled anti- PSA antibody, which is directed towards a different region of the antigen molecule) is then added. Following the incubation the wells are washed to remove unbound material; colour is developed using a chromogenic substrate (TMB). The colour reaction is stopped and the absorbance of the stopped reaction mixture is read in a microtitre plate reader. The intensity of colour developed is proportional to the antigen concentration in the sample. The measured ODs of the Calibrators are used to construct a calibration curve against which the unknown samples are calculated.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

1. PSA Calibrators (6 vials, 0.5 mL each)

CAL0 / Sample Diluent (10 mL)	REF DCE002/13706-0
CAL1	REF DCE002/13707-0
CAL2	REF DCE002/13708-0
CAL3	REF DCE002/13709-0
CAL4	REF DCE002/13710-0
CAL5	REF DCE002/13711-0

Contain non-mercury preservative.

2. Controls (2 vials, 0.5 mL each)

For concentration see the Certificate of Analysis.

Contains non-mercury preservative.

Low Control	REF DCE045/13701-0
High Control	REF DCE045/13702-0

3. Coated Microplate (1 breakable microplate)

Monoclonal antibody anti PSA adsorbed on microplate
REF DCE002/13703-0

4. Enzyme Conjugate (1 vial, 12 mL)

Monoclonal antibody anti PSA conjugated to horseradish peroxidase. Contains non-mercury preservative.

REF DCE002/13702-0

5. Substrate Solution (1 vial, 12 mL)

Contains TMB (tetramethylbenzidine) and H₂O₂

REF DCE004/13704-0

6. Stop Solution (1 vial, 14 mL)

Contains sulphuric acid (H₂SO₄) $\leq 5\%$ (exact concentration: 0.5M)

(avoid contact with the Stop Solution, it may cause skin irritations and burns)

REF DCE005/13705-0

3.2. Material required but not provided

- Distilled water.
- Precision micropipettes (volume: 25 μ L and 100 μ L) with disposable tips.
- ELISA photometer with 450 nm and 630 nm filters.
- Timer with 60 min. range or higher.
- Container for the proper handling of waste and samples after use.

3.3. Auxiliary materials and instrumentation

- Microplate washer.
- Vortex or similar mixing tools.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.

- Follow Good Laboratory Practice (GLP) for handling blood products.

- Use appropriate personal protective equipment while working with the reagents provided.

Serum and plasma samples should be treated as potentially infectious materials. Wear gloves and proper laboratory attire when handling sample materials. Do not eat, drink or smoke in areas where specimen or kit reagents are handled. Do not pipette with the mouth. In case of skin contact, wash with a germicidal soap and copious amounts of water. Seek medical advice if indicated.

- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the reagents should be handled in the same manner as potentially infectious material.

- Due to the potentially infectious character of samples and kit components all materials that have come in contact with these materials should be sterilized and disposed of according to local regulations. This also includes the liquid waste.

- The assay reagents contain preservatives, TMB, H₂O₂ or sulphuric acid and may be harmful if ingested. A direct skin or mucosa contact should be avoided. In case of skin contact, wash thoroughly with water and seek medical attention if required.

- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.

- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.

- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.

Furthermore, the instrumentation employed to dispense it should be thoroughly cleaned after use.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.

- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.

- Allow all kit components and specimens to reach room temperature (20-26°C) and mix well prior to use.

- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.

- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.

- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.

- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate

- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.

- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.

- Maximum precision is required for reconstitution and dispensation of the reagents.

- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.

- Plate readers measure vertically. Do not touch the bottom of the wells.

- Avoid reagent or sample carry-over by using fresh tips for solutions and samples.

- Do not use test kit if zip lock pouch or bottles have been damaged.

6. STORAGE AND STABILITY

Store the kit at 2-8°C; the kit is stable until the expiry date claimed on the kit label and in the Certificate of Analysis.

Do not use the kit or its components after the expiry date.

After opening, the kit is stable for 8 weeks at 2-8°C.

7. PROCEDURE

7.1. Preparation of Calibrators and Controls

The Calibrators are calibrated against the International Standard WHO 96/670 and have the following concentrations:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
ng/mL	0	1.56	3.12	6.25	12.50	25.00

Calibrators and controls are ready to use.

7.2. Preparation of the Conjugate

The component "Enzyme Conjugate" is ready to use.

7.3. Preparation of Samples

The kit can be used on serum and plasma samples.

Blood samples are collected by venipuncture. As different factors could influence the PSA level in blood, doctors should ensure that the patient has avoided the following conditions before taking the blood sample:

- conditions that may lead to an increase of PSA levels:
 - biking
 - sexual intercourse (ejaculation)
 - manipulation of the prostate during medical examinations like DRE, transrectal prostatic ultrasound etc.
 - prostatitis
 - liver dysfunction
- conditions that may lead to a decrease of PSA levels:
 - intake of 5-alpha-reductaseinhibitors, antiandrogens, or GnRH analogues

The preparation of serum or plasma samples is performed according to standard techniques. Serum or plasma should be prepared as soon as possible to avoid hemolysis and to improve the stability of PSA.

If not used immediately, the samples can be stored at 2-8°C or -20°C for 1 week.

5 freezing / thawing cycles of samples are permitted.

Samples with an expected PSA value higher than 25 ng/mL should be diluted with the "Cal 0 / Sample diluent" or reported as "> 25 ng/mL".

7.4. Procedure

- **Allow all reagents to reach room temperature (20-26°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C; avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve, two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample	Control
Calibrator C ₀ -C ₅	25 µL		
Control			25 µL
Sample		25 µL	
Incubate 5 minutes at room temperature (20-26°C).			
Enzyme Conjugate	100 µL	100 µL	100 µL

Mix by moving plate on the table (10 sec)
Incubate 1 hour at room temperature (20-26°C).
Remove solution from the wells by aspirating the liquid or by decanting it. If decanting, tap plate on adsorbent paper to remove residual liquid.

For washing fill plate with distilled water and wait 15 sec before removing the distilled water; wash 5x.

We recommend the following procedure: wash wells 5-times with 300 or 400 µL / well distilled water. Preferably use an automated washing procedure. If washing manually take care that the washing solution remains in each well for the same time. This is necessary to receive lowest possible CV-values.

Substrate Solution	100 µL	100 µL	100 µL
Incubate 20 minutes at room temperature (20-26°C) in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read absorbencies (OD) at 450 nm (blanking 630 nm)			

8. QUALITY CONTROL

- It is recommended that internal controls are used in every assay in duplicate. Control results should be within established ranges and should preferably represent low, medium, and high concentrations.
- The risk for the patient is mainly due to false negative results around the threshold (cut-off) value of PSA < 4.0 ng/mL. It is therefore highly recommended to validate the kit and laboratory via external trials (e.g. DGKC).

9. CALCULATION OF RESULTS

1. Calculate the average absorbance values for each set of Calibrators, controls and patient samples.
2. Using linear graph paper, construct a calibration curve by plotting the mean absorbance obtained from each Calibrator against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the calibration curve.
4. Automated method: The results in the IFU have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. 4 Parameter Logistics is the preferred method. Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this calibration curve. Samples with concentrations higher than that of the highest Calibrator have to be further diluted or reported as > 25 ng/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

9.1. Example of Typical Calibration Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Calibrator	Optical Units (450 nm)
Calibrator 0 (0 ng/mL)	0.05
Calibrator 1 (1,56 ng/mL)	0.24
Calibrator 2 (3,12 ng/mL)	0.39
Calibrator 3 (6,25 ng/mL)	0.74
Calibrator 4 (12.5 ng/mL)	1.27
Calibrator 5 (25,0 ng/mL)	2.01

10. REFERENCE VALUES

The generally recommended threshold for follow-up examinations is:

Cut-off value	4.0 ng PSA / mL
---------------	-----------------

Healthy men generally have a PSA concentration lower than 4.0 ng/mL. If the PSA concentration is equal or higher than 4.0 ng/ml follow-up examinations are highly recommended. This PSA concentration indicates an elevated risk for prostate cancer but might also be caused by BPH.

Please note that the 4 ng/ml threshold is only a guideline value. In the literature it is discussed that modifications according to age and ethnological background might be useful e.g. that for younger men the threshold should be lower than for older men. If possible, it is recommended for each laboratory to establish its own specific values that take into consideration a population indigenous to the area where the laboratory is located.

It is important to keep in mind that some prostate tumors do not cause elevated PSA levels so that PSA measurements should never replace DRE but should only be used in conjunction with DRE.

As elevated PSA levels might also be caused by non-cancerous conditions follow examinations might try to increase the diagnostic specificity of t-PSA values. In the literature PSA density, PSA velocity and the ratio of f-PSA to t-PSA are discussed to improve discrimination between cancerous and non-cancerous conditions and might be used to reduce unnecessary prostate biopsies. But only a prostate biopsy can finally show if a prostate carcinoma is present or not.

Note: PSA values can only be used to estimate the cancer risk. They should always be interpreted in conjunction with other clinical findings and should not be used as a sole basis for prostate cancer diagnosis.

11. PERFORMANCE CHARACTERISTICS

11.1. Sensitivity

Sensitivity of the assay has been investigated through the Limit of Blank (LoB), Limit of Detection (LoD), Limit of Quantitation (LoQ) and Analytical Sensitivity (AS). The following table shows the results obtained.

	Result (ng/mL)
LoB	0.045
LoD	0.216
LoQ	1.21
A.S.	0.054

11.2. Analytical Precision

Assay precision was established by analysing patient sera of different PSA concentrations. The results are shown below.

Intra-assay				
Patients	n° replicates	Mean ng/mL	SD ng/mL	CV %
1	24	12.5	0.65	6.0
2	24	3.4	0.13	3.9
3	32	0.8	0.07	8.8

Inter-assay				
Patients	n° replicates	n° lots	Mean ng/mL	CV %
1	75	3	12.3	6.7
2	75	3	3.3	8.0

Inter-lot			
Patients	n° replicates	Mean ng/mL	CV %
1	18	13.0	3.4
2	18	3.5	8.8
3	18	12.4	4.1
4	18	3.4	4.1

11.3. Analytical Specificity

The antibodies used in this kit are highly specific for total PSA (free PSA & PSA-ACT-complex), with a relatively low cross-reactivity to other proteins and polypeptides, lipids or chemotherapeutic agents that might be present in patient samples.

Antigens	Amount added	Cross reaction
Proteins		
AFP	10 µg/mL	No
CEA	10 µg/mL	No
HCG	10 µg/mL	No
Lactalbumin	10 µg/mL	No
PAP	1 µg/mL	No

	Amount added	Cross reaction
Interfering substances		
Bilirubin	0.2 mg/mL	No
Hemoglobin*	0.1 mg/ml	No
Triglyceride	15 mg/mL	No

* at higher concentration hemoglobin results in too high OD values, hemolytic samples should thus be avoided.

	Amount added	Cross reaction
Cytostatic drugs		
Carboplatin	700 µg/mL	No
Cisplatin	200 µg/mL	No
Calcium Folate	2.3 µg/mL	No
Cyclophosphamide	1000 µg/mL	No
5-Fluorouracil	500 µg/mL	No
Doxorubicin HCl	72 µg/mL	No
Dexamethasone	11 µg/mL	No
Diethylstilbestrol	2 µg/mL	No
Flutamide	10 µg/mL	No
Methotrexate	50 µg/mL	No
Hypertension drugs		
Simvastatin	0.1 µg/mL	No
Irbesartan	1.5 µg/mL	No
Sildenafil Citrate	5 µg/mL	No
Furosemide	200 µg/mL	No
Other substances		
Benzalkonium Chloride	0.5%	No
Plasmin-α2-Antiplasmin complex	5 µg/mL	No
Paclitaxel	5.5 µg/mL	No

Note: Diametra Total PSA kit does not discriminate between free PSA and complexed forms of PSA.

Until today no other substances (drugs) are known to us, which have an influence on the measurement of PSA in a sample.

11.4. Serum-plasma study

Plasma samples obtained using heparin, citrate and EDTA have been investigated.

Results from plasma samples have been compared to results from serum samples obtained from the same patient (recovery %).

A plasma sample concentration difference < 20% against serum sample concentration has been considered acceptable.

For each type of plasma, 10 different samples have been analyzed.

Results:

Plasma	Recovery (%)	Interference
Lithium-Heparin	104.3%	No
EDTA	97.8%	No
Citrate	100.1%	No

11.5. Accuracy: recovery

A known amount of PSA was added to three patient sera and the quantities recovered were measured. The results are shown in the following table.

Sample	Expected value (ng/mL)	Observed value (ng/mL)	Recovery %
1	6.4	6.4	99.8
2	4.6	4.7	97.6
3	10.1	10.9	92.6

11.6. Accuracy: dilution

4 samples containing different amounts of analyte were serially diluted with Cal 0, from 1:2 to 1:16. The percentage recovery was calculated by comparing the expected and measured values for the analyte; this percentage must be found within the acceptance range of 85-115%.

Sample	Conc. (ng/mL)	Mean recovery %
1	8.84	106.9
2	9.92	111.2
3	14.52	111.0
4	16.69	103.1

All linearity data are within the acceptance range.

11.7. Hook effect

The assay was tested for a high dose hook effect.

Up to a PSA concentration of 2000 ng/mL no hook effect was observed.

11.8. Correlation

Diametra Total PSA ELISA was compared with the Access Hybritech PSA assay.

57 samples have been assayed.

The results are the following:

$$y = 1.108x - 0.054$$

$$R^2 = 0.952$$

11.9. Clinical sensitivity and specificity

Clinical sensitivity and specificity have been studied on 74 samples; 11 samples have a PSA concentration > 4 ng/mL, 63 samples have a PSA concentration < 4 ng/mL.

The results have been compared between the Diametra ELISA kit and a commercially available ELISA kit (the "gold standard"). this study showed the following results:

The results are showed below:

Diagnostic sensitivity	81.2 %
Diagnostic specificity	> 99.9 %
Positive Predictive Value	> 99.9 %
Negative Predictive Value	96,9 %
Reproducibility	97,2 %

11.10. Normal Values

The normal values of the Diametra Total PSA kit were determined by measuring the values of apparently healthy individuals and men with suspected high PSA values.

The results are showed below:

	Healthy males	Suspicious males
Number	51	50
97.5th Percentile ng/mL	4.85	14.82
2.5th Percentile ng/mL	0.00	2.93
Mean (ng/mL)	1.20	6.56
Median (ng/mL)	0.86	5.51
Min Value (ng/mL)	0.00	1.33
Max Value (ng/mL)	5.35	16.94

The normal values corresponds to data available in the literature.

It is strongly recommended, that each laboratory should determine its own normal and abnormal values.

Please refer to the previous par. 10 for other information.

12. LIMITATIONS

12.1. Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact Diametra.

12.2. Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

12.3. Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

13. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

14. BIBLIOGRAPHY

1. Ornstein DK, Pruthi RS. Prostate-specific antigen. Expert Opin Pharmacother. 2000 Dec;1(7):1399-411.
2. Kostova MB, Brennen WN, Lopez D, Anthony L, Wang H, Platz E, Denmeade SR. PSA-alpha-2-macroglobulin complex is enzymatically active in the serum of patients with advanced prostate cancer and can degrade circulating peptide hormones. Prostate. 2018 Aug;78(11):819-829. doi: 10.1002/pros.23539. Epub 2018 Apr 16. PMID: 29659051.
3. Stamey TA, Yang N, Hay AR, McNeal JE, Freiha FS, Redwine E. Prostate-specific antigen as a serum marker for adenocarcinoma of the prostate. N Engl J Med. 1987 Oct 8;317(15):909-16.
4. Albertsen PC. Prostate cancer screening with prostate-specific antigen: Where are we going? Cancer. 2018 Feb 1;124(3):453-455. doi: 10.1002/cncr.31140. Epub 2017 Dec 12. PMID: 29231972.
5. Thompson IM, Pauler DK, Goodman PJ, Tangen CM, Lucia MS, Parnes HL, Minasian LM, Ford LG, Lippman SM, Crawford ED, Crowley JJ, Coltman CA Jr. Prevalence of prostate cancer among men with a prostate-specific antigen level < or =4.0 ng per milliliter. N Engl J Med. 2004 May 27;350(22):2239-46. doi: 10.1056/NEJMoa031918. Erratum in: N Engl J Med. 2004 Sep 30;351(14):1470. PMID: 15163773.
6. Schröder FH, Hugosson J, Roobol MJ, Tammela TL, Ciatto S, Nelen V, Kwiatkowski M, Lujan M, Lilja H, Zappa M, Denis LJ, Recker F, Páez A, Mänttinen L, Bangma CH, Aus G, Carlsson S, Villers A, Rebillard X, van der Kwast T, Kujala PM, Blijenberg BG, Stenman UH, Huber A, Taari K, Hakama M, Moss SM, de Koning HJ, Auvinen A; ESRPC Investigators. Prostate-cancer mortality at 11 years of follow-up. N Engl J Med. 2012 Mar 15;366(11):981-90. doi: 10.1056/NEJMoa1113135. Erratum in: N Engl J Med. 2012 May 31;366(22):2137. PMID: 22417251; PMCID: PMC6027585.
7. Moyer VA; U.S. Preventive Services Task Force. Screening for prostate cancer: U.S. Preventive Services Task Force recommendation statement. Ann Intern Med. 2012 Jul 17;157(2):120-34. doi: 10.7326/0003-4819-157-2-201207170-00459. PMID: 22801674.
8. Li B. Prostate cancer targeted therapy. In: Schwab M, ed. *Encyclopedia of Cancer*. Berlin, Germany: Springer; 2012:3062-3064.

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14. BIBLIOGRAFIA

- Ornstein DK, Pruthi RS. Prostate-specific antigen. Expert Opin Pharmacother. 2000 Dec;1(7):1399-411.
- Kostova MB, Brennen WN, Lopez D, Anthony L, *et al.* PSA-alpha-2-macroglobulin complex is enzymatically active in the serum of patients with advanced prostate cancer and can degrade circulating peptide hormones. Prostate. 2018 Aug;78(11):819-829. doi: 10.1002/pros.23539. Epub 2018 Apr 16. PMID: 29659051.
- Stamey TA, Yang N, Hay AR, McNeal JE, *et al.* Prostate-specific antigen as a serum marker for adenocarcinoma of the prostate. N Engl J Med. 1987 Oct 8;317(15):909-16.
- Alberts PC. Prostate cancer screening with prostate-specific antigen: Where are we going? Cancer. 2018 Feb 1;124(3):453-455. doi: 10.1002/cncr.31140. Epub 2017 Dec 12. PMID: 29231972.
- Thompson IM, Pauler DK, Goodman PJ, Tangen CM, *et al.* Prevalence of prostate cancer among men with a prostate-specific antigen level < or =4.0 ng per milliliter. N Engl J Med. 2004 May 27;350(22):2239-46. doi: 10.1056/NEJMoa031918. Erratum in: N Engl J Med. 2004 Sep 30;351(14):1470. PMID: 15163773.
- Schröder FH, Hugosson J, Roobol MJ, Tammela TL, *et al.*; ESRPC Investigators. Prostate-cancer mortality at 11 years of follow-up. N Engl J Med. 2012 Mar 15;366(11):981-90. doi: 10.1056/NEJMoa1113135. Erratum in: N Engl J Med. 2012 May 31;366(22):2137. PMID: 22417251; PMCID: PMC6027585.
- Moyer VA; U.S. Preventive Services Task Force. Screening for prostate cancer: U.S. Preventive Services Task Force recommendation statement. Ann Intern Med. 2012 Jul 17;157(2):120-34. doi: 10.7326/0003-4819-157-2-201207170-00459. PMID: 22801674.
- Li B. Prostate cancer targeted therapy. In: Schwab M, ed. *Encyclopedia of Cancer*. Berlin, Germany: Springer; 2012:3062-3064.

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ERROR POSSIBLE CAUSES / SUGGESTIONS**No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

Too low reaction (too low ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

Too high reaction (too high ODs)

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

Unexplainable outliers

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

DiaMetra	Packaging Information Sheet	Mod. PIS000-1
-----------------	------------------------------------	----------------------

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