DCM138-6  
Ed. 04/2020

## Free PSA

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### 14. BIBLIOGRAFIA

1. World Cancer Report. World Health Organization. 2014. Chapter 5.11 ISBN 978-9283204299
2. Henttu P and Vihko P. Prostate-specific antigen and human glandular kallikrein: two kallikreins of the human prostate. *Amm Med* 1994, 26(3):157-64.
3. Balk SP, Ko YL and Bubley GJ. Biology of prostate-specific antigen. *J Clin Oncol*. 2003, 15;21(2):383-91.
4. Zhou AM et al. Multiple forms of prostate-specific antigen in serum: differences in immunorecognition by monoclonal and polyclonal assays. *Clin Chem*. 1993, 39(12):2483-91.
5. Fritsche HA and Babalan RJ. Analytical performance goals for measuring prostate-specific antigen. *Clin Chem*, 1993, 39: 1529-1529.
6. Velonas VM, Woo HH, dos Remedios CG, Assinder SJ. Current status of biomarkers for prostate cancer. *Int J Mol Sci*. 2013, 24;14(6):11034-60.
7. Gion M et al. Percent free prostate-specific antigen in assessing the probability of prostate cancer under optimal analytical conditions. 1998, *Clin Chem* 44(12):2462-70.
8. Akdas et al. The role of free prostate specific antigen in the diagnosis of prostate cancer. *British J Urol*, 1997, 79: 920-923.
9. Chen YT et al. Using proportions of free to total prostate-specific antigen, age, and total prostate-specific antigen to predict the probability of prostate cancer. *Urology*, 1996, 47(4):518-24.
10. Catalona et al. Use of the percentage of free prostate-specific antigen to enhance differentiation of prostate cancer from benign prostatic disease: a prospective multicenter clinical trial. *JAMA*, 1998, 20;279(19):1542-7.
11. Liu J et al. Establishment of two new predictive models for prostate cancer to determine whether to require prostate biopsy when the PSA level is in the diagnostic gray zone (4-10 ng ml<sup>-1</sup>). *Asian J Androl*. 2019 21:1-4.
12. Huang Y, Li ZZ, Huang YL, Song HJ, Wang YJ. Value of free/total prostate-specific antigen (f/t PSA) ratios for prostate cancer detection in patients with total serum prostate-specific antigen between 4 and 10 ng/mL: A meta-analysis. *Medicine (Baltimore)*. 2018; 97(13).
13. Duffy MJ. Biomarkers for prostate cancer: prostate-specific antigen and beyond. *Clin Chem Lab Med*. 2019, 12. aop.
14. Carlsson SV and Roobol MJ. Improving the evaluation and diagnosis of clinically significant prostate cancer in 2017. *Curr Op in Urol*. 2017 27(3): 198–204.
15. Milford Ward A et al. Prostate specific antigen: biology, biochemistry and available commercial assays. *Ann Clin Biochem*, 2001, 38: 633-651.
16. Price CP et al., Pre-and post-analytical factors that may influence use of serum prostate specific antigen and its isoforms in a screening program for prostate cancer. *Ann Clin Biochem*, 2001, 38: 188-216.
17. Ferguson J et al. Continued provision of WHO International Standards for total and free PSA: Content and commutability of replacement preparations. *Clin Biochem*. 2019 71:58-66.

Ed. 04/2020

DCM138-6

#### Legal Manufacturer

Dia.Metra S.r.l.  
Via Pozzuolo 14  
06038 SPELLO (PG) Italy  
Tel. +39-0742-24851  
Fax +39-0742-316197  
E-mail: [info@diametra.com](mailto:info@diametra.com)

## Free PSA

for routine analysis

Enzyme immunoassay for the quantitative determination of free PSA (Prostate Specific Antigen) in human serum and plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 test

REF DKO138

### INTENDED USE

Diametra Free PSA ELISA is used for the quantitative determination of free Prostate Specific Antigen (f-PSA) in human serum or plasma samples. The determination of f-PSA levels is generally used in conjunction with a total PSA (t-PSA) measurement to determine the ratio between f-PSA and t-PSA. This ratio helps to estimate the risk for prostate cancer and to discriminate between elevated t-PSA levels caused by cancerous or non-cancerous conditions. F-PSA determinations are especially recommended for men with elevated t-PSA levels and negative results with digital rectal examination (DRE) in order to decide if a second prostate biopsy is indicated.

### 1. CLINICAL SIGNIFICANCE

Prostate cancer is the most frequent type of cancer found in man and is the second cause of death due to cancer in males. Until recently, digital rectal examination (DRE) was frequently used as only diagnostic modality for the detection of early stages of prostate cancer. In the recent years the determination of serum PSA levels has become the most accepted method to improve the diagnostic specificity of DRE. Although PSA is a tissue specific protein and is not solely tumor specific, it has become the most important marker for prostate carcinoma, showing a better specificity than other biochemical markers used in this context (PAP, total alkaline phosphatase, carcinoembryonic antigen, etc.)

In 1979, Wang et al. isolated a specific antigen for normal prostate tissue and called this protein PSA. PSA is a 33 kDa serine proteinase. Immunohistological studies have shown that PSA is localized in the cytoplasm of prostate acinar cells, ductal epithelium and in the secretion on the ductal lumina, present in normal, benign hyperplastic and malignant prostate tissues as well metastatic prostate cancer and in seminal plasma. If the structural integrity of the prostate is disturbed and/or the gland size is increased, the amount of PSA in the blood plasma may become elevated. In the blood plasma, most of the PSA forms complexes with various proteinase inhibitors. Only a small fraction of PSA circulates as free inactive PSA. Basically three major forms of PSA can be distinguished, only two of which are immunoreactive. The predominant form of PSA is a complex with  $\alpha$ 1-antichymotrypsin (ACT-PSA). Inactive free PSA (f-PSA) represents around 10-40% of the immunologically detectable PSA. The total amount of immunoreactive PSA is known as total PSA (t-PSA). PSA complexed with  $\alpha$ -2-macroglobulin cannot be detected by immunological assays and is therefore frequently called occult PSA (o-PSA).

Current methods of screening men for prostate cancer utilize the detection of t-PSA. Levels of 4.0 ng/ml or higher

are strong indicators of the possibility of prostatic cancer and are an indication for follow-up examinations of the patient. However, elevated serum PSA levels are frequently also attributed to benign prostatic hyperplasia, leading to a high percentage of false positive screening results. A potential solution to this problem involves the determination of free PSA levels. Studies have suggested that the percentage of free PSA is lower in patients with prostate cancer than those with benign prostatic hyperplasia. Thus, the measurement of free serum PSA in conjunction with total PSA, can improve specificity of prostate cancer screening in selected men with elevated total serum PSA levels, which would subsequently reduce unnecessary prostate biopsies with minimal effects on cancer detection rates.

### 2. PRINCIPLE

This f-PSA ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA). The microtiter wells are coated with an anti-f-PSA monoclonal antibody, directed towards an epitope of an antigen molecule. An aliquot of patient serum is incubated in the coated well with enzyme conjugated second antibody (E-Ab), directed towards a different region of the antigen molecule. After incubation the unbound E-Ab is washed off. The amount of bound E-Ab is proportional to the concentration of antigen in the sample. After adding the substrate solution, the intensity of colour developed is proportional to the antigen concentration in the sample. The measured ODs of the Calibrators are used to construct a calibration curve against which the unknown samples are calculated.

### 3. REAGENTS, MATERIALS AND INSTRUMENTATION

#### 3.1. Reagents and materials supplied in the kit

1. PSA Calibrators (6 vials, 0.5 mL each)  
CAL0 / Sample Diluent (10 mL) **REF DCE002/13806-0**  
CAL1 **REF DCE002/13807-0**  
CAL2 **REF DCE002/13808-0**  
CAL3 **REF DCE002/13809-0**  
CAL4 **REF DCE002/13810-0**  
CAL5 **REF DCE002/13811-0**

Contain non-mercury preservative.

2. Controls (2 vial, 0.5 mL each)  
Low Control **REF DCE045/13804-0**  
High Control **REF DCE045/13805-0**

For concentration see the Certificate of Analysis.

Contains non-mercury preservative.

3. Coated Microplate (1 breakable microplate)  
Monoclonal antibody anti PSA adsorbed on microplate  
**REF DCE002/13803-0**

4. Assay Reagent (1 vial, 6 mL)  
Contains non-mercury preservative  
**REF DCE043/13843-0**

5. Enzyme Conjugate (1 vial, 6 mL)  
Monoclonal antibody anti PSA conjugated to horseradish peroxidase. Contains non-mercury preservative.  
**REF DCE002/13802-0**

6. 40X Wash Buffer (1 vial, 30 mL)  
Dilute with distilled water before use.  
**REF DCE057/13857-0**

7. TMB Substrate (1 vial, 14 mL)  
Contains TMB (tetramethylbenzidine) and H<sub>2</sub>O<sub>2</sub>  
**REF DCE004/13804-0**

8. Stop Solution (1 vial, 14 mL)  
Contains sulphuric acid (H<sub>2</sub>SO<sub>4</sub>) 0.5M (avoid contact with the Stop Solution, it may cause skin irritations and burns)  
**REF DCE005/13805-0**

#### 3.2. Material required but not provided

- Precision micropipettes (volume: 25 µL and 100 µL) with disposable tips
- Distilled water
- ELISA photometer with 450 nm- and 630 nm-filters
- Timer with 60 min. range or higher
- Container for the proper handling of waste and samples after use

#### 3.3. Auxiliary materials and instrumentation

- Microplate washer.
- Vortex or similar mixing tools.

#### 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Use appropriate personal protective equipment while working with the reagents provided.  
Serum and plasma samples should be treated as potentially infectious materials. Wear gloves and proper laboratory attire when handling sample materials. Do not eat, drink or smoke in areas where specimen or kit reagents are handled. Do not pipette with the mouth. In case of skin contact, wash with a germicidal soap and

copious amounts of water. Seek medical advice if indicated.

- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the reagents should be handled in the same manner as potentially infectious material.
- Due to the potentially infectious character of samples and kit components all materials that have come in contact with these materials should be sterilized and disposed of according to local regulations. This also includes the liquid waste.
- The assay reagents contain preservatives, TMB, H<sub>2</sub>O<sub>2</sub> or sulphuric acid and may be harmful if ingested. A direct skin or mucosa contact should be avoided. In case of skin contact, wash thoroughly with water and seek medical attention if required.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested.  
To prevent chemical burns, avoid contact with skin and eyes.  
Furthermore, the instrumentation employed to dispense it should be thoroughly cleaned after use.

#### 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.  
Be sure to use the current valid revision of the Instruction for Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (20-26°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.

- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Avoid reagent or sample carry-over by using fresh tips for solutions and samples.
- Do not use test kit if zip lock pouch or bottles have been damaged.

## 6. STORAGE AND STABILITY

Store the kit at 2-8°C; the kit is stable until the expiry date claimed on the kit label and in the Certificate of Analysis. Do not use the kit or its components after the expiry date. After opening, the kit is stable 2 months at 2-8°C.

## 7. ASSAY PROCEDURE

### 7.1. Preparation of the Calibrators and Controls

The Calibrators are calibrated against WHO NIBSC 17/102 and have the following concentrations:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
ng/mL	0	0.75	1.50	3	6	12

Calibrators and controls are ready to use.

### 7.2. Preparation of the Conjugate

The component "Enzyme Conjugate" is ready to use.

### 7.3. Preparation of the "Assay Reagent"

The component "Assay Reagent" is ready to use.

### 7.4. Preparation of Samples

The kit can be used on serum and plasma (obtained with Lithium-Eparine, Citrate or EDTA) samples. Do not use samples containing Sodium Azide.

Blood samples are collected by venipuncture. As different factors could influence the PSA level in blood, doctors should ensure that the patient has avoided the following conditions before taking the blood sample:

- conditions that may lead to an increase of PSA levels:
  - biking
  - sexual intercourse (ejaculation)
  - manipulation of the prostate during medical examinations like DRE, transrectal prostatic ultrasound etc.
  - prostatitis
  - liver dysfunction
- conditions that may lead to a decrease of PSA levels:
  - intake of 5-alpha-reductaseinhibitors, antiandrogens, or GnRH analoga

The preparation of serum or plasma samples is performed according to standard techniques. Serum or plasma should be prepared as soon as possible to avoid hemolysis and to improve the stability of PSA.

Important note: samples with high concentration of rheumatoid factor and human anti-mouse antibodies (HAMA) could give erroneous results. Samples with an expected PSA value higher than 12 ng/mL should be diluted with the "Cal 0 / Sample diluent" or reported as "> 12 ng/mL".

## Samples storage

If the assay is not performed immediately, the samples can be stored:

- 3 days at 2-8°C
- 7 days at -20°C

It's recommended to perform only 1 freezing and thawing cycle.

## 7.5. Preparation of the Wash Solution

Before starting the assay the wash buffer must be diluted to the right concentration. Per well approximately 2 mL of diluted buffer are needed. Calculate the volume of buffer required for the testing. Take 1/40 of the volume wash buffer concentrate and dilute with 39/40 of the volume distilled water.

The diluted wash solution is stable up to 2 weeks at room temperature (20-26°C).

## 7.6. Test procedure

- **Allow all reagents to reach room temperature (20-26°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve, two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample	Controls
Calibrators C <sub>0</sub> -C <sub>5</sub>	25 µL		
Controls			25 µL
Sample		25 µL	
Assay Reagent	50 µL	50 µL	50 µL
Mix by moving plate on the table (10 sec); it's very important to mix completely the solution. Incubate 1 hour at room temperature (20-26°C).			
Enzyme Conjugate	50 µL	50 µL	50 µL
Mix by moving plate on the table (10 sec); it's very important to mix completely the solution. Incubate 1 hour at room temperature (20-26°C).			
Remove solution from the wells by aspirating the liquid or by decanting it. If decanting, tap plate on adsorbent paper to remove residual liquid.			
Washing step: fill the wells with 300 µL of diluted wash solution, wait 10sec and then remove the solution; repeat this washing procedure for a total of 5 times. In the case you use an automatic plate washer, use 400 µL of diluted wash solution.			
<u>Important note:</u> sensitivity and precision of the assay strictly depend on the correct execution of the washing step.			

Substrate Solution	100 µL	100 µL	100 µL
Incubate 15 minutes at room temperature (20-26°C) in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance at 450 nm against a reference wavelength of 620-630 nm. It's recommended to read the microplate within 10 minutes after adding the Stop Solution.			

## 8. QUALITY CONTROL

Good laboratory practice requires that controls have to be run with each calibration curve.

A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to the local active regulations to assure the day to day validity of results.

Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the Certificate of Analysis provided with the kit; these data always refer to the kit lot in use and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analyzing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

## 9. CALCULATION OF RESULTS


1. Calculate the mean absorbance for each duplicate.
2. Subtract the absorbance value of the Calibrator 0 from the mean absorbance values of Calibrators, control and samples.
3. Draw the calibration curve on lin-lin or log-log graph paper by plotting absorbance values of Calibrators against appropriate PSA concentrations.
4. Read off the f-PSA concentrations for the control and the samples.

Samples with an expected PSA value higher than 12 ng/mL should be diluted with the "Cal 0 / Sample diluent" or reported as "> 12 ng/mL".

## 9.1. Example of Calibration curve

The following data are for example use only and must not be used in any case.

Wells	Identity	A 450 nm		Conc. ng/mL
1-2	Cal 0 ng/mL	0.011	0,009	
3-4	Cal 0.75 ng/mL	0.109	0.096	
5-6	Cal 1.5 ng/mL	0.185	0.179	
7-8	Cal 3 ng/mL	0.393	0.392	
9-10	Cal 6 ng/mL	0.871	0.837	
11-12	Cal 12 ng/mL	1.901	1.844	
13-14	Control LOw	0,218	0.231	2,05
15-16	Control High	1,009	0,989	7,80

Note that the absolute OD values for the Calibrators might vary due to temperature influences or age of the conjugate. As long as the OD values form a calibration curve and remain within the specifications and the control shows the expected value, results for unknown PSA samples are valid.

## 10. EXPECTED VALUES LIMITATIONS OF THE PROCEDURE

In a study performed with apparently healthy men, using the Dia.Metra Srl Free PSA kit the following data were observed:

Population	Uomini, sani
n° samples	120
Mean (ng/mL)	0,2
Median (ng/mL)	0,1
Percentile (ng/mL)	0,0 - 1,2
Range (ng/mL)	0,0 - 4,6

When total PSA values are between 4-10 ng/mL the ratio of free PSA / total PSA might be employed to increase the diagnostic specificity of PSA testing.

A ratio  $\leq 10\%$  indicates 49% to 65% risk of prostate cancer depending on age; a free/total PSA ratio  $> 25\%$  indicates a 9% to 16% risk of prostate cancer, depending on age. According to the American Cancer Society and National Cancer Institute men with f-PSA below 7 % should undergo biopsy (see reference 10).

**Please note that an isolated f-PSA concentration is of no diagnostic value.**

The above values are only for user's guidance. It is recommended for each laboratory to establish its own specific values that take into consideration a population indigenous to the area where the laboratory is located.

**Note: PSA values and the free PSA / total PSA ratio can only be used to estimate the cancer risk. They should always be interpreted in conjunction with other clinical findings and should not be used as a sole basis for prostate cancer diagnosis.**

## 11. PERFORMANCE CHARACTERISTICS

### 11.1. Assay range

The range of the assay is 0,095 - 12,0 ng/mL.

### 11.2. Analytical sensitivity

Analytical sensitivity has been investigated through the LoB (Limit of Blank), LoD (Limit of Detection), LoQ (limit of Quantitation) and Analytical Sensitivity (A.S.). The table below shows the results obtained.

	Results (ng/mL)
LoB	0,001
LoD	0,014
LoQ	0,095
A.S.	0,01

### 11.3. Analytical Precision

- Intra-assay (within assay)**

Intra-assay precision was determined by measuring each sample 10 times.

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- Inter-assay (between assay)**

Inter-assay precision was determined by measuring each sample 10 times for 3 days.

Patients	Number of replicates	Mean ng/mL	CV %
1	30	0,7	10,0
2	30	4,6	8,0
3	30	6,5	6,8
4	30	8,3	6,3

- Inter-lot (between lot)**

Inter-lot precision was determined by measuring each sample 6 times with 3 different lots.

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### 11.4. Accuracy

- Recovery**

4 samples were prepared by adding aliquots of samples with known free PSA concentration. The recovery (%) was calculated as percentage ratio between measured values and expected values.

	Sample			
	1	2	3	4
Conc (ng/mL)	1,4	1,9	3,8	4,9
Average (%)	88,8	92,5	88,5	92,9
Recovery	from	87,8	89,5	86,9
	to	89,9	95,4	90,1

- Linearity**

4 samples were measured undiluted and diluted with Calibrator 0. Recovery (%) was calculated as percentage ratio between measured values and expected values.

	Sample			
	1	2	3	4
Conc (ng/mL)	1,43	2,65	4,30	10,92
Average (%)	100,5	94,1	104,1	108,2
Recovery	from	88,6	87,5	98,6
	to	107,9	98,1	107,4

### 11.5. Analytical Specificity

- Endogenous interfering substances**

Interferences by Triglycerides, Bilirubin and Hemoglobin have been investigated by adding the interfering substance (at different concentration) to a serum sample and by comparing its concentration to the unspiked sample. 4 samples were investigated.

The interference has been evaluated as "significant" if it causes a concentration bias >10% between spiked and unspiked sample.

The table below shows the results obtained:

Substance	Concentration	Interference
Trigliceridi	7,5 mg/mL	No
Bilirubina	0,5 mg/mL	No
Emoglobina	4 mg/mL	No

Conclusion: the study showed no interference in the assay by the aforementioned substances. It's anyway recommended to avoid using highly lipemic or hemolyzed samples.

- **Cross-reactivity**

The following substances have been tested for cross-reactivity of the assay:

Substance	Concentration (ng/mL)	Cross-reactivity (%)
AFP	193,6	0,0
CEA	100,0	0,0
b-HCG	21,8	0,0
Kallikrein	1200,0	0,0
Rheumatic Factor	1200,0	0,0
HAMA	1200,0	0,0

Conclusion: according to the stated specifications, no one of the tested substances have interference with the assay at the concentrations stated.

- **Drug interferences**

Interferences by the substances listed below have been investigated by adding the interfering substance (at different concentration) to a serum sample and by comparing its concentration to the unspiked sample. 4 samples were investigated.

The interference has been evaluated as "significant" if it causes a concentration bias >15% between spiked and unspiked sample.

The table below shows the results obtained:

Cytostatic drugs	Conc. tested (µg/mL)
Carboplatin	700
Cisplatin	200
Calcium Folate	2,3
Cyclophosphamide	550
Fluorouracil	520
Dexamethasone	11
Paclitaxel	5,3
Doxorubicin • HCl	72

Hypertension drugs	Conc. tested (µg/mL)
Simvastatin	0,1
Irbesartan	1,5
Sildenafil Citrate	5,0
Furosemide	0,2

Antimicrobial agent	Conc. tested (µg/mL)
Benzalkonium Chloride	0,5

Conclusion: according to the stated specifications, no one of the tested drugs have interference with the assay at the concentrations stated.

### 11.6. Hook effect

No hook effect has been observed up to 240 ng/mL of f-PSA added to the sample.

### 11.7. Study on plasma samples

Plasma samples obtained using Lithium-heparin, Citrate and EDTA were investigated.

The result of the resulting plasma assay obtained has been compared with the result of the assay obtained on the serum taken from the same patient from which the plasma has been obtained (% recovery).

A difference in serum - plasma concentration <20% has been considered acceptable.

For each type of plasma, 10 different samples have been analyzed.

Results:

Plasma	Interference
Lithium-eparin	No
Citrate	No
EDTA	No

Diametra Free PSA assay can be performed on serum and plasma (Lithium-heparin, Citrate and EDTA) samples.

### 11.8. Specificity and clinical sensitivity

The diagnostic performances for the Free PSA kit can not be calculated in relation to the "free PSA" molecule because there is no threshold value as in the case of total PSA.

For the estimate of the risk of prostate cancer, the ratio between free PSA and total PSA refers, therefore no analysis can be made in relation to positive and negative samples knowing only the value of free PSA without knowing the total PSA value.

### 11.9. Correlation

The present updated free PSA kit was compared to the previous version of the product.

The linear regression curve is the following:

$$Y = 0,928 \cdot X + 0,021$$

$$R^2 = 0,939$$

## 12. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

## 13. LEGAL ASPECTS

### 13.1. Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact Diametra.

### 13.2. Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under previous point 13.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived. The test result itself should never be the sole determinant for deriving any therapeutic consequences.

### 13.3. Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to previous point 13.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

### 14. BIBLIOGRAPHY

1. World Cancer Report. World Health Organization. 2014. Chapter 5.11 ISBN 978-9283204299
2. Henttu P and Vihko P. Prostate-specific antigen and human glandular kallikrein: two kallikreins of the human prostate. *Amm Med* 1994, 26(3):157-64.
3. Balk SP, Ko YL and Bublely GJ. Biology of prostate-specific antigen. *J Clin Oncol.* 2003, 15;21(2):383-91.
4. Zhou AM et al. Multiple forms of prostate-specific antigen in serum: differences in immunorecognition by monoclonal and polyclonal assays. *Clin Chem.* 1993, 39(12):2483-91.
5. Fritsche HA and Babalan RJ. Analytical performance goals for measuring prostate-specific antigen. *Clin Chem*, 1993, 39: 1529-1529.
6. Velonas VM, Woo HH, dos Remedios CG, Assinder SJ. Current status of biomarkers for prostate cancer. *Int J Mol Sci.* 2013, 24;14(6):11034-60.
7. Gion M et al. Percent free prostate-specific antigen in assessing the probability of prostate cancer under optimal analytical conditions. 1998, *Clin Chem* 44(12):2462-70.
8. Akdas et al. The role of free prostate specific antigen in the diagnosis of prostate cancer. *British J Urol*,1997, 79: 920-923.
9. Chen YT et al. Using proportions of free to total prostate-specific antigen, age, and total prostate-specific antigen to predict the probability of prostate cancer. *Urology*, 1996, 47(4):518-24.
10. Catalona et al. Use of the percentage of free prostate-specific antigen to enhance differentiation of prostate cancer from benign prostatic disease: a prospective multicenter clinical trial. *JAMA*, 1998, 20;279(19):1542-7.
11. Liu J et al. Establishment of two new predictive models for prostate cancer to determine whether to require prostate biopsy when the PSA level is in the diagnostic gray zone (4-10 ng ml<sup>-1</sup>). *Asian J Androl.* 2019 21:1-4.
12. Huang Y, Li ZZ, Huang YL, Song HJ, Wang YJ. Value of free/total prostate-specific antigen (f/t PSA) ratios for prostate cancer detection in patients with total serum prostate-specific antigen between 4 and 10 ng/mL: A meta-analysis. *Medicine (Baltimore).* 2018; 97(13).
13. Duffy MJ. Biomarkers for prostate cancer: prostate-specific antigen and beyond. *Clin Chem Lab Med.* 2019, 12. aop.
14. Carlsson SV and Roobol MJ. Improving the evaluation and diagnosis of clinically significant prostate cancer in 2017. *Curr Op in Urol.* 2017 27(3): 198–204.
15. Milford Ward A et al. Prostate specific antigen: biology, biochemistry and available commercial assays. *Ann Clin Biochem*, 2001, 38: 633-651.
16. Price CP et al., Pre-and post-analytical factors that may influence use of serum prostate specific antigen and its isoforms in a screening program for prostate cancer. *Ann Clin Biochem*, 2001, 38: 188-216.
17. Ferguson J et al. Continued provision of WHO International Standards for total and free PSA: Content and commutability of replacement preparations. *Clin Biochem.* 2019 71:58-66.

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**ERROR POSSIBLE CAUSES / SUGGESTIONS****No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

**Too low reaction (too low ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

**Too high reaction (too high ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

**Unexplainable outliers**

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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