Urine CrossLaps[®] (CTX-I) EIA

For the quantification of degradation products of C-terminal telopeptides of Type I collagen in human urine

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
Immunodiagnostic Systems Limited is not responsible for any other use of the kit or consequence hereof than the one specified above. Neither for misuse, e.g. use deviating from the procedure described in this manual. Furthermore, Immunodiagnostic Systems Limited is not to be made responsible for any diagnoses or conclusions made by the user or third party based on the results obtained with the Urine CrossLaps® (CTX-I) EIA kit nor for any consequences such interpretations may cause.

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INTRODUCTION

Intended use

The Urine CrossLaps® (CTX-I)EIA is an enzyme immunosorbent assay for quantitative determination of degradation products of C-terminal telopeptides of Type-I collagen in human urine.

The Urine CrossLaps® (CTX-I) EIA assay is intended for *in vitro* diagnostic use as an indication of human bone resorption and may be used as an aid in

A. Monitoring bone resorption changes of

- 1) Anti-resorptive therapies in postmenopausal women:
 - a) Hormone Replacement Therapies (HRT) with hormones and hormone like drugs.
 - b) Bisphosphonate therapies.
- 2) Anti-resorptive therapies in individuals diagnosed with osteopenia:
 - a) Hormone Replacement Therapies (HRT) with hormones and hormone like drugs.
 - b) Bisphosphonate therapies.

B. Predicting skeletal response (Bone Mineral Density) in postmenopausal women undergoing anti-resorptive therapies

- a) Hormone Replacement Therapies (HRT) with hormones and hormone like drugs.
- b) Bisphosphonate therapies.

Limitations

The use of the test has not been established to predict the development of osteoporosis or future fracture risk. The use of the test has not been established in hyperparathyroidism or hyperthyroidism. When using the test to monitor therapy, results may be confounded in patients afflicted with clinical conditions known to affect bone resorption e.g. bone metastases, hyperparathyroidism or hyperthyroidism.

Urine CrossLaps® (CTX-I)EIA results should be interpreted in conjunction with clinical findings and other diagnostic results and should not be used as a sole determinant in initiating or changing therapy. Do not interchange Serum CrossLaps® (CTX-I)EIA values with Urine CrossLaps® (CTX-I)EIA values.

Summary and explanation of the test

Type I collagen accounts for more than 90% of the organic matrix of bone and is synthesized primarily in bone (1). During renewal of the skeleton, Type-I collagen is degraded, and small peptide fragments are excreted into the urine. These fragments can be measured by Urine CrossLaps® EIA as has been reported as useful for follow up of anti resorptive treatment of patients with metabolic bone diseases (2-17).

Principle of the procedure

The Urine CrossLaps® (CTX-I)EIA assay is based on the competitive binding of the anti-CrossLaps antibodies to either soluble CrossLaps antigen or to CrossLaps antigen-coated micro-titre wells.

Standards, controls, or unknown samples are pipetted into the appropriate microtiter wells. Then the anti-CrossLaps antibodies are added and incubation takes place for 1 hour at 18-22°C. The wells are washed and peroxidase conjugated anti-rabbit immunoglobulin is added. The wells are then incubated again for 1 hour at 18-22°C. After a second washing step the wells are incubated for 15 min with a chromogenic substrate. The reaction is stopped and the absorbance is measured.

PRECAUTIONS

The following precautions should be observed in the laboratory:

- Do not eat, drink, smoke or apply cosmetics where immunodiagnostic materials are being handled
- Do not pipette by mouth
- Wear gloves when handling immunodiagnostic materials and wash hands thoroughly afterwards
- Cover working area with disposable absorbent paper

Warnings

For *in vitro* use only.

- All reagents and laboratory equipment should be handled and disposed of as if they were infectious
- Do not use kit components beyond the expiry date and do not mix reagents from different lots

Storage

Store the Urine CrossLaps® (CTX-I)EIA kit upon receipt at 2-8°C. Under these conditions the kit is stable up to the expiry date stated on the box.

MATERIALS

Specimen collection

For optimal results it is recommended to use urine from second morning void. Also for monitoring the individual patient, follow up samples should be collected under same conditions as the baseline sample. Levels of the peptide measured by the Urine CrossLaps® (CTX-I)EIA in urine do not vary with dietary intake, however for consistency, it is recommended to use second morning void. Also for monitoring the individual patient, follow-up samples should be collected under same conditions as the baseline sample (i.e. second morning void).

Please note that for optimal results it is recommendable to centrifuge the samples (e.g. 2000g; 10 min). Keep the urine sample refrigerated (2-8°C) for storage less than one week, or freeze the sample (<-18°C) for longer storage.

Specimen obviously contaminated with whole blood may interfere with assay performance. These specimens should be discarded and a specimen recollected.

Materials supplied

Before opening the kit, please read the section on **Precautions**. The kit contains reagents sufficient for 96 determinations.

Immunostrips MICROPLAT

Microwell strips (12 x 8 wells each) pre-coated with CrossLaps antigen. Supplied in a plastic frame.

CrossLaps Standard CAL 0

One vial (min. 11.0 mL) of ready-for-use buffer solution with protein stabilizer, detergent and preservative.

CrossLaps Standards CAL 1-5

Five vials (min. 0.4 mL/vial) of ready-for-use CrossLaps standard in buffer solution with protein stabilizer, detergent and preservative. The exact value of each Standard is printed on the QC Report.

Urine Controls CTRL 1-2

Two vials (min. 0.4mL/vial) of ready-for-use human urine in buffer solution with protein stabilizer, detergent and preservative. The exact value of each Control is printed on the QC Report.

Primary Antibody Solution AB

One vial (min. 12.0 mL) of ready-for-use anti- CrossLaps antibodies raised in rabbits. Prepared in buffer solution with protein stabilizer, detergent and preservative.

Peroxidase Conjugated Antibody ENZYMCNJ

One vial (min. 12.0 mL) of ready-for-use horseradish peroxidase-conjugated anti-rabbit antibodies in buffer solution with protein stabilizer, detergent and preservative.

Substrate Solution SUBS TMB

One vial (min. 12.0 mL) of ready-for-use tetramethylbenzidine (TMB) substrate in an acidic buffer. Please note that the chromogenic substrate might appear slightly bluish.

Stopping Solution H₂SO₄

One vial (min. 12.0 mL) of ready-for-use 0.18 mol/L sulfuric acid.

Washing Solution WASHBUF 50x

One vial (min. 20.0 mL) of concentrated washing buffer solution with detergent and preservative. Dilute 1+50 before use.

Sealing tape

Adhesive film for covering wells during incubation.

Materials required — not supplied

- Container for preparing the Washing Solution
- Precision micropipettes to deliver 15 µL
- Distilled water
- Precision 8- multipipette to deliver 100 µL
- Microwell mixing apparatus (300 rpm)
- Microtiter plate reader

ASSAY PROCEDURE

Assay Procedure

For optimal performance of the assay it is important to comply with the instructions given below. Prior to use, prepare and equilibrate all solutions to room temperature. **Perform the assay at 18-22°C.**

Determine the number of strips needed for the assay. It is recommended to test all samples in duplicate. In addition, for each run a total of 16 wells are needed for standards and control. Place the appropriate number of strips in the plastic frame. Store unused immunostrips in the tightly closed foil bag with desiccant capsules.

1 Pipette Standards or samples

Pipette 15 µL of **Standards** [CAL 0-5], **Controls** [CTRL 1-2] or unknown samples into appropriate wells.

2 Incubation in Immunostrips

Add 100 µL **Primary Antibody Solution** [AB] to each well. Cover the immunostrips with sealing tape and incubate for 60±5 minutes at 18-22°C on a microtiter plate mixing apparatus (300 rpm).

3 Washing

Wash the immunostrips 5 times manually with 300 µL diluted **Washing Solution** ([WASHBUF 50x diluted 1+50 in distilled water]).

Using an automatic plate washer, follow the instructions of the manufacturer or the guidelines of the laboratory. Usually, 5 washing cycles are adequate. Make sure that the wells are completely emptied after each manual or automatic washing cycle.

4 Incubation with Peroxidase Conjugated Antibody

Add 100 µL of **Peroxidase Conjugated Antibody** [ENZYMCONJ] to each well, cover with sealing tape, and incubate for 60±5 minutes at 18-22°C on the mixing apparatus (300 rpm).

5 Washing

Proceed as described in step 3.

6 Incubation with chromogenic substrate solution

Pipette 100 µL of the **Substrate Solution** [SUBS TMB] into each well, cover with sealing tape, and incubate for 15±2 minutes at 18-22°C in the dark on the mixing apparatus (300 rpm).

Do not pipette directly from the vial containing TMB substrate but transfer the needed volume to a clean reservoir. Remaining substrate in the reservoir should be discarded and not returned to vial.

7 Stopping of color reaction

Pipette 100 µL of **Stopping Solution** [H₂SO₄] into each well.

8 Measurement of absorbance

Measure the absorbance at 450 nm with 650 nm as reference within two hours.

Limitations of the procedure

If the absorbance of a sample exceeds that of **Standards 5**, the sample should be diluted in **Standard 0** and re-analysed.

QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analyzed with appropriate statistical methods.

RESULTS

Calculation of results

- Calculate the mean of the duplicate absorbance determinations. Construct a standard curve on log-linear graph paper by plotting the mean absorbencies of the six Standards (ordinate) against the corresponding CrossLaps concentrations (abscissa). Draw the best fitting curve. Alternatively, a four-parametric logistic curve fit can be used
- Determine the CrossLaps concentration of the Controls and each patient samples by interpolation

Example of results obtained:

Standards/ Controls/ Samples	CrossLaps Conc. (µg/L)	A ₄₅₀₋₆₅₀ (nm)	Mean A ₄₅₀₋₆₅₀ (nm)	Interpolated CrossLaps Conc. (µg/L)
Standard 0	0	1.950 / 2.027	1.989	
Standard 1	75	1.599 / 1.705	1.652	
Standard 2	250	1.236 / 1.275	1.256	
Standard 3	750	0.785 / 0.850	0.818	
Standard 4	2250	0.472 / 0.496	0.484	
Standard 5	6750	0.267 / 0.285	0.276	
Control 1		0.643 / 0.671	0.657	1193
Control 2		0.440 / 0.468	0.454	2462
Sample I		1.161 / 1.202	1.182	305
Sample II		0.729 / 0.765	0.747	919
Sample III		0.423 / 0.422	0.423	2829

Please note: The data above are for illustration only and should not be used to calculate the results of any run.

Calculation of corrected CrossLaps value

For each sample the CrossLaps concentration (µg/L) and the creatinine concentration (mM=mmol/L) should be determined. For determination of creatinine the method by Jaffe (3) or equivalent is recommended.

The following equation corrects the CrossLaps concentration for variation in urine concentration:

$$\text{Corr. CrossLaps Value (µg/mmol)} = \frac{\text{CrossLaps (µg/L)}}{\text{Creatinine (mM)}}$$

Performance characteristics

All performance data have been established using second morning void urine samples unless otherwise indicated.

Detection limit: 50 µg/L

This is the concentration corresponding to two standard deviations below the mean of 21 determinations of "CrossLaps Standard 0".

Precision

The following results were obtained when assaying five replicates of four samples each day in five consecutive days (25 determinations in total for each sample).

InterAssay Variation (n=25)

Mean (µg/L)	SD (µg/L)	CV (%)
240	13.68	5.7
480	13.92	2.9
1340	46.90	3.5
3790	20.47	5.4

IntraAssay Variation (n=25)

Mean (µg/L)	SD (µg/L)	CV (%)
240	22.56	9.4
480	25.92	4.7
1340	62.98	4.7
3790	32.59	8.6

Dilution/Linearity

The Urine CrossLaps® (CTX-I)EIA is linear in the range 50 µg/L to 6750 µg/L of CrossLaps.

Urine samples with the concentration of 3450-4340 µg/L CrossLaps were diluted with standard 0 and the concentration of CrossLaps were determined with Urine CrossLaps® (CTX-I)EIA. The urine neat sample is set to 100%.

Sample	Dilution	Expected µg/L	Measured µg/L	Corrected µg/L	Recovery %
1	neat	3450	3450	3450	100
	1:2	1730	1680	3360	97
	1:4	860	790	3160	92
	1:8	430	440	3520	102
2	neat	4060	4060	4060	100
	1:2	2030	2200	4400	108
	1:4	1015	1070	4280	105
	1:8	508	550	4400	108
3	neat	4340	4340	4340	100
	1:2	2170	2180	4380	101
	1:4	1085	1150	4600	106
	1:8	443	510	4080	94
Mean					101

Recovery

Spiked samples were prepared by adding varying amounts of the CrossLaps antigen to three different urine specimens.

Sample	Initial µg/L	Added µg/L	Expected µg/L	Measured µg/L	Recovery %
1	410	130	540	560	106
	410	380	790	870	110
	410	1130	1540	1670	108
2	500	130	630	560	89
	500	380	880	790	90
	500	1130	1630	1590	98
3	520	130	650	590	91
	520	380	900	790	88
	520	1130	1650	1390	84
Mean					96

Expected values

It is advisable for a laboratory to establish its own range of normal and pathological values. As an example, the mean values and standard deviations for various populations are given below; all values have been established using second morning void samples.

For further reading, please refer to the reference list.

Populations	Number of subjects	Mean Values Urine CrossLaps ($\mu\text{g}/\text{mmol Cr}$)	95% Confidence Interval ($\mu\text{g}/\text{mmol Cr}$)
Pre-menopausal women	175	191	67 – 544
Post-menopausal women	250	324	121 – 874
Males	249	174	54 - 559

Day to Day Intra-individual Variation:

The Day to Day Intra-individual Variation was assessed by analyzing urine samples (second morning void) from 14 healthy post menopausal women at five time points over 2 weeks.

Urine CrossLaps [®] EIA [$\mu\text{g}/\text{mmol creatinine}$]								
Woman no	Age	Visit 1	Visit 2	Visit 3	Visit 4	Visit 5	Mean	SD
001	75	368	331	299	290	307	319	31
002	75	239	261	240	267	247	251	13
003	72	200	193	182	183	185	188	8
004	74	220	218	203	236	203	216	14
005	65	76	78	83	81	72	78	4
006	63	71	52	65	57	59	61	7
007	72	101	135	153	125	106	124	21
008	65	110	104	129	148	141	126	19
009	74	128	134	134	98	131	125	15
010	73	107	131	140	141	117	127	15
011	60	109	124	107	109	114	113	7
012	54	98	112	107	104	113	107	6
013	53	362	433	424	392	479	418	44
014	73	117	142	136	111	99	121	18
Mean	68							16
SD	7.80							10.87

CLINICAL DATA

The Urine CrossLaps[®] (CTX-I) EIA has been used to monitor treatment in several clinical studies and the CrossLaps values have been compared to Bone Mineral Density ($\text{BMD}_{\text{spine}}$) measurements.

All the clinical studies presented below were performed according to the European Standard for good clinical practice (GCP and GLP). Most of the clinical studies presented here were conducted on white Danish women. However, several studies have been published showing that other demographic groups display similar CrossLaps decrease in response to anti-resorptive therapies (6-16). The Paget's disease study was performed on a French population.

For all the data presented below second morning void samples have been used.

The Bone mineral density was measured at the Lumbar spine (L1 - L4). The change in the bone mineral density is presented below as $\alpha\text{-BMD}$. $\alpha\text{-BMD}$ is defined as the slope of the linear regression line for $\text{BMD}_{\text{spine}}$ versus time (years) for the period of treatment. In most cases the calculation of $\alpha\text{-BMD}$ involves a minimum of 5 $\text{BMD}_{\text{spine}}$ measurements. The $\alpha\text{-BMD}$ thus represents the % change in $\text{BMD}_{\text{spine}}$ per year. Because there to this date is no universal agreement as to what constitutes positive BMD response we have calculated the sensitivity and specificity using two different cut-off values for $\alpha\text{-BMD}$; $\alpha\text{-BMD}>0$ and $\alpha\text{-BMD}>1$. The sensitivity is defined as the percent of the study population with a positive BMD response and who have a % change from baseline of Urine CrossLaps[®] EIA which is 40% or greater. The specificity is defined as the percent of the study population without a positive BMD response and who have a % change from baseline of Urine CrossLaps[®] EIA that is less than 40%.

Bisphosphonate studies

Below is shown the Urine CrossLaps® (CTX-I)EIA data [$\mu\text{g}/\text{mmol}$ creatinine] from two different bisphosphonate studies. For the mean values the observed ranges (95% confidence interval) are noted in parentheses.

Alendronate

- Women between age 40 and 59 years, 6 months to 3 years since menopause
- 12 participants on placebo (500 mg calcium)
- 40 participants on active treatment: (5 mg (n=14), 10 mg (n=13), 20 mg (n=13)) alendronate and 500 mg calcium
- Treatment period: 2 - 3 years

Urine CrossLaps® (CTX-I)EIA


	Placebo group		Alendronate group	
	Mean ($\mu\text{g}/\text{mmolL}$)	95% Confidence Interval ($\mu\text{g}/\text{mmolL}$)	Mean ($\mu\text{g}/\text{mmolL}$)	95% Confidence Interval ($\mu\text{g}/\text{mmolL}$)
Baseline	387	300 – 499	412	367 – 462
after 6 months treatment	250	176 – 357	55	43 – 71

Ibandronate

- Women less than 75 years, more than ten years after menopause and have a BMD forearm 1.5 SD or more below the standard for healthy pre-menopausal women
- 25 participants on placebo (1000 mg calcium)
- 41 participants on active treatment: (2.5 mg (n=23), 5 mg (n=18)) ibandronate and 1000 mg calcium
- Treatment period: 1 year

Urine CrossLaps® (CTX-I)EIA


	Placebo group		Ibandronate group	
	Mean ($\mu\text{g}/\text{mmolL}$)	95% Confidence Interval ($\mu\text{g}/\text{mmolL}$)	Mean ($\mu\text{g}/\text{mmolL}$)	95% Confidence Interval ($\mu\text{g}/\text{mmolL}$)
Baseline	377	319 – 445	337	278 – 290
after 6 months treatment	231	184 – 290	41	30 – 57

Using a cut-off for Urine CrossLaps® (CTX-I)EIA of 40% change from baseline the following sensitivities, and specificities and 95% confidence intervals are obtained. Below is indicated the actual fractions:

	Ibandronate		Alendronate	
	Sensitivity	Specificity	Sensitivity	Specificity
α -BMD>0	95% (88 - 102)	57% (31 - 83)	79% (68 - 90)	54% (27 - 81)
α -BMD>1	97% (91 - 103)	45% (23 - 67)	78% (66 -90)	38% (17 - 59)

Below are shown distribution plots for the combined bisphosphonate studies. The number over each bar indicates the number of participants in each class.

HRT studies

Below is shown Urine CrossLaps® (CTX-I)EIA [mg/mmol creatinine] data from three different HRT studies.

Tibolone

- Women less than 75 years and more than ten years after menopause
- 13 participants on placebo (400 mg calcium/day)
- 57 participants on active treatment: (1.25 mg (n=29), 2.5 mg (n=28)) Org OD14 and 400 mg calcium/day
- Treatment period 2 years

Urine CrossLaps® (CTX-I)EIA

	Placebo group		Tibolone group	
	Mean (µg/mmol)	95% Confidence Interval (µg/mmol)	Mean (µg/mmol)	95% Confidence Interval (µg/mmol)
Baseline	259	212 – 315	317	287 – 350
after 6 months treatment	214	174 – 262	99	82 – 119

HRT I

- Women more than 45 years, 1 to 6 years since menopause
 - 42 participants on placebo (400 mg calcium/day)
 - 125 participants on active treatment:

Days 1-16	Days 17-28
34; E 1mg	E 1mg + G 25 µg
29; E 2mg	E 2mg + G 25 µg
28; E 2mg	E 2mg + G 50 µg
34; E 1mg + G 25 µg	E 1mg + G 25 µg
- E = estradiol-15β, G = gestodene, active treatment also receive 400mg calcium/day
- Treatment period: 2 years

Urine CrossLaps® (CTX-I)EIA

	Placebo group		HRT I group	
	Mean (µg/mmol)	95% Confidence Interval (µg/mmol)	Mean (µg/mmol)	95% Confidence Interval (µg/mmol)
Baseline	352	308 – 402	332	312 – 353
after 6 months treatment	361	318 – 411	122	109 – 135


HRT II

- Women between 65 and 70 years and BMC forearm below 1SD of healthy pre-menopausal women
- 19 participants on placebo (1000 mg calcium/day)
- 17 participants on active treatment: 50 µg estradiol, 1 mg norethisterone and 1000 mg calcium/day
- Treatment period: 2 years

Urine CrossLaps® (CTX-I)EIA

	Placebo group		HRT II group	
	Mean (µg/mmoL)	95% Confidence Interval (µg/mmoL)	Mean (µg/mmoL)	95% Confidence Interval (µg/mmoL)
Baseline	251	206 – 305	217	163 – 291
after 6 months treatment	148	107 – 205	51	33 – 77




Urine CrossLaps® EIA versus α-BMD for patients treated with HRT

Using a cut-off for Urine CrossLaps® (CTX-I)EIA of 40% change from baseline the following sensitivities, specificities and 95% confidence intervals are obtained

	Tibolone		HRT I		HRT II	
	Sensitivity	Specificity	Sensitivity	Specificity	Sensitivity	Specificity
α-BMD>0	88% (80-96)	69% (44-96)	80% (73-87)	80% (68-92)	78% (92-94)	54% (24-88)
α-BMD>1	91% (83-99)	50% (30-70)	83% (76-90)	71% (59-83)	83% (68-98)	54% (27-81)

Below is a shown distribution plot for the combined HRT studies. The number over each bar indicates the number of participants in each class

Paget's disease

Below is shown the Urine CrossLaps® (CTX-I)EIA data [mg/mmol creatinine] from a Paget's disease study.


- Men and woman with Paget's disease documented by X-ray and bone scan
- 48 patients on active treatment: IV pamidronate, 3 injections of 60 mg
- Treatment period: 3 days

Urine CrossLaps® (CTX-I)EIA


	Pre treatment group		Post treatment group	
	Mean (µg/mmoL)	95% Confidence Interval (µg/mmoL)	Mean (µg/mmoL)	95% Confidence Interval (µg/mmoL)
Baseline	580	466 – 722	191	141 – 259

Below are shown distribution plots for the Paget's disease study. The number over each bar indicates the number of participants in each class.

Urine CrossLaps® EIA for Paget's patients

Urine CrossLaps® EIA change from baseline

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Example Version

REFERENCES

1. Burgeson RE. New collagens, new concepts. *Annu Rev Cell Biol* (1988); 4: 551-577.
2. Bonde M. et al., Immunoassay for quantifying type I collagen degradation products in urine evaluated. *Clin Chem* (1994); 40: 2022-2025.
3. Jaffé M. Über den niederschlag, welche pikrinsäure in normalem harn erzeugt und über eine neue reaktion des creatinins. *Hoppe Seylers Z Phys Chem* (1886); 10:391-400.
4. Garnero P. et al., Assessment of bone resorption with a new marker of collagen degradation in patients with metabolic bone disease. *J Clin Endocrinol Metab* (1994); 79: 780-785.
5. Bonde M. et al., Applications of an enzyme immunoassay for a new marker of bone resorption (CrossLaps) - follow-up on hormone replacement therapy and osteoporosis risk assessment. *J Clin Endocrinol Metab* (1995); 80: 864-868.
6. Delmas PD. et al., A combination of low doses of 17 beta-estradiol and norethisterone acetate prevents bone loss and normalizes bone turnover in postmenopausal women. *Osteoporos Int* (2000); 11: 177-187.
7. Ravn P. et al., Association between pharmacokinetics of oral ibandronate and clinical response in bone mass and bone turnover in women with postmenopausal osteoporosis. *Bone* (2002); 30: 320-324.
8. Buclin T. et al., Bioavailability and biological efficacy of a new oral formulation of salmon calcitonin in healthy volunteers. *J Bone Miner Res* (2002); 17: 1478-1485.
9. Ravn P. et al., Changes in biochemical markers and bone mass after withdrawal of ibandronate treatment: prediction of bone mass changes during treatment. *Bone* (1998); 22: 559-564.
10. Fledelius C. et al., Characterization of urinary degradation products derived from type I collagen. Identification of a beta-isomerized Asp-Gly sequence within the C-terminal telopeptide (alpha1) region. *J Biol Chem* (1997); 272: 9755-9763.
11. Garnero P. et al., Do markers of bone resorption add to bone mineral density and ultrasonographic heel measurement for the prediction of hip fracture in elderly women? The EPIDOS prospective study. *Osteoporos Int* (1998); 8: 563-569.
12. Bjarnason NH and Christiansen C. Early response in biochemical markers predicts long-term response in bone mass during hormone replacement therapy in early postmenopausal women. *Bone* (2000); 26: 561-569.
13. Delmas PD. et al., M. Monitoring individual response to hormone replacement therapy with bone markers. *Bone* (2000); 26: 553-560.
14. Tsai KS. et al., Raloxifene versus continuous combined estrogen/progestin therapy: densitometric and biochemical effects in healthy postmenopausal Taiwanese women. *Osteoporos Int* (2001); 12: 1020-1025.
15. Peichl P. et al., Serum osteocalcin and urinary crosslaps are suitable markers of bone turnover in response to short-term hormone replacement therapy. *Gynecol Endocrinol* (2000); 14: 374-381.
16. Hoshino H. et al., Urinary excretion of type I collagen degradation products in healthy women and osteoporotic patients with vertebral and hip fractures. *Calcif Tissue Int* (1998); 62: 36-39.

Doc: AC-03PL-A
Version: 08
Date: 28 Sep 2022

 EXP	GB Use By DE Verwendbar bis ES Fecha de caducidad IT Utilizzare entro FR Utiliser jusque NL Houdbaar tot DK Holdbar til CZ Použitelné do SK Použitelné do GR Ημερομηνία λήξης PT Prazo de validade HU Felhasználható SE Använd före PL Użyć przed	 LOT	GB Batch code DE Chargenbezeichnung ES Código de lote IT Codice del lotto Code du lot NL Lot nummer DK Lotnummer CZ Číslo šarže SK Číslo šarže GR Αριθμός Παρτίδας PT Código do lote HU Szrzszzám SE Lot nummer PL Kod partii
 REF	GB Catalogue number DE Bestellnummer ES Número de catálogo IT Numero di catalogo FR Référence du catalogue NL Catalogus nummer DK Katalognummer CZ Katalogové číslo SK Katalógové číslo GR Αριθμός καταλόγου PT Referência de catálogo HU Katalógusszám SE Katalognummer PL Numer katalogowy		GB Manufacturer DE Hersteller ES Fabricante IT Fabbicante FR Fabricant NL Fabrikant DK Producent CZ Výrobce SK Výrobca GR Κατασκευαστής PT Fabricante HU Gyártó SE Tillverkare PL Producent
 Σ	GB Contains sufficient for <n> tests DE Inhalt ausreichend für <n> Prüfungen ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi FR Contenu suffisant pour "n" tests NL Inhoud voldoende voor "n" testen DK Indeholder tilstrækkeligt til "n" test CZ Lze použít pro <n> testů SK Obsah postačuje na <n> stanovení GR Περιεχόμενο επαρκές για «n» εξετάσεις PT Conteúdo suficiente para "n" ensaios HU A doboz tartalma <n> vizsgálat elvégzéséhez elegendő SE Räckert till "n" antal tester PL Wystarczy na wykonanie <n> testów	 IVD	GB In-Vitro Diagnostic Medical Device DE In-Vitro-Diagnostikum ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro FR Dispositif médical de diagnostic in vitro NL Medisch hulpmiddel voor in-vitro diagnostiek DK Medicinsk udstyr til in vitro-diagnostik CZ In Vitro diagnostický zdravotnický prostředek SK Zdravotnícka pomocka in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro HU In vitro diagnosztikum SE Medicintekniska produkter för in vitro diagnostik PL Wyrób do diagnostyki In Vitro
	GB Temperature limitation DE Temperaturbegrenzung ES Límite de temperatura IT Limiti di temperatura FR Limites de température NL Temperatuurlimiet DK Temperaturbegrænsning CZ Teplotní rozmezí od do SK Teplotné rozmedzie od do GR Περιορισμοί θερμοκρασίας PT Limites de temperatura HU Hőmérsékletartomány SE Temperaturbegränsning PL Przestrzegać zakresu temperatury		GB Consult Instructions for Use DE Gebrauchsanweisung beachten ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso FR Consulter les instructions d'utilisation NL Raadpleeg de gebruiksaanwijzing DK Se brugsanvisning CZ Viz návod k použití SK Vid' návod na použitie GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização HU Nézze meg a Használati utasítást SE Se handhavandebeskrivningen PL Sprawdź w instrukcji obsługi

	GB EU Importer DE EU-Importeur ES UE Importador IT Importatore UE FR Importateur UE NL EU Importeur DK EU Importør CZ EU Dvozcce SK Dovozca z EÚ GR Εισαγωγέας ΕΕ PT Importador UE HU EU importőr SE EU Importör PL Importer w UE
---	--

Immunodiagnostic Systems Limited,
10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK
Tel.: +44 191 519 0660
Fax: +44 191 519 0760
E-Mail: info.uk@idsplc.com
www.idsplc.com

IDS Immunodiagnostic Systems GmbH
Herriotstrasse 1
60528 Frankfurt am Main
Germany