

fastGEN QC Recommendation

1. Threshold default filter setting table:
 - Minimal read depth: 500x
 - Minimal variant allele frequency (VAF): 1%
 - Synonymous variants: No

2. Provide a sheet with recommendations like min. raw reads per sample, min. coverage (1000), etc.
 - Cluster density: 800–1 200K clusters/mm² (MiSeq)
 - Total reads per sample in Genovesa:
 - Solid/Solid II Cancer kit: 100 000 reads per sample
 - Lung Cancer kit: 50 000 reads per sample
 - Brain Cancer kit: 25 000 reads per sample
 - POLE Cancer kit: 50 000 reads per sample
 - POLE/CTNNB1 Cancer kit: 70 000 reads per sample
 - PIK3CA Cancer kit: 125 000 reads per sample
 - TP53 Cancer kit: 200 000 reads per sample
 - Minimal variant read depth: 500x
 - Minimal Coverage/Depth
 - Wild type samples > 1 000 for each exon/region
 - Mutated samples
 - Minimal Total Depth > 500 for mutated variant VAF > 5 % (thus minimal variant count = 25)
 - Minimal Total Depth > 1000 for mutated variant VAF 2–5 % (thus minimal variant count = 20–50)
 - Minimal Total Depth > 2000 for mutated variant VAF 1–2 % (thus minimal variant count = 20–40)

3. Support e-mail address: support@biovendor-mdx.com