

ENG

Instructions for Use:
HUMAN MMP-9 ELISA

Catalogue number:
RBL002R

For research use only!

 **BioVendor**
R&D[®]



BioVendor – Laboratorní medicína a.s.

Karásek 1767/1, 621 00 Brno, Czech Republic

+420 549 124 185

info@biovendor.com

sales@biovendor.com

www.biovendor.com

1. INTENDED USE	3
2. STORAGE, EXPIRATION	3
3. INTRODUCTION	3
4. TEST PRINCIPLE	3
5. PRECAUTIONS	4
6. REAGENT SUPPLIED	4
7. MATERIAL REQUIRED BUT NOT SUPPLIED	4
8. PREPARATION OF REAGENTS	5
9. PREPARATION OF SAMPLES	5
10. ASSAY PROCEDURE	6
11. PERFORMANCE CHARACTERISTICS	6
12. CALCULATION	7
13. REFERENCES	8
14. EXPLANATION OF SYMBOLS	9
15. ASSAY PROCEDURE - SUMMARY	10

HISTORY OF CHANGES

Previous version	Current version
	ENG.001.A
New edition	

1. INTENDED USE

Enzyme Immunoassay for the quantitative determination of Matrix metalloproteinase 9 (MMP-9) in human plasma and serum.

2. STORAGE, EXPIRATION

- The kit must be stored at 2 – 8°C.
- The opened components can be stored for one week at 2 – 8°C.

3. INTRODUCTION

Matrix metalloproteinases (MMPs) are zinc-dependent endopeptidases, responsible for the integrity of the basement membrane (BM) via degradation of extracellular matrix and BM components.¹ MMPs are presented in various types of cells, including cardiomyocytes, cerebral neurons, hepatocytes, and many others, thus influencing various processes.

Increased matrix metalloproteinase 9 (MMP9), also called Gelatinase B, after myocardial infarction (MI) exacerbates ischemia-induced chronic heart failure (CHF).² It is also upregulated in the diabetic heart, and ablation of MMP9 decreases infarct size in the non-diabetic myocardial infarction heart.³ There is also strong evidence that the increased levels of MMP-9 are associated with many inflammatory conditions, such as coronary artery disease, COPD, arthritis, and metabolic syndrome.⁴ It was also proved that MMP-9 is associated with progression of atherosclerosis and higher risk of cardiovascular events. It was also suggested that elevated levels of MMP-9 (combined with MMP-2) are associated with highly tumorigenic cancers.

4. TEST PRINCIPLE

The microtiter plate is coated with the antibody specifically binding the Matrix metalloproteinase 9. The human serum or plasma is incubated in the plate with the capture antibody.

The specimen is washed out and the specifically bound protein is incubated with biotin-labelled detection antibody. Following another washing step, Streptavidin-HRP conjugate is added into the well.

Unbound reagent is then washed out. Horseradish peroxidase (HRP) bound in the complex reacts with the chromogenic substrate (TMB) creating the blue colour. The reaction is stopped by addition of STOP solution (H₂SO₄).

The absorbance values are measured at 450 nm (optionally 450/630 nm) and are proportional to the concentration of MMP-9 in the specimen. The concentration of MMP-9 in unknown samples is determined from the calibration curve which is created by plotting the absorbance values against the standard concentration values.

5. PRECAUTIONS

- For research use only
- For professional laboratory use
- The reagents with different lot numbers should not be mixed
- To prevent cross sample contamination, use disposable labware and pipette tips
- To protect laboratory stuff, wear protective gloves and protective clothing
- The substrate solution should remain colourless, keep it protected from light
- The test should be performed at standard laboratory conditions (temperature $25^{\circ}\text{C} \pm 2^{\circ}\text{C}$).

6. REAGENT SUPPLIED

Item	Qty.
Antibody Coated Microtiter Plate	96 wells
Biotin-labelled Antibody	13 mL
Streptavidin-HRP Conjugate	13 mL
Master Standard (lyophilized)	1 vial
Quality Control A (human serum, lyophilized)	1 vial
Quality Control B (human serum, lyophilized)	1 vial
Dilution Buffer	2x13 mL
Wash Buffer 15x conc.	50 mL
Substrate Solution	13 mL
STOP Solution	13 mL

7. MATERIAL REQUIRED BUT NOT SUPPLIED

- Glassware and test tubes
- Microtiter plate washer
- Precision pipettes (various volumes) with tips
- Orbital shaker
- Microtiter plate reader capable of measuring absorbance at 450 nm or 450/630 nm with software for data generation

8. PREPARATION OF REAGENTS

- Use new pipette tip for pipetting different reagents and samples to prevent cross-contamination.
- All reagents and samples should be allowed to reach the temperature $25^{\circ}\text{C} \pm 2^{\circ}\text{C}$.

8.1 Preparation of Standards

Reconstitute lyophilized Human MMP-9 Standard in Dilution Buffer, for the volume information see the Certificate of Analysis. Let it rehydrate for 15 min prior to use. The concentration of human MMP-9 in Master Standard is 20 ng/mL. Use the Master Standard for serial dilution (as below). Mix each tube thoroughly before the next transfer. The Dilution Buffer serves as Blank. Prepare set of Standard solution as follows:

	Volume of Standard	Dilution Buffer	Concentration
Std1	Standard 20 ng/mL (lyophilized)	See CofA	20 ng/mL
Std2	250 μL of Std1	250 μL	10 ng/mL
Std3	250 μL of Std2	250 μL	5 ng/mL
Std4	250 μL of Std3	250 μL	2.5 ng/mL
Std5	250 μL of Std4	250 μL	1.25 ng/mL
Std6	250 μL of Std5	250 μL	0.625 ng/mL
Blank	-	200 μL	0 ng/mL

8.2 Preparation of Quality Control A and B

Reconstitute the lyophilized human serum Quality Controls in deionized/distilled water, for the volume information see the Certificate of Analysis. Let the QCs rehydrate for 15 min and dilute them 1:50 in Dilution Buffer, prior to use, see Preparation of samples.

8.3 Preparation of Wash Buffer 1x

Prepare a working solution of Wash Buffer by adding 50 mL of Wash Buffer 15x conc. to 700 mL of deionized/ distilled water (dH_2O). Mix well. Store at 4°C for two weeks or at -20°C for long term storage.

9. PREPARATION OF SAMPLES

Human plasma or serum may be used with this assay. For long-term storage the samples should be frozen at minimum -70°C . Lipemic or haemolytic samples may cause false results. **Pay attention to a possibly elevated serum level of human MMP-9 due to MMP-9 release by platelets during sampling process. This may cause variable and irreproducible results.** Use a silica-based activator with polymer gel tubes for serum separation (SST/ BD Vacutainer, Serum gel (coagulation activator) / Sarstedt S-Monovette), or platelet-poor plasma. Platelet-poor plasma can be prepared in two-step procedure: 1/ centrifugation 20 minutes at $1500 \times g$, and 2/ additional centrifugation step of the plasma at $10\,000 \times g$ for 10 minutes for platelet removal. Recommended dilution of samples is 1:50, i.e., 5 μL of sample + 245 μL of Dilution Buffer, for duplicates and for singlets. Do not store the diluted samples.

10. ASSAY PROCEDURE

1. Prepare the reagents as described in the previous chapter.
2. Pipette 100 μL of set of Standards, Quality Controls, diluted Samples and Dilution Buffer = Blank into each well. Incubate for **1 hour** at $25^{\circ}\text{C} \pm 2^{\circ}\text{C}$, shaking at 300 rpm.
3. Wash the wells 3-times with 1x Wash Buffer (350 μL /well). When finished, tap the plate against the paper towel to remove the liquid completely.
4. Pipette 100 μL of Biotin-labelled Antibody into each well. Incubate for **1 hour** at $25^{\circ}\text{C} \pm 2^{\circ}\text{C}$, shaking at 300 rpm.
5. Wash the wells as described in point 3.
6. Pipette 100 μL of Streptavidin-HRP into each well. Incubate for **30 min** at $25^{\circ}\text{C} \pm 2^{\circ}\text{C}$, shaking at 300 rpm.
7. Wash the wells as described in point 3.
8. Pipette 100 μL Substrate solution, incubate for **10 min**, at $25^{\circ}\text{C} \pm 2^{\circ}\text{C}$. Avoid exposure to the light during this step.
9. Pipette 100 μL of STOP solution.
10. Read the signal at 450 or 450/630 nm within 15 min.

11. PERFORMANCE CHARACTERISTICS

Samples used in the tests were diluted 1:50 as recommended and assayed. The results are multiplied by the dilution factor.

11.1 Sensitivity

The limit of detection, defined as a concentration of human MMP-9 giving absorbance higher than absorbance of blank + 3 standard deviations, is better than 0.16 ng/mL of sample.

11.2 Precision

11.2.1 Intra-assay

Sample	Mean (ng/mL)	SD	CV (%)
1	527	9.7	2
2	224	6.7	3

11.2.2 Inter-assay (Run – to – run)

Sample	Mean (ng/mL)	SD	CV (%)
1	389	31	8
2	255	14	6

11.3 Accuracy

11.3.1 Dilution linearity

Sample	Dilution	Measured concentration (ng/mL)	Expected concentration (ng/mL)	Yield (%)
1		519	-	-
	2x	247	260	95


	4x	123	130	95
	8x	59	65	91
2		582	-	-
	2x	289	291	99
	4x	151	145	104
	8x	69	73	95

11.3.2 Spiking Recovery

Sample	Spike (ng/ml)	Measured concentration (ng/mL)	Expected concentration (ng/mL)	Yield (%)
1	-	223	-	-
	250.0	480	473	102
	125.0	357	348	103
	62.5	290	285	102

12. CALCULATION

The standard curve needs to be measured in every test. Most of the microplate reader can automatically calculate the analyte concentration using 4-parameter algorithm or alternative functions to fit the standard points properly. The concentrations need to be multiplied by the dilution factor, either automatically by reader or manually.












13. REFERENCES

- ¹ Łukaszewicz-Zajac M, Mroczko B, Słowik A. Matrix metalloproteinases (MMPs) and their tissue inhibitors (TIMPs) in amyotrophic lateral sclerosis (ALS). *J Neural Transm (Vienna)*. 2014 Nov;121(11):1387-97. doi: 10.1007/s00702-014-1205-3. Epub 2014 Jul 22. PMID: 25047909; PMCID: PMC4210652.
- ² (Nandi SS, Katsurada K, Sharma NM, Anderson DR, Mahata SK, Patel KP. MMP9 inhibition increases autophagic flux in chronic heart failure. *Am J Physiol Heart Circ Physiol*. 2020 Dec 1;319(6):H1414-H1437. doi: 10.1152/ajpheart.00032.2020. Epub 2020 Oct 16. PMID: 33064567; PMCID: PMC7792705.)
- ³ Yadav SK, Kambis TN, Kar S, Park SY, Mishra PK. MMP9 mediates acute hyperglycemia-induced human cardiac stem cell death by upregulating apoptosis and pyroptosis in vitro. *Cell Death Dis*. 2020 Mar 13;11(3):186. doi: 10.1038/s41419-020-2367-6. PMID: 32170070; PMCID: PMC7070071.
- ⁴ Snitker S. Correlation of Circulating MMP-9 with White Blood Cell Count in Humans: Effect of Smoking. *PLoS ONE*. 2013; 8(6): e66277

Example Version

14. EXPLANATION OF SYMBOLS

	Catalogue number
	Batch code
	Caution
	Use by date
	Temperature limit
	Manufacturer
 <p data-bbox="256 1182 464 1205">www.biovendor.com</p>	Read electronic instructions for use - eIFU
	The content is sufficient for 96 tests
	Biological risks

15. ASSAY PROCEDURE - SUMMARY

Add 100 μ L of Standards, diluted QCs and Samples to the wells



Incubate for 1 hour at 25°C, shaking at 300 rpm

3-times wash the wells (350 μ L/well)



Add 100 μ L of Biotin-labelled Antibody to the wells



Incubate for 1 hour at 25°C, shaking at 300 rpm

3-times wash the wells (350 μ L/well)



Add 100 μ L of SAV-HRP to the wells



Incubate for 30 min at 25°C, shaking at 300 rpm

3-times wash the wells (350 μ L/well)



Add 100 μ L of Substrate Solution to the wells



Incubate for 10 min in the dark at 25°C, NO shaking

Add 100 μ L of Stop Solution to the wells



Read the signal at 450 nm (450/630 nm) within 15 min



BioVendor – Laboratorní medicína a.s.

Karásek 1767/1, 621 00 Brno, Czech Republic

+420 549 124 185

info@biovendor.com

sales@biovendor.com

www.biovendor.com



FOLLOW US