IGFBP-2 ELISA

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ENGLISH - IGFBP-2 ELISA E05 INSTRUCTIONS FOR USE

TABLE OF CONTENTS

ENGLISH - IGFBP-2 ELISA E05 INSTRUCTIONS FOR USE	20
1 INTENDED USE	20
2 INTRODUCTION	20
3 PRINCIPLE	21
4 WARNINGS AND PRECAUTIONS	22
5 SAMPLES	23
6 MATERIALS	24
7 TECHNICAL NOTES	25
8 ASSAY PROCEDURE	26
9 QUALITY CONTROL	27
10 EVALUATION OF RESULTS	27
11 INTERPRETATION OF RESULTS	28
12 LIMITATIONS OF PROCEDURE	28
13 REFERENCE VALUES	29
14 PERFORMANCE CHARACTERISTICS	31
15 SCIENTIFIC APPLICATION	33
16 LITERATUR/ LITERATURE	34
17 INTERNATIONAL ASSAY DESCRIPTION	35
18 ASSAY PROCEDURE	36

ENGLISH - IGFBP-2 ELISA E05 INSTRUCTIONS FOR USE

IGFBP-2 ELISA E05	96 Determinations
CE	DE/CA40/00809/13/1
Principle of the test	Enzyme Immunoassay
Duration (incubation period)	1,75 h
Antibodies	specific, monoclonal antibody and high-affinity polyclonal antiserum, ready for use
Buffer	Ready for use and 20fold conc.
Calibrator	Single Calibrators with 2; 10; 20; 40 und 80 ng/mL IGFBP-2 rec. hIGFBP-2 , lyophilised
Calibration	Calibrated against recombinant IGFBP-2
Assay Range	0.2 ng/mL – 1680 ng/mL
Control	2 Controls, lyophilised
Sample	Serum und EDTA-Plasma
Required sample volume	15 µL
Sample dilution	1:21
Analytical sensitivity	ø 0.2 ng/mL
Intra- / Interassay Variance	ø ≤ 10 %
Reference values	Normal values for healthy individuals (1.5 to > 70 years) were evaluated for this assay.

1 INTENDED USE

Measurement of human IGFBP-2 in human serum, EDTA-plasma.

2 INTRODUCTION

Insulin-like growth factors (IGFs) regulate the proliferation, differentiation, apoptosis, cell adhesion and metabolism in various tissues and cell types. The IGF-I, which is produced mainly in liver under the influence of growth hormone (GH), regulates as hormone the linear growth of the bones and the process of sexual maturity, while IGF-II is mainly a growth factor of fetal tissue (11-13). The biological actions of IGF over the IGF-Type-I receptor are modulated variably through the IGF binding proteins (IGFBP-1 to-6) (14). IGFBP-2 is, after IGFBP-3, the second most frequent IGFBP in the human blood. IGFs, especially tumor typical pro-IGF-forms and hormones regulate the expression of IGFBP-2, GH effect is thereby inhibiting. At cellular level IGFBP-2 seems to stimulate the proliferation and dissemination of solid tumors via an IGF-independent mechanism (15,16).

IGFBP-2 is a unglycosylated polypeptide of 31.3 kDa, which forms binary IGF-complexes and shows no circadian rhythm in the circulation. The serum concentration of IGFBP-2 increases in fasting, after major surgery and after trauma, but the increasing of the concentration is most intensive in malignant diseases. The correlation of the IGFBP-2 level to the degree of progression is a striking feature in various tumor types as is the normalization of the IGFBP-serum levels after remission (5-10). During the GH-therapy, e.g. in short stature and in GH-abuse (doping) the IGFBP-2 level decreases. In Trisomy 18 IGFBP-2 in maternal serum is decreased and IGFBP-1 is increased; therefore the ratio IGFBP-2 /IGFBP-1 is a marker for this chromosome abnormality (17).

Clinical Implications

The IGFBP-2 concentration is age-dependent in blood (3).

Supplementary parameter to IGFBP-3 in the diagnosis of growth disorders (IGFBP-2/IGFBP-3 ratio), IGFBP-2 is an inhibitor of growth hormone action (3,4).

Progression-dependent tumor marker in leukaemia (5), astrocytic CNS tumors (6,7), prostate- (8), suprarenal cortex- (9)-, hepatocellular (10) and other carcinomas.

Anti-aging parameter: IGFBP-2 as a marker of physiological functionality (20).

3 PRINCIPLE

The Mediagnost ELISA for IGFBP-2 E05 is a so-called Sandwich-Assay using two specific and high affinity antibodies. The IGFBP-2 in the samples binds to the first antibody coated on the microtiter plate. In the following step the second specific anti-IGFBP-2-Antibody binds in turn to the immobilised IGFBP-2. The second antibody is biotinylated and will be applied in a mixture with a Streptavidin-Peroxidase-Enzyme Conjugate. In the closing substrate reaction the turn of the colour will be catalysed quantitatively depending on the IGFBP-2-level of the samples.

It recognizes IGFBP-2 quantitatively and is not influenced by increased IGF-I or IGF-II values. Related molecules such as IGFBP-3 do not show any cross-reactions in the test.

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4 WARNINGS AND PRECAUTIONS

For In Vitro Diagnostic Use only. For Professional use only.

The Mediagnost kit is suitable only for in vitro diagnostics and not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. Mediagnost will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided.

Do not use obviously damaged or microbial contaminated or spilled material.

Caution: This kit contains material of human and/or animal origin. Therefore all components and patient's specimens should be treated as potentially infectious.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations. Material Safety Data Sheet is available on request.

Human Serum

Following components contain human serum: **Control CTR1 and CTR 2**

Source human serum for the control sera provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore all components and patient's specimens should be treated as potentially infectious.

Reagents CAL A-E, DET, EC, DIL, WB

Contain as preservative a mixture of **5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one** (<0.015%)

H317	May cause an allergic skin reaction.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P501	Dispose of contents/ container in accordance with local/ regional/ national/ international regulations.

Substrate S

The TMB-Substrate (S) contains 3,3',5,5' Tetramethylbencidine (<0.05%)

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.

Stop Solution STP

The Stop Solution contains 0.2 M acid sulphuric acid (H₂SO₄)

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P301+P330+	IF SWALLOWED: rinse mouth.
P331	Do NOT induce vomiting.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.
P309+P310	IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician.

4.1 General first aid procedures:

Skin contact: Wash affected area rinse immediately with plenty of water at least 15 minutes. Remove contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing spread the eyelids.

Ingestion: After swallowing the product, if the affected person is conscious, rinse out the mouth with plenty of water: seek medical advice immediately.

5 SAMPLES

5.1 Sample type

Serum and EDTA-Plasma

5.2 Specimen collection

The blood sample for serum preparation should be gained according to standardized venipuncture procedure. Hemolytic reactions have to be avoided. The blood has to be allowed to clot and after complete clotting, serum is separated by centrifugation.

Whole blood should ideally be separated after about 2 hours and processed after 12 hours at the latest and stored at -20 ° C until the measurement.

Blood samples may be taken at any time of the day. Haemolytic reactions are to be avoided.

5.3 Required sample volume: 15 µL

5.4 Sample stability

In firmly closed sample vials:

- Storage at RT max. 24 hours
- Storage at -20°C min. 2 years

Samples are not allowed to have more than 10 freeze/thaw cycles.

5.5 Interference

Triglycerides bilirubin and hemoglobin in the sample do not interfere up to a concentration of 100 mg/mL, 200 µg/mL and 1 mg/mL respectively. The use of hemolytic, lipemic or icteric samples should nevertheless be validated beforehand by the user.

5.6 Sample dilution

- Dilution: **1:21** with Dilution Buffer **DIL**
- Example: 15 µL sample are added to 300 µL Dilution Buffer **DIL** (dilution factor 21). After mixing this solution, 2 x 100 µL are used in the assay.
- The serum samples must be diluted 1:10 – 1:30 with Dilution Buffer **DIL** prior to measurement, depending on the expected IGFBP-2 values. In general, a dilution of 1:21 is suitable (and thus the minimum required sample volume for a duplicate determination is 15 µL serum). An extraction step is not required.

6 MATERIALS

6.1 Materials provided

The reagents listed below are sufficient for 96 wells including the calibration curve.

MTP	Microtiter plate , ready for use, coated with mouse-anti-hlGFBP-2-antibody. Wells are separately breakable.	(8x12) wells
CAL A-E	Calibrators , lyophilized, (rec. hlGFBP-2), concentrations are given on vial labels and on quality certificate.	5 x 750 µL
CTR1	Control 1 , lyophilized, (human serum), concentration is given on quality certificate in ng/mL.	1 x 100 µL
CTR2	Control 2 , lyophilized, (human serum), concentration is given on quality certificate in ng/mL.	1 x 100 µL
DET	Antibody-HRP-Conjugate , rabbit-anti-hlGFBP-2-antibody biotinylated + streptavidin horseradish peroxidase conjugate ready for use.	1 x 12 mL
DIL	Dilution Buffer , ready for use	1 x 50 mL
WB	Washing Buffer , 20-fold concentrated solution	1 x 50 mL
S	Substrate , ready for use, horseradish-peroxidase-(HRP) substrate, stabilised Tetramethylbenzidine.	1 x 12 mL
STP	Stop Solution , ready for use, 0.2 M sulphuric acid.	1 x 12 mL
-	Sealing Tape , for covering the microtiter plate .	2 x
	Instructions for use	1 x
-	Quality Certificate	1 x

6.2 Materials required, but not provided

- Distilled (Aqua destillata) or deionized water for dilution of the Washing Buffer **WB (A. dest.)**, **950 mL**.
- Precision pipettes and multichannel pipettes with disposable plastic tips
- Polyethylene PE/Polypropylene PP tubes for dilution of samples
- Vortex-mixer
- Microtiter plate shaker (350 rpm)
- Microtiter plate washer (recommended)
- Micro plate reader ("ELISA-Reader") with filter for 450 and ≥ 590 nm

7 TECHNICAL NOTES

Storage Conditions

Store the kit at 2-8°C after receipt until its expiry date. The lyophilized reagents should be stored at -20 °C after reconstitution. Avoid repeated thawing and freezing.

Storage Life

The shelf life of the components **after initial opening** is warranted for **4 weeks**, store the unused strips and microtiter wells **airtight** together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The **reconstituted components** Calibrators **A-E** and Control **CTR1** and **CTR 2** must be stored at -20°C (max. 4 weeks). For further use, thaw quickly but gently (avoid temperature increase above room temperature and avoid excessive vortexing). Up to 3 of the freeze-thaw cycles did not influence the assay. The 1:20 diluted Washing Buffer **WB** is 4 weeks stable at 2-8°C

Preparation of reagents

Bring all reagents to room temperature (20 - 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

Reconstitution

The Calibrators **A – E** and Control **CTR1** and **CTR2** are reconstituted with the Dilution Buffer **DIL**. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

Dilution

After reconstitution dilute the Control **CTR1** and **CTR2** with the Dilution Buffer **DIL** in the same ratio (1:21) as the sample. The required volume of Washing Buffer **WB** is prepared by 1:20 dilution of the provided 20-fold concentrate with Aqua dest.

Assay Procedure

When performing the assay, Blank, Calibrators **A-E**, Controls **CTR1** and **CTR2** and the samples should be pipette as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times, Antibody-POD-Conjugate **DET** as well as the succeeding Substrate **S** should be added to the plate in the same order and in the same time interval as the samples. Stop Solution **STP** should be added to the plate in the same order as Substrate **S**. All determinations (Blank, Calibrators **A-E**, Control **CTR1** and **CTR2** and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

Incubation

Incubation at room temperature means: Incubation at 20 - 25°C. The Substrate **S**, stabilised Tetramethylbencidine, is photosensitive—store and incubation in the dark.

Shaking

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtiter plate shaker. We are recommending approx. 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must be adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/ or false values, excessive shaking may result in high optical densities and/ or false values.

Washing

Proper washing is of basic **importance** for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer **WB** diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

The danger of handling with potentially infectious material must be taken into account.

Manual washing should be used. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. Decant contents into a biohazard bin, then blot plate on absorbent tissue. Wash the plate by adding 300µL Washing Buffer WB/well, then decant and blot on absorbent tissue. Repeat this step 4 more times for total of 5 washes.

8 ASSAY PROCEDURE

Preparation of reagents		Reconstitution:	Dilution
CAL A-E	Calibrators	in 750 µL Dilution Buffer DIL	-
CTR1	Control 1	in 100 µL Dilution Buffer DIL	1:21 with Dilution Buffer DIL
CTR2	Control 2	in 100 µL Dilution Buffer DIL	1:21 with Dilution Buffer DIL
WB	Washing Buffer Conc.	-	1:20 with Aqua dest. → WB 1:20
Dilute samples 1:21 in Dilution Buffer DIL, mix immediately, Use 100 µL for each well in the assay.			
Before assay procedure bring all reagents to room temperature +20-25°C.			
Assay procedure in double determination			
Pipette	Reagents		Position
100 µL	Dilution Buffer DIL (Blank)		A1/A2
100 µL	Calibrator A (2 ng/mL)		B1/B2
100 µL	Calibrator B (10 ng/mL)		C1/C2
100 µL	Calibrator C (20 ng/mL)		D1/D2
100 µL	Calibrator D (40 ng/mL)		E1/E2
100 µL	Calibrator E (80 ng/mL)		F1/F2
100 µL	Control CTR1 (1:21 diluted)		G1/G2
100 µL	Control CTR2 (1:21 diluted)		H1/H2
100 µL	Sample SPE (1:21 diluted)		in the rest of the wells according the requirements
Cover the wells with the sealing tape.			
Sample-Incubation: 1 h at +20-25°C, 350 rpm			
5x 300 µL	Decant the contents of the wells and wash 5 x with 300 µL each Washing Buffer WB 1:20 / well		In each well
100 µL	Antibody-HRP-Conjugate DET		In each well
Cover the wells with the sealing tape.			
Incubation: 30 Minutes at +20-25°C, 350 rpm			
5x 300 µL	Decant the contents of the wells and wash 5 x with 300 µL each Washing Buffer WB 1:20 / well		In each well
100 µL	Substrate S		In each well
Incubation: 15 Minutes in the Dark at +20-25°C			
100 µL	Stop Solution STP		In each well
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.			

9 QUALITY CONTROL

Good laboratory practice requires that controls are included in each assay. The test results are only valid if the test has been performed following the instructions. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state or local standards/laws.

9.1 Quality criteria

For the evaluation of the assay it is required that the absorbance values of the blank should be below 0.25, and the absorbance of **Calibrator E** should be above 1.00. Samples, which yield higher absorbance values than **Calibrator E**, should be re-tested with a higher dilution.

10 EVALUATION OF RESULTS

10.1 Establishing of the calibration curve

Calibrators are provided in the following IGFBP-2 concentrations

Calibrator	A	B	C	D	E
ng/mL	2	10	20	40	80

- 1) Calculate the **mean absorbance** value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all samples, controls and calibrators.
- 3) Plot the calibrator concentrations on the x-axis versus the mean value of the absorbance of the calibrators on the y-axis.
- 4) Calculation of the calibration curve should be done by using a computer program, because the curve is in general (without respective transformation) not ideally described by linear regression. A **four parametric logistic (4-PL) curve fit** should be used for recalculation of IGFBP-2 concentrations.
- 5) The IGFBP-2 concentrations of the diluted samples or the diluted control CTR1 & CTR2 in ng/mL is calculated in this way, the IGFBP-2 concentration of the **undiluted sample** and of CTR1 & CTR2 is calculated **by multiplication** with the respective dilution factor.

10.2 Example of a typical calibration curve

The exemplary shown calibration curve in Fig.1 **cannot be used** for calculation of your test results. You have to establish a calibration curve for each test you conduct!


Fig 1. Exemplary Calibration Curve with a polynomial 3rd degree as curve fit

Exemplary calculation of the IGFBP-2 concentration of undiluted sample:

Measured extinction of your sample	0.37
Measured extinction of the blank	0.06

Your **measurement programm** will calculate the IGFBP-2 concentration of the diluted sample automatically by using the difference (0.31) of sample and blank for the calculation. You only have to determine the most suitable curve fit (here: polynomial 3rd degree).

In this exemplary case the following equation is solved by the programm to calculate the IGFBP-2 concentration in the sample:

$$0.31 = -0.0012048 + 0.039581x + 5.1788 \cdot 10^{-5} \cdot x^2 - 1.8929 \cdot 10^{-6} \cdot x^3$$

$$7.93 = x$$

if the dilution factor (**1:21**) is taken into account the IGFBP-2 concentration of the undiluted sample is $7.93 \cdot 21 = 166.55 \text{ ng/mL}$

11 INTERPRETATION OF RESULTS

The test results should not be the only base for therapeutic decisions. The results should be interpreted in regard to anamnesis, further clinical observations and results of other diagnostic investigations. Further, it is recommended to establish reference and cut-off values corresponding to the relevant group of patients for each laboratory.

12 LIMITATIONS OF PROCEDURE

Deviation from the reference range can be expected especially in hypothyroidism, after major surgery, in polytrauma, in Diabetes mellitus (due to insulin therapy), in fasting and in malignant diseases.

The Mediagnost IGFBP-2 ELISA is based on mono- and polyclonal antibodies. Generally the result of any immunological test system can be influenced by heterophilic antibodies, anti-species antibodies or rheumatic factors. The assay design reduces these potential influences to a minimum but they cannot be excluded in any case.

13 REFERENCE VALUES

The IGFBP-2 concentration in serum is depended on Body Mass Index (BMI; Table 1) and on age (Table 2). For data collection of these reference values IGFBP-2 levels were determined in serum of over 400 normal children and adults (see table 2 figure 2); (3).

Table 1: BMI dependent reference values, adults between 20 and 80 years

BMI		mean		percentiles		
[kg/m ²]	n	IGFBP-2 [ng/mL]	SD	5.	50.	95.
15	12	612	110	431	612	793
17.5	14	568	126	361	568	775
20	76	509	144	271	509	746
22.5	124	449	162	182	449	716
25	101	398	165	127	398	670
27.5	52	348	147	106	348	590
30	25	315	118	120	315	510
32.5	15	282	90	135	282	430
35	4	251	80	119	251	383
37.5	4	220	71	104	220	336

SD = Standard Deviation, n=number

Table 2: IGFBP-2 serum levels (in ng/mL) of > 400 healthy individuals. The normal range is given by the 5., 50. and 95. percentile for age classes.

Age-dependent normal range of serum IGFBP-2

age (years)	5. percentile (ng/mL)	50. percentile (ng/mL)	95. percentile (ng/mL)
1	408	545	728
2	359	500	696
3	317	460	668
4	277	421	640
5	243	388	617
6	217	361	602
7	194	336	583
8	173	312	562
9	154	289	542
10	138	268	522
11	123	249	503
12	111	232	486
13	101	219	477
14	94	212	470
15	89	207	465
16	86	207	460
17	84	214	466
18	84	223	483
19	84	232	500
25	99	280	580
35	110	381	686
45	130	403	702
55	140	410	715
65	151	418	727
75	153	427	740
80	156	430	744


Fig. 2: IGFBP-2 serum levels (in ng/mL) of > 400 healthy individuals. The normal range is given by the 5th, 50th and 95th percentile.

14 PERFORMANCE CHARACTERISTICS

14.1 Sensitivity

Sensitivity was assed by measuring the blank and calculating the theoretical concentration of the blank + 2SD. The analytical sensitivity of the E05 is 0.2 ng/mL.

14.2 Specificity

The specificity of this ELISA was investigated. Neither 100 ng/mL IGFBP-1 nor 20 µg/mL IGFBP-3 containing samples resulted in a signal significantly different of the blank.

14.3 Precision

Intra-Assay Variance

The inter- and intra assay coefficients of variability are below 10%. Exemplary determinations are shown in table 3 and 4.

Table 3: Intra-Assay-Variation

Sample 1 ng/mL	322	375	298	305	318	311	320	325	302	301	305	317	6.48% CV
Sample 2 ng/mL	612	609	616	648	594	597	620	613	617	611	636	698	4.49% CV

CV: Coefficient of Variation

14.4 Inter-Assay

Serum samples were measured in independent assays. Exemplary results are shown in Table 4.

Table 4: Inter-Assay-Variation

Sample 1 (ng/mL)	137	159	152	8% CV
Sample 2 (ng/mL)	672	697	688	2% CV
Sample 3 (ng/mL)	928	929	956	2% CV

CV: Coefficient of Variation

14.5 Linearity

The Mediagnost IGFBP-2 ELISA E05 can be diluted over a very wide range. The linearity of serum dilutions is shown in table 5.

Table 5: Linearity of the sample dilution

Dilution	Serum 1 (ng/mL)	Serum 2 (ng/mL)
1:10	938	582
1:20	1061	673
1:40	1055	719
1:80	1004	691
1:160	952	668

14.6 Recovery

Recombinant IGFBP-2 was added in three different concentrations to human serum. The IGFBP-2 concentration was measured and the mean relative recovery in comparison to buffer was 108%. Some exemplary data are shown in table 6.

Table 6: Recovery of recombinant human IGFBP-2 in Serum

IGFBP-2 [ng/mL]	+1000 ng/mL	+500 ng/mL	+100 ng/mL	Mean [ng/mL]
% Recovery	100,00	112,00	114,00	108,67

14.7 Interference

Interference of bilirubin, triglycerides and hemoglobin was tested by adding different amounts of these substances to human serum containing IGFBP-2. For comparison the same amount of buffer without any substance was also added to the serum. Table 7 demonstrates that neither bilirubin nor triglycerides exert any influence on the measurement of IGFBP-2 in human serum.

Table 7: Interference

Bilirubin		Triglycerides		Hemoglobin	
[µg/mL]	%	[mg/mL]	%	[mg/mL]	%
25	95.07	12.5	100.79	0.125	99
50	92.80	25	101.01	0.25	105
100	93.83	50	103.65	0.5	100
200	88.15	100	101.34	1	100

Instructions for use for scientific application

15 SCIENTIFIC APPLICATION

IGFBP-2 is present in different concentrations in various body fluids and in conditioned cell culture media of many cell lines.

15.1 Samples suitable for scientific application

Cerebrospinal fluid, breast milk, amniotic fluid, saliva and in cell culture medium.

Table 8: Linearity of the sample dilution.

Dilution	Cerebrospinal fluid (ng/mL)	Amniotic fluid (ng/mL)
1:10	426	Not determined
1:20	428	460
1:40	379	483
1:80	318	431
1:160	426	415

Table 9: Expected values of IGFBP-2 in body fluids of human origin and in cell culture supernatants:

Sample	Expected Value [ng/mL]
Serum	[100 - 1000]
Cerebrospinal fluid	[100 - 300]
Amniotic fluid	[200 - 10000]
Seminal plasma	[5000 - 15000]
Breast milk	[1500 - 3000]
Cell culture supernatants	[5 - 300]

Table 10: Results of saliva-sample matrix tests:.

Recombinant IGFBP-2: 38,2 [ng/mL] % Recovery		
Saliva + DIL	2.02	n.a.
Saliva + rec. IGFBP-2	35.7	89
Saliva + rec. IGFBP-2	37.36	93
Saliva + rec. IGFBP-2	35.67	89
Saliva + rec. IGFBP-2	36.17	90

15.2 Cross reactions with animal samples

This assay is specific for human IGFBP-2, low degree of cross reactions was found with commercial dog, horse, donkey, cat and goat sera. No cross-reactivity was found with pig, bovine, rabbit, mouse, chicken, rat, guinea pig, sheep sera.

16 LITERATUR/ LITERATURE

Literature

1. Elmlinger MW, Wimmer K, Biemer E, Blum WF, Ranke MB, Dannecker GE, (1996) Insulin-like growth factor binding protein 2 is differentially expressed in leukaemic T-and B-cell lines. *Growth Regulation* 6, 152-157
2. Hoeflich A, Nedbal S, Blum WF, Erhard M, Lahm H, Brem G, Kolb HJ, Wanke R, Wolf E. (2001) Growth inhibition in giant growth hormone transgenic mice by overexpression of insulin-like growth factor-binding protein-2. *Endocrinology* 142:1889-98
3. Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2000) Significance of basal IGF-I, IGFBP-3 and IGFBP-2 measurements in the diagnostics of short stature in children. *Hormone Research* 54: 60-68
4. Ranke MB, Schweizer R, Elmlinger MW, Weber K, Binder G, Schwarze CP, Wollmann HA (2001). Relevance of IGF-I, IGFBP-3, and IGFBP-2 measurements during GH treatment of GH-deficient and non-GH-deficient children and adolescents. *Hormone Research* 55: 115-124
5. Crofton PM, Ahmed SF, Wade JC, Elmlinger MW, Ranke MB, Kelnar CJH, Wallace WHB (1999) Effects of a third intensification block of chemotherapy on bone and collagen turnover, insulin-like growth factor I, its binding proteins and short-term growth in children with acute lymphoblastic leukemia. *European Journal of Cancer* 35: 960-967
6. Muller HL, Oh Y, Lehnbecher T, Blum WF, Rosenfeld RG.(1994) Insulin-like growth factor-binding protein-2 concentrations in cerebrospinal fluid and serum of children with malignant solid tumors or acute leukemia. *J Clin Endocrinol Metab.* 79: 428-34
7. Martin W Elmlinger, Martin H Deininger, Burkhardt S Schuett, Richard Meyermann, Frank Duffner, Ernst H Grote, and Michael B Ranke (2001) In vivo expression of the insulin-like growth factor binding protein-2 in human gliomas increases with the tumor grade. *Endocrinology* 142: 1652-1658
8. Cohen P, Peehl DM, Stamey TA, Wilson KF, Clemmons DR, Rosenfeld RG (1993) Elevated levels of insulin-like growth factor-binding protein-2 in the serum of prostate cancer patients. *J Clin Endocrinol Metab.* 76: 1031-1035
9. Boulle N, Baudin E, Gicquel C, Logie A, Bertherat J, Penfornis A, Bertagna X, Luton JP, Schlumberger M, Le Bouc Y. (2001) Evaluation of plasma insulin-like growth factor binding protein-2 as a marker for adrenocortical tumors. *Eur J Endocrinol.* 144: 29-36
10. Ranke MB, Maier KP, Schweizer R, Stadler B, Schleicher s, Elmlinger MW, Flehmig B (2003), Pilot study of elevated levels of insulin-like growth factor-binding protein-2 as indicators of hepatocellular carcinoma. *Horm Res* 60:174-180
11. Rosenfeld RG, Roberts CT Jr. (eds.) (1999) *The IGF system: Contemporary Endocrinology Series*; Humana Press
12. Jones JI, Clemmons DR. (1995) Insulin-like growth factors and their binding proteins: biological actions. *Endocr Rev.* 16: 3-34
13. Chard T (1994) Insulin-like growth factors and their binding proteins in normal and abnormal fetal growth. *Growth Reg* 4: 91-100
14. Ranke MB, Elmlinger MW (1997) Functional role of insulin-like growth factor binding proteins. *Hormone Research* 48 (suppl 4): 9-15
15. Elmlinger MW, Bell M, Schütt, BS, Langkamp M, Kutoh E, Ranke MB (2001) Transactivation of the IGFBP-2 promoter in human tumor cell lines. *Molecular and Cellular Endocrinology* 175: 211-218
16. Hoeflich A, Reisinger R, Lahm H, Kiess W, Blum WF, Kolb HJ, Weber MM, Wolf E. (2001) Insulin-like growth factor-binding protein 2 in tumorigenesis: protector or promoter? *Cancer Res.*15: 8601-8610
17. Miell JP, Langford KS, Jones JS, Noble P, Westwood M, White A, Nicolaides KH. (1997) The maternal insulin-like growth factor (IGF) and IGF-binding protein response to trisomic pregnancy during the first trimester: a possible diagnostic tool for trisomy 18 pregnancies. *J Clin Endocrinol Metab.*82: 287-292
18. Elmlinger MW., Regine Grund, Michael Buck, Hartmut A Wollmann, Norman Feist, Matthias M. Weber, Christian C. Speer and Michael B. Ranke (1999). Limited proteolysis of the insulin-like growth factor binding protein (IGFBP)-2 by a specific serine protease activity in early breast milk. *Pediatric Research* 46: 76-81
19. Elmlinger MW, Michael S Sanatani, Mathias Bell, Günther E Dannecker, Michael B Ranke (1998) Elevated insulin-like growth factor (IGF) binding protein (IGFBP)-2 and IGFBP-4 expression of leukemic T-cells is affected by autocrine/paracrine IGF-II action but not by IGF type I receptor expression. *European Journal of Endocrinology* 138: 337-343
20. van den Beld, Blum W, Pols H, Grobbee D, Lamberts S (2003) Serum insulin-like growth factor binding protein-2 levels as an indicator of functional ability in elderly men. *European Journal of Endocrinology* 148: 627-634

17 INTERNATIONAL ASSAY DESCRIPTION

International Assay Description

CAL A-E	Rec in 750 µL DIL	-
CTR1	Rec in 100 µL DIL	1:21 DIL
CTR2	Rec in 100 µL DIL	1:21 DIL
WB 20x	-	1:20 A. dest. → WB 1:20
SPE	-	1:21 DIL
°C +20-25°C ; max 350 rpm		
100 µL	DIL	A1/A2
100 µL	CAL A (2 ng/mL)	B1/B2
100 µL	CAL B (10 ng/mL)	C1/C2
100 µL	CAL C (20 ng/mL)	D1/D2
100 µL	CAL D (40 ng/mL)	E1/E2
100 µL	CAL E (80 ng/mL)	F1/F2
100 µL	CTR1 1:21 DIL	G1/H1
100 µL	CTR1 1:21 DIL	H1/H2
100 µL	SPE 1:21 DIL	
TAPE		
1 h °C +20-25°C 350 rpm		
5x 300 µL	5x WB 1:20	
100 µL	DET	
TAPE		
0.5 h °C +20-25°C 350 rpm		
5x 300 µL	5x WB 1:20	
100 µL	S	
15 min °C +20-25°C		
STP		
MEASURE		

18 ASSAY PROCEDURE

Preparation of reagents		Reconstitution:	Dilution
CAL A-E	Calibrators	in 750 µL Dilution Buffer DIL	-
CTR1	Control 1	in 100 µL Dilution Buffer DIL	1:21 with Dilution Buffer DIL
CTR2	Control 2	in 100 µL Dilution Buffer DIL	1:21 with Dilution Buffer DIL
WB	Washing Buffer Conc.	-	1:20 with Aqua dest. → WB 1:20
Dilute samples 1:21 in Dilution Buffer DIL, mix immediately, Use 100 µL for each well in the assay.			
Before assay procedure bring all reagents to room temperature +20-25°C.			
Assay procedure in double determination			
Pipette	Reagents		Position
100 µL	Dilution Buffer DIL (Blank)		A1/A2
100 µL	Calibrator A (2 ng/mL)		B1/B2
100 µL	Calibrator B (10 ng/mL)		C1/C2
100 µL	Calibrator C (20 ng/mL)		D1/D2
100 µL	Calibrator D (40 ng/mL)		E1/E2
100 µL	Calibrator E (80 ng/mL)		F1/F2
100 µL	Control CTR1 (1:21 diluted)		G1/G2
100 µL	Control CTR2 (1:21 diluted)		H1/H2
100 µL	Sample SPE (1:21 diluted)		in the rest of the wells according the requirements
Cover the wells with the sealing tape.			
Sample-Incubation: 1 h at +20-25°C, 350 rpm			
5x 300 µL	Decant the contents of the wells and wash 5 x with 300 µL each Washing Buffer WB 1:20 / well		In each well
100 µL	Antibody-POD-Conjugate DET		In each well
Cover the wells with the sealing tape.			
Incubation: 30 Minutes at +20-25°C, 350 rpm			
5x 300 µL	Decant the contents of the wells and wash 5 x with 300 µL each Washing Buffer WB 1:20 / well		In each well
100 µL	Substrate S		In each well
Incubation: 15 Minutes in the Dark at +20-25°C			
100 µL	Stop Solution STP		In each well
Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.			

Copy of the page 26