IGF-II ELISA

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ENGLISH**Instructions for use**

IGF-II ELISA	96 Determinations
CE	DE/CA40/00809/26
Principle of the test	Enzyme Immunoassay
Duration (incubation period)	3 h
Antibody Conjugate	ready for use
Enzyme Conjugate	ready for use
Buffer	ready for use and 20fold concentrate
Substrate	ready for use
Calibrators	5 single Calibrators: 0.45 - 9 ng/mL, lyophilized, human recombinant IGF-II
Reference material	Calibrated against the International Standard WHO/NIBSC 96/538
Assay Range	0.06 – 3636 ng/mL
Control	2 Controls, lyophilised
Sample	human serum / plasma
Required sample volume	10 µL
Sample dilution	1:404
Analytical sensitivity	Ø 0.06 ng/mL
average Intra- / Inter-Assay Variance	Ø < 10%
Reference Values	Blum W.F et.al.1990 Insulin-Like Growth Factors and Their Binding Proteins. In: Ranke MB, Mullins P.E.(ed): Diagnostics of Endocrine Function in Children and Adolescents. Basel, Karger, 2011, pp.157-181

1 INTENDED USE

This enzyme immunoassay kit is suited for measuring IGF-II, human Insulin-like Growth Factor-II, in human serum and plasma for diagnostic and scientific purposes.

2 INTRODUCTION

The insulin-like growth factors (IGF)-I and -II play a pivotal role in the regulation of proliferation and differentiation of several tissue types (1-3). IGF-I also called Somatomedin C (4) has a molecular weight of 7.469 kDa (5). Its expression is mainly regulated by Growth Hormone and nutrition (6). But several hormones and peptide factors are known to influence IGF-II synthesis in different tissues. Bioavailability of the IGFs is regulated by specific binding proteins (IGFBP). Beside the high affinity Insulin-like Growth Factor Binding Proteins 1-6, IGFs are also bound by IGFBP-related Proteins (7, 8, 22). These binding proteins bind IGF-I and IGF-II with the same affinity or prefer IGF-II (9, 10). Direct measurement of IGFs in serum samples without pretreatment results in false values because of the extremely slow dissociation of the IGF/IGFBP complexes during the assay incubation only a part of the IGF-II in the specimen can bind to the antibodies and be detected.

Therefore, various techniques were applied to physically separate IGF-II from its binding proteins before measurement, including (a) size exclusion chromatography under acidic conditions, (b) solid-phase extraction and (c) acid-ethanol extraction (2, 12, 13) These

techniques, however, are either inconvenient or time-consuming or give incomplete and not-reproducible recoveries.

This assay is easy, fast and results do not depend on the binding protein concentration of the sample. It is based on the high specificity of the employed antibodies for IGF-II. There is virtually no cross-reactivity with IGF-I. This allows the separation of IGF-II from the binding proteins by acidification and blocking of the free binding proteins with IGF-I. Thus, the endogenous IGF-II is free in solution (15, 21, 23).

Indication

Scientific studies in the context of neonatal hypertrophy or hypotrophy (IGF-II is a fetal growth factor) and malignancy (IGF-II is an oncogenic growth factor). Age-related reference values are shown in Table 1. IGF-II appears to be suitable for differential diagnosis in various malignant diseases. For example, IGF-II can be used to differentiate between adrenocortical tumors and adenomas. In prostate cancer, too, tumor staging and the differentiation between carcinoma and hyperplasia can be improved by measuring IGF-II in serum (24, 25). Recent results from neurology show that the IGF system is also important in the development of Alzheimer's and Parkinson's disease (26).

3 ASSAY PRINCIPLE

The Mediagnost ELISA for IGF-II E30 is a so-called Sandwich-Assay. It utilizes two specific and high affinity antibodies for this protein. The first antibody, immobilized on the microtiter plate, and the added second biotinylated antibody are binding the IGF-II in the sample. The Streptavidin-Peroxidase Enzyme Conjugate subsequently binds to the complex. In the closing substrate reaction, the turn of the colour will be high specific catalysed, quantitatively depending on the IGF-II-level of the samples.

IGF-II-IGFBP complex is dissociated by dilution in an acidic buffer. IGFBPs are blocked by IGF-I excess, thus allowing the measurement of free IGF-II. With this method, the IGFBPs are not removed, but their function and therefore their interference in the assay is neutralized.

Due to the low cross-reactivity of the IGF-II antibody with IGF-I, excess IGF-I does not disturb the interaction of the first antibody with IGF-II.

4 WARNINGS AND PRECAUTIONS

For In Vitro Diagnostic Use only. For Professional use only.

The Mediagnost kit is suitable only for in vitro diagnostics and not for internal use in humans and animals. Follow strictly the test protocol. Use the valid version of the package insert provided with the kit. Be sure that everything has been understood. Mediagnost will not be held responsible for any loss or damage (except as required by statute) howsoever caused, arising out of noncompliance with the instructions provided.

Do not use obvious damaged or microbial contaminated or spilled material.

Caution: This kit contains material of human and/or animal origin. Therefore, all components and patient's specimens should be treated as potentially infectious.

Appropriate precautions and good laboratory practices must be used in the storage, handling and disposal of the kit reagents. The disposal of the kit components must be made according to the local regulations.

Human Serum

Following components contain human serum: **Controls CTR1 / CTR2**

Source human serum for the control sera provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV). No known methods can offer total security of the absence of infectious agents; therefore, all components and patient's specimens should be treated as potentially infectious.

Reagents DET, EC, SB, A-E, WB

Contain as preservative **5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one** (<0.015%)

H317	May cause an allergic skin reaction.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P501	Dispose of contents/ container in accordance with local/ regional/ national/ international regulations.

Substrate S

The TMB-Substrate (S) contains 3,3',5,5' Tetramethylbenzidine (<0.05%)

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapours/ spray.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.

Stop Solution STP

The Stop Solution contains 0.2 M acid sulphuric acid (H₂SO₄)

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
P280	Wear protective gloves/ protective clothing/ eye protection face protection.
P301+P330+	IF SWALLOWED: rinse mouth.
P331	Do NOT induce vomiting.
P305+P351+	IF IN EYES: Rinse cautiously with water for several minutes.
P338	Remove contact lenses, if present and easy to do. Continue rinsing.
P309+P310	If exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician.

4.1 General first aid procedures:

Skin contact: Wash affected area rinse immediately with plenty of water at least 15 minutes. Remove contaminated cloths and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water at least 15 minutes. In order to assure an effectual rinsing, spread the eyelids.

Ingestion: After swallowing the product, if the affected person is conscious, rinse out the mouth with plenty of water: seek medical advice immediately.

5 SAMPLES

5.1 Sample type

Serum and Plasma

Serum and Heparin/ EDTA/ Citrate Plasma yield comparable values.
Possible dilution of the sample by the anticoagulant must be considered.

5.2 Specimen collection

Use standard venipuncture for the blood sampling.
Haemolytic reactions have to be avoided.

5.3 Required sample volume: 10 µL

5.4 Sample stability

Samples should be handled as recommended in general: as fast as possible and chilled as soon as possible. In case there will be a longer period between the sample withdrawal and determination store the undiluted samples frozen -20°C or below in tightly closable plastic tubes. Avoid on principal repeated freeze-thaw cycles of serum/plasma (if required, please subaliquote) although IGF-II levels were found to be unaffected by few cycles in our experiments.

5.5 Interference

Triglyceride, bilirubin and hemoglobin in the sample do not interfere to a concentration of 100 mg/mL, 200 µg/mL or 1 mg/mL, respectively. However, the use of haemolytic, lipemic or icteric samples should be validated by the user.

5.6 Sample dilution

Suggestion for dilution protocol:

- Dilution: **1:404** with Sample Buffer **SB**.
- **1-step dilution:**
Please pipette **2015 µL** Sample Buffer **SB** in PE/PP-Tubes (application of a multi-stepper is recommended in larger series); subsequently **add 5 µL** sample (dilution 1:404). Incubate at least for 15 Minutes, max. 2 h, use 50 µL in Assay
- **2-step dilution:**
Because the pipetting accuracy can rise by the use of 10 µL sample, a 2-step dilution is alternatively possible, for this place **1000 µL** Sample Buffer **SB** in in PE/PP-Tubes, add **10 µL sample** (samples are 1:101 diluted), mix from this solution **50 µL** with **150 µL** Sample Buffer **SB** (samples are 1:404 diluted). Incubate at least for 15 Minutes, max. 2 h, use 50 µL in Assay.

6 MATERIALS

6.1 Materials provided

The reagents listed below are sufficient for 96 wells including the calibration curve.

MTP	Microtiter plate , ready for use, coated with mouse-anti-hIGF-II antibody. Wells are separately breakable.	(8x12) wells
CAL A-E	Calibrators , lyophilized, (recombinant human IGF-II), concentrations are given on vial labels and on the QC-certificate.	5 x 500 µL
CTR1	Control 1 , lyophilised, (human serum), concentration is given on the QC-certificate.	1 x 250 µL
CTR2	Control 2 , lyophilised, (human serum), concentration is given on the QC-certificate.	1 x 250 µL
DET	Antibody Conjugate , ready for use, contains coat anti-hIGF-II antibody.	1 x 6 mL
EC	Enzyme Conjugate , ready for use, contains horseradish-peroxidase conjugate to streptavidin,	1 x 12 mL
SB	Sample Buffer , ready for use, Please shake before use!	1 x 125 mL
WB	Washing Buffer , 20-fold concentrated solution	1 x 50 mL
S	Substrate , ready for use, horseradish-peroxidase-substrate, stabilised Tetramethylbencidine.	1 x 12 mL
STP	Stop Solution , ready for use, 0.2 M sulphuric acid.	1 x 12 mL
-	Sealing Tape , for covering the microtiter plate .	2 x
	Instructions for use	1 x
--	Quality Control Certificate	1 x

6.2 Material is required, but not provided

- Distilled (Aqua destillata) or deionized water for dilution of the Washing Buffer **WB (A. dest.)**, 950 mL.
- Precision pipettes and multichannel pipettes with disposable plastic tips
- Polyethylene PE/Polypropylene PP tubes for dilution of samples
- Vortex-mixer
- Microtiter plate shaker (350 rpm)
- Microtiter plate washer (recommended)
- Micro plate reader ("ELISA-Reader") with filter for 450 and ≥ 590 nm

7 TECHNICAL NOTES

Storage Conditions

Store the kit at 2-8°C after receipt until its expiry date. The lyophilized reagents should be stored at -20 °C after reconstitution. Avoid repeated thawing and freezing.

Storage Life

The shelf life of the components **after initial opening** is warranted for **4 weeks**, store the unused strips and microtiter wells **airtight** together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The **reconstituted components** Calibrators **A-E** and Controls **CTR1** and **CTR2** must be stored at -20°C (max. 4 weeks). For further use, thaw quickly but gently (avoid temperature increase above room temperature and avoid excessive vortexing). Up to 3 of the freeze-thaw cycles did not influence the assay. The 1:20 diluted Washing Buffer **WB** is 4 weeks stable at 2-8°C

Preparation of reagents

Bring all reagents to room temperature (20 - 25°C) before use. Possible precipitations in the buffers have to be resolved before usage by mixing and / or warming. Reagents with different lot numbers cannot be mixed.

Reconstitution

The Calibrators **A – E** and Controls **CTR1** and **CTR2** are reconstituted with the Sample Buffer **SB**. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a Vortex mixer.

Dilution

After reconstitution dilute the Controls **CTR1** and **CTR2** with the Sample Buffer **SB** in the same ratio (1:404) as the sample. The required volume of Washing Buffer **WB** is prepared by 1:20 dilution of the provided 20fold concentrate with Aqua dest.

Assay Procedure

When performing the assay, Blank, Calibrators **A-E**, Controls **CTR1** and **CTR2** and the samples should be pipette as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times, Antibody-Conjugate **DET**, **Enzyme Conjugate**, as well as the succeeding Substrate **S** should be added to the plate in the same order and in the same time interval as the samples. Stop Solution **STP** should be added to the plate in the same order as Substrate **S**.

All determinations (Blank, Calibrators **A-E**, Controls **CTR1** and **CTR2** and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

Incubation

Incubation at room temperature means: Incubation at 20 - 25°C. The Substrate **S**, stabilised Tetramethylbencidine, is photosensitive—store and incubation in the dark.

Shaking

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtiter plate shaker. We are recommending 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must be adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/ or false values, excessive shaking may result in high optical densities and/ or false values.

Washing

Proper washing is of basic **importance** for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results calculations of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer **WB** diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

The danger of handling with potentially infectious material must be taken into account.

When using an **automatic microtiter** plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

Manual washing is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamical swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

8 ASSAY PROCEDURE

Preparation of reagents		Reconstitution:	Dilution
CAL A-E	Calibrators	in 500 µL Sample Buffer SB	-
CTR1	Control 1	in 250 µL Sample Buffer SB	1:404 with Sample Buffer SB
CTR2	Control 2	in 250 µL Sample Buffer SB	1:404 with Sample Buffer SB
WB	Washing Buffer conc.	-	1:20 with Aqua dest. → WB 1:20

Sample **SPE and Control (CTR1, CTR2) dilution:** each with Sample Buffer **SB 1:404**, mix immediately, incubate at least for 15 Minutes, max. 2 h. Use **50 µL** for each well in the assay.

Before assay procedure bring all reagents to room temperature **20-25°C**.

Assay Procedure in Double Determination:		
Pipette	Reagents	Position
50 µL	Antibody Conjugate DET	Pipette in all required number of wells
50 µL	Sample Buffer SB as Blank	A1/A2
50 µL	Calibrator A (0.45 ng/mL)	B1/B2
50 µL	Calibrator B (1.5 ng/mL)	C1/C2
50 µL	Calibrator C (3 ng/mL)	D1/D2
50 µL	Calibrator D (5.63 ng/mL)	E1/E2
50 µL	Calibrator E (9 ng/mL)	F1/F2
50 µL	Control CTR1 (1:404 diluted)	G1/G2
50 µL	Control CTR2 (1:404 diluted)	H1/G2
50 µL	Sample SPE (1:404 diluted)	in the rest of the wells according the requirements
Cover the wells with the sealing tape.		
Sample Incubation: 2 h at 20-25°C, 350 rpm		
5 x 300 µL	Aspirate the contents of the wells and wash 5 x with 300 µL each Washing Buffer WB 1:20 / well	In each well
100 µL	Enzyme Conjugate EC	In each well
Cover the wells with the sealing tape.		
Incubation: 30 Minutes at 20-25°C, 350 rpm		
5 x 300 µL	Aspirate the contents of the wells and wash 5 x with 300 µL each Washing Buffer WB 1:20 /well	In each well
100 µL	Substrate S	In each well
Substrate S Incubation: 30 Minutes in the Dark at 20-25°C		
100 µL	Stop Solution STP	In each well
Measure the absorbance within 30 Minutes at 450 nm with ≥ 590 nm as reference wavelength.		

9 QUALITY CONTROL

Good laboratory practice requires that controls are included in each assay. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance. The test results are only valid if the test has been performed following the instructions. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state or local standards/laws. All Calibrators and kit controls must be found within the acceptable ranges as stated on the QC Certificate. If the criteria are not met, the run is not valid and should be repeated. Each laboratory should use known samples as further controls.

9.1 Quality criteria

For the evaluation of the assay it is required that the absorbance values of the blank should be below **0.25**, and the absorbance of Calibrator **E** should be above **1.00**. Samples, which yield higher absorbance values than Calibrator **E**, should be re-tested with a higher dilution.

10 EVALUATION OF RESULTS

10.1 Establishing of the standard curve

The Calibrators provided contain the following concentrations of recombinant human IGF-II

Calibrator	CAL A	CAL B	CAL C	CAL D	CAL E
ng/mL	0.45	1.5	3	5.63	9

- 1) Calculate the **mean absorbance** value for the blank from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the blank from the mean absorbances of all other samples, controls and calibrators.
- 3) Plot the calibrator concentrations on the x-axis versus the mean value of the absorbance of the Calibrators on the y-axis.
- 4) Recommendation: Calculation of the calibration curve should be done by using a computer program, because the curve is in general (without respective transformation) not ideally described by linear regression. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).
- 5) The **IGF-II** concentration in ng/mL of the samples can be calculated by **multiplication** with the respective **dilution factor**.

10.2 Example of a typical calibration curve

The following data is for demonstration only and cannot be used in place of data generation at the time of assay.

	Blank	CAL A	CAL B	CAL C	CAL D	CAL E
ng/mL	0	0.45	1.5	3	5.63	9
OD (450-620 nm)	0.0329	0.127	0.4195	0.8595	1.3695	1.855


Figure 1 Exemplary standard curve

The exemplary shown standard curve in Figure 1 **cannot** be used for calculation of your test results. You have to establish a standard curve for each test you conduct!

Exemplary calculation of IGF-II concentrations

Sample dilution: 1:404

Measured extinction of your sample	0.449
Measured extinction of the blank	0.0329

Your measurement program will calculate the IGF-II concentration of the diluted sample automatically by using the difference of sample and blank for the calculation. You only have to determine the most suitable curve fit (here: polynomial 2nd degree).

In this exemplary case the following equation is solved by the program to calculate the IGF-II concentration in the sample:

$$y = -0.0115x^2 + 0.3136x - 0.0433$$

$$0.4161 = -0.0115x^2 + 0.3136x - 0.0433$$

$$x = 1.49 \text{ ng/mL}$$

If the dilution factor (**1:404**) is taken into account the IGF-II concentration of the undiluted sample is

$$1.49 \text{ ng/mL} \times 404 = 601.96 \text{ ng/mL} = 0.602 \text{ mg/L}$$

11 INTERPRETATION OF RESULTS

The test results should not be the only base for therapeutic decisions. The results should be interpreted in regard to anamnesis, further clinical observations and results of other diagnostic investigations. Further, it is recommended to establish reference and cut-off values corresponding to the relevant group of patients for each laboratory.

12 LIMITATION OF PROCEDURE

The Mediagnost IGF-II ELISA, E30, is based on specific antibodies. Generally, this technique is sensible to heterophilic antibodies in the sample. The influence of heterophilic antibodies is reduced by assay design, but cannot be excluded completely.

13 REFERENCE VALUES

Table 1 Serum levels of IGF-II in ng/mL in healthy persons at various ages*

Percentile			
Age group	5th	50th	95th
Newborns	158	284	516
1-4 weeks	350	486	673
1-6 months	348	551	871
6-12 months	388	582	876
1-3 years	384	596	926
3-5 years	397	617	920
5-7 years	419	638	973
7-9 years	433	656	997
9-11 years	442	662	994
11-13 years	448	671	1006
13-15 years	455	679	1014
15-17 years	452	686	1042
20-30 years	436	679	1058
30-40 years	442	680	1049
40-50 years	407	650	1039
50-60 years	396	644	1049
60-70 years	373	611	1000

*Measurement was performed after acid-ethanol extraction, and values were corrected for recovery (correction factor 1.2). Blum W., Schweizer R.: Insulin-like growth factors and their binding proteins; in Ranke MB (ed): Diagnostics of endocrine function in children and adolescents. Basel, Karger, 2003, pp 166-199 (22).

14 PERFORMANCE CHARACTERISTICS

14.1 Sensitivity

Sensitivity was assessed by measuring the blank and calculating the theoretical concentration of the 2fold standard deviation of the blank. The mean analytical sensitivity of the E30 is 0.06 ng/mL (Range 0.01 to 0.223, n = 7)

14.2 Specificity

IGF-I, a homologous molecule to IGF-II, was diluted in three different concentrations (250, 750, 1250 ng/mL) and tested with the IGF-II ELISA. None of the samples gave a significant signal.

14.3 Reproducibility and Precision

The Inter- and Intra-Assay variation coefficients are in mean <10%. Exemplary determinations are represented in the Table 2 and 3.

Table 2 Inter-Assay- Variation within one lot measured within 10 months.

	Number [n]	Mean value [ng/mL]	Standard deviation [ng/mL]	Variation Coefficient [%]
Sample 1	18	1019	86.4	8.49
Sample 2	18	692	55.4	8.00
Sample 3	18	570	41.2	7.28

Table 3 Intra-Assay Variation

	Number [n]	Mean value [ng/mL]	Standard deviation [ng/mL]	Variation Coefficient [%]
Sample 1	16	666	20	3.07
Sample 2	16	875	58	6.61
Sample 3	20	621	26	4.17
Sample 4	19	670	41	6.13
Sample 5	20	795	42	5.29

14.4 Linearity

Three serum samples were measured in different dilutions. Standard assay dilution is 1:404. Here sample dilutions of 1:100 – 1:800 were used. Figure 2 demonstrates the linearity of sample dilution. These dilutions cover a concentration range of 4.5 to 0.6 ng/mL.


Figure 2 Linearity. Shown are the measured concentrations in different dilutions of three serum samples.

14.5 Interference

Interference of haemoglobin, bilirubin and triglycerides was tested by adding different amounts of these substances to human serum containing IGF-II. For comparison the same amount of buffer without any substance was also added to the serum. Table 4 demonstrates that there is no significant influence of these substances on the measurement of IGF-II in human serum.

Table 4 Interference. Serum samples were enriched with the respective amount of a potentially interfering substance. The recovery (%) in comparison to the non-enriched sample was calculated to assess the influence of haemoglobin, bilirubin and triglycerides on IGF-II measurement.

	Triglyceride 100 mg/mL	Bilirubin 100 µg/mL	Hemoglobin 1 mg/mL
Recovery	%	%	%
Sample 1	97	103	102
Sample 2	95	95	91
Sample 3	90	111	119

15 COMPARISON WITH OTHER ASSAYS

The test system was compared with an independent in-house assay, used for clinical routine by the University Children's Hospital Tübingen. A good coefficient of determination was measured ($R^2=0.8$). Passing Bablok regression analysis revealed a regression line described by $y = 18 + 1.3x$ (RSD 87). The systematic error is negligible and the proportional bias of 30% is probably resulting from usage of different antibodies and calibration material.


Figure 3 Assay Comparison with an in-house assay of the University Children's Hospital Tübingen (n = 18)

Example Version

Information for use for scientific application

16 SCIENTIFIC APPLICATION

IGF-II is present in low concentrations in various body fluids and in conditioned cell culture media of many cell lines.

16.1 Samples

Serum samples as well as Heparin-, EDTA- and Citrate-Plasma samples are suited. Possible dilution of the sample by the anticoagulant must be considered.

Furthermore, suitable samples are: urine, saliva (low concentration, at least 1:10 dilution), cerebrospinal fluid (at least 1:10 dilution) and cell culture medium (incl. 5% FCS, at least 1:5 dilution). IGF-II concentration in other body fluids or cell culture supernatants can deviate strongly from the serum values.

Table 5 Results of sample matrix tests 4 ng/mL IGF-II was added to the respectively diluted samples. Enriched samples were measured without further dilution. Shown is the relative recovery of added IGF-II of the value measured in enriched assay buffer.

	IGF-II in the native sample [ng/mL]	IGF-II in the enriched sample [ng/mL]	Recovery [%]
Cerebrospinal fluid			
undiluted	2.86	4.39	31.1
1:2	3.2	5.51	57.8
1:5	3.84	6.78	73.5
1:10	3.03	6.48	86.3
Urine			
undiluted	0.36	4.5	103.5
1:2	0.35	4.92	92.9
1:5	0.35	5.17	98.0
1:10	0.35	5.06	95.7
Saliva			
undiluted	0.38	3.06	67.0
1:2	0.35	3.99	74.0
1:5	0.24	4.1	78.5
1:10	0.34	4.32	80.9
Cell Culture Medium			
undiluted	0.38	3.06	67.0
1:2	0.35	3.99	91.0
1:5	0.24	4.1	96.5
1:10	0.34	4.32	99.5
Cell Culture Medium			
undiluted	0.4	3.62	80.5
1:2	0.4	3.53	78.3
1:5	0.7	5.55	121.3
1:10	0.5	5.74	131.0

16.2 Cross reactions with animal samples




It has been shown that the test can be used as a heterologous assay for IGF-II measurement in serum samples from bovine and pigs. Species-specific calibration must be carried out by the user.

17 LITERATUR / REFERENCES

- 1) Baxter RC. 1986 The somatomedins:insulin-like growth factors. *Adv Clin Chem.*25:49-115
- 2) Daughaday WH, Rotwein P. 1989 Insulin-like growth factors I and II. Peptide, messenger ribonucleic acid and gene structures, serum, and tissue concentrations. *Endocr Rev.* 10:68-91
- 3) Spencer EM (Ed.) 1991 *Modern Concepts of Insulin-Like Growth Factors.* New York: Elsevier.
- 4) Klapper DG, Svoboda ME, Van Wyk JJ. 1983 Sequence analysis of somatomedin-C: confirmation of identity with insulin-like growth factor-I. *Endocrinology* 112:2215-2217.
- 5) Rinderknecht E, Humbel RE. 1978 The amino acid sequence of human insulin-like growth factor I and its structural homology with proinsulin. *J Biol Chem.* 253:2769-2276.
- 6) Clemmons DR, Van Wyk JJ. 1984 Factors controlling blood concentration of somatomedin C. *Clin Endocrinol Metab.* 13:113-143.
- 7) Ballard J, Baxter R, Binoux M, et al. 1989 On the nomenclature of the IGF binding proteins. *Acta Endocrinol (Copenh).* 121:751-752.
- 8) Drop SLS. 1992 Report on the nomenclature of the IGF binding proteins. *J. Clin Endocrinol Metab.* 74:1215-1216.
- 9) Martin JL, Baxter RC. 1986 Insulin-like growth factor binding protein from human plasma. Purification and characterization. *J Biol Chem.* 261:8754-8760.
- 10) Binkert C, Landwehr J, Mary JL, Schwander J, Heinrich G. 1989 Cloning, sequence analysis and expression of a cDNA encoding a novel insulin-like growth factor binding protein (IGFBP-2). *EMBO J.* 8:2497-2502.
- 11) Furlanetto RW, Underwood LE, Van Wyk JJ, D'Ercole AJ. 1977 Estimation of somatomedin-C levels in normals and patients with pituitary disease by radioimmunoassay. *J. Clin Invest.* 60:648-657.
- 12) Daughaday WH, Kapadia M, Mariz I. 1987 Serum somatomedin binding proteins: physiologic significance and interference in radioligand assay. *J Lab Clin Med.* 109:355-363.
- 13) Breier BH, Gallaher BW, Gluckman PD. 1991 Radioimmunoassay for insulin-like growth factor-I: solutions to some potential problems and pitfalls. *J Endocrinol.* 128:347-357.
- 14) Daughaday WH, Mariz IK, Blethen SL. 1980 Inhibition of access of bound somatomedin to membrane receptor and immunobinding sites: a comparison of radioreceptor and radioimmunoassay of somatomedin in native and acid-ethanol-extracted serum. *J Clin Endocrinol Metab.* 51:781-788.
- 15) Blum WF, Gallaher B, Ranke MB. 1992 An IGFBP-blocked IGF-I RIA that measures what it pretends to measure: IGF-I. 74th Annual Meeting of the American Endocrine Society. 293.
- 16) Rosenfeld RG, Wilson DM, Lee PDK, Hirtz RL. 1986 Insulin-like growth factors I and II in evaluation of growth retardation. *J Pediatr.* 109:428-433.
- 17) Clemmons DR, Van-Wyk JJ, Ridgway EC, Kliman B, Kjellberg RN, Underwood LE. 1979 Evaluation of acromegaly by radioimmunoassay of somatomedin-C. *N.Engl J Med.* 301:1138-1142
- 18) Zapf J, Walter H, Froesch ER. 1981 Radioimmunological determination of insulin-like growth factors I and II in normal subjects and in patients with growth disorders and extrapancreatic tumor hypoglycemia. *J Clin Invest.* 68:1321-1330.
- 19) Blum WF. 1992 Insulin-like growth factors and their binding proteins. In: Ranke MB, ed. *Functional Endocrinologic Diagnostics in Children and Adolescence.* Mannheim: J + J Verlag; 102-117.
- 20) Rieu M, Girard F, Bricaire H, Binoux M. 1982 The importance of insulin-like growth factor (somatomedin) measurements in the diagnosis and surveillance of acromegaly. *J Clin Endocrinol Metab.* 55:147-153.
- 21) Blum WF, Ranke MB, Bierich JR. 1988 A specific radioimmunoassay for insulin-like growth factor II: the interference of IGF binding proteins can be blocked by excess IGF-I. *Acta Endocrinol (Copenh).*118:374-380.
- 22) Blum W., Schweizer R.: Insulin-like growth factors and their binding proteins; in Ranke MB (ed): *Diagnostics of endocrine function in children and adolescents.* Basel, Karger, 2003, pp 166-199.
- 24) Blum WH, Breier BH (1994) Radioimmunoassays for IGFs and IGFBPs. *Growth Regulation* 4 (Suppl. 1):11-19
- 25) Trojan L, Bode C, Weiss C, Mayer D, Grobholz R, Alken P, Michel MS: IGF-II serum levels increase discrimination between benign prostatic hyperplasia and prostate cancer and improve the predictive value of PSA in clinical staging. *Eur Urol* (2006), 49(2):286-292
- 26) Soon PS, Gill AJ, Benn DE, Clarkson A, Robinson BG, McDonald KC, Sidhu SB: Microarray gene expression and immunohistochemistry analyses of adrenocortical tumors identify IGF2 and Ki-67 as useful in differentiating carcinomas from adenomas. *Endocr Relat Cancer* (2009), 16(2):573-583
- 27) Freude S et al: Neuronal IGF-1 resistance reduces Abeta accumulation and protects against premature death in a model of Alzheimer's disease. *FASEB J* (2009) Jun1 Epub ahead of print.

18 INTERNATIONAL ASSAY DESCRIPTION

International Test Description

CAL A-E	Rec in 500 µL SB	-
CTR1	Rec in 250 µL SB	1:404 SB
CTR2	Rec in 250 µL SB	1:404 SB
WB 20x	-	1:20 A. dest. → WB 1:20
SPE	-	1:404 SB ↔
 20-25°C  15 min		
50 µL	DET	A1 - End
50 µL	SB	A1/A2
50 µL	CAL A (0.45 ng/mL)	B1/B2
50 µL	CAL B (1.5 ng/mL)	C1/C2
50 µL	CAL C (3 ng/mL)	D1/D2
50 µL	CAL D (5.63 ng/mL)	E1/E2
50 µL	CAL E (9 ng/mL)	F1/F2
50 µL	CTR1 1:404 SB	G1/G2
50 µL	CTR2 1:404 SB	H1/H2
50 µL	SPE 1:404 SB	
TAPE		
 2 h  20-25°C  350 rpm		
5x 300 µL	5x WB 1:20	
100 µL	EC	
TAPE		
 0.5 h  20-25°C  350 rpm		
5x 300 µL	5x WB 1:20	
100 µL	S	
 0.5 h  20-25°C 		
STP		
MEASURE		

19 ASSAY PROCEDURE

Preparation of reagents		Reconstitution:	Dilution
CAL A-E	Calibrators	in 500 µL Sample Buffer SB	-
CTR1	Control 1	in 250 µL Sample Buffer SB	1:404 with Sample Buffer SB
CTR2	Control 2	in 250 µL Sample Buffer SB	1:404 with Sample Buffer SB
WB	Washing Buffer conc.	-	1:20 with Aqua dest. → WB 1:20

Sample and Control dilution: each with Sample Buffer **SB 1:404**, mix immediately, incubate at least for 15 Minutes, max. 2 h. Use **50 µL** for each well in the assay.

Before assay procedure bring all reagents to room temperature **20-25°C**.

Assay Procedure in Double Determination:		
Pipette	Reagents	Position
50 µL	Antibody Conjugate DET	Pipette in all required number of wells
50 µL	Sample Buffer SB as Blank	A1/A2
50 µL	Calibrator A (0.45 ng/mL)	B1/B2
50 µL	Calibrator B (1.5 ng/mL)	C1/C2
50 µL	Calibrator C (3 ng/mL)	D1/D2
50 µL	Calibrator D (5.63 ng/mL)	E1/E2
50 µL	Calibrator E (9 ng/mL)	F1/F2
50 µL	Control CTR1 (1:404 diluted)	G1/G2
50 µL	Control CTR2 (1:404 diluted)	H1/G2
50 µL	Sample SPE (1:404 diluted)	in the rest of the wells according the requirements
Cover the wells with the sealing tape.		
Sample Incubation: 2 h at 20-25°C, 350 rpm		
5 x 300 µL	Aspirate the contents of the wells and wash 5 x with 300 µL each Washing Buffer WB 1:20 /well	In each well
100 µL	Enzyme Conjugate EC	In each well
Cover the wells with the sealing tape.		
Incubation: 30 Minutes at 20-25°C, 350 rpm		
5 x 300 µL	Aspirate the contents of the wells and wash 5 x with 300 µL each Washing Buffer WB 1:20 /well	In each well
100 µL	Substrate S	In each well
Substrate S Incubation: 30 Minutes in the Dark at 20-25°C		
100 µL	Stop Solution STP	In each well
	Measure the absorbance within 30 Minutes at 450 nm with ≥ 590 nm as reference wavelength.	

(Copy of Page 26)