ALPHA CrossLaps[®] (CTX-I) EIA

For the quantification of degradation products of C-terminal telopeptides of Type I collagen in human urine

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Immunodiagnostic Systems Limited is not responsible for any other use of the kit or consequence hereof than the one specified above. Neither for misuse, e.g., use deviating from the procedure described in this manual.

Furthermore, Immunodiagnostic Systems Limited is not to be made responsible for any diagnoses or conclusions made by the user or third party based on the results obtained with the ALPHA CrossLaps® (CTX-I) EIA kit nor for any consequences such interpretations may cause.

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INTRODUCTION

Intended use

The ALPHA CrossLaps® (CTX-I) EIA is an enzyme-linked immunosorbent assay for the quantification of non-isomerized fragments of C-terminal telopeptides of Type-I collagen (α CTX-I) in human urine.

The ALPHA CrossLaps® (CTX-I) EIA is intended for use as an indicator of degradation of non-isomerized bone collagen and may be used as an aid in the identification of skeletal metastases in patients with breast and prostate cancer.

In addition, the ALPHA CrossLaps® (CTX-I) EIA may be applied to humans as an aid in the diagnosis of Paget's disease and the monitoring of anti-resorptive therapy in patients with Paget's disease

Limitations

The use of the test has not been established for determination of the level of bone resorption.

Summary and explanation of the test

Type I collagen accounts for more than 90% of the organic matrix of bone and is synthesized primarily in bone. During renewal of the skeleton, Type I collagen is degraded, and small peptide fragments are excreted into the urine. These fragments can be measured by ALPHA CrossLaps® (CTX-I) EIA.

The measurements of the specific degradation products of Type I collagen (ALPHA CrossLaps) in human urine have been reported as useful assessment of bone resorption in Paget's disease (1) and for detection of bone metastases in prostate (2, 5) and breast cancer (3, 4, 5, 6).

Principle of the procedure

The ALPHA CrossLaps® (CTX-I) EIA is based on one highly specific monoclonal antibodies against the amino acid sequence of EKAHDGGR. In order to obtain a specific signal in the ALPHA CrossLaps® (CTX-I) EIA, two chains of EKAHDGGR must be cross linked.

Standards, control, or unknown urine samples are pipetted into the appropriate microtitre wells coated with streptavidin, followed by application of a mixture of biotinylated antibody and peroxidase-conjugated antibody. Then, a complex between ALPHA CrossLaps antigens, biotinylated antibody and peroxidase-conjugated antibody is generated, and this complex binds to the streptavidin surface via the biotinylated antibody. Following the one step incubation at 2-8°C, the wells are emptied and washed. A chromogenic substrate is added, and the colour reaction is stopped with sulfuric acid. Finally, the absorbance is measured.

PRECAUTIONS

The following precautions should be observed in the laboratory:

- Do not eat, drink, smoke or apply cosmetics where immunodiagnostic materials are being handled
- Do not pipette by mouth.
- Wear gloves when handling immunodiagnostic materials and wash hands thoroughly afterwards
- Cover working area with disposable absorbent paper

Warnings

For in vitro use only.

- All reagents and laboratory equipment should be handled and disposed of as if they were infectious.
- Do not use kit components beyond the expiry date and do not mix reagents from different lots.

Storage

Store the ALPHA CrossLaps® (CTX-I) EIA kit upon receipt at 2-8°C. Under these conditions the kit is stable up to the expiry date stated on the box.

MATERIALS

Specimen collection

It is recommended to use second morning void urine specimens from fasting individuals.

Urine samples are stable for 7 days at 4°C. For longer storage, the urine samples should be stored frozen (< -18°C).

Prior to use, urine specimens should be shaken and sedimentation allowed for a minimum of 30 minutes.

Materials supplied

Before opening the kit, read the section on **Precautions**. The kit contains reagents sufficient for 96 determinations.

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Microwell strips (12x8 wells) pre-coated with streptavidin. Supplied in a plastic frame.

Standard CAL 0

One vial (min. 12.0 mL) of ready-for-use solution with protein stabiliser, detergent and preservative.

Standards CAL 1 - 5

Five vials (min. 0.4 mL) of ready-for-use, standard in a buffered solution with protein stabiliser, detergent and preservative. The exact value of each Standard is printed on the QC Report.

Control CTRL 1 - 2

Two vials (min. 0.4 mL) of ready-for-use, di-peptide in a buffered solution with protein stabiliser, detergent and preservative. Please refer to enclosed QC Report for control range.

Biotinylated Antibody AB BIOTIN

One vial (min. 0.2 mL) of a concentrated solution of a biotinylated monoclonal murine antibody specific for degradation products of C-terminal telopeptides of Type I collagen. Prepared in a buffered solution with protein stabiliser, detergent and preservative.

Peroxidase Conjugated Antibody ENZYMCONE

One vial (min. 0.25 mL) of a concentrated solution of a peroxidase conjugated murine monoclonal antibody specific for degradation products of C-terminal telopeptides of Type I collagen. Prepared in a buffered solution with protein stabiliser, detergent and preservative.

Incubation Buffer BUF

One vial (min. 15.0 mL) of a ready-for-use buffered solution with protein stabiliser, detergent and preservative.

Substrate Solution SUBS TMB

One vial (min. 12.0 mL) of a ready-for-use tetramethylbenzidine (TMB) substrate in an acidic buffer. Please note that the chromogenic substrate might appear slightly blueish.

Stopping Solution H2O4

One vial (min. 12.0 mL) of ready-for-use 0.18 mol/L sulfuric acid.

Washing Buffer WASHBUF 50x

One vial (min. 20 mL) of a concentrated washing buffer with detergent and preservative.

Sealing tape

Adhesive film for covering wells during incubation.

Materials required — not supplied

- Containers for preparing the Antibody Solution and the Washing Solution
- Precision micropipettes to deliver 25-200 µL
- Distilled or deionized water
- Precision 8- or 12-channel multipipette to deliver 100 µL
- ELISA plate reader with 450 nm, and 650 nm as reference wavelength

ASSAY PROCEDURE

Prior to use, prepare and equilibrate all solutions to room temperature (18-22°C). Mix all reagents and samples before use (avoid foam). **Perform the assay at 18-22°C.**

Determine the number of strips needed for the assay. It is recommended to test all samples in duplicate. In addition, for each run a total of 16 wells are needed for the standards and controls of the kit. Place the appropriate number of strips in the plastic frame. Store unused immunostrips in the tightly closed foil bag with desiccant capsules.

1 Pre-dilution of urine samples

Dilute kit controls **CTRL 1 - 2** supplied with the kit and urine samples 1+7 in **Standard 0 CAL 0** prior to testing (e.g. 25+175 µL).

2 Preparation of the Antibody Solution:

ATTENTION: Prepare the **Antibody Solution immediately before use**. The **Antibody Solution** is prepared by mixing the Biotinylated Antibody **AB BIOTIN**, Peroxidase Conjugated Antibody **ENZYMCNJ** and Incubation Buffer **BUF** in the volumetric ratio 1+1+100 in an empty container. Mix carefully and avoid formation of foam. **Prepare a fresh solution before each run of the assay**.

3 One Step incubation

Pipette 25 µL of either **Standards CAL 0 - 5**, Control or unknown samples into appropriate wells followed by 100 µL of the **Antibody Solution**. Cover the immunostrips with sealing tape and incubate for **60±5 minutes at 2-8°C** without shaking.

4 Washing

Wash the immunostrips 5 times manually with 300 µL **Washing Buffer (WASHBUF)** diluted 1+50 in distilled or deionized water). Using an automated plate washer, follow the instructions of the manufacturer or the guidelines of the laboratory. Usually, 5 washing cycles are adequate. Make sure that the wells are **completely emptied** after each manual or automatic washing cycle.

5 Incubation with chromogenic substrate solution

Pipette 100 µL of the **Substrate Solution SUBS TMB** into each well and incubate for 15±2 minutes at room temperature (18-22°C) in the dark without shaking. Use sealing tape. Do not pipette directly from the vial containing TMB substrate but transfer the needed volume to a clean reservoir. Remaining substrate in the reservoir should be discarded and not returned to the vial with TMB.

6 Stopping of colour reaction

Pipette 100 µL of the **Stopping Solution H₂SO₄** into each well.

7 Measurement of absorbance

Measure the absorbance at 450 nm with 650 nm as reference within two hours.

Limitations of the procedure

If the absorbance of a (pre-diluted) sample is above 5.0 ng/mL, the sample should be additionally diluted in **Standard 0** and re-analysed.

QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

RESULTS

Calculation of results

It is recommended to use a **4-parameter** curve fit.

Example of results obtained:

Standards/ Controls/ Samples	ALPHA CrossLaps conc. (ng/mL)	A ₄₅₀₋₆₅₀ (Abs) Obs 1 / Obs 2	Mean A ₄₅₀₋₆₅₀ (Abs)	ALPHA CrossLaps Interpolated conc. (ng/mL)	ALPHA CrossLaps Conc. corrected for dilution x8 (ng/mL)
Standard 0	0.00	0.008 / 0.009	0.009		
Standard 1	0.50	0.031 / 0.029	0.030		
Standard 2	1.12	0.145 / 0.135	0.140		
Standard 3	3.22	0.929 / 0.873	0.901		
Standard 4	4.59	1.534 / 1.350	1.442		
Standard 5	5.92	1.975 / 1.888	1.932		
Control 1		0.655 / 0.535	0.595	2.49	19.9
Control 2		1.170 / 1.146	1.158	3.86	30.9
Sample I		0.569 / 0.600	0.585	2.03	16.2
Sample II		0.293 / 0.285	0.289	1.23	9.84
Sample III		0.159 / 0.155	0.157	0.66	5.26

Please note:

The data above are for illustration only and should not be used to calculate the results of any run.

Calculation of corrected ALPHA CrossLaps value

For each urine sample the ALPHA CrossLaps concentration (ng/mL) and the creatinine concentration (mM= mmol/L) should be determined using an enzymatic colometric method for clinical chemistry analyzers (e.g. CREA plus for Roche/Hitachi analyzers) or equivalent.

The following equation corrects the ALPHA CrossLaps concentration for variation in urine concentration:

$$\text{ALPHA CrossLaps (ng/mL)}$$

$$\text{Corr. ALPHA CrossLaps Value (\mu g/mmol)} =$$

$$\frac{\text{ALPHA CrossLaps (ng/mL)}}{\text{Creatinine (mmol/L)}}$$

Performance characteristics

Detection limit: 0.80 ng/mL ALPHA CrossLaps

This is the concentration corresponding to three standard deviations above the mean of 21 determinations of the blank (Standard 0) multiplied by the dilution factor 8.

Precision

The precision of the ALPHA CrossLaps® (CTX-I) EIA was evaluated for three urine samples. The results summarised in the table below (the results have not been corrected by dilutionfactor).

InterAssay Variation (n=10)

IntraAssay Variation (n=10)

Mean (ng/mL)	SD (ng/mL)	CV (%)
5,21	0,30	5,8
3,19	0,17	5,3
1,15	0,11	9,4

Mean (ng/mL)	SD (ng/mL)	CV (%)
5,21	0,07	1,4
3,19	0,07	2,3
1,15	0,02	2,0

Dilution/Linearity

The ALPHA CrossLaps® (CTX-I) EIA is linear in the range 0.10 ng/mL to 5.00 ng/mL (before correction of dilution) of ALPHA CrossLaps.

Urine samples with the concentration of 6.32-13.36 (corrected for dilution) were diluted in **Standard 0** and the concentration of ALPHA CrossLaps were determined with ALPHA CrossLaps® (CTX-I) EIA. The urine neat sample is set to 100%

Dilution	Sample 1			Sample 2			Sample 3		
Sample vs. diluent	Obs. (ng/mL)	Exp. (ng/mL)	RC (%)	Obs. (ng/mL)	Exp. (ng/mL)	RC (%)	Obs. (ng/mL)	Exp. (ng/mL)	RC (%)
Neat	0.79	0.79	100	1.67	1.67	100	1.16	1.16	100
1+1	0.38	0.39	97	1.00	0.83	120	0.59	0.58	102
1+3	0.19	0.20	97	0.49	0.42	117	0.28	0.29	98
1+8	0.10	0.10	104	0.21	0.21	100	0.14	0.14	96

Obs.: Observed, Exp.: Expected, RC: Recovery.

Interference:

To be determined.

Expected values

It is advisable for a laboratory to establish its own range of normal and pathological values. As an example, the mean values and 95% CI for various populations are given below. For further reading, please refer to the reference list. All samples were morning fasting samples from healthy individuals.

Populations	Number of subjects (n)	Age (Years)	Geometric Mean Values (µg/mmol)	95% CI
Pre-menopausal women	81	32-49	0.32	0.10 – 0.99
Post-menopausal women	320	43-75	0.62	0.17 – 2.26
Males	209	30-87	0.38	0.13 – 1.13

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REFERENCES

1. Alexandersen et al. Non-Isomerized C-Telopeptide fragments are highly sensitive markers for monitoring disease activity and treatmentefficacy in Padget´s disease of bone.
2. Cloos et al. Investigation of bone diseases using isomerised and racemised fragments of type I collagen. Calcif Tissue Int (2003);72:8-17.
3. Cloos et al. Breast cancer patients with bone metastases are characterised by increasing levels of nonisomerised type I collagen fragments. Breast Cancer Res (2003);5:R103-9.
4. Cloos et al. An immunoassay for measuring fragments of newly synthesized collagen type I produced during metastatic invasion of bone. Cli. Lab. (2004); 50:279-289.
5. Leeming DJ, Koizumi M, Byrjalsen I, Li B, Qvist P, Tanko LB: The relative use of eight collagenous and noncollagenous markers for diagnosis of skeletal metastases in breast, prostate, or lung cancer patients. Cancer Epidemiol Biomarkers Prev 15: 32-8, 2006
6. Leeming DJ, Delling G, Koizumi M, Henriksen K, Karsdal MA, Li B, Qvist P, Tanko LB, Byrjalsen I: Alpha CTX as a biomarker of skeletal invasion of breast cancer: immunolocalization and the load dependency of urinary excretion. Cancer Epidemiol Biomarkers Prev 15: 1392-5, 2006

	GB Use By DE Verwendbar bis ES Fecha de caducidad IT Utilizzare entro FR Utiliser jusque NL Houdbaar tot DK Holdbar til CZ Použitelné do SK Použiteľné do GR Ημερομνία λήξης PT Prazo de validade HU Felhasználható SE Använd före PL Użyć przed		GB Batch code DE Chargenbezeichnung ES Código de lote IT Codice del lotto FR Code du lot NL Lot nummer DK Lotnummer CZ Číslo šarže SK Číslo šarže GR Αριθμός Παρτίδας PT Código do lote HU Sarzszám SE Lot nummer PL Kod partii
	GB Catalogue number DE Bestellnummer ES Número de catálogo IT Numero di catalogo FR Référence du catalogue NL Catalogus nummer DK Katalognummer CZ Katalogové číslo SK Katalógové číslo GR Αριθμός καταλόγου PT Referência de catálogo HU Katalógusszám SE Katalognummer PL Numer katalogowy		GB Manufacturer DE Hersteller ES Fabricante IT Fabbricante FR Fabricant NL Fabrikant DK Producent CZ Výrobce SK Výrobca GR Κατασκευαστής PT Fabricante HU Gyártó SE Tillverkare PL Producent
	GB Contains sufficient for <n> tests DE Inhalt ausreichend für <n> Prüfungen ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi FR Contenu suffisant pour "n" tests NL Inhoud voldoende voor "n" testen DK Indeholder tilstrækkeligt til "n" test CZ Lze použít pro <n> testů SK Obsah postačuje na <n> stanovení GR Περιεχόμενο επαρκές για «n» εξετάσεις PT Conteúdo suficiente para "n" ensaios HU A doboz tartalma <n> vizsgálat elvégzéséhez elegendő SE Räcker till "n" antal tester PL Wystarczy na wykonanie <n> testów		GB In Vitro Diagnostic Medical Device DE In-Vitro-Diagnostikum ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro FR Dispositif médical de diagnostic in vitro NL Medisch hulpmiddel voor in-vitro diagnostiek DK Medicinsk udstyr til in vitro-diagnostik CZ In Vitro diagnostický zdravotnický prostředek SK Zdravotnícka pomôcka in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro HU In vitro diagnosztikum SE Medicintekniska produkter för in vitro diagnostik PL Wyrob do dignistyki In Vitro
	GB Temperature limitation DE Temperaturbegrenzung ES Límite de temperatura IT Limiti di temperatura FR Limites de température NL Temperatuurlimiet DK Temperaturbegrensning CZ Teplotní rozmezí od do SK Teplotné rozmedzie od do GR Τεμπερατοριού θερμοκρασίας PT Limites de temperatura HU Hőmérsékletartamány SE Temperaturbegränsning PL Przestrzegać zakresu temperatury		GB Consult Instructions for Use DE Gebrauchsanweisung beachten ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso FR Consulter les instructions d'utilisation NL Raadpleeg de gebruiksaanwijzing DK Se brugsanvisning CZ Viz návod k použití SK Viď návod na pužtie GR Συμβουλεύετε τις οδηγίες χρήσης PT Consulte as instruções de utilização HU Nézze meg a Használati utasítást SE Se handhavandebeskrivningen PL Sprawdź w instrukcji obsługi

	GB EU Importer DE EU-Importeur DK EU Importør PT Importador UE SE EU Importör

Immunodiagnostic Systems Limited,
10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK
Tel.: +44 191 519 0660
Fax: +44 191 519 0760
E-Mail: info.uk@idspcl.com • www.idspcl.com

IDS Immunodiagnostic Systems GmbH
Herriotstrasse 1
60528 Frankfurt am Main
Germany