Urine BETA CrossLaps[®] (CTX-I)ELISA

For the quantification of degradation products of C-terminal telopeptides of Type I collagen in human urine

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Immunodiagnostic Systems Limited is not responsible for any other use of the kit or consequence hereof than the one specified above. Neither for misuse, e.g. use deviating from the procedure described in this manual. Furthermore, Immunodiagnostic Systems Limited is not to be made responsible for any diagnoses or conclusions made by the user or third party based on the results obtained with the Urine BETA CrossLaps® (CTX-I)ELISA kit nor for any consequences such interpretations may cause.

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INTRODUCTION

Intended use

The Urine BETA CrossLaps® (CTX-I)ELISA is an enzyme immunological test for the quantification of degradation products of C-terminal telopeptides of Type-I collagen in human urine.

The Urine BETA CrossLaps® (CTX-I)ELISA assay is intended for *in vitro* diagnostic use as an indication of human bone resorption as an aid in

A. Monitoring bone resorption changes of

- 1) Anti-resorptive therapies in postmenopausal women:
 - a) Hormone Replacement Therapies (HRT) with hormones and hormone like drugs
 - b) Bisphosphonate therapies

B. Predicting skeletal response (Bone Mineral Density) in postmenopausal women undergoing anti-resorptive therapies

- a) Hormone Replacement Therapies (HRT) with hormones and hormone like drugs
- b) Bisphosphonate therapies

Limitations

The use of the test has not been established to predict the development of osteoporosis or future fracture risk.

Summary and explanation of the test

Type I collagen accounts for more than 90% of the organic matrix of bone and is synthesized primarily in bone (1). During renewal of the skeleton, type I collagen is degraded, and small peptide fragments are excreted into the urine. These fragments can be measured by Urine BETA CrossLaps® (CTX-I)ELISA. Urine BETA CrossLaps® (CTX-I)ELISA demonstrate a high correlation to corresponding measurement of serum samples in Serum CrossLaps® (CTX-I)ELISA. Serum CrossLaps® (CTX-I)ELISA has been reported to be useful for follow-up of anti-resorptive treatment of patients with metabolic bone diseases (3 17).

Principle of the procedure

The Urine BETA CrossLaps® (CTX-I)ELISA is based on two highly specific monoclonal antibodies against the amino acid sequence of EKAHD- β -GGR, where the aspartic acid residue (D) is β -isomerized. In order to obtain a specific signal in the Urine BETA CrossLaps® (CTX-I)ELISA, two chains of EKAHD- β -GGR must be cross-linked.

Standards, controls, or unknown urine samples are pipetted into the appropriate microtitre wells coated with streptavidin, followed by application of a mixture of a biotinylated antibody and a peroxidase-conjugated antibody. Then, a complex between CrossLaps antigens, biotinylated antibody and peroxidase-conjugated antibody is generated, and this complex binds to the streptavidin surface via the biotinylated antibody. Following the one-step incubation at room temperature, the wells are emptied and washed. A chromogenic substrate is added and the colour reaction is stopped with sulfuric acid. Finally, the absorbance is measured.

PRECAUTIONS

The following precautions should be observed in the laboratory:

- Do not eat, drink, or smoke where immunodiagnostic materials are being handled
- Do not pipette by mouth.
- Wear gloves when handling immunodiagnostic materials and wash hands thoroughly afterwards
- Cover working area with disposable absorbent paper

Warnings

For *in vitro* use only.

- All reagents and laboratory equipment should be handled and disposed of as if they were infectious.
- Do not use kit components beyond the expiry date and do not mix reagents from different lots.

Storage

Store the Urine BETA CrossLaps® (CTX-I)ELISA kit upon receipt at 2-8°C. Under these conditions the kit is stable up to the expiry date stated on the box.

MATERIALS

Specimen collection

As bone resorption has a marked circadian variation primarily regulated by food intake (18), it is recommended to use second morning void urines obtained after an overnight fast.

Keep the urine sample refrigerated (2-8°C) for storage less than one week, or freeze the sample

(<-18°C) for longer storage. Please note that for optimal results it is recommended to mix the samples and then centrifuge the samples (e.g. 2000g; 10 min) before testing.

Materials supplied

Before opening the kit, read the section on **Precautions**. The kit contains reagents sufficient for 96 determinations.

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Microwell strips (12x8 wells) pre-coated with streptavidin. Supplied in a plastic frame.

CrossLaps Standard **CAL - 0**

One vial (min. 10.0 mL/vial) of ready-for-use PBS buffered solution with protein stabiliser and preservative.

CrossLaps Standards **CAL 1 - 5**

Five vials (min. 0.4 mL/vial) of ready-for-use, CrossLaps standard in a PBS-buffered solution with protein stabiliser and preservative. The exact value of each Standard is printed on the QC Report.

Controls **CTRL 1 - 2**

Two vials (min. 0.4mL/vial) CrossLaps peptide in a PBS-buffered solution with protein stabiliser and preservative. Please refer to enclosed QC Report for control range.

Biotinylated Antibody **Ab BIOTIN**

One vial (min. 0.25 mL) of a concentrated solution of a biotinylated monoclonal murine antibody specific for degradation products of C-terminal telopeptides of type I collagen. Prepared in a buffered solution with protein stabiliser and preservative.

Peroxidase Conjugated Antibody **ENZYMCONJ**

One vial (min. 0.25 mL) of a concentrated solution of a peroxidase conjugated murine monoclonal antibody specific for degradation products of C-terminal telopeptides of type I collagen. Prepared in a buffered solution with protein stabiliser and preservative.

Incubation Buffer **BUF**

One vial (min. 19 mL) of a ready-for-use buffered solution with protein stabiliser, detergent, and preservative.

Substrate Solution **SUBS TMB**

One vial (min. 12 mL) of a ready-for-use tetramethylbenzidine (TMB) substrate in an acidic buffer. Please note that the chromogenic substrate might appear slightly bluish.

Stopping Solution **H2S04**

One vial (min. 12 mL) of ready-for-use 0.18 mol/L sulfuric acid.

Washing Buffer **WASHBUF 50x**

One vial (min. 20 mL) of a concentrated washing buffer with detergent and preservative.

Sealing tape

Adhesive film for covering wells during incubation.

Materials required — not supplied

- Containers for preparing the Antibody Solution and the Washing Solution
- Precision micropipettes to deliver 10-200 µL
- Distilled water
- Precision 8- or 12-channel multipipette to deliver 100 µL and 150 µL
- Microwell mixing apparatus
- Microtiter plate reader

ASSAY PROCEDURE

Assay Procedure

For optimal performance of the assay, it is important to comply with the instructions given below. Equilibrate all reagents to room temperature (18-22°C) prior to use. Determine the number of strips needed for the assay. It is recommended to test all samples in duplicate. In addition, for each run a total of 16 wells are needed for standards and controls. Place the appropriate number of strips in the plastic frame. Store unused immunostrips in the tightly closed foil bag with desiccant capsules. Mix all reagents and samples before use (avoid foam).

1 Preparation of the Antibody Solution:

ATTENTION: Prepare the **Antibody Solution** maximum 30 minutes before starting the assay. Mix the Biotinylated Antibody **Ab BIOTIN** Peroxidase Conjugated Antibody **ENZYMCONJ** and Incubation Buffer **BUF** in the volumetric ratio 1+1+100 in an empty container. Mix carefully and avoid formation of foam. **Prepare a fresh solution before each run of the assay.**

2 Pre-dilution of test specimens

All specimens (unknown samples and **CTRL 1 - 2**), except the standards delivered with the kit must be pre-diluted 1+3 in standard 0 prior to testing (e.g. 20 µL+60 µL).

3 One Step incubation

Pipette 10 µL of un-diluted **Standards CAL 0 -5**, pre-diluted **Controls**, and pre-diluted unknown samples into appropriate wells followed by 150 µL of the **Antibody Solution**. Cover the immunostrips with sealing tape and incubate for 120±5 minutes at room temperature (18-22°C) on a microtitre plate mixing apparatus (300 rpm).

4 Washing

Wash the immunostrips 5 times manually with 300 µL **Washing Buffer (WASHBUF 50x)** diluted 1+50 in distilled water). Using an automated plate washer, follow the instructions of the manufacturer or the guidelines of the laboratory. Usually 5 washing cycles are adequate. Make sure that the wells are completely emptied after each manual or automatic washing cycle.

5 Incubation with chromogenic substrate solution

Pipette 100 µL of the Substrate Solution **SUBS TMB** into each well and incubate for 15±2 minutes at room temperature (18-22°C) in the dark on the mixing apparatus (300 rpm). Use sealing tape.

Do not pipette directly from the vial containing TMB substrate but transfer the needed volume to a clean reservoir. Remaining substrate in the reservoir should be discarded and not returned to vial.

6 Stopping of colour reaction

Pipette 100 µL of the **Stopping Solution** into each well.

7 Measurement of absorbance

Measure the absorbance at 450 nm with 650 nm as reference within two hours.

Limitations of the procedure

If the absorbance of a sample exceeds that of **Standard 5**, the sample should be further diluted in **Standard 0** and re-analysed.

QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

RESULTS

Calculation of results

A quadratic curve fit should be used.

Alternatively, calculate the mean of the duplicate absorbance determinations. Construct a standard curve on graph paper by plotting the mean absorbances of the six standards (ordinate) against the corresponding CrossLaps concentrations (abscissa). Determine the CrossLaps concentration of the controls and each patient sample by interpolation.

Example of results obtained

Standards/ Controls/ Samples	CrossLaps conc. (µg/L)	A ₄₅₀₋₆₅₀ (nm) Obs 1/ Obs 2	Mean A ₄₅₀₋₆₅₀ (nm)	Interpolated CrossLaps conc. (µg/L)	Conc. Corrected for 4X dilution (µg/L)
Standard 0	0.00	0.046 / 0.045	0.046		
Standard 1	1.43	0.164 / 0.160	0.162		
Standard 2	4.24	0.402 / 0.383	0.393		
Standard 3	8.27	0.742 / 0.700	0.721		
Standard 4	16.6	1.461 / 1.366	1.414		
Standard 5	23.7	2.027 / 1.923	1.975		
Control 1		0.327 / 0.319	0.323	3.38	13.5
Control 2		1.216 / 1.173	1.195	14.0	55.9
Sample I		0.162 / 0.150	0.156	1.36	5.44
Sample II		0.507 / 0.500	0.504	5.56	22.2
Sample III		1.173/ 1.082	1.128	13.2	52.6

Please note:

The data above are for illustration only and should not be used to calculate the results of any run.

Calculation of corrected CrossLaps value

For each sample the CrossLaps concentration (ng/mL) and the creatinine concentration should be determined. Determine the concentration of creatinine (mmol/L) in the sample using a routine enzymatic colorimetric method for clinical chemistry analysers, and perform the correction using the equation:

$$\text{Corr. CrossLaps Value } (\mu\text{g}/\text{mmol}) = \frac{\text{CrossLaps (ng/mL)}}{\text{Creatinine (mM)}}$$

Performance characteristics

Detection limit: 0.80 µg/L CrossLaps

This is the concentration corresponding to three standard deviations above the mean of 21 determinations of the blank ("CrossLaps Standard 0") multiplied by the dilution factor 4.

Precision ≤ 6.9%

The precision of the Urine BETA CrossLaps® (CTX-I)ELISA was evaluated for three urine samples. The results are summarised in the table below (the results have not been corrected by dilutionfactor).

Intraassay ≤ 3.9% (n=10)			Interassay ≤ 6.9% (n=10)		
Mean (µg/L)	SD (µg/L)	CV (%)	Mean (µg/L)	SD (µg/L)	CV (%)
5.38	0.05	3.9	5.38	0.05	3.6
21.6	0.15	2.9	21.6	0.37	6.9
56.0	0.42	3.0	56.0	0.66	4.7

Dilution/Linearity

The Urine BETA CrossLaps® (CTX-I)ELISA is linear in the range 0.20 ng/mL to 40.7 ng/mL of CrossLaps.

Urine samples with the concentration of 16.7 – 52.8 ng/mL CrossLaps were diluted with standard 0 and the

concentration of CrossLaps were determined with Urine BETA CrossLaps® (CTX-I)ELISA.

The data below is calculated from 3 different urine samples

Dilution Procedure								
Urine [%]	Std. 0 [%]	Sample 1		Sample 2		Sample 3		Overall recovery (%)
		Exp. (µg/L)	RC (%)	Exp. (µg/L)	RC (%)	Exp. (µg/L)	RC (%)	
100	0	16,68	100	23,36	100	52,80	100	100
90	10	15,00	103	21,00	105	47,52	110	106
80	20	13,32	100	18,68	104	42,24	109	104
70	30	11,68	95	16,36	99	36,96	110	101
60	40	10,00	96	14,00	98	31,68	104	99
50	50	8,32	97	11,68	88	26,40	99	95
40	60	6,68	96	9,32	87	21,12	94	93
30	70	5,00	96	7,00	88	15,84	84	90
20	80	3,32	96	4,68	86	10,56	84	89
10	90	1,68	87	2,32	90	5,28	76	85
Mean								96

RC: Recovery, Exp: Expected

Expected values

It is advisable for a laboratory to establish its own range of normal and pathological values.

Populations	Number of subjects (n)	Age (Years)	Geometric Mean Values (µg/mmol)	95% CI
Pre-menopausal women	25	40-49	1.66	0.83-3.32
Post-menopausal women	272	52-75	2.27	0.73-7.07
Males	25	50-78	1.46	0.53-4.04


Day to Day Individual Variation

The Day to Day Intra-individual Variation was assessed by analyzing urine samples (morning fasting) from 11 healthy post menopausal women at five time points over 2 weeks

Subject No	Urine BETA CrossLaps® (CTX-I) ELISA (µg/mmol)					Mean (µg/mmol)	SD (µg/mmol)	CV (%)
	Visit 1	Visit 2	Visit 3	Visit 4	Visit 5			
1	3.04	3.17	3.88	3.50	4.06	3.53	0.44	12
2	0.96	1.21	0.99	1.17	1.17	1.04	0.14	13
3	1.15	1.04	1.56	1.78	1.78	1.34	0.31	23
4	1.33	1.51	1.65	0.99	0.99	1.42	0.27	19
5	1.35	1.60	1.54	1.37	1.37	1.43	0.14	10
6	2.40	1.91	1.88	1.90	1.90	2.11	0.29	14
7	1.29	1.04	1.18	1.32	1.32	1.21	0.11	9
8	1.18	1.21	1.20	1.07	1.07	1.22	0.14	11
9	5.77	5.67	4.92	4.78	4.78	5.59	0.81	15
10	1.42	1.40	1.28	1.22	1.22	1.26	0.18	14
11	4.33	4.47	5.61	5.14	5.14	5.07	0.64	13

CLINICAL DATA

Correlation between serum samples measured in Serum CrossLaps® (CTX-I)ELISA and corresponding urine samples measured in Urine BETA CrossLaps® (CTX-I)ELISA (n=159). The 159 samples are from pre- and postmenopausal women before and after anti-resorptive therapy, as well as males. The women were treated with Ibandronate, Raloxifene, Drospirone, or Aerodin

Anti-resorptive therapies

Below is the Urine BETA CrossLaps® (CTX-I)ELISA data [$\mu\text{g}/\text{mmol}$ creatinine] from two anti-resorptive studies. The geometric mean value and the 95% confidence interval of mean are shown.

Bisphosphonate

- Postmenopausal women between age 55 and 75 years, more than 5 years since menopause, and a BMD of spine or femoral neck at osteoporotic level
- 20 participants received placebo (500 calcium and 400 mg vitamin D daily)
- 20 participants received active treatment (2.5 mg Ibandronate, 500 mg calcium, and 400 mg vitamin D daily)
- Treatment period 1 year

	Placebo group Mean [95% CI] ($\mu\text{g}/\text{mmol}$)	Bisphosphonate group Mean [95% CI] ($\mu\text{g}/\text{mmol}$)
Baseline	2.69 [2.11;3.42]	2.17 [1.49;3.14]
After 6 months of treatment	1.97 [1.36;2.85]	0.67 [0.43;1.04]
After 12 months of treatment	1.92 [1.39;2.65]	0.53 [0.29;0.95]

HRT

- Postmenopausal women between age 45 and 65 years, more than 2 years since menopause, and a BMD of hip within the reference range of premenopausal women
- 57 participants received placebo (500 mg calcium)
- 175 participants received active treatment (1 mg of estradiol continuously combined with 1 mg of drospirenone (n=59), 2 mg of drospirenone (n=58), or 3 mg of drospirenone (n=58))
- Treatment period 2 years

	Placebo group Mean [95% CI] ($\mu\text{g}/\text{mmol}$)	HRT group Mean [95% CI] ($\mu\text{g}/\text{mmol}$)
Baseline	2.11 [1.84;2.42]	2.29 [2.11;2.49]
After 6 months of treatment	1.86 [1.46;2.36]	0.46 [0.40;0.53]
After 12 months of treatment	1.70 [1.39;2.09]	0.44 [0.38;0.51]
After 24 months of treatment	2.09 [1.75;2.51]	0.66 [0.57;0.76]

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REFERENCES

1. Rosenquist C. et al., Serum CrossLaps® ELISA. First application of monoclonal antibodies for measurement in serum of bone-related degradation products from C-terminal telopeptides of type I collagen. *Clin. Chem* (1998); 44; 2281-2289.
2. Christgau S. et al., Clinical evaluation of the Serum CrossLaps® ELISA, a new assay measuring the serum concentration of bone-derived degradation products of type I collagen C telopeptides. *Clin Chem* (1998); 44; 2290-2300.
3. Bjarnason NH. et al., Early response in biochemical markers predicts long-term response in bone mass during hormone replacement therapy in early postmenopausal women. *Bone* (2000); 26; 561-569.
4. Buclin T. et al., Bioavailability and biological efficacy of a new oral formulation of salmon calcitonin in healthy volunteers. *J Bone Miner Res* (2002); 17: 1478-1485.
5. Ravn P. et al., Association between pharmacokinetics of oral ibandronate and clinical response in bone mass and bone turnover in women with postmenopausal osteoporosis. *Bone* (2002); 30: 320-324.
6. Ravn P, Clemmesen B, Christiansen C. Biochemical markers can predict the response in bone mass during alendronate treatment in early postmenopausal women. *Bone* (1999); 24: 237-244.
7. Rosen HN. et al., Serum CTX: A new marker of bone resorption that shows treatment effect more often than other markers because of low coefficient of variability and large changes with bisphosphonate therapy. *Calcif Tissue Int* (2000); 66: 100-103.
8. Tanko LB. et al., Safety and efficacy of a novel salmon Calcitonin (sCT) technology-based oral formulation in healthy postmenopausal women: acute and 3-month effects on biomarkers of bone turnover. *J Bone Miner Res* (2004); 19:1531-1538.
9. Warming L. et al., Safety and efficacy of drospirenone used in a continuous combination with 17beta-estradiol for prevention of postmenopausal osteoporosis. *Climacteric* (2004); 7: 103-111.

Doc: AC-05PL-A

Issue: 08

Date: 11 October 2022

Example Version

	<p>GB <i>Use By</i> DE <i>Verwendbar bis</i> ES <i>Fecha de caducidad</i> IT <i>Utilizzare entro</i> FR <i>Utiliser jusque</i> NL <i>Houdbaar tot</i> DK <i>Holdbar til</i> CZ <i>Použitelné do</i> SK <i>Použitelné do</i> GR <i>Ημερομηνία λήξης</i> PT <i>Prazo de validade</i> HU <i>Felhasználható</i> SE <i>Använd före</i> PL <i>Użyć przed</i></p>	<p>LOT</p>	<p>GB <i>Batch code</i> DE <i>Chargenbezeichnung</i> ES <i>Código de lote</i> IT <i>Codice del lotto</i> FR <i>Code du lot</i> NL <i>Lot nummer</i> DK <i>Lotnummer</i> CZ <i>Číslo šarže</i> SK <i>Číslo šarže</i> GR <i>Αριθμός Παρτίδας</i> PT <i>Código do lote</i> HU <i>Sarzszzám</i> SE <i>Lot nummer</i> PL <i>Kod partii</i></p>
<p>REF</p>	<p>GB <i>Catalogue number</i> DE <i>Bestellnummer</i> ES <i>Número de catálogo</i> IT <i>Numero di catalogo</i> FR <i>Référence du catalogue</i> NL <i>Catalogus nummer</i> DK <i>Katalognummer</i> CZ <i>Katalogové číslo</i> SK <i>Katalógové číslo</i> GR <i>Αριθμός καταλόγου</i> PT <i>Referência de catálogo</i> HU <i>Katalógusszám</i> SE <i>Katalognummer</i> PL <i>Numer katalogowy</i></p>		<p>GB <i>Manufacturer</i> DE <i>Hersteller</i> ES <i>Fabricante</i> IT <i>Fabbricante</i> FR <i>Fabricant</i> NL <i>Fabrikant</i> DK <i>Producent</i> CZ <i>Výrobce</i> SK <i>Výrobca</i> GR <i>Κατασκευαστής</i> PT <i>Fabricante</i> HU <i>Gyártó</i> SE <i>Tillverkare</i> PL <i>Producent</i></p>
	<p>GB <i>Contains sufficient for <n> tests</i> DE <i>Inhalt ausreichend für <n> Prüfungen</i> ES <i>Contenido suficiente para <n> ensayos</i> IT <i>Contenuto sufficiente per "n" saggi</i> FR <i>Contenu suffisant pour "n" tests</i> NL <i>Inhoud voldoende voor "n" testen</i> DK <i>Indeholder tilstrækkeligt til "n" test</i> CZ <i>Lze použít pro <n> testů</i> SK <i>Obsah postačuje na <n> stanovení</i> GR <i>Περιεχόμενο επαρκές για «n» εξετάσεις</i> PT <i>Conteúdo suficiente para "n" ensaios</i> HU <i>A doboz tartalma <n> vizsgálat elvégzéséhez elegendő</i> SE <i>Räcker till "n" antal tester</i> PL <i>Wystarczy na wykonanie <n> testów</i></p>	<p>IVD</p>	<p>GB <i>In Vitro Diagnostic Medical Device</i> DE <i>In-Vitro-Diagnostikum</i> ES <i>Producto sanitario para diagnóstico in vitro</i> IT <i>Dispositivo medico-diagnostico in vitro</i> FR <i>Dispositif médical de diagnostic in vitro</i> NL <i>Medisch hulpmiddel voor in-vitro diagnostiek</i> DK <i>Medicinsk udstyr til in vitro-diagnostik</i> CZ <i>In Vitro diagnostický zdravotnický prostředek</i> SK <i>Zdravotnícka pomocka in vitro</i> GR <i>Ιn Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν</i> PT <i>Dispositivo médico para diagnóstico in vitro</i> HU <i>In vitro diagnosztikum</i> SE <i>Medicintekniska produkter för in vitro diagnostik</i> PL <i>Wyrób do diagnostyki In Vitro</i></p>
	<p>GB <i>Temperature limitation</i> DE <i>Temperaturbegrenzung</i> ES <i>Límite de temperatura</i> IT <i>Limiti di temperatura</i> FR <i>Limites de température</i> NL <i>Temperatuurlimiet</i> DK <i>Temperaturbegrænsning</i> CZ <i>Teplotní rozmezí od do</i> SK <i>Teplotné rozmedzie od do</i> GR <i>Περιορισμοί θερμοκρασίας</i> PT <i>Limites de temperatura</i> HU <i>Hőmérséklettartomány</i> SE <i>Temperaturbegränsning</i> PL <i>Przestrzegać zakresu temperatury</i></p>		<p>GB <i>Consult Instructions for Use</i> DE <i>Gebrauchsanweisung beachten</i> ES <i>Consulte las instrucciones de uso</i> IT <i>Consultare le istruzioni per l'uso</i> FR <i>Consulter les instructions d'utilisation</i> NL <i>Raadpleeg de gebruiksaanwijzing</i> DK <i>Se brugsanvisning</i> CZ <i>Viz návod k použití</i> SK <i>Vid' návod na použitie</i> GR <i>Συμβουλευτείτε τις οδηγίες χρήσης</i> PT <i>Consulte as instruções de utilização</i> HU <i>Nézze meg a Használati utasítást</i> SE <i>Se handhavandebeskrivningen</i> PL <i>Sprawdź w instrukcji obsługi</i></p>

	<p>GB <i>EU Importer</i> DE <i>EU-Importeur</i> IT <i>Importatore UE</i> DK <i>EU Importør</i> CZ <i>EU Dovoze</i> PT <i>Importador UE</i> SE <i>EU Importör</i></p>
---	--

Immunodiagnostic Systems Limited,
10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK
Tel.: +44 191 519 0660
Fax: +44 191 519 0760
E-Mail: info.uk@idsplc.com
www.idsplc.com

IDS Immunodiagnostic Systems GmbH
Herriotstrasse 1
60528 Frankfurt am Main
Germany