Urine CartiLaps® (CTX-II) EIA

Instructions For Use

Product Name	Urine CartiLaps® (CTX-II) EIA	REF	AC-10F1
Abbreviated Product Name	Urine CTX-II EIA		

1. Intended Use

For In Vitro Diagnostic Use.

The Urine CartiLaps® (CTX-II) EIA is an enzyme immunological test for the quantification of degradation products of C-terminal telopeptides of Type II collagen in human urine. Results are to be used in conjunction with other clinical and laboratory data for the indication of degradation of cartilage and may be used as an aid for:

- quantitative assessment of disease activity (structural damage of articular cartilage) in patients with RA and OA
- prognosis of disease activity in patients with RA and OA, and
- early assessment of long-term effect of therapy in patients with RA

2. Summary and Explanation

Disruption of the structural integrity of cartilage is the major histological finding in osteoarthritis and rheumatoid arthritis. Type II collagen is the major organic constituent of cartilage and fragments of type II collagen (CTX-II) are being released into circulation and subsequently secreted into urine following degradation of cartilage. In urine, the CTX-II fragments can be quantified by Urine CartiLaps® (CTX-II) EIA.

The Urine CartiLaps® (CTX-II) EIA has been reported to be useful in prediction of progression of osteoarthritis¹⁻² and in other clinical and pre-clinical investigations³⁻⁹.

3. Method Description

Urine CartiLaps® (CTX-II) EIA is based on the competitive binding of a monoclonal antibody to urinary fragments of type II collagen⁴ or to biotinylated, synthetic peptides bound to the surface of microtitre plates coated with streptavidin.

Initially, biotinylated, synthetic peptides are bound to the surface of streptavidin-coated wells of the microtitre plate. After washing, calibrators, controls, and urine samples are pipetted into the wells followed by addition of the monoclonal antibody. The wells are washed, and a solution of peroxidase-conjugated anti-mouse immunoglobulin (rabbit) is added to the wells. Following the second washing step, the wells are incubated with a chromogenic substrate. The reaction is stopped and the absorbance is measured.

4. Warnings and Precautions

The Urine CartiLaps® (CTX-II) EIA is for *in vitro* diagnostic use only and is not for internal use in humans or animals. This product must be used strictly in accordance with the instructions set out in these Instructions For Use (IFU). Immunodiagnostic Systems Limited (IDS) will not be held responsible for any loss or damage (except as required by statute), howsoever caused, arising out of non-compliance with the instructions provided.

CAUTION: This kit contains material of animal origin. Handle kit reagents as if capable of transmitting an infectious agent. Appropriate precautions and good laboratory practice must be used in the storage, handling and disposal of the kit reagents. Disposal of kit reagents should be in accordance with local regulations.

Human materials

Human material used in the preparation of this product has been tested by FDA recommended assays for the presence of antibody to Human Immunodeficiency Virus (HIV I and II), Hepatitis B surface antigen, antibody to Hepatitis C, and found negative. As no test can offer complete assurance that infectious agents are absent, the reagents should be handled according to Biosafety Level 2.

Reagents containing Sodium Azide

Some reagents in this kit contain sodium azide (NaN₃) <0.1 % (w/w) which may react with lead, copper or brass plumbing to form highly explosive metal azides. On disposal, flush with large volumes of water to prevent azide build up.

Classification under CLP:

EUH208

Hazard statements:

EUH208 Contains a mixture of: 5-chloro-2-methyl-2h-isothiazol-3-one [ec no 247-500-7] and 2-methyl-2h-isothiazol-3-one [ec no 220-239-6]. May produce an allergic reaction.

Precautionary statements:

N/A

5. Shelf Life And Storage Of Reagents

Store the Urine CartiLaps® (CTX-II) EIA kit upon receipt at 2-8°C. Under these conditions the kit is stable up to the expiry date stated on the box.

Do not use kit components beyond the expiry date and do not mix reagents from different lots.

Indications of possible deterioration of kit reagents include:

- The presence of abnormal particulate matter in any of the reagents.
- A decrease in the maximum binding.
- A high non-specific binding.
- A shift in the slope of the curve from its normal position.

6. Sample Collection and Storage

- For optimal results it is recommended to use urine from second morning void. Spot urine sample may be used.
- For monitoring purpose, follow up samples should be collected under same conditions as the baseline sample.
- Prior to use, urine specimens should be shaken and sedimentation allowed for a minimum of 30 minutes.
- The urine sample can be stored at 2-8°C for a maximum of 24 hours; store the samples at -20°C or lower for longer storage.
- Specimens obviously contaminated with whole blood may interfere with assay performance. These specimens should be discarded and a new specimen collected

Note:

- Samples containing particulate matter must be centrifuged before performing the assay.
- Samples displaying microbial contamination should not be assayed with the kit.
- Before performing assays, make sure that samples, calibrators and controls are at room temperature (18 - 22 °C).
- Each laboratory should follow the guidelines or requirements of local, state, and/or federal regulations or accrediting organizations to establish its own specimens handling and storage stability. For guidance on appropriate practices, please refer to the CLSI GP44-A4, "Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests".

7. Materials

Materials Provided

MICROPLAT	Streptavidin coated microtiter plate Microwell strips (12x8 wells) pre-coated with streptavidin. Supplied in a plastic frame.
CAL0	Calibrator Ready to use buffered solution with protein stabiliser, detergent and preservative; 1 vial, 3.0 mL
CAL 1 -5	Calibrators Ready to use buffered solution containing a synthetic peptide with protein stabilisers, detergent and preservative; 1 each of 5 concentration levels, 0.4 mL per vial The exact value of each standard is printed on the QC report.
CTRL 1 - 2	Controls Ready to use buffered solution containing a synthetic peptide with protein stabilisers, detergent and preservative; 1 each of 2 concentration levels, 0.4 mL per vial The established ranges for the controls are printed on the QC report.
Ag BIOTIN	Biotinylated Antigen Biotinylated synthetic peptide prepared in a buffered solution with protein stabiliser, detergent and preservative; 1 vial, 12.0 mL
Ab	Primary Antibody Ready to use buffered solution containing monoclonal antibody with protein stabiliser, detergent, preservative and red dye; 1 vial, 12.0 mL
ENZYMCONJ	Peroxidase conjugated antibody Ready to use peroxidase conjugated anti mouse immunoglobulins (rabbit) in a buffered solution with protein stabiliser, detergent, preservative and blue dye; 1 vial, 12.0 mL.
SUBS TMB	Substrate Solution Ready to use tetramethylbenzidine (TMB) substrate in an acidic buffer; 1 vial, 12.0 mL Please note that the chromogenic substrate might appear slightly blueish.
H₂SO₄	Stopping Solution Ready to use solution of 0.18 mol/L sulphuric acid; 1 vial, 12.0 mL
WASHBUF 50x	Washing Buffer Concentrated washing buffer with detergent and preservative; 1 vial, 20.0 mL

Adhesive Plate Sealer Adhesive film for covering wells during incubation.

Documentation Instructions for Use and QC report.

Materials required — not supplied

- Containers for preparing the Washing Solution
- Precision pipetting devices to deliver 40 µL
- Distilled water
- Precision 8 or 12 channel multipipette to deliver 100 µL to 300 µL
- Vortex mixer
- Automatic microplate washer (optional)
- Photometric microplate reader and data analysis equipment

8. Preparation Of Reagents

Allow all reagents to come to room temperature (18 - 22°C) for a minimum of 60 minutes before use. Do not interchange kit components from different lots.

Wash buffer preparation

Prepare by adding 1-part Wash Concentrate **WASHBUF 50x** to 50-parts distilled water.

All other reagents are supplied ready for use.

N.B. To avoid potential microbial and / or chemical contamination, unused reagents should never be returned into the original vials.

9. Assay Procedure

Prepare reagents as described in § 8. Preparation of Reagents. Mix all reagents and samples before use (avoid formation of foam).

NOTE: To ensure consistent results between runs, between operators, and to minimise any drift effect; strictly adhere to the following procedure:

- a. Bring all reagents to room temperature (18 – 22 °C) prior to use – this will take approximately 60 minutes.
- b. Seal the plate during incubations using the plate sealers which are supplied with the assay kit.
- c. Do not stack plates during incubation in order to ensure a consistent temperature for all plates.
- d. Do not under or over-fill the assay wells during the washing steps.
- e. Add reagents in the same sequence each time to reduce time deviation between reactions

Do not pipette directly from the vial containing TMB substrate. The required volume should first be transferred to a clean container. Solution remaining in the container should be discarded following use and NOT returned to the stock vial **SUBS TMB**

Determine the number of strips needed for the assay; it is recommended to test all samples in duplicate. In addition, for each run a total of 16 wells are needed for the standards and controls. Place the appropriate number of strips in the plastic frame. Store any unused strips in the tightly closed foil bag with desiccant capsules.

1. Add 100 µL of Biotinylated Antigen **Ag BIOTIN** to each well, cover with sealing tape, and incubate for 30±5 minutes at room temperature (18-22°C) without shaking.
2. Wash all wells 5 times with wash buffer

Automatic plate wash	Set plate washer to dispense 300 µL of wash solution per well
	Fill and aspirate for 5 cycles
Manual wash	Decant the contents of the wells by inverting sharply
	Pipette 300 µL of wash solution into each well, decant and repeat 5 times
	Remove excess wash buffer by tapping firmly on absorbent tissue before proceeding
	Make sure wells are completely emptied after each manual or automatic washing cycle.
	Proceed immediately to the next step at the end of washing.
3. Pipette 40 µL of Calibrators **CAL 0 – 5**, Controls **CTRL 1 - 2** or unknown samples into appropriate wells followed by 100 µL of the Antibody Solution **Ab**.
4. Cover the immunostrips with sealing tape and incubate for 21±3 hrs. at 2-8°C without shaking.
5. Wash all wells as step 2
6. Add 100 µL of the Peroxidase-Conjugated Antibody **ENZYMCONJ** to each well.
7. Cover the immunostrips with sealing tape and incubate for 60±5 minutes at room temperature (18-22°C) without shaking.
8. Wash all wells as step 2
9. Pipette 100 µL of the Substrate Solution **SUBS TMB** into each well and incubate for 15±2 minutes at room temperature (18-22°C) in the dark. Use sealing tape.

NOTE: do not pipette directly from the vial containing TMB substrate. The required volume should first be transferred to a clean container. Solution remaining in the container should be discarded following use and NOT returned to the stock vial **SUBS TMB**
10. Pipette 100 µL of the Stopping Solution **H₂SO₄** into each well.
11. Measure the absorbance at 450 nm with 650 nm as reference within two hours.

N.B. Microplate readers measure vertically; when pipetting, do not touch the bottom of the wells

Automated Platforms

The Urine CartiLaps® (CTX-II) EIA kit was designed and developed to be performed manually using the protocol described above. The protocol is not necessarily applicable to automated platforms.

If automated platforms are used it is the responsibility of the user to ensure the kit has been appropriately tested. To improve the performance of the kit on automated platforms, it is recommended to increase the number of wash cycles at each wash step.

10. Calculation of Results

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic curve fit should be used.

If the absorbance of a pre-diluted sample is above **Calibrator 5**, the sample should be diluted in **Calibrator 0** and re-analysed.

For each urine sample the Urine CTX-II EIA concentration ($\mu\text{g/L}$) and the creatinine concentration ($\text{mM} = \text{mmol/L}$) should be determined using an enzymatic colorimetric method for clinical chemistry analyser.

The following equation corrects the CTX-II concentration for variation in urine concentration:

$$\text{Corrected CTX-II Value (ng/mmol)} = [1000 \times \text{CTX-II } (\mu\text{g/L})] \div \text{Creatinine (mmol/L)}$$

11. Quality Control

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples and the results analysed with appropriate statistical methods.

The two kit controls provided in the kit are intended to assist in assessing the validity of results obtained with each assay plate.

IDS recommend the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories¹⁸.

12. Measurement Range

Detection limit 0.20 $\mu\text{g/L}$

The detection limit was determined to 0.20 $\mu\text{g/L}$, which is the concentration corresponding to three standard deviations below the mean of 21 determinations of the absorbance of the Urine CartiLaps Standard 0.

13. Limitations of Use

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays¹⁹. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.
- The following substances do not interfere in the Urine CartiLaps[®] (CTX-II) EIA when the concentrations presented in the following table are below the stated threshold.

Potentially Interfering Agent	Threshold Concentration
Urea	30 g/L
Creatinine	10 mg/L
Glucose	5 mg/L
Ascorbic acid	5 mg/L
Albumin	50 mg/L
Ibuprofen	50 g/L
Acetyl-salicylic acid	50 g/L
Paracetamol	50 g/L

14. Expected Values

It is advisable for each laboratory to establish its own range of healthy and pathological CTX-II values. As an example, the geometric mean values and 95% confidence interval (CI) for various populations are given below. For further information, please refer to the reference list at the end of these instructions.

Population	Number of subjects	Age (years)	Mean CTX-II (ng/mmol)	95% CI (ng/mmol)
All women	459	20 – 85	299	79 – 1137
Premenopausal	38	20 – 30	464	103 – 2086
	165	30 - 59	200	65 – 618
	28	48 – 53	164	66 – 410
Postmenopausal	38	48 – 53	318	89 – 1132
	256	46 – 85	363	112 – 1172
All male	247	22 – 87	278	87 – 895
Male	27	20 – 30	501	214 – 1171
	141	30 – 60	236	89 – 628
	79	> 60	305	85 – 1096

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

15. Performance Data

Representative performance data are shown. Results obtained at individual laboratories may vary.

15.1 Precision

The precision was determined using ten analytical runs, each with duplicate determinations of urine samples.

Sample	Mean Conc. (µg/L)	Intra-assay		Inter-assay	
		SD (µg/L)	CV%	SD (µg/L)	CV%
Low	0.52	0.04	7.8	0.06	12.2
Medium	1.84	0.08	4.6	0.20	10.8
High	5.50	0.28	5.2	0.38	6.9

15.2 Dilution /Linearity


The dilution recovery of the Urine CartiLaps® (CTX-II) EIA was determined as 96%. Four urine samples were diluted in Urine CartiLaps Standard 0. The concentration of CTX-II was determined in the UrineCartiLaps® (CTX-II) EIA and the recovery calculated after correction with the dilution factor.

Dilution Factor	Recovery (%)				Mean
	1	2	3	4	
1 x	100	100	100	100	96
2 x	92	104	92	98	
4 x	88	105	86	93	
8 x	84	104	89	108	

15.3 Specificity

The epitope detected in the Urine CartiLaps® (CTX-II) EIA is highly conserved and therefore the test can be applied to urine samples from many other species, including non-human primates, bovines, horses, pigs, rabbits, rats and mice.

16. Symbols used

	Catalogue Number
	<i>In Vitro</i> Diagnostic Device
	Manufacturer
	Applied in accordance with directive 98/79/EC
	Caution Federal Law restricts this device to sale by or on the order of a (licensed healthcare practitioner) (Only for US)
	EU Importer

17. REFERENCES


1. Ceunick F De. et al., Urinary collagen type II C-telopeptide fragments are sensitive markers of matrix metallo-proteinase dependent cartilage degradation in rat adjuvant induced arthritis. *J Rheumatol* (2003); 30: 1561-1564.
2. Christgau S. et al., Suppression of Elevated Cartilage Turnover in Postmenopausal Women and in Ovariectomized Rats by Estrogen and a Selective Estrogen-Receptor Modulator (SERM). *Menopause* (2004); 11: 508-518.
3. Christgau S. et al., Cartilage Degradation In Glucosamine Sulphate Treated Knee Osteoarthritis Patients With Elevated Levels Of Urinary Collagen Type II C-Telopeptide Fragments. *Clin Exp Rheumatol.* (2004) 22(1):36-42.
4. Christgau S. et al., Collagen type II degradation products in urine as an index of cartilage degradation. *Bone* (2001); 29: 209-215.
5. Forsblad d'Elia H. et al., Hormone replacement therapy, calcium and vitamin D3 versus calcium and vitamin D3 alone decreases markers of cartilage and bone metabolism in rheumatoid arthritis: a randomized controlled trial. *Arthritis Res Ther* (2004); 6: R457-68.
6. Garnero P. et al., Association of 10 molecular markers of bone, cartilage and synovium with disease activity and joint damage in hip osteoarthritis patients: the ECHODIAH cohort. *ACR* 2003.
7. Garnero P. et al., Association of baseline levels of markers of bone and cartilage degradation with long-term progression of joint damage in patients with early rheumatoid arthritis: the COBRA Study. *Arthritis & Rheum* (2002); 46:2847-2856.
8. Garnero P. et al., Uncoupling of type II collagen synthesis and degradation predicts progression of joint damage in patients with knee osteoarthritis. *Arthritis & Rheum* (2002); 46:2613-2624.
9. Garnero P. et al., Association of baseline levels of urinary glucosyl-galactosyl pyridinoline and type II collagen C-telopeptide with progression of joint destruction in patients with early rheumatoid arthritis. *Arthritis & Rheum* (2002); 46: 21-30.
10. Garnero P. et al., The bisphosphonate Zoledronate decreases type II collagen breakdown in patients with Paget's disease of bone. *Bone* (2001); 28: 461-464.
11. Garnero P. et al., Cross sectional evaluation of biochemical markers of bone, cartilage, and synovial tissue metabolism in patients with knee osteoarthritis: relations with disease activity and joint damage. *Ann Rheum Dis* (2001); 60: 619-26.
12. Høegh-Andersen P. et al., Ovariectomized Rats as a Model of Postmenopausal Osteoarthritis. Validation and Application. *Annals of Rheum Dis.* (2004); 6(2): R169-80.
13. Jung M. et al., Elevated urinary concentration of collagen type II C-telopeptide fragments in patients with osteoarthritis. *Pathobiology.* (2004); 71(2): 70-6
14. Lehmann HJ. et al., The effects of bisphosphonates on CartiLaps: A new marker for cartilage degradation. *Ann Rheum Dis* (2002); 61:530-533.
15. Mazières B. et al., Molecular markers of cartilage breakdown and synovitis are strong independent predictors of structural progression of hip osteoarthritis (OA). the ECHODIAH cohort. *ACR* 2003.
16. Mouritzen U. et al., CartiLaps: A novel marker of Cartilage Degradation. The influence of age, gender, menopause, hormone replacement therapy and bone mass index. *Annals Rheum Dis.* (2003); 62: 332-336.
17. Roy-Beaudry M. et al., Entothelin-1 promotes osteoarthritic cartilage degradation via mmp-1 and mmp-13 induction. *Arthritis & Rheum* (2003); 48:2855-2864.
18. Basic QC Practices On-line Course; <http://www.Westgard.com>.
19. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

Immunodiagnostic Systems Limited
 10 Didcot Way, Boldon Business Park
 Boldon, Tyne & Wear, NE35 9PD, UK
 Tel.: +44 191 519 0660
 Fax: +44 191 519 0760
 e-mail: info.uk@idsplc.com • www.idsplc.com
 IDS Immunodiagnostic Systems GmbH
 Herriotstrasse 1
 60528 Frankfurt am Main
 Germany

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EU Dovoze

17. LITERATURA

1. Ceunick F De. et al., Urinary collagen type II C-telopeptide fragments are sensitive markers of matrix metallo-proteinase dependent cartilage degradation in rat adjuvant induced arthritis. *J Rheumatol* (2003); 30: 1561-1564.
2. Christgau S. et al., Suppression of Elevated Cartilage Turnover in Postmenopausal Women and in Ovariectomized Rats by Estrogen and a Selective Estrogen-Receptor Modulator (SERM). *Menopause* (2004); 11: 508-518.
3. Christgau S. et al., Cartilage Degradation In Glucosamine Sulphate Treated Knee Osteoarthritis Patients With Elevated Levels Of Urinary Collagen Type II C-Telopeptide Fragments. *Clin Exp Rheumatol.* (2004) 22(1):36-42.
4. Christgau S. et al., Collagen type II degradation products in urine as an index of cartilage degradation. *Bone* (2001); 29: 209-215.
5. Forsblad d'Elia H. et al., Hormone replacement therapy, calcium and vitamin D3 versus calcium and vitamin D3 alone decreases markers of cartilage and bone metabolism in rheumatoid arthritis: long-term progression of joint damage in patients with early rheumatoid arthritis: the COBRA Study. *Arthritis & Rheum* (2002); 46:2847-2856.
6. Garnero P. et al., Uncoupling of type II collagen synthesis and degradation predicts progression of joint damage in patients with knee osteoarthritis. *Arthritis & Rheum* (2002); 46:2613-2624.
7. Garnero P. et al., Association of baseline levels of urinary glucosyl-galactosyl pyridinoline and type II collagen C-telopeptide with progression of joint destruction in patients with early rheumatoid arthritis. *Arthritis & Rheum* (2002); 46: 21-30.
8. Garnero P. et al., The bisphosphonate Zoledronate decreases type II collagen breakdown in patients with Paget's disease of bone. *Bone* (2001); 28: 461-464.
9. Garnero P. et al., Cross sectional evaluation of biochemical markers of bone, cartilage, and synovial tissue metabolism in patients with knee osteoarthritis: relations with disease activity and joint damage. *Ann Rheum Dis* (2001); 60: 619-26.
10. Høegh-Andersen P. et al., Ovariectomized Rats as a Model of Postmenopausal Osteoarthritis. Validation and Application. *Annals of Rheum Dis.* (2004); 6(2): R169-80.
11. Jung M. et al., Elevated urinary concentration of collagen type II C-telopeptide fragments in patients with osteoarthritis. *Pathobiology.* (2004); 71(2): 70-6
12. Lehmann HJ. et al., The effects of bisphosphonates on CartiLaps: A new marker for cartilage degradation. *Ann Rheum Dis* (2002); 61:530-533.
13. Mazières B. et al., Molecular markers of cartilage breakdown and synovitis are strong independent predictors of structural progression of hip osteoarthritis (OA). the ECHODIAH cohort. *ACR* 2003.
14. Mouritzen U. et al., CartiLaps: A novel marker of Cartilage Degradation. The influence of age, gender, menopause, hormone replacement therapy and bone mass index. *Annals Rheum Dis.* (2003); 62: 332-336.
15. Roy-Beaudry M. et al., Entothelin-1 promotes osteoarthritic cartilage degradation via mmp-1 and mmp-13 induction. *Arthritis & Rheum* (2003); 48:2855-2864.

16. Basic QC Practices On-line Course; <http://www.Westgard.com>.
17. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

Immunodiagnostic Systems Limited

10 Didcot Way, Boldon Business Park
Boldon, Tyne & Wear, NE35 9PD, UK
Tel.: +44 191 519 0660
Fax: +44 191 519 0760

e-mail: info.uk@idsplc.com • www.idsplc.com

IDS Immunodiagnostic Systems GmbH
Herriotstrasse 1
60528 Frankfurt am Main
Germany

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