Ostase® BAP EIA

Immunoenzymetric Assay

For the quantitative measurement of
Bone-specific Alkaline Phosphatase (BAP) in
human serum

INTENDED USE

The Ostase® BAP EIA Immunoenzymetric Assay is an in vitro device indicated for the quantitative measurement of bone-specific alkaline phosphatase (BAP), an indicator of osteoblastic activity, in human serum. This device is intended to be used as an aid in the management of postmenopausal osteoporosis and Paget's disease.

SUMMARY AND EXPLANATION

Bone is a dynamic tissue in which bone formation and bone removal (also referred to as resorption) continue throughout life in a process called remodeling. The remodeling process is a function of complex interactions between two types of bone cells: osteoblasts for the formation of bone, and osteoclasts for the resorption of bone (1-3). Bone formation and resorption are interdependent processes that are, under normal circumstances, tightly coupled (2,4). This coupled relationship is integral to maintaining the biochemical competence of the skeleton, thereby preserving the organization of bone structure, form, and strength (2,3,5). Serum levels of BAP are believed to reflect the metabolic status of osteoblasts (6,7). An accurate assessment of bone metabolism is critical for determining the severity of metabolic bone disease and responses to therapy. Measurement of serum levels of BAP has been shown to be useful in evaluating patients with Paget's disease, osteomalacia, primary hyperparathyroidism, renal osteodystrophy, osteoporosis and metastases to bone (6-10). Total alkaline phosphatase determinations have been the accepted method for the diagnosis and monitoring of patients with Paget's disease.

Paget's disease of bone is a common skeletal disorder in which there is a focal proliferation of the normal cellular components of bone. Paget's disease is more prevalent than once thought with the incidence rate in certain populations at 3%-4% in middle-aged patients and 10%-15% in the elderly (11). This disease does not affect young individuals. The majority of patients with Paget's disease have no symptoms and often go undiagnosed unless an abnormal X-ray or serum alkaline phosphatase level is found in the course of a medical evaluation for unrelated reasons. The most common complaints in symptomatic patients are pain and deformity.

The risk of osteoporosis, another bone remodeling disorder, depends in part upon skeletal development, the attainment of peak bone mass, and in later life, the amount of bone lost. In healthy children, bone formation is favored over bone resorption, which results in bone development and normal skeletal growth (3). In healthy young adults, bone formation and bone resorption are balanced, resulting in no net increase or decrease in skeletal mass. With advancing age, men and women experience an imbalance in bone remodeling in which resorption is slightly greater than formation, resulting in a continuous net loss of bone mass with time (1,2,4,12). If this imbalance persists, bone mass may decline until the skeleton is insufficient to withstand normal mechanical stresses, and it becomes abnormally susceptible to fractures. The excessive loss of bone mass with an increased susceptibility to fractures is a disorder known as osteoporosis (5).

The most common form of osteoporosis occurs in postmenopausal women and is the result of estrogen deficiency (2,12,13). Rapid bone loss accompanies the decline of estrogen levels at the onset of menopause or as a result of surgical removal of the ovaries. Rapid bone loss occurs as a result of the combined effects of an imbalance in bone remodeling and an increase in bone turnover (5,14-16). In the United States, osteoporosis affects some 25 million postmenopausal women and is the cause of approximately

1.5 million fractures annually, including approximately 500,000 vertebral crush fractures, 250,000 hip fractures, and 200,000 distal radius fractures (2,5,17).

Hormone replacement therapy is currently the most widely prescribed therapy for the prevention of osteoporotic fractures in postmenopausal women (4,5,18-20). However many women cannot, or will not, avail themselves of hormone replacement therapy because of the potential for the increased risk of cancer and the resumption of menstrual bleeding. For this reason, other compounds such as bisphosphonates, a standard treatment for Paget's disease of bone, have been developed to treat osteoporosis. The anti-resorptive properties of bisphosphonates decrease bone remodeling and, consequently, decrease the overall loss of bone.

Biochemical markers are useful in monitoring metabolic bone disease. Urinary hydroxyproline and total serum alkaline phosphatase have been used for monitoring the treatment of Paget's disease. Osteoporosis, however, represents a more subtle modification of the bone remodeling process; therefore, more specific and sensitive markers are needed.

The Ostase® BAP assay is an in vitro device for the quantitative measurement of bone-specific alkaline phosphatase (BAP) in human serum. Changes in BAP have been shown to be useful in patients undergoing therapy for metabolic bone disorders (6,7,10,21,22).

PRINCIPLES OF THE PROCEDURE

The Ostase® EIA assay is a solid phase, monoclonal antibody immunoenzymetric assay. Samples containing BAP are reacted with a solution containing a biotin-labeled, BAP-specific monoclonal antibody. The reaction takes place in plastic well strips (solid phase) coated with streptavidin and enclosed in a plastic frame. Following the formation of a solid phase/capture antibody/BAP complex, the microplate is washed to remove unbound BAP and is then incubated with an enzyme substrate. The amount of substrate turnover is determined colorimetrically by measuring the absorbance of the quenched reaction at 405 nm in a

microplate reader. The absorbance is proportional to the concentration of BAP present in the test sample. The calculation of BAP concentration in the sample is based on concurrent testing of BAP calibrators and Zero Calibrator/Diluent.

PRODUCT INFORMATION

Components	Cat No.	AC-20F1
Conjugate CONJ		96 Tests 1 x 14 mL
Anti-BAP (mouse monoclonal IgG) with biotin in a bovine/horse protein matrix with 0.09% sodium azide.		
Microplates MICROPLAT		12 x 8 wells
Streptavidin coated plastic well strips in a plastic tray. Desiccant: silica gel.		
Zero Calibrator/Diluent (0) CAL 0		1 x 14 mL
A bovine protein matrix containing no detectable concentration of BAP (0 µg BAP/L) and 0.09% sodium azide.		
Calibrators (1-5) CAL 1 - 5		5 x 1 mL
A bovine protein matrix containing approximately 7, 15, 30, 60, and 90 µg human BAP/L and 0.09% sodium azide.		
Low Control (1) CTRL 1		1 x 1 mL
A bovine protein matrix containing approximately 11 µg human BAP/L and 0.09% sodium azide. Refer to Quality Control Report for assigned range.		
High Control (2) CTRL 2		1 x 1 mL
A bovine protein matrix containing approximately 45 µg human BAP/L and 0.09% sodium azide. Refer to Quality Control Report for assigned range.		
Wash Concentrate WASHBUF 20X		1 x 50mL
Phosphate buffered saline containing Tween.		
Substrate SUBS pNPP		1 x 20 mL
p-nitrophenyl phosphate in a stabilizing buffer containing preservatives		
Quench Reagent NaOH		1 x 14 mL
1 N sodium hydroxide.		
Quality Control Report, Package Insert		1 each

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only.
- Do not pipette by mouth.
- Do not eat, drink or smoke in designated work areas.
- Wash hands thoroughly after handling specimens and kit reagents.
- Some components in this kit contain material of human origin which has been tested using FDA-approved methods and has been found negative for antibody to human immunodeficiency virus (HIV-1 and HIV-2), antibody to Hepatitis C virus and for Hepatitis B surface antigen (HBsAg). No known test method can offer total assurance that HIV-1, HIV-2, Hepatitis B virus, Hepatitis C virus or other infectious agents are absent. Handle these reagents as if they were potentially infectious (23).
- Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up (24).
- The Quench Reagent for this kit contains 1 N NaOH, a corrosive which causes burns. Avoid contact with skin and eyes. Do not ingest or inhale. Safety glasses, gloves and lab coat should be worn when working in the laboratory.
- Avoid microbial contamination of the reagents when dispensing aliquots.

STORAGE AND STABILITY

- Ostase® BAP EIA assay reagents are to be stored between 2 and 8°C, at which temperatures they are stable until the expiration date printed on the box label.
- Wash and Quench reagents are stable at 2-30°C until the expiration dates printed on the bottle labels.
- All reagents must be brought to room temperature (18-25°C) prior to use. After use all reagents except the Wash solution should be stored at 2-8°C.
- Recovery of the kit control concentrations should fall within the stated ranges.
- Unused well strips should be returned to the plastic storage tray with the desiccant pouch

provided and kept at 2–8°C.

- Do not use kit components beyond the expiration date.
- The specimens storage and stability information stated above are general recommendations for use in a variety of settings of laboratories. Each laboratory should follow the guidelines or requirements of local, state, and/or federal regulations or accrediting organizations to establish its own specimens handling and storage stability. For guidance on appropriate practices, please refer to the CLSI GP44-A4, Procedures for Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline – Fourth Edition.

SPECIMEN COLLECTION AND PREPARATION

- No special preparation of the patient is necessary.
- A whole blood specimen should be obtained by acceptable medical technique.
- Allow the blood to clot and separate the serum by centrifugation.
- Specimens should be collected in such a way as to avoid hemolysis.
- Serum is required for the Ostase® BAP EIA assay. Performance with plasma samples has not been established.
- If the serum sample is to be assayed in 24 to 48 hours after collection, the specimen should be stored in a refrigerator at 2–8°C.
- Specimens held for longer times (up to 39 months) should be frozen at –80°C.
- Turbid serum samples or samples containing particulate matter should be centrifuged prior to assay.

MATERIALS REQUIRED BUT NOT PROVIDED

- Disposable tip precision pipettors: fixed volume or adjustable for 50, 100, and 150 µL (± 1%). Optionally, a multichannel precision pipette can be used together with disposable V-shaped troughs for addition of anti-BAP assay conjugate, Substrate and Quench reagents.
- Disposable pipette tips for 50, 100, 150 µL
- Test tubes for sample dilutions
- Microplate washer
- Aspiration device
- Distilled water
- Timer
- Container for storage of Wash solution
- Microplate reader (405 nm primary wavelength, 600–650 nm background subtraction) and data reduction software
- Horizontal microplate rotator (500 rpm to 900 rpm range)

For information on a complete reagent handling system, data reduction system, and liquid handling devices, contact your local sales representative.

PREPARATION OF REAGENTS

- Bring all reagents to room temperature (18–25°C) prior to use.
- Thoroughly mix reagents before each use by gentle agitation or swirling.
- Use a clean pipette tip for each specimen, calibrator or control to avoid contamination.
- Wash solution: To prepare Wash solution, add Wash Concentrate to 950 mL of distilled water and mix.

PROCEDURE

- The procedure for the Ostase® BAP EIA assay is performed at room temperature. Bring all serum specimens and kit components to room temperature (18–25°C) and mix well before use.
- The Zero Calibrator/Diluent, calibrators and controls should be tested in duplicate in the same plate holder.
- Because the termination of each incubation stops a reaction that is in progress (i.e., antibody binding or substrate turnover), reliable calibration of the assay depends on ensuring that the incubation times are essentially the same for all of the wells.

The assay procedure is as follows:

1. Label test well strips/holders appropriately and load strips into holder frame. Run in duplicate each of 6 calibrators, low and high controls for each assay.
2. Pipette 50 µL of Zero Calibrator/calibrators B-F/controls/patient specimens into the bottom of each assigned well.
3. Pipette 100 µL of conjugate into each well.
4. Incubate for 1 hour at 18–25°C using a horizontal rotator set at 500 rpm to 900 rpm.
5. Wash the microplate wells 3 times by:
 - a) Aspirating the liquid from the first strip.
 - b) Pipetting 300 µL of Wash solution into the first strip.
 - c) Repeating steps a and b on all subsequent strips.
 - d) Repeating steps a and c twice.
6. Dispense 150 µL of Substrate reagent into each well. Proceed to next step without delay.
7. Incubate for 13–15 minutes at 18–25°C using a horizontal rotator set at 500 rpm to 900 rpm. Proceed to next step without delay.
8. Dispense 100 µL of Quench Reagent into each well.
9. Read the absorbance of each well at 405 nm in a microplate reader, subtracting a blank reading at 600–650 nm for each well. Read the plate within 1 hour of quenching the reaction.
10. Calculate the results as described in "Instrumentation and Calculation of Results."

PROCEDURAL COMMENTS

1. Note on plate washing: Immunoassays require efficient washing to remove the unbound, biotin-labeled antibody. Therefore, it is very important to wash each well efficiently, removing the last droplets of Wash solution to achieve optimal results.
2. If a specimen was found to contain BAP at a concentration greater than the highest

calibrator, the specimen should be diluted with the Zero Calibrator/Diluent and assayed according to the assay Procedure. The dilution factor must be incorporated in the calculation of results. Each diluted specimen should be mixed thoroughly prior to testing. The recommended dilutions for specimens containing BAP greater than the highest calibrator are 1:3, 1:5, or 1:10.

However, it is desirable to dilute serum specimens with BAP greater than the highest calibrator so that the diluted sample reads greater than 10 µg BAP/L.

3. Because absorbance is a function of temperature and duration of the Substrate reagent incubation, it is very important that this incubation be the same for all wells/plates. This can be accomplished by ensuring that the elapsed time for pipetting reagents from beginning to end – without interruption – is exactly the same for both the Substrate reagent addition step and the Quench Reagent addition step. To ensure best results, the addition of these reagents should not exceed 90 seconds and total substrate incubation should not exceed 15 minutes.
4. For convenience, repeating or multichannel pipettors may be used for dispensing assay conjugate, Wash solution, Substrate and Quench reagents. Pipettors with disposable tips are recommended for pipetting calibrators, controls and specimens. The pipette tips should be changed after each sample is pipetted to avoid potential sample carryover and contamination of the reagents or specimens.
5. Do not mix materials from different kit lots.

INSTRUMENTATION AND CALCULATION OF RESULTS

Results may be calculated by using computer-assisted methods or manually on linear graph paper.

Computer-Assisted Method

A point-to-point curve fit is recommended. Point-to-point software that connects a straight line

between the means of calibrator replicates – including the 0 µg BAP/L calibrator – provides good results with the procedure and calibration method described. For additional information on computer-assisted data reduction, consult your local sales representative.

Manual Method


The Ostase® BAP EIA calibration curve may be constructed manually on linear graph paper by plotting the average absorbance for each calibrator on the y-axis versus the concentration of BAP in the Calibrator on the x-axis. A point-to-point curve should be drawn through the calibration points. Do not force the curve to a straight line.

To determine the concentration of BAP in the controls and patient specimens, extend a horizontal line from the absorbance value for the test sample to the calibration curve. At the point of intersection of the horizontal line and the curve, drop a vertical line to the x-axis and read the concentration of BAP.

If the absorbance for any replicate of the sample is greater than the absorbance of the highest Calibrator the specimen must be diluted and re-assayed. The observed concentration of the diluted sample must be multiplied by the dilution factor.

Example Data

Well #	Description	Abs. 405 nm	Mean Abs.	BAP µg/L
1	Calibrator (0)	0.193		
2		0.187	0.190	0.0
3	Calibrator (1)	0.401		
4		0.402	0.402	7.2
5	Calibrator (2)	0.595		
6		0.572	0.584	15.5
7	Calibrator (3)	0.908		
8		0.866	0.887	31.9
9	Calibrator (4)	1.499		
10		1.460	1.479	58.4
11	Calibrator (5)	2.000		
12		2.051	2.026	85.8
13	Patient Sample	0.911		33.0

QUALITY CONTROL AND ACCEPTABILITY OF RESULTS

Good laboratory practices include the use of control specimens within an assay run to ensure that all reagents and protocols are performing properly. The Ostase® BAP EIA antibody set contains controls which can be used to verify assay performance.

1. Recovery of the control concentrations should fall within the stated ranges.
2. The coefficient of variation (%CV) of the 405 nm absorbance readings for each calibrator and control sample should be less than 10%.

LIMITATIONS OF THE PROCEDURE

HAMA Interference: Some individuals have antibodies to mouse protein (HAMA), which can cause interference in immunoassays that employ antibodies derived from mice. In particular, it has been reported that serum samples from patients who have undergone therapy or diagnostic procedures that include infusion of mouse monoclonal antibody may produce erroneous results in such assays. Therefore, Ostase® BAP results for such patients should be used only in conjunction with results from some other diagnostic procedure and with information available from the clinical evaluation of the patient.

The immunoreactivity of liver ALP has been determined in the Ostase® BAP EIA assay: 100 U/L of liver ALP activity gives a result of 2.8 to 6.2 µg/L in the Ostase® BAP EIA assay. Serum samples with significant elevations of liver ALP activity may yield elevated results in the Ostase® BAP EIA assay. Patients with metabolic bone disorders who have low levels of disease activity may have bone-specific ALP levels that fall within the Ostase® BAP EIA assay expected values.

Ostase® BAP EIA results should be used only in conjunction with information available from the clinical evaluation of the patient and other diagnostic procedures. Therefore, the Ostase® BAP EIA assay is not recommended for use as a screening procedure to detect the presence of osteoporosis in the general population. Furthermore, the Ostase® BAP EIA assay cannot be used to assess the rate of bone formation or bone remodeling.

EXPECTED VALUES

The BAP results in the Expected Values and Clinical Studies sections were generated with the Tandem-R Ostase assay. However, a correlation study of the Ostase® BAP EIA and Tandem-R Ostase assays demonstrated good agreement between the assays ($y = 1.02x + 0.28$, $r = 0.9700$, $n = 136$).

The Tandem-R Ostase assay was evaluated in studies involving apparently healthy adults (20 to 89 years of age) at 6 test sites. The mean BAP concentration, standard deviation (SD), median, and 95th percentiles for males (N = 217), premenopausal females (N = 228), and postmenopausal females (N = 529) are presented in the following table.


Summary of BAP Concentrations in Apparently Healthy Adults*

	N	BAP Mean µg/L	SD	BAP Median µg/L	BAP 95th Percentile µg/L
Males	217	12.3	4.3	11.6	20.1
Pre-menopausal Females	228	8.7	2.9	8.5	14.3
Post-menopausal Females	529	13.2	4.7	12.5	22.4

* Results generated with the Tandem-R Ostase Immunoradiometric Assay

The results presented above show that mean BAP concentrations in a population of postmenopausal women are elevated over mean BAP concentrations in premenopausal women ($p \leq 0.0001$). This increase in mean BAP reflects the increase in bone remodeling associated with estrogen deficiency in a postmenopausal population compared to a premenopausal population (2,8,12,13). However, there is considerable overlap in BAP concentrations in the populations shown in the distribution below.

Distribution of BAP Concentrations in Apparently Healthy Adults*

* Results generated with the Tandem-R Ostase Immunoradiometric assay.

For each of the populations shown in the plot: 1) The horizontal lines in each box represent median concentrations. 2) The upper and lower margins of each box represent interquartile ranges – i.e., the upper margin delimits 25% of the values above the median while the lower margin delimits 25% of the values below the median. 3) The tip of each lower whisker represents the 5th percentile concentration, while the tip of each upper whisker represents the 95th percentile concentration.

The ranges of observations reported above are representative of this study only, and do not necessarily reflect the ranges that will be observed in a particular clinical laboratory. Each laboratory should establish its own reference interval.

CLINICAL STUDIES


Paget's Disease

Correlation studies were performed to compare serum BAP concentrations obtained using the Tandem-R Ostase assay to those results obtained using two currently marketed electrophoretic methods and two enzyme activity assays for total alkaline phosphatase (TAP). These studies were conducted at clinical investigation sites using 100 samples from patients with Paget's disease. The correlation coefficient (r) obtained between the Tandem-R Ostase assay and the electrophoretic methods is 0.9418. The correlation coefficient (r) obtained for the correlation between Tandem-R Ostase and the total alkaline phosphatase enzyme activity assays is 0.9459.


Postmenopausal Osteoporosis – Bisphosphonate Treatment

In order to demonstrate the ability of the Tandem-R Ostase assay to reflect alterations in bone remodeling in response to therapy in postmenopausal osteoporosis, a study was conducted in patients with clinically-defined osteoporosis who were treated with a bisphosphonate (alendronate sodium). This study was a multi-site, prospective, double-blind, placebo-controlled study (25). Data are available for subjects who were randomized to placebo (N = 148) or to alendronate 10 mg (N = 74). All subjects took 500 mg of supplemental calcium daily. Bone mineral density (BMD) measurements were made at the lumbar spine using dual-energy x-ray absorptiometry. Serial BAP determinations were made using the Tandem Ostase assay.

Mean percent changes from baseline and standard error bars for BMD in the placebo and treated groups are plotted at each timepoint in the figure below. A 7.8% increase in BMD over baseline is observed in the alendronate group at 24 months, indicating that the bisphosphonate therapy resulted in a positive effect on bone. As expected, a 0.8% decrease in BMD is observed in the calcium-supplemented placebo group at 24 months.

A summary of two-year data consisting of serial BAP concentrations as measured by the Tandem Ostase assay for the placebo and the alendronate 10 mg groups is presented below. Mean percent decreases from baseline and standard error bars for BAP for the placebo (baseline BAP = 17.7 µg/L ± 6.0) and alendronate-treated (baseline BAP = 17.0 µg/L ± 4.6) groups are plotted at each timepoint. The data show a time-dependent decrease in BAP concentrations for the alendronate-treated group, reaching a nadir at 6 months. As expected a small, transient change in marker levels is observed in the calcium-supplemented placebo group (26,27). The decrease and subsequent stabilization of the BAP changes in the alendronate-treated group suggest that a new steady state in bone remodeling occurred as a result of the bisphosphonate treatment (24,28,29).

The observed decreases in BAP concentrations for the alendronate-treated group, as measured by the Tandem Ostase assay, are significantly different from baseline ($p < 0.0001$) as early as 3 months and at all subsequent timepoints (6, 12, and 24 months). In contrast, BAP concentrations decreased by 11% at 3 and 6 months in the placebo group and returned to baseline by 24 months. The BAP changes and BMD changes in this study are consistent with current knowledge of bone remodeling and the mechanism of action of anti-osteoporotic therapies (24). From these results, it may be concluded that the decrease in BAP concentrations as measured by the Tandem Ostase assay in alendronate-treated subjects reflects changes in bone remodeling that result from anti-osteoporotic therapy.

The percent change between two serum marker values must exceed both the biological (within-subject) variation of the marker and the analytical variation of the assay in order to have clinical significance (30). This threshold is referred to as the minimum percent change. Percent changes in a marker that exceed the minimum percent change may be attributed to a change in the clinical status of the subject, such as response to therapy (30). The minimum percent change for a biochemical marker has been described by Soletormos (30) and others (31) with the formula

$$\text{Minimum Percent Change} = \sqrt{2} \times Z \times \sqrt{CV_p^2 + CV_a^2}$$

in which CV_p is within-subject variability, CV_a is assay imprecision, and Z is the Z-statistic, which depends on the probability selected for significance and on whether the change expected is unidirectional ($Z = 1.645$) or bidirectional ($Z = 1.96$).

The within-subject variability of the Tandem Ostase assay was determined from serum samples from 17 healthy postmenopausal women collected every day over a five-day period. The average biological coefficient of variation was calculated as 4.2%.

Analytical variation was based on between-run precision data found in the SPECIFIC PERFORMANCE CHARACTERISTICS section of this insert.


The minimum percent change for BAP as measured by the Tandem Ostase assay was calculated as 25%. Therefore, changes in BAP concentrations in postmenopausal women that exceed 25% may be attributed to changes in bone remodeling.

Percent changes of BAP relative to baseline for individual subjects are presented below. Data are presented for both the placebo and the alendronate 10 mg subjects at the 3-month and 6-month time-points. The horizontal line at zero represents baseline, and the horizontal line at -25% represents the minimum percent change for BAP as measured by the Tandem Ostase assay. For the alendronate 10 mg group, 77.0% (57/74) of the subjects had an BAP decrease from baseline (time 0) of 25% or greater at 3 months, and 85.1% (63/74) of the subjects had an BAP decrease from baseline of 25% or greater at the nadir (6 months). By 24 months, 90.3% (65/72) had a decrease in BAP of 25% or more. For the calcium-supplemented placebo group, 75.7% (112/148) of the subjects had an BAP decrease from baseline (time 0) that was less than 25% (at 6 months). By 24 months, 88.1% (126/143) of the placebo group had BAP decreases of less than 25%. These data are further summarized in the table provided below. This table shows the number of treated and placebo subjects (and the percent of subjects) who had BAP decreases from baseline < 25% or > 25% (at 3 and 6 months, respectively).

Summary of Individual BAP Percent Changes from Baseline for Alendronate 10 mg and Placebo Subjects Number of Subjects/Total Number of Subjects

BAP Percent Change Cutoff	3 Months		6 Months	
	Placebo N=148	Alendronate 10mg N=74	Placebo N=148	Alendronate 10mg N=74
25% or Greater Decrease	33/148 (22.3%)	57/74 (77.0%)	36/148 (24.3%)	63/74 (85.1%)
Less than 25% Decrease	115/148 (77.7%)	17/74 (23.0%)	112/148 (75.7%)	11/74 (14.9%)

In the alendronate 10 mg group, a mean BAP decrease of 45.7% from baseline to 6 months was followed by a mean increase in BMD of 7.8% at 24 months. Conversely, the small change from baseline for BAP at 6 months for the calcium-supplemented placebo group (-11.3%) was accompanied by a small decrease in BMD of -0.8% at 24 months.

In summary, it has been demonstrated that:

1. BAP concentrations reflect the increase in bone remodeling associated with estrogen deficiency in apparently healthy postmenopausal females.
2. BAP concentrations reflect changes in bone remodeling that occur as a result of anti-osteoporotic therapy with alendronate in postmenopausal women.
3. Clinically significant decreases in BAP concentrations at 3 and 6 months are indicators of changes in bone remodeling. Effective therapy with alendronate 10 mg in postmenopausal osteoporotic women is indicated by increases in BMD at 24 months.
4. A decrease in BAP of 25% or more from baseline, as measured by the Tandem Ostase assay, may indicate changes in bone remodeling resulting from anti-osteoporotic therapy.


Subjects undergoing anti-osteoporotic therapy with alendronate whose BAP levels do not decrease by 25% from baseline at early timepoints should be re-tested at later timepoints. Of the 11 subjects who demonstrated less than 25% decrease in BAP from baseline at 6 months, 9 of them eventually demonstrated greater than 25% decrease in BAP from baseline by 24 months. All of these subjects responded to alendronate 10 mg as determined by increases in BMD. Treated subjects whose BAP levels do not decrease by 25% from baseline should be evaluated by other clinical means to determine treatment efficacy.

Postmenopausal Women – Hormone Replacement Therapy (HRT)

In addition to the bisphosphonate study described above, other studies were initiated in order to demonstrate the ability of the Tandem Ostase assay to reflect alterations in bone remodeling in response to estrogen therapy (29,32). Data are presented below from one of the studies. In this study, healthy non-hysterectomized postmenopausal women (n = 12), aged 65-75 years, received daily estrogen/progestin therapy (Premarin 0.625 mg and Provera 2.5 mg) for two years. In addition, 6 healthy hysterectomized women received estrogen alone for a total study population of 18.

Skeletal ALP was measured by the Tandem Ostase assay at baseline and at 6, 12, and 24 months. BMD at the lumbar spine was evaluated at baseline and at 24 months.

Mean percent changes from baseline and standard error bars for BAP (left y-axis) and for BMD (right y-axis) are plotted at each timepoint in the figure below. A 6.6% increase in BMD at the lumbar spine over baseline is observed for the estrogen-treated subjects, indicating that the estrogen therapy resulted in a positive effect on bone. The data presented in the figure below also show a time-dependent decrease in BAP concentrations for the treated group, reaching a nadir at 12 months. The observed decreases in BAP concentrations for the estrogen treated group, as measured by the Tandem Ostase assay, are significantly different from baseline ($p \leq 0.0001$) at 6 months and at all subsequent timepoints (12 and 24 months).

The BAP and BMD changes in this study are consistent with current knowledge of bone remodeling and the mechanism of action of anti-osteoporotic therapies (5,33). From these results, it may be concluded that the decrease in BAP concentrations as measured by the Tandem Ostase assay in estrogen-treated subjects reflects changes in bone turnover that result from anti-osteoporotic therapy.

Tandem Ostase results should be used only in conjunction with information available from the clinical evaluation of the patient and other diagnostic procedures.

SPECIFIC PERFORMANCE CHARACTERISTICS

Within-Run Precision

Within-run precision was determined by assaying four serum pools containing various concentrations of BAP. Each serum pool was assayed in replicates of 20. The data are presented below.

Serum Pool	1	2	3	4
Number of Replicates	20	20	20	20
Mean µg/L	7.4	27.3	52.7	79.5
Standard Deviation	0.48	0.78	1.35	3.54
%CV	6.5	2.9	2.6	4.5

Between-Run Precision

Between-run precision was determined by duplicate measurements of four serum pools over a series of 20 individually calibrated runs:

Serum Pool	A	B	C	D
Number of Assays	20	20	20	20
Mean µg/L	8.4	29.2	55.6	81.1
Standard Deviation	0.47	1.88	2.03	4.9
%CV	5.8	6.4	3.7	2

Recovery and Dilution

Various quantities of a serum sample containing elevated levels of BAP were added to human sera containing endogenous BAP and the samples were assayed in triplicate.

Spike µg/L	Expected Concentration µg/L	Observed Concentration µg/L	n	% Recovery*
14.0	24.6	24.8	24	100.9
39.9	50.5	48.8	24	96.6
59.2	69.8	66.0	24	94.6
92.2	102.8	92.6	24	90.1

*% recovery equals the observed concentration divided by the expected concentration

A serum sample containing an elevated BAP concentration was diluted with the Zero Calibrator/Diluent and assayed at multiple dilutions.

Dilution	Expected Concentration µg/L	Observed Concentration µg/L	% Recovery*
Neat	N/A	54.6	N/A
1:4	27.3	27.9	102.3
1:6	18.2	17.4	95.4
1:8	13.7	12.7	92.8

*% recovery equals the observed concentration divided by the expected concentration times 100.

Interfering Substances

- Hemoglobin, which was tested at concentrations up to 500 mg/dL, does not interfere with the Ostase® BAP EIA assay.
- Unconjugated and conjugated bilirubin, which were tested at concentrations up to 40 mg/dL and 20 mg/dL, respectively, do not interfere with the Ostase® BAP EIA assay.
- Triglycerides, which were tested at concentrations up to 2000 mg/dL, do not interfere with the Ostase® BAP EIA assay.
- Total protein, which was tested at concentrations between 3 and 14 g/dL, does not interfere with the Ostase® BAP EIA assay.
- 100 U/L Intestinal ALP produced a result of 1.0 µg/L in the Ostase® BAP EIA assay.
- 100 U/L Placental ALP did not produce a detectable result in the Ostase® BAP EIA assay.

Liver ALP Reactivity

The liver ALP reactivity in the Ostase® BAP EIA assay was determined using serum samples from liver disease patients and serum samples from patients with Paget's disease. Samples were screened by electrophoresis and shown to contain > 95% liver ALP or BAP.

Two methods were used to evaluate liver ALP reactivity. The first procedure, described by Moss and Whitby (35), used a heat inactivation method that minimizes the contribution of endogenous BAP in the samples with elevated liver ALP. Using this method, 100 U/L liver ALP activity produced a result of 2.8 to 3.4 µg/L in the Ostase® BAP EIA assay.

The second procedure, described by Price, et al. (36), used the slope method Ostase® BAP EIA assay vs. total ALP), where no pretreatment of the samples occurred. In this study, Ostase® BAP EIA quantitation of the liver ALP and BAP samples (y-axis) was plotted against the total ALP activity (x-axis) in each sample. From the slope values of the liver ALP and BAP samples, it was determined that:

- 100 U/L of liver ALP activity produced a result of 6.2 µg/L in the Ostase® BAP EIA assay; and
- 100 U/L of BAP activity produced a result of 36.9 µg/L in the Ostase® BAP EIA assay

Interference by Drugs

Various concentrations of drugs were added to three separate serum pools containing BAP and assayed in quadruplicate. The drugs and highest concentrations tested are listed below. The testing was based on NCCLS Guideline EP-7P (Interference Testing in Clinical Testing).

acetaminophen	20 mg/dL
alendronate	5 mg/dL
aspirin	50 mg/dL
calcitonin-salmon	112 IU/dL
calcium	20 mg/dL
estrogen	400 ng/dL
etidronate	105 mg/dL
ibuprofen	40 mg/dL
pamidronate	18 mg/dL
progesterone	25 mg/dL
vitamin D	80,500 IU/dL

These drugs did not interfere with the recovery of BAP from the serum pools in the Ostase® BAP EIA assay.

Minimum Detectable Concentration

The minimum detectable concentration of BAP in the Ostase® BAP EIA assay is estimated to be 0.7 µg BAP/L. The minimum detectable concentration is defined as the concentration of BAP that corresponds to the absorbance that is two standard deviations greater than the mean absorbance of 20 replicate determinations of the Zero Calibrator/Diluent.

Purchase of this kit licenses its use under U.S. Patent Nos. 4,376,110 and 4,486,530.

Tandem® and Ostase® are registered trademarks of Hybritech Incorporated, a subsidiary of Beckman Coulter, Inc.

*Tween is a Trademark of ICI Americas, Inc.

REFERENCES/BIBLIOGRAPHIE/LITERATUR

- Riggs L. Causes of age-related bone loss and fractures, *Osteoporosis: Physiological Basis, Assessment, and Treatment*. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 7-16, 1990.
- Status of Research in Osteoporosis. *Osteoporosis Research, Education and Health Promotion*. U.S. Department of Health and Human Services, Public Health Service, National Institutes of Health.
- Canalis E. Regulation of bone remodeling, *Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism*. Favus MJ, ed., Raven Press (New York), pp. 33-37, 1993.
- Kleerekoper M and Avioli LV. Evaluation and treatment of postmenopausal osteoporosis, *Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism*. Favus MJ, ed., Raven Press (New York), pp. 223-229, 1993.
- Lindsay R. *Osteoporosis. A Guide to Diagnosis, Prevention, and Treatment*. Raven Press (New York), pp. 1-39, 1992.
- Epstein S. Serum and urinary markers of bone remodeling: assessment of bone turnover. *Endocrine Reviews* 9:437-448, 1988.
- Whyte MP. Alkaline phosphatase and the measurement of bone formation, *Clinical Disorders of Bone and Mineral Metabolism*. Potts JT and Frame B, eds., pp. 120-125, 1983.
- Garnero P and Delmas P. Assessment of the serum levels of bone alkaline phosphatase with a new immunoradiometric assay in patients with metabolic bone disease. *J Clin Endocrinol Metab* 77:1046-1053, 1993.
- Cooper EH, Whelan P, and Puryes D. Bone alkaline phosphatase and prostate-specific antigen in the monitoring of prostate cancer. *The Prostate* 25:236-242, 1994.
- Van Straalen JP, Sanders E, Prummel MF, and Sanders GTB. Bone alkaline phosphatase as indicator of bone formation. *Clin Chim Acta* 201:27-33, 1991.
- Merkow RL and Lane JM. Metabolic bone disease and Paget's disease in the elderly: II. Paget's Disease. *Clin Rheum Dis* 12:70, 1986.
- Lindsay R. Sex steroids in the pathogenesis and prevention of osteoporosis, *Osteoporosis: Etiology, Diagnosis, and Management*. Riggs L and Melton LJ, eds., Raven Press (New York), pp. 333-358, 1988.
- Peck W. The pathogenesis of postmenopausal osteoporosis, *Osteoporosis: Physiological Basis, Assessment and Treatment*. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 3-6, 1990.
- Christiansen C. Osteoporosis: Diagnosis and management today. *Bone* 17:513S-516S, 1995.
- Turner R, Riggs B and Spelsberg T. Skeletal effects of estrogen. *Endocrine Reviews* 15:275-300, 1994.
- Kanis J. Pathogenesis of osteoporosis and fracture, *Osteoporosis*. London: Blackwell Science Ltd.; pp. 22-55, 1994.
- Mack TM and Ross RK. A current perception of HRT risks and benefits, *Osteoporosis: Physiological Basis, Assessment and Treatment*. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 161-178, 1990.
- Riggs L and Melton LJ. The prevention and treatment of osteoporosis. *N Engl J Med* 327:620-627, 1992.
- Guidelines for Preclinical and Clinical Evaluation of Agents Used in the Prevention or Treatment of Postmenopausal Osteoporosis. Division of Metabolism and Endocrine Drug Products, Food and Drug Administration. April, 1994.
- Ribot C, Tremollieres F and Pouilles JM. Prevention of postmenopausal bone loss by long-term parenteral administration of 17B estradiol: comparison of percutaneous and transdermal route, *Osteoporosis: Physiological Basis, Assessment and Treatment*. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 137-145, 1990.
- Delmas P. Biochemical markers of bone turnover. Methodology and clinical use in osteoporosis. *Am J Med* 91:59S-63S, 1991.
- Duda RJ, O'Brien JF, Katzman JA, Peterson JM, Mann KG, and Riggs BL. Concurrent assays of circulating bone gla-protein and bone alkaline phosphatase: effects of sex, age, and metabolic bone disease. *J Clin Endocrinol Metab* 66:951-957, 1988.
- HHS Publication No. (CDC) 93-8395. 3rd edition. May 1993. Biosafety in Microbiological and Biomedical Laboratories. Washington, DC: US Government Printing Office, 1993.
- Manual Guide – Safety Management No. CDC-22. Decontamination of Laboratory Sink Drains to Remove Azide Salts. April 30, 1976. Atlanta, GA: Centers for Disease Control.
- Hill CS and Wolfert RL. The preparation of monoclonal antibodies which react preferentially with human bone alkaline phosphatase and not liver alkaline phosphatase. *Clin Chim Acta* 186:315-320, 1989.
- Garnero P, Shih WJ, Gineyts E, Karpf DB, and Delmas PD. Comparison of new biochemical markers of bone turnover in late postmenopausal osteoporotic women in response to alendronate treatment. *J Clin Endocrinol Metab* 79:1693-1700, 1994.
- Lieberman UA, Weiss SR, Broll J, Minne HW, Quan H, Bell NH, Rodriguez-Portales J, Downs RW, Dequeker J, Favus M, Seeman E, Recker RR, Capizzi T, Santora AC, Lombardi A, Shah RV, Hirsch LJ and Karpf DB. For the Alendronate Phase III Osteoporosis Treatment Study Group. Effect of oral alendronate on bone mineral density and the incidence of fractures in postmenopausal osteoporosis. *N Engl J Med* 333:1437-1443, 1995.
- Dawson-Hughes B, Dallal GE, Krall EA, Sadowski L, Sahyoun N and Tannenbaum S. A controlled trial of the effect of calcium supplementation on bone density in postmenopausal women. *N Engl J Med* 323:878-883, 1990.
- Reid IR, Ames RW, Evans MC, Gamble GD, and Sharpe SJ. Long-term effects of calcium supplementation on bone loss and fractures in postmenopausal women: a randomized controlled trial. *Am J Med* 98:331-335, 1995.
- Devogelaer JP, Broll H, Correa-Rotter R, Cumming DC, De Deuxchaisnes CN, Geusens P, Hosking D, Jaeger P, Kaufman JM, Leite M, Leon J, Lieberman U, Menkes CJ, Meunier PJ, Reid I, Rodriguez J, Romanowicz A, Seeman E, Vermeulen A, Hirsch LJ, Lombardi A, Plezia K, Santora AC, Yates AJ and Yuan W. Oral alendronate induces progressive increases in bone mass of the spine, hip, and total body over 3 years in postmenopausal women with osteoporosis. *Bone* 18:141-150, 1996.
- Overgaard K, Alexandersen P, Riis BJ and Christiansen C. Evaluation of a new commercial IRMA for bone-specific alkaline phosphatase during treatment with hormone replacement therapy and calcitonin. *Clin Chem* 42:973-974, 1996.
- Soletormos G, Schioler V, Nielsen D, Skovsgaard T, and Dombrowsky P. Interpretation of results for tumor markers on the basis of analytical imprecision and biological variation. *Clin Chem* 10:2077-2083, 1993.
- Blumsohn A, Hannon RA, Al-Dehaimi AW, and Eastell R. Short-term intraindividual variability of markers of bone turnover in healthy adults. *J Bone Min Res* 9:S153, 1994.
- Raisz L, Wiita B, Artis A, Bowen A, Schwartz S, Trahiotis M, Shoukri K and Smith J. Comparison of the effects of estrogen alone and estrogen plus androgen on biochemical markers of bone formation and resorption in postmenopausal women. *J Clin Endocrinol Metab* 81:37-43, 1996.
- Kanis J. Treatment of generalized osteoporosis with inhibitors of bone resorption, *Osteoporosis*. London: Blackwell Science Ltd.; pp. 168-195, 1994.
- Kennedy J, Carey R, Coolen R, Garber C, Hartman A, Lee H, Leitz V, Levine J, McLean M, Osberg I, Steindel S, and Sylvester E. User Evaluation of Precision Performance of Clinical Chemistry Devices. Tentative Guideline, NCCLS 4:8 EP5-T2, 1984.
- Moss D and Whitby L. A simplified heat-inactivation method for investigating alkaline phosphatase isoenzymes in serum. *Clin Chim Acta* 61:63, 1975.
- Price C, Mitchell C, Moriarty J, Gray M and Noonan K. Mass versus activity of an immunometric assay for bone alkaline phosphatase in serum. *Ann Clin Biochem* 32:405, 1995.

39.

 EXP	<p>GB Use By DE Verwendbar bis ES Fecha de caducidad IT Utilizzare entro FR Utiliser jusque NL Houdbaar tot DK Holdbar til CZ Použitelné do SK Použitelné do GR Ημερομηνία λήξης PT Prazo de validade HU Felhasználható SE Använd före PL Użyty przed</p>	 LOT	<p>GB Batch code DE Chargenbezeichnung ES Código de lote IT Codice del lotto FR Code du lot NL Lot nummer DK Lotnummer CZ Číslo šarže SK Číslo šarže GR Αριθμός Παρτίδας PT Código do lote HU Szarzsám SE Lot nummer PL Kod partii</p>
 REF	<p>GB Catalogue number DE Bestellnummer ES Número de catálogo IT Numero di catalogo FR Référence du catalogue NL Catalogus nummer DK Katalognummer CZ Katalogové číslo SK Katalogové číslo GR Αριθμός καταλόγου PT Referência de catálogo HU Katalógusszám SE Katalognummer PL Numer katalogowy</p>		<p>GB Manufacturer DE Hersteller ES Fabricante IT Fabbicante FR Fabricant NL Fabrikant DK Producent CZ Výrobce SK Výrobca GR Κατασκευαστής PT Fabricante HU Gyártó SE Tillverkare PL Producent</p>
	<p>GB Contains sufficient for <n> tests DE Inhalt ausreichend für <n> Prüfungen ES Contenido suficiente para <n> ensayos IT Contenuto sufficiente per "n" saggi FR Contenu suffisant pour "n" tests NL Inhoud voldoende voor "n" testen DK Indeholder tilstrækkeligt til "n" test CZ Lze použít pro <n> testů SK Obsah postačuje na <n> stanovení GR Περιεχόμενο επαρκές για «n» εξετάσεις PT Conteúdo suficiente para "n" ensaios HU A doboz tartalma <n> vizsgálat elvégzéséhez elegendő SE Räckert till "n" antal tester PL Wystarczy na wykonanie <n> testów</p>	 IVD	<p>GB In Vitro Diagnostic Medical Device DE In-Vitro-Diagnostikum ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro FR Dispositif médical de diagnostic in vitro NL Medisch hulpmiddel voor in-vitro diagnostiek DK Medicinsk udstyr til in vitro-diagnostik CZ In Vitro diagnostický zdravotnický prostředek SK Zdravotnícka pomôcka in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro HU In vitro diagnosztikum SE Medicintekniska produkter för in vitro diagnostik PL Wzrób do diagnostyki In Vitro</p>
	<p>GB Temperature limitation DE Temperaturbegrenzung ES Límite de temperatura IT Limiti di temperatura FR Limites de température NL Temperatuurlimiet DK Temperaturbegrænsning CZ Teplotní rozmezí od do SK Teplotné rozmedzie od do GR Περιορισμοί θερμοκρασίας PT Limites de temperatura HU Hőmérsékletkorlát SE Temperaturbegränsning PL Przestrzegać zakresu temperatury</p>		<p>GB Consult Instructions for Use DE Gebrauchsanweisung beachten ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso FR Consulter les instructions d'utilisation NL Raadpleeg de gebruiksaanwijzing DK Se brugsanvisning CZ Viz návod k použití SK Viď návod na použitie GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização HU Nézz meg a Használati utasítást SE Se handhavandebeskrivningen PL Sprawdź w instrukcji obsługi</p>

EU Importer

Immunodiagnostic Systems Limited,
10 Didcot Way, Boldon Business Park,
Boldon, Tyne & Wear, NE35 9PD, UK
Tel.: +44 191 519 0660
Fax: +44 191 519 0760
E-Mail: info.uk@idsplc.com
www.idsplc.com

IDS Immunodiagnostic Systems GmbH
Herriotstrasse 1
60528 Frankfurt am Main
Germany

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REFERENCES/BIBLIOGRAPHIE/LITERATUR

1. Riggs L. Causes of age-related bone loss and fractures, Osteoporosis: Physiological Basis, Assessment, and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 7-16, 1990.
2. Status of Research in Osteoporosis. Osteoporosis Research, Education and Health Promotion. U.S. Department of Health and Human Services, Public Health Service, National Institutes of Health.
3. Canalis E. Regulation of bone remodeling. Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. Favus MJ, ed., Raven Press (New York), pp. 33-37, 1993.
4. Kleerekoper M and Avioli LV. Evaluation and treatment of postmenopausal osteoporosis, Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. Favus MJ, ed., Raven Press (New York), pp. 223-229, 1993.
5. Lindsay R. Osteoporosis. A Guide to Diagnosis, Prevention, and Treatment. Raven Press (New York), pp. 1-39, 1992.
6. Epstein S. Serum and urinary markers of bone remodeling: assessment of bone turnover. Endocrine Reviews 9:437-448, 1988.
7. Whyte MP. Alkaline phosphatase and the measurement of bone formation, Clinical Disorders of Bone and Mineral Metabolism. Potts JT and Frame B, eds., pp. 120-125, 1983.
8. Garnero P and Delmas P. Assessment of the serum levels of bone alkaline phosphatase with a new immunoradiometric assay in patients with metabolic bone disease. J Clin Endocrin Metab 77:1046- 1053, 1993.
9. Cooper EH, Whelan P, and Purves D. Bone alkaline phosphatase and prostate-specific antigen in the monitoring of prostate cancer. The Prostate 25:236-242, 1994.
10. Van Straalen JP, Sanders E, Prummel MF, and Sanders GTB. Bone alkaline phosphatase as indicator of bone formation. Clin Chim Acta 201:27-33, 1991.
11. Merkow RL and Lane JM. Metabolic bone disease and Paget's disease in the elderly: II. Paget's Disease. Clin Rheum Dis 12:70, 1986.
12. Lindsay R. Sex steroids in the pathogenesis and prevention of osteoporosis, Osteoporosis: Etiology, Diagnosis, and Management. Riggs L and Melton LJ, eds., Raven Press (New York), pp. 333-358, 1988.
13. Peck W. The pathogenesis of postmenopausal osteoporosis, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 3-6, 1990.
14. Christiansen C. Osteoporosis: Diagnosis and management today. Bone 17:513S-516S, 1995.
15. Turner R, Riggs B and Spelsberg T. Skeletal effects of estrogen. Endocrine Reviews 15:275-300, 1994.
16. Kanis J. Pathogenesis of osteoporosis and fracture, Osteoporosis. London: Blackwell Science Ltd.; pp. 22-55, 1994.
17. Mack TM and Ross RK. A current perception of HRT risks and benefits, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 161-178, 1990.
18. Riggs L and Melton LJ. The prevention and treatment of osteoporosis. N Engl J Med 327:620-627, 1992.
19. Guidelines for Preclinical and Clinical Evaluation of Agents Used in the Prevention or Treatment of Postmenopausal Osteoporosis. Division of Metabolism and Endocrine Drug Products, Food and Drug Administration. April, 1994.
20. Ribot C, Tremolieres F and Pouilles JM. Prevention of postmenopausal bone loss by long-term parenteral administration of 17B estradiol: comparison of percutaneous and transdermal route, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 137-145, 1990.
21. Delmas P. Biochemical markers of bone turnover. Methodology and clinical use in osteoporosis. Am J Med 91:59S-63S, 1991.
22. Duda RJ, O'Brien JF, Katzman JA, Peterson JM, Mann KG, and Riggs BL. Concurrent assays of circulating bone gla-protein and bone alkaline phosphatase: effects of sex, age, and metabolic bone disease. J Clin Endocrin Metab 66:951-957, 1988.
23. HHS Publication No. (CDC) 93-8395. 3rd edition. May 1993. Biosafety in Microbiological and Biomedical Laboratories. Washington, DC: US Government Printing Office, 1993.
24. Manual Guide – Safety Management No. CDC-22. Decontamination of Laboratory Sink Drains to Remove Azide Salts. April 30, 1976. Atlanta, GA: Centers for Disease Control.
25. Hill CS and Wolfert RL. The preparation of monoclonal antibodies which react preferentially with human bone alkaline phosphatase and not liver alkaline phosphatase. Clin Chim Acta 186:315-320, 1989.
26. Garnero P, Ship WJ, Gineyts E, Karpf DB, and Delmas PD. Comparison of new biochemical markers of bone turnover in late postmenopausal osteoporotic women in response to alendronate treatment. J Clin Endocrin Metab 79:1693-1700, 1994.
27. Liberman UA, Weiss SR, Broll J, Minne HW, Quan H, Bell NH, Rodriguez-Portales J, Downs RW, Dequeker J, Favus M, Seeman E, Recker RR, Capizzi T, Santora AC, Lombardi A, Shah RV, Hirsch LJ and Karpf DB. For the Alendronate Phase III Osteoporosis Treatment Study Group. Effect of oral alendronate on bone mineral density and the incidence of fractures in postmenopausal osteoporosis. N Engl J Med 333:1437-1443, 1995.
28. Dawson-Hughes B, Dallal GE, Krall EA, Sadowski L, Sahyoun N and Tannenbaum S. A controlled trial of the effect of calcium supplementation on bone density in postmenopausal women. N Engl J Med 323:878-883, 1990.
29. Reid IR, Ames RW, Evans MC, Gamble GD, and Sharpe SJ. Long-term effects of calcium supplementation on bone loss and fractures in postmenopausal women: a randomized controlled trial. Am J Med 98:331-335, 1995.
30. Devogelaer JP, Broll H, Correa-Rotter R, Cumming DC, De Deuxchaisnes CN, Geusens P, Hosking D, Jaeger P, Kaufman JM, Leite M, Leon J, Liberman U, Menkes CJ, Meunier PJ, Reid I, Rodriguez J, Romanowicz A, Seeman E, Vermeulen A, Hirsch LJ, Lombardi A, Plezia K, Santora AC, Yates AJ and Yuan W. Oral alendronate induces progressive increases in bone mass of the spine, hip, and total body over 3 years in postmenopausal women with osteoporosis. Bone 18:141-150, 1996.
31. Overgaard K, Alexandersen P, Riis BJ and Christiansen C. Evaluation of a new commercial IRMA for bone-specific alkaline phosphatase during treatment with hormone replacement therapy and calcitonin. Clin Chem 42:973-974, 1996.
32. Soletormos G, Schioler V, Nielsen D, Skovsgaard T, and Dombrowsky P. Interpretation of results for tumor markers on the basis of analytical imprecision and biological variation. Clin Chem 10:2077- 2083, 1993.
33. Blumsohn A, Hannon RA, Al-Dehaimi AW, and Eastell R. Short-term intraindividual variability of markers of bone turnover in healthy adults. J Bone Min Res 9:S153, 1994.
34. Raisz L, Wiita B, Artis A, Bowen A, Schwartz S, Trahiotis M, Shoukri K and Smith J. Comparison of the effects of estrogen alone and estrogen plus androgen on biochemical markers of bone formation and resorption in postmenopausal women. J Clin Endocrin Metab 81:37-43, 1996.
35. Kanis J. Treatment of generalized osteoporosis with inhibitors of bone resorption, Osteoporosis. London: Blackwell Science Ltd.; pp. 168-195, 1994.
36. Kennedy J, Carey R, Coolen R, Garber C, Hartman A, Lee H, Leitz V, Levine J, McLean M, Osberg I, Steindel S, and Sylvester E. User Evaluation of Precision Performance of Clinical Chemistry Devices. Tentative Guideline, NCCLS 4:8 EP5-T2, 1984.
37. Moss D and Whitby L. A simplified heat-inactivation method for investigating alkaline phosphatase isoenzymes in serum. Clin Chem Acta 61:63, 1975.
38. Price C, Mitchell C, Moriarty J, Gray M and Noonan K. Mass versus activity of an immunometric assay for bone alkaline phosphatase in serum. Ann Clin Biochem 32:405, 1995.

 EXP	<p>GB <i>Use By</i> DE <i>Verwendbar bis</i> ES <i>Fecha de caducidad</i> IT <i>Utilizzare entro</i> FR <i>Utiliser jusque</i> NL <i>Houdbaar tot</i> DK <i>Holdbar til</i> CZ <i>Použitelné do</i> SK <i>Použitelné do</i> GR <i>Ημερομηνία λήξης</i> PT <i>Prazo de validade</i> HU <i>Felhasználható</i> SE <i>Använd före</i> PL <i>Użyć przed</i></p>		<p>GB Batch code DE Chargenbezeichnung ES Código de lote IT Codice del lotto FR Code du lot NL Lot nummer DK Lotnummer CZ Číslo šarže SK Číslo šarže GR Αριθμός Παρτίδας PT Código do lote HU Szarzsorszám SE Lot nummer PL Kod partii</p>
	<p>GB <i>Catalogue number</i> DE <i>Bestellnummer</i> ES <i>Número de catálogo</i> IT <i>Numero di catalogo</i> FR <i>Référence du catalogue</i> NL <i>Catalogus nummer</i> DK <i>Katalognummer</i> CZ <i>Katalogové číslo</i> SK <i>Katalogové číslo</i> GR <i>Αριθμός καταλόγου</i> PT <i>Referência de catálogo</i> HU <i>Katalógusszám</i> SE <i>Katalognummer</i> PL <i>Numer katalogowy</i></p>		<p>GB Manufacturer DE Hersteller ES Fabricante IT Fabbricante FR Fabricant NL Fabrikant DK Producent CZ Výrobce SK Výrobca GR Κατασκευαστής PT Fabricante HU Gyártó SE Tillverkare PL Producent</p>
	<p>GB <i>Contains sufficient for <n> tests</i> DE <i>Inhalt ausreichend für <n> Prüfungen</i> ES <i>Contenido suficiente para <n> ensayos</i> IT <i>Contenuto sufficiente per "n" saggi</i> FR <i>Contenu suffisant pour "n" tests</i> NL <i>Inhoud voldoende voor "n" testen</i> DK <i>Indeholder tilstrækkeligt til "n" test</i> CZ <i>Lze použít pro <n> testů</i> SK <i>Obsah postačuje na <n> stanovení</i> GR <i>Περιεχόμενο επαρκές για «n» εξετάσεις</i> PT <i>Conteúdo suficiente para "n" ensaios</i> HU <i>A doboz tartalma <n> vizsgálat elvégzéséhez elegendő</i></p> <p>SE <i>Räcker till "n" antal tester</i> PL <i>Wystarczy na wykonanie <n> testów</i></p>		<p>GB In Vitro Diagnostic Medical Device DE In-Vitro-Diagnostikum ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro FR Dispositif médical de diagnostic in vitro NL Medisch hulpmiddel voor in-vitro diagnostiek DK Medicinsk udstyr til in vitro-diagnostik CZ In Vitro diagnostický zdravotnický prostředek SK Zdravotnícka pomôcka in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro HU In vitro diagnosztikum SE Medicintekniska produkter för in vitro diagnostik PL Wyrób do diagnostyki In Vitro</p>
	<p>GB <i>Temperature limitation</i> DE <i>Temperaturbegrenzung</i> ES <i>Límite de temperatura</i> IT <i>Limiti di temperatura</i> FR <i>Limites de température</i> NL <i>Temperatuurlimiet</i> DK <i>Temperaturbegrænsning</i> CZ <i>Teplotní rozmezí od do</i> SK <i>Teplotné rozmedzie od do</i> GR <i>Περιορισμοί θερμοκρασίας</i> PT <i>Limites de temperatura</i> HU <i>Hőmérsékletartomány</i> SE <i>Temperaturbegränsning</i> PL <i>Przestrzeżać zakresu temperatury</i></p>		<p>GB Consult Instructions for Use DE Gebrauchsanweisung beachten ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso FR Consulter les instructions d'utilisation NL Raadpleeg de gebruiksaanwijzing DK Se brugsanvisning CZ Viz návod k použití SK Vid' návod na použitie GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização HU Nézze meg a Használati utasítást SE Se handhavandebeskrivningen PL Sprawdź w instrukcji obsługi</p>

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
REFERENCES/BIBLIOGRAPHIE/LITERATUR

1. Riggs L. Causes of age-related bone loss and fractures. Osteoporosis: Physiological Basis, Assessment, and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 7-16, 1990.
2. Status of Research in Osteoporosis. Osteoporosis Research, Education and Health Promotion. U.S. Department of Health and Human Services, Public Health Service, National Institutes of Health.
3. Canalis E. Regulation of bone remodeling, Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. Favus MJ, ed., Raven Press (New York), pp. 33-37, 1993.
4. Kleerekoper M and Avioli LV. Evaluation and treatment of postmenopausal osteoporosis, Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. Favus MJ, ed., Raven Press (New York), pp. 223-229, 1993.
5. Lindsay R. Osteoporosis. A Guide to Diagnosis, Prevention, and Treatment. Raven Press (New York), pp. 1-39, 1992.
6. Epstein S. Serum and urinary markers of bone remodeling: assessment of bone turnover. Endocrine Reviews 9:437-448, 1988.
7. Whyte MP. Alkaline phosphatase and the measurement of bone formation. Clinical Disorders of Bone and Mineral Metabolism. Potts JT and Frame B, eds., pp. 120-125, 1983.
8. Gamaro P and Delmas P. Assessment of the serum levels of bone alkaline phosphatase with a new immunoradiometric assay in patients with metabolic bone disease. J Clin Endocrinol Metab 77:1046-1053, 1993.
9. Cooper EH, Whelan P, and Purves D. Bone alkaline phosphatase and prostate-specific antigen in the monitoring of prostate cancer. The Prostate 25:236-242, 1994.
10. Van Straalen JP, Sanders E, Prummel MF, and Sanders GTB. Bone alkaline phosphatase as indicator of bone formation. Clin Chim Acta 201:27-33, 1991.
11. Merkow RL and Lane JM. Metabolic bone disease and Paget's disease in the elderly: II. Paget's Disease. Clin Rheum Dis 12:70, 1986.
12. Lindsay R. Sex steroids in the pathogenesis and prevention of osteoporosis, Osteoporosis: Etiology, Diagnosis, and Management. Riggs L and Melton LJ, eds., Raven Press (New York), pp. 333-358, 1988.
13. Peck W. The pathogenesis of postmenopausal osteoporosis, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 3-6, 1990.
14. Christiansen C. Osteoporosis: Diagnosis and management today. Bone 17:513S-516S, 1995.
15. Turner R, Riggs B and Spelsberg T. Skeletal effects of estrogen. Endocrine Reviews 15:275-300, 1994.
16. Kanis J. Pathogenesis of osteoporosis and fracture, Osteoporosis. London: Blackwell Science Ltd.; pp. 22-55, 1994.
17. Mack TM and Ross RK. A current perception of HRT risks and benefits, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 161-178, 1990.
18. Riggs L and Melton LJ. The prevention and treatment of osteoporosis. N Engl J Med 327:620-627, 1992.
19. Guidelines for Preclinical and Clinical Evaluation of Agents Used in the Prevention or Treatment of Postmenopausal Osteoporosis. Division of Metabolism and Endocrine Drug Products, Food and Drug Administration. April, 1994.
20. Ribot C, Tremolieres F and Pouilles JM. Prevention of postmenopausal bone loss by long-term parenteral administration of 17β estradiol: comparison of percutaneous and transdermal route, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 137-145, 1990.
21. Delmas P. Biochemical markers of bone turnover. Methodology and clinical use in osteoporosis. Am J Med 91:59S-63S, 1991.
22. Duda RJ, O'Brien JF, Katzman JA, Peterson JM, Mann KG, and Riggs BL. Concurrent assays of circulating bone gla-protein and bone alkaline phosphatase: effects of sex, age, and metabolic bone disease. J Clin Endocrinol Metab 66:951-957, 1988.
23. HHS Publication No. (CDC) 93-8395. 3rd edition. May 1993. Biosafety in Microbiological and Biomedical Laboratories. Washington, DC: US Government Printing Office, 1993.
24. Manual Guide – Safety Management No. CDC-22. Decontamination of Laboratory Sink Drains to Remove Azide Salts. April 30, 1976. Atlanta, GA: Centers for Disease Control.
25. Hill CS and Wolfert RL. The preparation of monoclonal antibodies which react preferentially with human bone alkaline phosphatase and not liver alkaline phosphatase. Clin Chim Acta 186:315-320, 1989.
26. Gamaro P, Shih WJ, Gineyts E, Karpf DB, and Delmas PD. Comparison of new biochemical markers of bone turnover in late postmenopausal osteoporotic women in response to alendronate treatment. J Clin Endocrinol Metab 79:1693-1700, 1994.
27. Liberman JA, Weiss SR, Broll J, Minne HW, Quan H, Bell NH, Rodriguez-Portales J, Downs RW, Dequeker J, Favus M, Seeman E, Recker RR, Capizzini T, Santora AC, Lombardi A, Shah RV, Hirsch LJ and Karpf DB. For the Alendronate Phase III Osteoporosis Treatment Study Group. Effect of oral alendronate on bone mineral density and the incidence of fractures in postmenopausal osteoporosis. N Engl J Med 333:1437-1443, 1995.
28. Dawson-Hughes B, Dallal GE, Krall EA, Sadowski L, Sahyoun N and Tannenbaum S. A controlled trial of the effect of calcium supplementation on bone density in postmenopausal women. N Engl J Med 323:878-883, 1990.
29. Reid IR, Ames RW, Evans MC, Gamble GD, and Sharpe SJ. Long-term effects of calcium supplementation on bone loss and fractures in postmenopausal women: a randomized controlled trial. Am J Med 98:331-335, 1995.

30. Devogelaer JP, Broll H, Correa-Rotter R, Cumming DC, De Deuxchaisnes CN, Geusens P, Hosking D, Jaeger P, Kaufman JM, Leite M, Leon J, Liberman U, Menkes CJ, Meunier PJ, Reid I, Rodriguez J, Romanowicz A, Seeman E, Vermeulen A, Hirsch LJ, Lombardi A, Plezia K, Santora AC, Yates AJ and Yuan W. Oral alendronate induces progressive increases in bone mass of the spine, hip, and total body over 3 years in postmenopausal women with osteoporosis. *Bone* 18:141-150, 1996.
31. Overgaard K, Alexandersen P, Riis BJ and Christiansen C. Evaluation of a new commercial IRMA for bone-specific alkaline phosphatase during treatment with hormone replacement therapy and calcitonin. *Clin Chem* 42:973-974, 1996.
32. Soletormos G, Schioler V, Nielsen D, Skovsgaard T, and Dombrowsky P. Interpretation of results for tumor markers on the basis of analytical imprecision and biological variation. *Clin Chem* 10:2077- 2083, 1993.
33. Blumsohn A, Hannon RA, Al-Dehaimi AW, and Eastell R. Short-term intraindividual variability of markers of bone turnover in healthy adults. *J Bone Min Res* 9:S153, 1994.
34. Raisz L, Wiita B, Artis A, Bowen A, Schwartz S, Trahiotis M, Shoukri K and Smith J. Comparison of the effects of estrogen alone and estrogen plus androgen on biochemical markers of bone formation and resorption in postmenopausal women. *J Clin Endocrin Metab* 81:37-43, 1996.
35. Kanis J. Treatment of generalized osteoporosis with inhibitors of bone resorption, Osteoporosis. London: Blackwell Science Ltd.; pp. 168-195, 1994.
36. Kennedy J, Carey R, Coolen R, Garber C, Hartman A, Lee H, Leitz V, Levine J, McLean M, Osberg I, Steindel S, and Sylvester E. User Evaluation of Precision Performance of Clinical Chemistry Devices. Tentative Guideline, NCCLS 4:8 EP5-T2, 1984.
37. Moss D and Whitby L. A simplified heat-inactivation method for investigating alkaline phosphatase isoenzymes in serum. *Clin Chem Acta* 61:63, 1975.
38. Price C, Mitchell C, Moriarty J, Gray M and Noonan K. Mass versus activity of an immunometric assay for bone alkaline phosphatase in serum. *Ann Clin Biochem* 32:405, 1995.

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REFERENCES/BIBLIOGRAPHIE/LITERATUR

- Riggs L. Causes of age-related bone loss and fractures, Osteoporosis: Physiological Basis, Assessment, and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 7-16, 1990.
- Status of Research in Osteoporosis. Osteoporosis Research, Education and Health Promotion. U.S. Department of Health and Human Services, Public Health Service, National Institutes of Health.
- Canalis E. Regulation of bone remodeling, Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. Favus MJ, ed., Raven Press (New York), pp. 33-37, 1993.
- Kleerekoper M and Avioli LV. Evaluation and treatment of postmenopausal osteoporosis, Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. Favus MJ, ed., Raven Press (New York), pp. 223-229, 1993.
- Lindsay R. Osteoporosis. A Guide to Diagnosis, Prevention, and Treatment. Raven Press (New York), pp. 1-39, 1992.
- Epstein S. Serum and urinary markers of bone remodeling: assessment of bone turnover. Endocrine Reviews 9:437-448, 1988.
- Whyte MP. Alkaline phosphatase and the measurement of bone formation, Clinical Disorders of Bone and Mineral Metabolism. Potts JT and Frame B, eds., pp. 120-125, 1983.
- Garnero P and Delmas P. Assessment of the serum levels of bone alkaline phosphatase with a new immunoradiometric assay in patients with metabolic bone disease. J Clin Endocrin Metab 77:1046-1053, 1993.
- Cooper EH, Whelan P, and Purves D. Bone alkaline phosphatase and prostate-specific antigen in the monitoring of prostate cancer. The Prostate 25:236-242, 1994.
- Van Straalen JP, Sanders E, Prummel MF, and Sanders GTB. Bone alkaline phosphatase as indicator of bone formation. Clin Chim Acta 201:27-33, 1991.
- Merkow RL and Lane JM. Metabolic bone disease and Paget's disease in the elderly: II. Paget's Disease. Clin Rheum Dis 12:70, 1986.
- Lindsay R. Sex steroids in the pathogenesis and prevention of osteoporosis, Osteoporosis: Etiology, Diagnosis, and Management. Riggs L and Melton LJ, eds., Raven Press (New York), pp. 333-358, 1988.
- Peck W. The pathogenesis of postmenopausal osteoporosis, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 3-6, 1990.
- Christiansen C. Osteoporosis: Diagnosis and management today. Bone 17:513S-516S, 1995.
- Turner R, Riggs B and Spelsberg T. Skeletal effects of estrogen. Endocrine Reviews 15:275-300, 1994.
- Kanis J. Pathogenesis of osteoporosis and fracture, Osteoporosis. London: Blackwell Science Ltd.; pp. 22-55, 1994.
- Mack TM and Ross RK. A current perception of HRT risks and benefits, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 161-178, 1990.
- Riggs L and Melton LJ. The prevention and treatment of osteoporosis. N Engl J Med 327:620-627, 1992.
- Guidelines for Preclinical and Clinical Evaluation of Agents Used in the Prevention or Treatment of Postmenopausal Osteoporosis. Division of Metabolism and Endocrine Drug Products, Food and Drug Administration. April, 1994.
- Ribot C, Tremollieres F and Pouilles JM. Prevention of postmenopausal bone loss by long-term parenteral administration of 17B estradiol: comparison of percutaneous and transdermal route, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 137-145, 1990.
- Delmas P. Biochemical markers of bone turnover. Methodology and clinical use in osteoporosis. Am J Med 91:59S-63S, 1991.
- Duda RJ, O'Brien JF, Katzman JA, Peterson JM, Mann KG, and Riggs BL. Concurrent assays of circulating bone gla-protein and bone alkaline phosphatase: effects of sex, age, and metabolic bone disease. J Clin Endocrin Metab 66:951-957, 1988.
- HHS Publication No. (CDC) 93-8395. 3rd edition. May 1993. Biosafety in Microbiological and Biomedical Laboratories. Washington, DC: US Government Printing Office, 1993.
- Manual Guide - Safety Management No. CDC-22. Decontamination of Laboratory Sink Drains to Remove Azide Salts. April 30, 1976. Atlanta, GA: Centers for Disease Control.
- Hill CS and Wolfert RL. The preparation of monoclonal antibodies which react preferentially with human bone alkaline phosphatase and not liver alkaline phosphatase. Clin Chim Acta 186:315-320, 1989.
- Garnero P, Shih WJ, Gineyts E, Karpf DB, and Delmas PD. Comparison of new biochemical markers of bone turnover in late postmenopausal osteoporotic women in response to alendronate treatment. J Clin Endocrin Metab 79:1693-1700, 1994.
- Lieberman UA, Weiss SR, Broll J, Minne HW, Quan H, Bell NH, Rodriguez-Portales J, Downs RW, Dequeker J, Favus M, Seeman E, Recker RR, Capizzi T, Santora AC, Lombardi A, Shah RV, Hirsch LJ and Karpf DB. For the Alendronate Phase III Osteoporosis Treatment Study Group. Effect of oral alendronate on bone mineral density and the incidence of fractures in postmenopausal osteoporosis. N Engl J Med 333:1437-1443, 1995.
- Dawson-Hughes B, Dallal GE, Krall EA, Sadowski L, Sahyoun N and Tannenbaum S. A controlled trial of the effect of calcium supplementation on bone density in postmenopausal women. N Engl J Med 323:878-883, 1990.
- Reid IR, Ames RW, Evans MC, Gamble GD, and Sharpe SJ. Long-term effects of calcium supplementation on bone loss and fractures in postmenopausal women: a randomized controlled trial. Am J Med 98:331-335, 1995.
- Devogelaer JP, Broll H, Correa-Rotter R, Cumming DC, De Deuxchaisnes CN, Geusens P, Hosking D, Jaeger P, Kaufman JM, Leite M, Leon J, Lieberman U, Menkes CJ, Meunier PJ, Reid I, Rodriguez J, Romanowicz A, Seeman E, Vermeulen A, Hirsch LJ, Lombardi A, Plezia K, Santora AC, Yates AJ and Yuan W. Oral alendronate induces progressive increases in bone mass of the spine, hip, and total body over 3 years in postmenopausal women with osteoporosis. Bone 18:141-150, 1996.
- Overgaard K, Alexandersen P, Riis BJ and Christiansen C. Evaluation of a new commercial IRMA for bone-specific alkaline phosphatase during treatment with hormone replacement therapy and calcitonin. Clin Chem 42:973-974, 1996.
- Soletormos G, Schioler V, Nielsen D, Skovsgaard T, and Dombrowsky P. Interpretation of results for tumor markers on the basis of analytical imprecision and biological variation. Clin Chem 10:2077-2083, 1993.
- Blumsohn A, Hannon RA, Al-Dehaimi AW, and Eastell R. Short-term intraindividual variability of markers of bone turnover in healthy adults. J Bone Min Res 9:S153, 1994.
- Raisz L, Wiita B, Artis A, Bowen A, Schwartz S, Trahiotis M, Shoukri K and Smith J. Comparison of the effects of estrogen alone and estrogen plus androgen on biochemical markers of bone formation and resorption in postmenopausal women. J Clin Endocrin Metab 81:37-43, 1996.
- Kanis J. Treatment of generalized osteoporosis with inhibitors of bone resorption, Osteoporosis. London: Blackwell Science Ltd.; pp. 168-195, 1994.
- Kennedy J, Carey R, Coolen R, Garber C, Hartman A, Lee H, Leitz V, Levine J, McLean M, Osberg I, Steindel S, and Sylvester E. User Evaluation of Precision Performance of Clinical Chemistry Devices. Tentative Guideline, NCCLS 4:8 EP5-T2, 1984.
- Moss D and Whitby L. A simplified heat-inactivation method for investigating alkaline phosphatase isoenzymes in serum. Clin Chem Acta 61:63, 1975.
- Price C, Mitchell C, Moriarty J, Gray M and Noonan K. Mass versus activity of an immunometric assay for bone alkaline phosphatase in serum. Ann Clin Biochem 32:405, 1995.

 EXP	GB <i>Use By</i> DE <i>Verwendbar bis</i> ES <i>Fecha de caducidad</i> IT <i>Utilizzare entro</i> FR <i>Utiliser jusque</i> NL <i>Houdbaar tot</i> DK <i>Holdbar til</i> CZ <i>Použitelné do</i> SK <i>Použitelné do</i> GR <i>Ημερομηνία λήξης</i> PT <i>Prazo de validade</i> HU <i>Felhasználható</i> SE <i>Använd före</i> PL <i>Użyty przed</i>		GB Batch code DE Chargenbezeichnung ES Código de lote IT Codice del lotto FR Code du lot NL Lot nummer DK Lotnummer CZ Číslo šarže SK Číslo šarže GR Αριθμός Παρτίδας PT Código do lote HU Szarzsorszám SE Lot nummer PL Kod partii
 REF	GB <i>Catalogue number</i> DE <i>Bestellnummer</i> ES <i>Número de catálogo</i> IT <i>Numero di catalogo</i> FR <i>Référence du catalogue</i> NL <i>Catalogus nummer</i> DK <i>Katalognummer</i> CZ <i>Katalogové číslo</i> SK <i>Katalogové číslo</i> GR <i>Αριθμός καταλόγου</i> PT <i>Referência de catálogo</i> HU <i>Katalógusszám</i> SE <i>Katalognummer</i> PL <i>Numer katalogowy</i>		GB Manufacturer DE Hersteller ES Fabricante IT Fabbicante FR Fabricant NL Fabrikant DK Producent CZ Výrobce SK Výrobca GR Κατασκευαστής PT Fabricante HU Gyártó SE Tillverkare PL Producent
	GB <i>Contains sufficient for <n> tests</i> DE <i>Inhalt ausreichend für <n> Prüfungen</i> ES <i>Contenido suficiente para <n> ensayos</i> IT <i>Contenuto sufficiente per "n" saggi</i> FR <i>Contenu suffisant pour "n" tests</i> NL <i>Inhoud voldoende voor "n" testen</i> DK <i>Indeholder tilstrækkeligt til "n" test</i> CZ <i>Lze použít pro <n> testů</i> SK <i>Obsah postačuje na <n> stanovení</i> GR <i>Περιεχόμενο επαρκές για «n» εξετάσεις</i> PT <i>Conteúdo suficiente para "n" ensaios</i> HU <i>A doboz tartalma <n> vizsgálat elvégzéséhez elegendő</i> SE <i>Räcker till "n" antal tester</i> PL <i>Wystarczy na wykonanie <n> testów</i>		GB <i>In Vitro Diagnostic Medical Device</i> DE <i>In-Vitro-Diagnostikum</i> ES <i>Producto sanitario para diagnóstico in vitro</i> IT <i>Dispositivo medico-diagnostico in vitro</i> FR <i>Dispositif médical de diagnostic in vitro</i> NL <i>Medisch hulpmiddel voor in-vitro diagnostiek</i> DK <i>Medicinsk udstyr til in vitro-diagnostik</i> CZ <i>In Vitro diagnostický zdravotnický prostředek</i> SK <i>Zdravotnícka pomocka in vitro</i> GR <i>In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν</i> PT <i>Dispositivo médico para diagnóstico in vitro</i> HU <i>In vitro diagnosztikum</i> SE <i>Medicintekniska produkter för in vitro diagnostik</i> PL <i>Wyrob do diagnostyki In Vitro</i>
	GB <i>Temperature limitation</i> DE <i>Temperaturbegrenzung</i> ES <i>Límite de temperatura</i> IT <i>Limiti di temperatura</i> FR <i>Limites de temperature</i> NL <i>Temperatuurlimiet</i> DK <i>Temperaturbegrænsning</i> CZ <i>Teplotní rozmezí od do</i> SK <i>Teplotné rozmedzie od do</i> GR <i>Περιορισμοί θερμοκρασίας</i> PT <i>Limites de temperatura</i> HU <i>Hőmérsékletkorlát</i> SE <i>Temperaturbegränsning</i> PL <i>Przestrzeżenie zakresu temperatury</i>		GB Consult Instructions for Use DE Gebrauchsanweisung beachten ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso FR Consulter les instructions d'utilisation NL Raadpleeg de gebruiksaanwijzing DK Se brugsanvisning CZ Viz návod k použití SK Viz návod na použitie GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização HU Nézze meg a Használati utasítást SE Se handhavandebeskrivningen PL Sprawdź w instrukcji obsługi

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 Immunodiagnostic Systems Limited,
 10 Didcot Way, Boldon Business Park,
 Boldon, Tyne & Wear, NE35 9PD, UK
 Tel.: +44 191 519 0660
 Fax: +44 191 519 0760
 E-Mail: info.uk@idsplc.com
www.idsplc.com

IDS Immunodiagnostic Systems GmbH
 Herriotstrasse 1
 60528 Frankfurt am Main
 Germany

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REFERENCES/BIBLIOGRAPHIE/LITERATUR

1. Riggs L. Causes of age-related bone loss and fractures. Osteoporosis: Physiological Basis, Assessment, and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 7-16, 1990.
2. Status of Research in Osteoporosis. Osteoporosis Research, Education and Health Promotion. U.S. Department of Health and Human Services, Public Health Service, National Institutes of Health.
3. Canalis E. Regulation of bone remodeling. Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. Favus MJ, ed., Raven Press (New York), pp. 33-37, 1993.
4. Kleerekoper M and Avioli LV. Evaluation and treatment of postmenopausal osteoporosis. Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. Favus MJ, ed., Raven Press (New York), pp. 223-229, 1993.
5. Lindsay R. Osteoporosis. A Guide to Diagnosis, Prevention, and Treatment. Raven Press (New York), pp. 1-39, 1992.
6. Epstein S. Serum and urinary markers of bone remodeling: assessment of bone turnover. Endocrine Reviews 9:437-448, 1988.
7. Whyte MP. Alkaline phosphatase and the measurement of bone formation. Clinical Disorders of Bone and Mineral Metabolism. Potts JT and Frame B, eds., pp. 120-125, 1983.
8. Gamero P and Delmas P. Assessment of the serum levels of bone alkaline phosphatase with a new immunoradiometric assay in patients with metabolic bone disease. J Clin Endocrin Metab 77:1046- 1053, 1993.
9. Cooper EH, Whelan P, and Purves D. Bone alkaline phosphatase and prostate-specific antigen in the monitoring of prostate cancer. The Prostate 25:236-242, 1994.
10. Van Straalen JP, Sanders E, Prummel MF, and Sanders GTB. Bone alkaline phosphatase as indicator of bone formation. Clin Chim Acta 201:27-33, 1991.
11. Merkow RL and Lane JM. Metabolic bone disease and Paget's disease in the elderly: II. Paget's Disease. Clin Rheum Dis 12:70, 1986.
12. Lindsay R. Sex steroids in the pathogenesis and prevention of osteoporosis. Osteoporosis: Etiology, Diagnosis, and Management. Riggs L and Melton LJ, eds., Raven Press (New York), pp. 333-358, 1988.
13. Peck W. The pathogenesis of postmenopausal osteoporosis. Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 3-6, 1990.
14. Christiansen C. Osteoporosis: Diagnosis and management today. Bone 17:513S-516S, 1995.
15. Turner R, Riggs B and Spelsberg T. Skeletal effects of estrogen. Endocrine Reviews 15:275-300, 1994.
16. Kanis J. Pathogenesis of osteoporosis and fracture. Osteoporosis. London: Blackwell Science Ltd.; pp. 22-55, 1994.
17. Mack TM and Ross RK. A current perception of HRT risks and benefits. Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 161-178, 1990.

18. Riggs L and Melton LJ. The prevention and treatment of osteoporosis. *N Engl J Med* 327:620-627, 1992.
19. Guidelines for Preclinical and Clinical Evaluation of Agents Used in the Prevention or Treatment of Postmenopausal Osteoporosis. Division of Metabolism and Endocrine Drug Products, Food and Drug Administration. April, 1994.
20. Ribot C, Tremolieres F and Pouilles JM. Prevention of postmenopausal bone loss by long-term parenteral administration of 17 β estradiol: comparison of percutaneous and transdermal route. *Osteoporosis: Physiological Basis, Assessment and Treatment*. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 137-145, 1990.
21. Delmas P. Biochemical markers of bone turnover. Methodology and clinical use in osteoporosis. *Am J Med* 91:59S-63S, 1991.
22. Duda RJ, O'Brien JF, Katzman JA, Peterson JM, Mann KG, and Riggs BL. Concurrent assays of circulating bone gla-protein and bone alkaline phosphatase: effects of sex, age, and metabolic bone disease. *J Clin Endocrin Metab* 66:951-957, 1988.
23. HHS Publication No. (CDC) 93-8395. 3rd edition. May 1993. Biosafety in Microbiological and Biomedical Laboratories. Washington, DC: US Government Printing Office, 1993.
24. Manual Guide – Safety Management No. CDC-22. Decontamination of Laboratory Sink Drains to Remove Azide Salts. April 30, 1976. Atlanta, GA: Centers for Disease Control.
25. Hill CS and Wolfert RL. The preparation of monoclonal antibodies which react preferentially with human bone alkaline phosphatase and not liver alkaline phosphatase. *Clin Chim Acta* 186:315-320, 1989.
26. Gamero P, Shih WJ, Gineyts E, Karpf DB, and Delmas PD. Comparison of new biochemical markers of bone turnover in late postmenopausal osteoporotic women in response to alendronate treatment. *J Clin Endocrin Metab* 79:1693-1700, 1994.
27. Liberman UA, Weiss SR, Broll J, Minne HW, Quan H, Bell NH, Rodriguez-Portales J, Downs RW, Dequeker J, Favus M, Seeman E, Recker RR, Capizzi T, Santora AC, Lombardi A, Shah RV, Hirsch LJ and Karpf DB. For the Alendronate Phase III Osteoporosis Treatment Study Group. Effect of oral alendronate on bone mineral density and the incidence of fractures in postmenopausal osteoporosis. *N Engl J Med* 333:1437-1443, 1995.
28. Dawson-Hughes B, Dallal GE, Krall EA, Sadowski L, Sahyoun N and Tannenbaum S. A controlled trial of the effect of calcium supplementation on bone density in postmenopausal women. *N Engl J Med* 323:878-883, 1990.
29. Reid IR, Ames RW, Evans MC, Gamble GD, and Sharpe SJ. Long-term effects of calcium supplementation on bone loss and fractures in postmenopausal women: a randomized controlled trial. *Am J Med* 98:331-335, 1995.
30. Devogelaer JP, Broll H, Correa-Rotter R, Cumming DC, De Deuxchaisnes CN, Geusens P, Hosking D, Jaeger P, Kaufman JM, Leite M, Leon J, Liberman U, Menkes CJ, Meunier PJ, Reid I, Rodriguez J, Romanowicz A, Seeman E, Vermeulen A, Hirsch LJ, Lombardi A, Plezia K, Santora AC, Yates AJ and Yuan W. Oral alendronate induces progressive increases in bone mass of the spine, hip, and total body over 3 years in postmenopausal women with osteoporosis. *Bone* 18:141-150, 1996.
31. Overgaard K, Alexandersen P, Riis BJ and Christiansen C. Evaluation of a new commercial IRMA for bone-specific alkaline phosphatase during treatment with hormone replacement therapy and calcitonin. *Clin Chem* 42:973-974, 1996.
32. Soletormos G, Schioler V, Nielsen D, Skovsgaard T, and Dombrowsky P. Interpretation of results for tumor markers on the basis of analytical imprecision and biological variation. *Clin Chem* 10:2077- 2083, 1993.
33. Blumsohn A, Hannon RA, Al-Dehaimi AW, and Eastell R. Short-term intraindividual variability of markers of bone turnover in healthy adults. *J Bone Min Res* 9:S153, 1994.
34. Raisz L, Wiita B, Artis A, Bowen A, Schwartz S, Trahiotis M, Shoukri K and Smith J. Comparison of the effects of estrogen alone and estrogen plus androgen on biochemical markers of bone formation and resorption in postmenopausal women. *J Clin Endocrin Metab* 81:37-43, 1996.
35. Kanis J. Treatment of generalized osteoporosis with inhibitors of bone resorption. *Osteoporosis*. London: Blackwell Science Ltd., pp. 168-195, 1994.
36. Kennedy J, Carey R, Coolen R, Garber C, Hartman A, Lee H, Leitz V, Levine J, McLean M, Osberg I, Stendel S, and Sylvester E. User Evaluation of Precision Performance of Clinical Chemistry Devices. Tentative Guideline, NCCLS 4:8 EP5-T2, 1984.
37. Moss D and Whitby L. A simplified heat-inactivation method for investigating alkaline phosphatase isoenzymes in serum. *Clin Chem Acta* 61:63, 1975.
38. Price C, Mitchell C, Moriarty J, Gray M and Noonan K. Mass versus activity of an immunometric assay for bone alkaline phosphatase in serum. *Ann Clin Biochem* 32:405, 1995.

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- Riggs L. Causes of age-related bone loss and fractures, Osteoporosis: Physiological Basis, Assessment, and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 7-16, 1990.
- Status of Research in Osteoporosis. Osteoporosis Research, Education and Health Promotion. U.S. Department of Health and Human Services, Public Health Service, National Institutes of Health.
- Canalis E. Regulation of bone remodeling, Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. Favus MJ, ed., Raven Press (New York), pp. 33-37, 1993.
- Kleerekopper M and Avioli LV. Evaluation and treatment of postmenopausal osteoporosis, Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. Favus MJ, ed., Raven Press (New York), pp. 223-229, 1993.
- Lindsay R. Osteoporosis. A Guide to Diagnosis, Prevention, and Treatment. Raven Press (New York), pp. 1-39, 1992.
- Epstein S. Serum and urinary markers of bone remodeling: assessment of bone turnover. Endocrine Reviews 9:437-448, 1988.
- Whyte MP. Alkaline phosphatase and the measurement of bone formation, Clinical Disorders of Bone and Mineral Metabolism. Potts JT and Frame B, eds., pp. 120-125, 1983.
- Garnero P and Delmas P. Assessment of the serum levels of bone alkaline phosphatase with a new immunoradiometric assay in patients with metabolic bone disease. J Clin Endocrin Metab 77:1046-1053, 1993.
- Cooper EH, Whelan R, and Purves D. Bone alkaline phosphatase and prostate-specific antigen in the monitoring of prostate cancer. The Prostate 25:236-242, 1994.
- Van Straalen JP, Sanders E, Prummel MF, and Sanders GTB. Bone alkaline phosphatase as indicator of bone formation. Clin Chim Acta 201:27-33, 1991.
- Merkow RL and Lane JM. Metabolic bone disease and Paget's disease in the elderly: II. Paget's Disease. Clin Rheum Dis 12:70, 1986.
- Lindsay R. Sex steroids in the pathogenesis and prevention of osteoporosis, Osteoporosis: Etiology, Diagnosis, and Management. Riggs L and Melton LJ, eds., Raven Press (New York), pp. 333-358, 1988.
- Peck W. The pathogenesis of postmenopausal osteoporosis, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 3-6, 1990.
- Christiansen C. Osteoporosis: Diagnosis and management today. Bone 17:513S-516S, 1995.
- Turner R, Riggs B and Spelsberg T. Skeletal effects of estrogen. Endocrine Reviews 15:275-300, 1994.
- Kanis J. Pathogenesis of osteoporosis and fracture, Osteoporosis. London: Blackwell Science Ltd.; pp. 22-55, 1994.
- Mack TM and Ross RK. A current perception of HRT risks and benefits, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 161-178, 1990.
- Riggs L and Melton LJ. The prevention and treatment of osteoporosis. N Engl J Med 327:620-627, 1992.
- Guidelines for Preclinical and Clinical Evaluation of Agents Used in the Prevention or Treatment of Postmenopausal Osteoporosis. Division of Metabolism and Endocrine Drug Products, Food and Drug Administration. April, 1994.
- Ribot C, Tremollieres F and Pouilles JM. Prevention of postmenopausal bone loss by long-term parenteral administration of 17B estradiol: comparison of percutaneous and transdermal route, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 137-145, 1990.
- Delmas P. Biochemical markers of bone turnover. Methodology and clinical use in osteoporosis. Am J Med 91:59S-63S, 1991.
- Duda RJ, O'Brien JF, Katzman JA, Peterson JM, Mann KG, and Riggs BL. Concurrent assays of circulating bone gla-protein and bone alkaline phosphatase: effects of sex, age, and metabolic bone disease. J Clin Endocrin Metab 66:951-957, 1988.
- HHS Publication No. (CDC) 93-8395. 3rd edition. May 1993. Biosafety in Microbiological and Biomedical Laboratories. Washington, DC: US Government Printing Office, 1993.
- Manual Guide - Safety Management No. CDC-22. Decontamination of Laboratory Sink Drains to Remove Azide Salts. April 30, 1976. Atlanta, GA: Centers for Disease Control.
- Hill CS and Wolfert RL. The preparation of monoclonal antibodies which react preferentially with human bone alkaline phosphatase and not liver alkaline phosphatase. Clin Chim Acta 186:315-320, 1989.
- Garnero P, Shih WJ, Gineyts E, Karpf DB, and Delmas PD. Comparison of new biochemical markers of bone turnover in late postmenopausal osteoporotic women in response to alendronate treatment. J Clin Endocrin Metab 79:1693-1700, 1994.
- Liberman JA, Weiss SR, Broll J, Minne HW, Quan H, Bell NH, Rodriguez-Portales J, Downs RW, Dequeker J, Favus M, Seeman E, Recker RR, Capizzi T, Santora AC, Lombardi A, Shah RV, Hirsch LJ and Karpf DB. For the Alendronate Phase III Osteoporosis Treatment Study Group. Effect of oral alendronate on bone mineral density and the incidence of fractures in postmenopausal osteoporosis. N Engl J Med 333:1437-1443, 1995.
- Dawson-Hughes B, Dallal GE, Krall EA, Sadowski L, Sahyoun N and Tannenbaum S. A controlled trial of the effect of calcium supplementation on bone density in postmenopausal women. N Engl J Med 323:878-883, 1990.
- Reid IR, Ames RW, Evans MC, Gamble GD, and Sharpe SJ. Long-term effects of calcium supplementation on bone loss and fractures in postmenopausal women: a randomized controlled trial. Am J Med 98:331-335, 1995.
- Devogelaer JP, Broll H, Correa-Rotter R, Cumming DC, De Deuxchaisnes CN, Geusens P, Hosking D, Jaeger P, Kaufman JM, Leite M, Leon J, Liberman U, Menkes CJ, Meunier PJ, Reid I, Rodriguez J, Romanowicz A, Seeman E, Vermeulen A, Hirsch LJ, Lombardi A, Plezia K, Santora AC, Yates AJ and Yuan W. Oral alendronate induces progressive increases in bone mass of the spine, hip, and total body over 3 years in postmenopausal women with osteoporosis. Bone 18:141-150, 1996.
- Overgaard K, Alexandersen P, Riis BJ and Christiansen C. Evaluation of a new commercial IRMA for bone-specific alkaline phosphatase during treatment with hormone replacement therapy and calcitonin. Clin Chem 42:973-974, 1996.
- Soletormos G, Schioler V, Nielsen D, Skovsgaard T, and Dombernowsky P. Interpretation of results for tumor markers on the basis of analytical imprecision and biological variation. Clin Chem 10:2077-2083, 1993.
- Blumsohn A, Hannon RA, Al-Dehaimi AW, and Eastell R. Short-term intraindividual variability of markers of bone turnover in healthy adults. J Bone Min Res 9:S153, 1994.
- Raisz L, Wiita B, Artis A, Bowen A, Schwartz S, Trahiotis M, Shoukri K and Smith J. Comparison of the effects of estrogen alone and estrogen plus androgen on biochemical markers of bone formation and resorption in postmenopausal women. J Clin Endocrin Metab 81:37-43, 1996.
- Kanis J. Treatment of generalized osteoporosis with inhibitors of bone resorption, Osteoporosis. Lon- don: Blackwell Science Ltd.; pp. 168-195, 1994.
- Kennedy J, Carey R, Coolen R, Garber C, Hartman A, Lee H, Leitz V, Levine J, McLean M, Osberg I, Steindel S, and Sylvester E. User Evaluation of Precision Performance of Clinical Chemistry Devices. Tentative Guideline, NCCLS 4:8 EP5-T2, 1984.
- Moss D and Whitby L. A simplified heat-inactivation method for investigating alkaline phosphatase isoenzymes in serum. Clin Chem Acta 61:63, 1975.
- Price C, Mitchell C, Moriarty J, Gray M and Noonan K. Mass versus activity of an immunometric assay for bone alkaline phosphatase in serum. Ann Clin Biochem 32:405, 1995.

 EXP	<p>GB <i>Use By</i> DE <i>Verwendbar bis</i> ES <i>Fecha de caducidad</i> IT <i>Utilizzare entro</i> FR <i>Utiliser jusque</i> NL <i>Houdbaar tot</i> DK <i>Holdbar til</i> CZ <i>Použitelné do</i> SK <i>Použitelné do</i> GR <i>Ημερομηνία λήξης</i> PT <i>Prazo de validade</i> HU <i>Felhasználható</i> SE <i>Använd före</i> PL <i>Użyty przed</i></p>		<p>GB <i>Batch code</i> DE <i>Chargenbezeichnung</i> ES <i>Código de lote</i> IT <i>Codice del lotto</i> FR <i>Code du lot</i> NL <i>Lot nummer</i> DK <i>Lotnummer</i> CZ <i>Číslo šarže</i> SK <i>Číslo šarže</i> GR <i>Αριθμός Παρτίδας</i> PT <i>Código do lote</i> HU <i>Sarzszzám</i> SE <i>Lot nummer</i> PL <i>Kod partii</i></p>
	<p>GB <i>Catalogue number</i> DE <i>Bestellnummer</i> ES <i>Número de catálogo</i> IT <i>Numero di catalogo</i> FR <i>Référence du catalogue</i> NL <i>Catalogus nummer</i> DK <i>Katalognummer</i> CZ <i>Katalogové číslo</i> SK <i>Katalógové číslo</i> GR <i>Αριθμός καταλόγου</i> PT <i>Referência de catálogo</i> HU <i>Katalógusszám</i> SE <i>Katalognummer</i> PL <i>Numer katalogowy</i></p>		<p>GB <i>Manufacturer</i> DE <i>Hersteller</i> ES <i>Fabricante</i> IT <i>Fabbricante</i> FR <i>Fabricant</i> NL <i>Fabrikant</i> DK <i>Producent</i> CZ <i>Výrobce</i> SK <i>Výrobca</i> GR <i>Κατασκευαστής</i> PT <i>Fabricante</i> HU <i>Gyártó</i> SE <i>Tillverkare</i> PL <i>Producent</i></p>
	<p>GB <i>Contains sufficient for <n> tests</i> DE <i>Inhalt ausreichend für <n> Prüfungen</i> ES <i>Contenido suficiente para <n> ensayos</i> IT <i>Contenuto sufficiente per "n" saggi</i> FR <i>Contenu suffisant pour "n" tests</i> NL <i>Inhoud voldoende voor "n" testen</i> DK <i>Indeholder tilstrækkeligt til "n" test</i> CZ <i>Lze použít pro <n> testů</i> SK <i>Obsah postačuje na <n> stanovení</i> GR <i>Περιεχόμενο επαρκές για «n» εξετάσεις</i> PT <i>Conteúdo suficiente para "n" ensaios</i> HU <i>A doboz tartalma <n> vizsgálat elvégzéséhez elegendő</i> SE <i>Räcker till "n" antal tester</i> PL <i>Wystarczy na wykonanie <n> testów</i></p>		<p>GB <i>In Vitro Diagnostic Medical Device</i> DE <i>In-Vitro-Diagnostikum</i> ES <i>Producto sanitario para diagnóstico in vitro</i> IT <i>Dispositivo medico-diagnostico in vitro</i> FR <i>Dispositif médical de diagnostic in vitro</i> NL <i>Medisch hulpmiddel voor in-vitro diagnostiek</i> DK <i>Medicinsk udstyr til in vitro-diagnostik</i> CZ <i>In Vitro diagnostický zdravotnický prostředek</i> SK <i>Zdravotnícka pomocka in vitro</i> GR <i>In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν</i> PT <i>Dispositivo médico para diagnóstico in vitro</i> HU <i>In vitro diagnosztikum</i> SE <i>Medicintekniska produkter för in vitro diagnostik</i> PL <i>Wyrób do diagnostyki In Vitro</i></p>
	<p>GB <i>Temperature limitation</i> DE <i>Temperaturbegrenzung</i> ES <i>Límite de temperatura</i> IT <i>Limiti di temperatura</i> FR <i>Limites de température</i> NL <i>Temperatuurlimiet</i> DK <i>Temperaturbegrænsning</i> CZ <i>Teplotní rozmezí od do</i> SK <i>Teplotné rozmedzie od do</i> GR <i>Περιορισμοί θερμοκρασίας</i> PT <i>Limites de temperatura</i> HU <i>Hőmérsékletartomány</i> SE <i>Temperaturbegränsning</i> PL <i>Przestrzeżać zakresu temperatury</i></p>		<p>GB <i>Consult Instructions for Use</i> DE <i>Gebrauchsanweisung beachten</i> ES <i>Consulte las instrucciones de uso</i> IT <i>Consultare le istruzioni per l'uso</i> FR <i>Consulter les instructions d'utilisation</i> NL <i>Raadpleeg de gebruiksaanwijzing</i> DK <i>Se brugsanvisning</i> CZ <i>Viz návod k použití</i> SK <i>Vid' návod na použitie</i> GR <i>Συμβουλευτείτε τις οδηγίες χρήσης</i> PT <i>Consulte as instruções de utilização</i> HU <i>Nézze meg a Használati utasítást</i> SE <i>Se handhavandebeskrivningen</i> PL <i>Sprawdź w instrukcji obsługi</i></p>

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REFERENCES/BIBLIOGRAPHIE/LITERATUR

1. Riggs L. Causes of age-related bone loss and fractures, Osteoporosis: Physiological Basis, Assessment, and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 7-16, 1990.
2. Status of Research in Osteoporosis. Osteoporosis Research, Education and Health Promotion. U.S. Department of Health and Human Services, Public Health Service, National Institutes of Health.
3. Canalis E. Regulation of bone remodeling, Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. Favus MJ, ed., Raven Press (New York), pp. 33-37, 1993.
4. Kleerekoper M and Avioli LV. Evaluation and treatment of postmenopausal osteoporosis, Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. Favus MJ, ed., Raven Press (New York), pp. 223-229, 1993.
5. Lindsay R. Osteoporosis. A Guide to Diagnosis, Prevention, and Treatment. Raven Press (New York), pp. 1-39, 1992.
6. Epstein S. Serum and urinary markers of bone remodeling: assessment of bone turnover. Endocrine Reviews 9:437-448, 1988.
7. Whyte MP. Alkaline phosphatase and the measurement of bone formation. Clinical Disorders of Bone and Mineral Metabolism. Potts JT and Frame B, eds., pp. 120-125, 1983.
8. Garnero P and Delmas P. Assessment of the serum levels of bone alkaline phosphatase with a new immunoradiometric assay in patients with metabolic bone disease. J Clin Endocrin Metab 77:1046-1053, 1993.
9. Cooper EH, Whelan P, and Purves D. Bone alkaline phosphatase and prostate-specific antigen in the monitoring of prostate cancer. The Prostate 25:236-242, 1994.
10. Van Straalen JP, Sanders E, Prummel MF, and Sanders GTB. Bone alkaline phosphatase as indicator of bone formation. Clin Chim Acta 201:27-33, 1991.
11. Merkow RL and Lane JM. Metabolic bone disease and Paget's disease in the elderly: II. Paget's Disease. Clin Rheum Dis 12:70, 1986.
12. Lindsay R. Sex steroids in the pathogenesis and prevention of osteoporosis, Osteoporosis: Etiology, Diagnosis, and Management. Riggs L and Melton LJ, eds., Raven Press (New York), pp. 333-358, 1988.
13. Peck W. The pathogenesis of postmenopausal osteoporosis, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 3-6, 1990.
14. Christiansen C. Osteoporosis: Diagnosis and management today. Bone 17:513S-516S, 1995.
15. Turner R, Riggs B and Spelsberg T. Skeletal effects of estrogen. Endocrine Reviews 15:275-300, 1994. Kanis J. Pathogenesis of osteoporosis and fracture, Osteoporosis. London: Blackwell Science Ltd.; pp. 22-55, 1994.
16. Mack TM and Ross RK. A current perception of HRT risks and benefits, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 161-178, 1990.
17. Riggs L and Melton LJ. The prevention and treatment of osteoporosis. N Engl J Med 327:620-627, 1992.
18. Guidelines for Preclinical and Clinical Evaluation of Agents Used in the Prevention or Treatment of Postmenopausal Osteoporosis. Division of Metabolism and Endocrine Drug Products, Food and Drug Administration. April, 1994.
19. Ribot C, Tremolieres F and Pouilles JM. Prevention of postmenopausal bone loss by long-term parenteral administration of 17B estradiol: comparison of percutaneous and transdermal route, Osteoporosis: Physiological Basis, Assessment and Treatment. Deluca HF and Mazess R, eds., Elsevier (New York), pp. 137-145, 1990.
20. Delmas P. Biochemical markers of bone turnover. Methodology and clinical use in osteoporosis. Am J Med 91:59S-63S, 1991.
21. Duda RJ, O'Brien JF, Katzman JA, Peterson JM, Mann KG, and Riggs BL. Concurrent assays of circulating bone gla-protein and bone alkaline phosphatase: effects of sex, age, and metabolic bone disease. J Clin Endocrin Metab 66:951-957, 1988.
22. HHS Publication No. (CDC) 93-8395. 3rd edition. May 1993. Biosafety in Microbiological and Biomedical Laboratories. Washington, DC: US Government Printing Office, 1993.
23. Manual Guide – Safety Management No. CDC-22. Decontamination of Laboratory Sink Drains to Remove Azide Salts. April 30, 1976. Atlanta, GA: Centers for Disease Control.
24. Hill CS and Wolfert RL. The preparation of monoclonal antibodies which react preferentially with human bone alkaline phosphatase and not liver alkaline phosphatase. Clin Chim Acta 186:315-320, 1989.
25. Garnero P, Shih WJ, Gineyts E, Karpf DB, and Delmas PD. Comparison of new biochemical markers of bone turnover in late postmenopausal osteoporotic women in response to alendronate treatment. J Clin Endocrin Metab 79:1693-1700, 1994.
26. Liberman UA, Weiss SR, Broll J, Minne HW, Quan H, Bell NH, Rodriguez-Portales J, Downs RW, Dequeker J, Favus M, Seeman E, Recker RR, Capizzi T, Santora AC, Lombardi A, Shah RV, Hirsch LJ and Karpf DB. For the Alendronate Phase III Osteoporosis Treatment Study Group. Effect of oral alendronate on bone mineral density and the incidence of fractures in postmenopausal osteoporosis. N Engl J Med 333:1437-1443, 1995.
27. Dawson-Hughes B, Dallal GE, Krall EA, Sadowski L, Sahyoun N and Tannenbaum S. A controlled trial of the effect of calcium supplementation on bone density in postmenopausal women. N Engl J Med 323:878-883, 1990.
28. Reid IR, Ames RW, Evans MC, Gamble GD, and Sharpe SJ. Long-term effects of calcium supplementation on bone loss and fractures in postmenopausal women: a randomized controlled trial. Am J Med 98:331-335, 1995.
29. Devogelaer JP, Broll H, Correa-Rotter R, Cumming DC, De Deuxchaisnes CN, Geusens P, Hosking D, Jaeger P, Kaufman JM, Leite M, Leon J, Liberman U, Menkes CJ, Meunier PJ, Reid I, Rodriguez J, Romanowicz A, Seeman E, Vermeulen A, Hirsch LJ, Lombardi A, Plezia K, Santora AC, Yates AJ and Yuan W. Oral alendronate induces progressive increases in bone mass of the spine, hip, and total body over 3 years in postmenopausal women with osteoporosis. Bone 18:141-150, 1996.
30. Overgaard K, Alexandersen P, Riis BJ and Christiansen C. Evaluation of a new commercial IRMA for bone-specific alkaline phosphatase during treatment with hormone replacement therapy and calcitonin. Clin Chem 42:973-974, 1996.
31. Soletormos G, Schioler V, Nielsen D, Skovsgaard T, and Dombrowsky P. Interpretation of results for tumor markers on the basis of analytical imprecision and biological variation. Clin Chem 10:2077-2083, 1993.
32. Blumsohn A, Hannon RA, Al-Dehaimi AW, and Eastell R. Short-term intraindividual variability of markers of bone turnover in healthy adults. J Bone Min Res 9:S153, 1994.
33. Raisz L, Wiita B, Artis A, Bowen A, Schwartz S, Trahiotis M, Shoukri K and Smith J. Comparison of the effects of estrogen alone and estrogen plus androgen on biochemical markers of bone formation and resorption in postmenopausal women. J Clin Endocrin Metab 81:37-43, 1996.
34. Kanis J. Treatment of generalized osteoporosis with inhibitors of bone resorption, Osteoporosis. London: Blackwell Science Ltd.; pp. 168-195, 1994.
35. Kennedy J, Carey R, Coolen R, Garber C, Hartman A, Lee H, Leitz V, Levine J, McLean M, Osberg I, Steindel S and Sylvester E. User Evaluation of Precision Performance of Clinical Chemistry Devices. Tentative Guideline, NCCLS 4:8 EP5-T2, 1984.
36. Moss D and Whitby L. A simplified heat-inactivation method for investigating alkaline phosphatase isoenzymes in serum. Clin Chem Acta 61:63, 1975.
37. Price C, Mitchell C, Moriarty J, Gray M and Noonan K. Mass versus activity of an immunometric assay for bone alkaline phosphatase in serum. Ann Clin Biochem 32:405, 1995.

<p>bar bis e caducidad e entro isque r tot til é do é do via λήξης validade álható öre ed</p>		<p>GB Batch code DE Chargenbezeichnung ES Código de lote IT Codice del lotto FR Code du lot NL Lot nummer DK Lotnummer CZ Číslo šarže SK Číslo šarže GR Αριθμός Παρτίδας PT Código do lote HU Sarzszám SE Lot nummer PL Kod partii</p>
<p>ie number immer de catálogo di catalogo e du catalogue is nummer ummer vé číslo vé číslo καταλόγου ia de catálogo sszám ummer atalogowy</p>		<p>GB Manufacturer DE Hersteller ES Fabricante IT Fabbricante FR Fabricant NL Fabrikant DK Producent CZ Výrobce SK Výrobca GR Κατασκευαστής PT Fabricante HU Gyártó SE Tillverkare PL Producent</p>
<p>sufficient for <n> tests sreichend für <n> Prüfungen o suficiente para <n> ensayos o sufficiente per "n" saggi suffisant pour "n" tests oldoende voor "n" testen er tilstrækkeligt til "n" test it pro <n> testú ostačuje na <n> stanovení eno επαρκές για «n» εξετάσεις o suficiente para "n" ensaios tartalma <n> vizsgálat elvégzéséhez elegendő ll "n" antal tester y na wykonanie <n> testów</p>		<p>GB In Vitro Diagnostic Medical Device DE In-Vitro-Diagnostikum ES Producto sanitario para diagnóstico in vitro IT Dispositivo medico-diagnostico in vitro FR Dispositif médical de diagnostic in vitro NL Medisch hulpmiddel voor in-vitro diagnostiek DK Medicinsk udstyr til in vitro-diagnostik CZ In Vitro diagnostický zdravotnický prostředek SK Zdravotnícka pomôcka in vitro GR In Vitro Διαγνωστικό Ιατροτεχνολογικό προϊόν PT Dispositivo médico para diagnóstico in vitro HU In vitro diagnosztikum SE Medicintekniska produkter för in vitro diagnostik PL Wyrób do diagnostyki In Vitro</p>
<p>ture limitation turbegrenzung temperature temperature e temperature tuurimiet turbegrænsning rozmezi od do rozmedzie od do μοι θερμοκρασίας e temperature klettartomány turbegrænsning egac zakresu temperatury</p>		<p>GB Consult Instructions for Use DE Gebrauchsanweisung beachten ES Consulte las instrucciones de uso IT Consultare le istruzioni per l'uso FR Consulter les instructions d'utilisation NL Raadpleeg de gebruiksaanwijzing DK Se brugsanvisning CZ Viz návod k použití SK Vid' návod na použitie GR Συμβουλευτείτε τις οδηγίες χρήσης PT Consulte as instruções de utilização HU Nézze meg a Használati utasítást SE Se handhavandebeskrivningen PL Sprawdź w instrukcji obsługi</p>

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