

biocolor

life science assays

Internet Manual Downloaded from www.biocolor-assays.com

Hyaluronan Assay Protocol

Detection Limit: 0.2 ug

Hyaluronan abbreviated to HA

All test samples will require protein digestion.

Protein removal

[1] Prepare tissue sample, typically 50 mg \pm 10 mg for skin and cartilage and 250 mg \pm 50 mg for soft tissues.

Chop the test samples into small fragments using a surgical scalpel and record the weights.

Transfer samples to 2.0 ml screw capped round bottom microcentrifuge tubes. Add 400ul TRIS-HCl buffer (pH 7.6) containing 20 ul Proteinase K. Digest tissue protein overnight at 55°C. Centrifuge extracts at 12,000 r.p.m.for 10 min. **DISCARD RESIDUES**

Recovery of GAG

- [2] Transfer supernatants to a set of 1.5 ml conical microcentrifuge tubes. Add 1.0 ml **GAG Precipitation Reagent** to each tube. Mix and then leave undisturbed for 15 minutes. Centrifuge at 12,000 r.p.m for 10 min. **RETAIN RESIDUES**
- [3] To the residues add 360 ul water to each tube.

 Extract for 15 minutes by mixing with a vortex mixer at intervals to fully dissolve the GAG. The pellet must be completely dislodged from tube wall and fully dispersed. It may be necessary to manually dislodge the pellet using a fine tipped pipette tip.

Add 40 ul concentrated **NaCl** and mix before adding 95 ul **CPC**. (Cetylpyridinium chloride). Leave for 30 minutes. Centrifuge at 12,000 r.p.m for 10 min. **DISCARD RESIDUES**

- All Ensure the complete transfer of the supernatants to a new set of 1.5 ml conical centrifuge tubes.

 Add 1.0 ml GAG Precipitation Reagent to each tube.

 Mix and leave undisturbed for 15 minutes.

 Centrifuge at 12,000 r.p.m for 10 min. RETAIN RESIDUES
- [5] Add 300 ul water to each residue, extract for 15 min, with mechanical mixing or occasional vortexing. Add 33 ul concentrated NaCl and mix before adding 77 ul **CPC.** Leave 30 minutes before centrifuging at 12,000 r.p.m for 10 min. **DISCARD RESIDUES**

ASSAY PROTOCOL CONTINUED ON INSIDE BACK COVER

ASSAY PROTOCOL CONTINUED FROM INSIDE FRONT COVER

- HA [6] Transfer supernatants to a new set of 1.5 ml conical centrifuge tubes.

 Isolation Add 1.0 ml GAG Precipitation Reagent to each tube.

 Mix and leave undisturbed for 15 min before centrifuging at 12,000 r.p.m. for 10 min. RETAIN RESIDUES
 - [7] To the residues add 500 ul +98% ethanol without sodium acetate. Do not mix. Centrifuge at 12,000 r.p.m. for 5 min then drain and discard supernatant. Leave tubes inverted on absorbent paper to drain dry.
 - [8] Add 100 ul water to each residue and use intermittent mixing to complete the solubility the HA pellets. Allow 30 min for HA to fully hydrate.

Measurement of HA

[9] Transfer 20 ul aliquots of Standards / Reagent Blanks / Test Samples to the flat bottom wells of a 96 microwell plate.

Add 200 ul of **Dye Reagent.** Mix by gentle rotation of the plate on the lab bench. Full colour development occurs within 10 min and is stable for one hour when stored in the dark.

IMPORTANT: Avoid exposure to light, hold microplate within a light proof drawer or box until transfer to the microplate reader.

[10] Using a microplate spectrophotometer read the Absorbance values at 655 nm and print data obtained.

Data Analysis

With the aid of a Standard Curve (see page 8) convert the Abs @ 655 nm into ug HA. Multiply this value by five (20 ul x 5 = 100 ul).

The 100 ul of purified HA [Step 8] was recovered from x mg of wet tissue (see [1] for the recorded wet weights of each tissue sample).

Convert the ug HA extracted from x mg wet tissue to ug HA contained in one gram of wet tissue.

If a test sample volume of 20 ul produces an Absorbance value greater than 2.0 then use a smaller aliquot of [Step 8] made up to 20 ul of water before adding 200 ul Dye Reagent.

Biocolor's

Extracellular Matrix Assays

Sircol Soluble Collagen Assay

A dye-binding method for the analysis of cold pepsin-soluble collagen fibres The Assay can assess the rate of newly synthesised collagen produced during periods of rapid growth, development, tissue repair, remodelling and wound healing.

Sircol Insoluble Collagen Assay

A dye-binding method for the analysis of residual, insoluble collagen fibres.

The Assay replaces hydroxyproline analysis and can be completed within 1 working day.

The ratio of soluble and insoluble collagens within a sample can be measured using the same Sircol dye.

Blyscan Sulfated Glycosaminoglycan Assay

A quantitative dye-binding method for the analysis of sulfated proteoglycans and glycosaminoglycans.

Measures the total sGAG content and can determine the O and N-sulfated glycosaminoglycan content.

Test material can be assayed directly when present in soluble form, or following papain extraction from biological material

Fastin Elastin Assay

A quantitative dye-biding method for the analysis of elastins extracted from mammalian sources.

The Assay can measure soluble tropoelastins, and insoluble elastins following hot acid solubilisation

The dye reagent binds to 'basic' and 'non-polar' amino acid sequences found in mammalian elastins

Purple-Jelley Assay

Hyaluronan / Hyaluronic acid

Assay manuals are available to download free of charge from our website www.biocolor-assays.com

Purple-Jelley HYALURONAN [Hyaluronic acid] Assay

Assay designed for research work with animal tissues. The assay is not suitable for blood and urine analysis


Handle the Assay using Good Laboratory Practice

CONTENTS

ASSAV PROTOCOL

ASSAT TROTOCOL	mside none and back covers
The nature of Hyaluronan	2
Storage of tissue	3
Sampling of tissue	4
Hyaluronan Assay components	5
The nature of the dye label for H	IA 6
Supplement Notes to the ASSAY	PROTOCOL
The need to isolate HA prior to	dye labelling 7
HA Standard Curve	8

Published by
Biocolor Ltd
8 Meadowbank Road, Carrickfergus
BT38 8YF, Northern Ireland, U.K



Nature of HYALURONAN [HA]

A hydrated molecule of HA is a gentle giant of a carbohydrate polymer. It is composed of a long flexible non-branching chain formed with a repeating disaccharide pattern. The disaccharide is made up of alternative uronic acid and aminosugar units.

In the human umbilical cord and synovial fluid, the size of HA is reported to be about 3,000,000 Da. (a flexible chain of 4000 disaccharide units).

Synthesis:

HA is formed by Hyaluronan synthases embedded within the cell membrane. HA is attached to the outside of the cell membrane and can be released from intact cells grown in culture by brief treatment with trypsin. The newly formed HA interacts with the ECM, the cell membrane and with water to form and maintain the formation of a pericellular substrata around the cell, the HA generating a territorial gel, or in cartilage to produce a collagen free lagoon for the cell.

Degradation:

HA is degraded principally by a group of seven hyaluronidases and may account for the variations in the sizes of HA reported. The activity of some hyaluronidases may remain during post mortem acid glycolysis when the animal body is still warm.

Fig. 1 Carbohydrate polymer of HA composed of a repeating disaccharide of an uronic acid and an aminosugar neither of which are sulphated.

The role of Hyaluronan

HA has the ability to exert control of the water balance of tissues. This property allows selected cells to migrate during growth and is essential for animal development from fertilization to birth.

HA however can have adverse effects in that it will also provide migratory tumour cells the ability to relocate at extended distances from the source of a primary tumour. The wandering tumour cells on relocation may induce metastasis to generate multiple secondary tumours.

Hyaluronan Reviews

Chemistry and Biology of Hyaluronan, Edited by H.G. Garg and C.A. Hales, (2004).

Publisher: Elsevier Oxford, UK

Hyaluronan in Cancer Biology, Edited by R.Stern, (2009) Publisher: Academic Press, San Diego, U.S.A

Storage of tissues

Until tissue is required for analysis it should be kept frozen to prevent the action of hyaluronidases.

Enzyme degradation of HA produces a mix of tetrasaccharide and disaccharide fragments that do not react with the assay dye reagent.

Blood and urine contains HA fragments that have been reduced to tetrasaccharides and disaccharides and these fluids are unsuitable for this assay.

In vitro cell culture, in monolayer or in suspension, do not generate sufficient HA for dye analysis.

Purple-Jelley HYALURONAN Assay

Pack Sizes and Storage Conditions

Standard Assay Kit: Product Code H1000 (100 assays)
Economy Assay Kit: Product Code H2000 (400 assays)

All components are stable for one year when stored unopened at 10 to 25°C. Once opened the Hyaluronan Standard should be stored at +5°C. The Dye Reagent must be stored in the dark.

Sampling of tissue, wet weight

The scheme is described on the inside front cover of this manual.

The actual weight of each sample should be recorded and the sample placed into 2.0 ml screw capped tubes suitably labelled for subsequent protein removal by the action of Proteinase K.

The distribution of HA in tissues

4

as produced in this assay procedure is listed in Table 1 and are given only as a outline guide.

In adult animals the $H\Lambda$ percentage gradually increase during ageing as muscle and fat mass decreases

Table 1. Distribution of HA in tissue

Species	Tissue	Concentration
		ug/g
Mouse	heart	48
	kidney	11
	lung	14
	leg muscle	20
	skin	307
Rat	heart	83
	kidney	28
	lung	17
	leg muscle	67
	skin	510
Rabbit	heart	138
	kidney	27
	lung	45
	leg muscle	50
	skin	428

Hyaluronan Assay Components

PROTEIN REMOVAL

Protein Digestion Reagent: A buffer containing 50 mM TRIS-HCl, pH 7.6 supplied in a tablet format, (one tablet dissolved in 15 ml dH₂O); use within 2 days.

For each tissue sample add 20 ul Proteinase K per 400 ul of buffer. A suitable Proteinase K suspension is available from Sigma Aldrich (Code: P4850

Screw Capped Round Bottomed Centrifuge Tubes, 2 ml capacity for the protein digestion process.

HYALURONAN / HA ASSAY

GAG Precipitation Reagent is a sodium acetate saturated ethanol solution. The kit contains 2×170 ml bottles part filled with water saturated suspension of Na acetate.

Add 136 ml Ethanol (or Industrial Methylated Spirits containing 99% ethanol) to the Na acetate in each bottle. Ensure complete mixing by inverting the bottle and gently shaking the contents..

Leave for 30 min to allow excess Na acetate to settle out at the bottom of the bottle. Do not shake contents again so as to avoid solid Na acetate being resuspended.

Concentrated sodium chloride supplied as a sterile solution in a 20 ml vial (23.3 gm % w/v)

CPC, [Cetylpyridinium chloride] supplied as a sterile solution in a 20 ml vial Warm before use, to ensure complete solution, (2.0 gm % w/v)

HA,[Hyaluronan / Hyaluronic acid Standard] as a sterile solution containing 200 ug /ml (4 ug HA/20 ul) in a 10 ml vial, store at 50 C.

Purple Dye Reagent; 3,3'-Diethyl-9-methyl-4,5,4',5'-dibenzothia carbocyanine bromide supplied in 55% solution of dimethyl sulfoxide (DMSO) in a 20 ml vial.

THE DYE REAGENT MUST BE PROTECTED FROM LIGHT.

REQUIRED FOR THE ASSAY, BUT NOT SUPPLIED

Proteinase K product is available from Sigma Aldrich (Code: P4850)

Ethanol (including Industrial Methylated Spirits) is subject to Government Alcohol Regulations and cannot be supplied by Biocolor.

The Dye used to Assay HA

In 1936 Edwin Jelley sent a 'Letter to the Editor' of Nature, (138 1009 -1010), regarding the unusual behaviour of some cyanine dyes.

When dissolved in 5 M NaCl these dyes produce a third absorbance peak at a longer wave-length (650nm) whereas in deionised water only a double peak occurs at \sim 540 nm and \sim 570 nm.

The 650 nm peak in concentration dye solution induces aggregation of the dye molecules and has been described as a '*J-aggregate*' named after Jelley.

The cyanine dye used in the assay is not a planar shaped molecule (Fig 2):

$$\begin{array}{c|c} & & & \\ &$$

Fig.2 Structure of 3,3'-Diethyl-9-methyl-4,5,4',5'- dibenzothiacarbo-cyanine bromide.

Individual dye molecules may form a stacking format to produce a supramolecular complex as seen in peak 3, (Image on front cover of the Manual)

Further studies in the 1960s notably by Kay et.al. (J. Physical Chem. **68** 1896 – 1906) found many biological polymers including proteins, DNA, polar lipids and glycosaminoglycans could induce the third peak when using high dye concentrations.

The dye was renamed 'Stains-all' by Dhlberg, Dingman and Peacock in 1969 (J. Mol. Biol. 41, 132).




Fig. 3 Absorbance Spectrum of Stains-all when mixed with glycosaminoglycans. Sulfated glycosaminoglycans need to be removed before measuring the isolated HA, from a two step critical electrolyte salting out process (CEC).

RED line spectrum: HA only (2 ug).

BLUE line pectrum chondroitin 4-sulfate only (2 ug).

GREEN line spectrum: equal quantities of HA and chondroitin 4-sulfate.

GREY line spectrum HA plus twice the concentration of chondroitin 4-sulfate.

Assay protocol - Supplementary Notes

- [1] Removal of tissue protein is essential for isolating HA from dye binding proteins.

 Proteinase K (EC 3.4.21.64), a broad serine Proteinase was selected in preference to Papain.
 - Proteinase K was used at 55°C at pH 7.6 from 3 hours to overnight for soft tissue. For hard tissues, including cartilage, it was used overnight.
- Recovery of GAG from the protein digest. This requires a sequence of precipitation steps using (a) ethanol saturated with Na acetate and (b) CPC containing Na Cl. Industrial Methylated Spirits containing 99% ethanol is suitable but may requires a Government Alcohol Licence for your laboratory.


Assay Protocol - Supplementary Notes

For comfort breaks,(coffee / lunch /etc) select the end of section [4] and/or [6] and/or [8]. The same breakpoints applies to overnight stops where the precipitate should be stored at 50 C.

Assay Protocol: [9] Color Stability:

No decrease of Absorbance readings at 655 nm occurred over one hour from mixing the dye and HA, providing the microwell plates were wrapped in aluminium foil between Absorbance readings.





Linear regression fit

Fig.4. Typical straight line Standard Curve

Microplate Readers vary in their design and performance so this Figure should be considered as a guide only.

The HA Standard aliquots sampled were adjusted to 20ul with water before mixing with 200ul Dye Reagent.