


# BioVendor

Research  
and Diagnostic Products



## HUMAN CYCLOPHILIN A ELISA

Product Data Sheet

Cat. No.: RD191329200R

For Research Use Only

# CONTENTS

---

1.	INTENDED USE	3
2.	STORAGE, EXPIRATION	3
3.	INTRODUCTION	4
4.	TEST PRINCIPLE	6
5.	PRECAUTIONS	6
6.	TECHNICAL HINTS	6
7.	REAGENT SUPPLIED	7
8.	MATERIAL REQUIRED BUT NOT SUPPLIED	7
9.	PREPARATION OF REAGENTS	8
10.	PREPARATION OF SAMPLES	9
11.	ASSAY PROCEDURE	9
12.	CALCULATIONS	12
13.	PERFORMANCE CHARACTERISTICS	13
14.	DEFINITION OF THE STANDARD	16
15.	PRELIMINARY POPULATION AND CLINICAL DATA	16
16.	METHOD COMPARISON	17
17.	TROUBLESHOOTING AND FAQs	18
18.	REFERENCES	19
19.	EXPLANATION OF SYMBOLS	21

**»» This kit is manufactured by:  
BioVendor-Laboratorní medicína a.s.**

**»» Use only the current version of Product Data Sheet enclosed with the kit!**

## 1. INTENDED USE

---

The RD191329200R Human Cyclophilin A ELISA is a sandwich enzyme immunoassay for the quantitative measurement of human cyclophilin A (CYPA).

### »» Features

- **It is intended for research use only**
- The total assay time is less than 3.5 hours
- The kit measures CYPA in serum and plasma (EDTA, citrate, heparin)
- Assay format is 96 wells
- Standard is recombinant protein based
- Components of the kit are provided ready to use, concentrated or lyophilized

## 2. STORAGE, EXPIRATION

---

Store the complete kit at 2-8°C. Under these conditions, the kit is stable until the expiration date (see label on the box).

For stability of opened reagents see Chapter 9.

### 3. INTRODUCTION

---

Cyclophilins belong to a group of proteins collectively known as immunophilins that have peptidyl-prolyl cis-trans isomerase activity which catalyze the isomerization of peptides and may play crucial roles in many biological conditions including protein folding, trafficking, assembly, T cell activation and cell signaling [4]. Cyclophilins also have varying degrees of affinity for the immunosuppressive drug Cyclosporine A (CsA), a cyclic 11-amino-acid peptide produced by the fungus *Tolypocladium infatum*.

Cyclophilin A (CYPA), an 18 kDa protein consisting of 165 amino acid residues, is the major intracellular receptor for Cyclosporine A [14].

In mammals, the CsA-CYPA complex binds to and inhibits calcineurin, a calcium-calmodulin-activated serine/threonine-specific protein phosphatase. The inhibition of calcineurin blocks the translocation of nuclear factor of activated T cells from the cytosol to the nucleus, thus preventing the transcription of genes encoding cytokines such as interleukin-2 [14].

Although CYPA is primarily located intracellularly, it can be secreted into the extracellular environment in various cell types [4].

Secreted CYPA activates cardiovascular cells resulting in a variety of cardiovascular diseases, including vascular remodeling, abdominal aortic aneurysm formation, atherosclerosis, cardiac hypertrophy and myocardial ischemic reperfusion injury [9].

CYPA is also an important mediator of angiotensin II-induced cardiac hypertrophy. Patients with acute coronary syndrome have high plasma concentrations of CYPA, and CYPA is strongly expressed in the atheromatous plaques of patients with acute myocardial infarction (AMI) [2].

Intracellular CYPA is secreted from cells in response to inflammatory stimuli such as hypoxia, infection and oxidative stress [10]. Its expression increases in inflammatory conditions including rheumatoid arthritis, autoimmune disease and cancer [9]. CYPA levels are higher in the serum and synovial fluids of rheumatoid arthritis patients and are also elevated in tumors including non-small cell lung cancer, pancreatic adenocarcinoma, hepatocellular carcinoma, oral cancer, buccal squamous cell carcinomas and breast cancer [6]. It has been shown that CYPA and its vascular receptor basigin (BSG) are crucial for hypoxia-induced pulmonary hypertension by inducing growth factor secretion, inflammatory cell recruitment and vascular smooth muscle cells proliferation [16].

Extracellular CYPA has a potent chemotactic effect on leukocytes, monocytes, and lymphocytes [7]. Studies with CYPA protein demonstrate that CYPA can induce chemotaxis of leukocytes and signaling via two distinct pathways: PPlase activity and extracellular binding to CD147 receptor. At present, the precise molecular mechanisms of CYPA/CD147 interaction have still not been elucidated in detail [4].

A growing body of research suggested CYPA and CD147 involvement in key processes of kidney disease pathologies (renal cell carcinoma, acute kidney injury, nephritis and renal fibrosis). Based on various studies reviewed, it is clear that in kidney diseases, both CD147 and CYPA have multifunctional properties, both independently and as an interacting complex [4].

CYPA is also known to play critical roles in the proliferation of a number of viruses, including human immunodeficiency virus type 1 (HIV-1), influenza virus, hepatitis virus C (HCV), vesicular stomatitis virus (VSV), vaccinia virus, severe acute respiratory syndrome coronavirus (SARS-CoV), rotavirus (RV) and human papillomavirus (HPV) [3].

Areas of investigation:

Cardiovascular Disease

Immune response, Infection and Inflammation

Oncology

## 4. TEST PRINCIPLE

---

In the BioVendor Human Cyclophilin A ELISA, standards and samples are incubated in microplate wells pre-coated with polyclonal anti-human CYPA antibody. After 60 minutes incubation and washing, biotin labelled polyclonal anti-human CYPA antibody is added and incubated for 60 minutes with captured human CYPA. After another washing, streptavidin-HRP conjugate is added. After 30 minutes incubation and the last washing step, the remaining conjugate is allowed to react with the substrate solution (TMB). The reaction is stopped by addition of acidic solution and absorbance of the resulting yellow product is measured. The absorbance is proportional to the concentration of CYPA. A standard curve is constructed by plotting absorbance values against concentrations of standards, and concentrations of unknown samples are determined using this standard curve.

## 5. PRECAUTIONS

---

- **For professional use only**
- Wear gloves and laboratory coats when handling immunodiagnostic materials
- Do not drink, eat or smoke in the areas where immunodiagnostic materials are being handled
- This kit may contain components of human or animal origin. These materials should be handled as potentially infectious
- Avoid contact with the acidic Stop Solution and Substrate Solution, which contains hydrogen peroxide and tetramethylbenzidine (TMB). Wear gloves and eye and clothing protection when handling these reagents. Stop and/or Substrate Solutions may cause skin/eyes irritation. In case of contact with the Stop Solution and the Substrate Solution wash skin/eyes thoroughly with water and seek medical attention, when necessary
- The materials must not be pipetted by mouth

## 6. TECHNICAL HINTS

---

- Reagents with different lot numbers should not be mixed
- Use thoroughly clean glassware
- Use deionized (distilled) water, stored in clean containers
- Avoid any contamination among samples and reagents. For this purpose, disposable tips should be used for each sample and reagent
- Substrate Solution should remain colourless until added to the plate. Keep Substrate Solution protected from light

- Stop Solution should remain colourless until added to the plate. The colour developed in the wells will turn from blue to yellow immediately after the addition of the Stop Solution. Wells that are green in colour indicate that the Stop Solution has not mixed thoroughly with the Substrate Solution
- Dispose of consumable materials and unused contents in accordance with applicable national regulatory requirements

## 7. REAGENT SUPPLIED

<i>Kit Components</i>	<i>State</i>	<i>Quantity</i>
Antibody Coated Microtiter Strips	ready to use	96 wells
Biotin Labelled Antibody Conc. (100x)	concentrated	0.13 ml
Streptavidin-HRP Conjugate	ready to use	13 ml
Master Standard	lyophilized	2 vials
Biotin-Ab Diluent	ready to use	13 ml
Dilution Buffer	ready to use	20 ml
Wash Solution Conc. (10x)	concentrated	100 ml
Substrate Solution	ready to use	13 ml
Stop Solution	ready to use	13 ml
Product Data Sheet + Certificate of Analysis	-	1 pc

## 8. MATERIAL REQUIRED BUT NOT SUPPLIED

- Deionized (distilled) water
- Test tubes for diluting samples
- Glassware (graduated cylinder and bottle) for Wash Solution (Dilution Buffer)
- Precision pipettes to deliver 5-1000  $\mu$ l with disposable tips
- Multichannel pipette to deliver 100  $\mu$ l with disposable tips
- Absorbent material (e.g. paper towels) for blotting the microtitre plate after washing
- Vortex mixer
- Orbital microplate shaker capable of approximately 300 rpm
- Microplate washer (optional). [Manual washing is possible but not preferable.]
- Microplate reader with  $450 \pm 10$  nm filter, preferably with reference wavelength 630 nm (alternatively another one from the interval 550-650 nm)
- Software package facilitating data generation and analysis (optional)

## 9. PREPARATION OF REAGENTS

- All reagents need to be brought to room temperature prior to use
  - Always prepare only the appropriate quantity of reagents for your test
  - Do not use components after the expiration date marked on their label
- Assay reagents supplied ready to use:

### Antibody Coated Microtiter Strips

#### Stability and storage:

Return the unused strips to the provided aluminium zip-sealed bag with desiccant and seal carefully. Remaining Microtiter Strips are stable 3 months stored at 2-8°C and protected from the moisture.

### Streptavidin-HRP Conjugate

#### Biotin-Ab Diluent

#### Dilution Buffer

#### Substrate Solution

#### Stop Solution

#### Stability and storage:

Opened reagents are stable 3 months when stored at 2-8°C.

- Assay reagents supplied concentrated or lyophilized:

### Human CYPA Master Standard

**Refer to the Certificate of Analysis for current volume of Dilution Buffer needed for reconstitution of standard!!!**

Reconstitute the lyophilized Master Standard with Dilution Buffer just prior to the assay. Let it dissolve at least 15 minutes with occasional gentle shaking (not to foam).

The resulting concentration of CYPA in the stock solution is **25 ng/ml**.

Prepare set of standards using Dilution Buffer as follows:

<i>Volume of Standard</i>	<i>Dilution Buffer</i>	<i>Concentration</i>
Stock	-	25 ng/ml
250 µl of stock	250 µl	12.5 ng/ml
250 µl of 12.5 ng/ml	250 µl	6.25 ng/ml
250 µl of 6.25 ng/ml	250 µl	3.13 ng/ml
250 µl of 3.13 ng/ml	250 µl	1.56 ng/ml
250 µl of 1.56 ng/ml	250 µl	0.78 ng/ml
250 µl of 0.78 ng/ml	250 µl	0.39 ng/ml

**Prepared Standards are ready to use, do not dilute them.**

Stability and storage:

**Do not store the reconstituted Master Standard and/or diluted standard solutions.**

#### **Biotin Labelled Antibody Conc. (100x)**

Prepare the working Biotin Labelled Antibody solution by adding 1 part Biotin Labelled Antibody Concentrate (100x) with 99 parts Biotin-Ab Diluent.

Example: 10 µl of Biotin Labelled Antibody Concentrate (100x) + 990 µl of Conjugate Diluent for 1 strip (8 wells). **Mix well** (not to foam).

Stability and storage:

Opened Biotin Labelled Antibody Concentrate (100x) is stable 3 months when stored at 2–8°C.

**Do not store the diluted Biotin Labelled Antibody Solution.**

#### **Wash Solution Conc. (10x)**

Dilute Wash Solution Concentrate (10x) ten-fold in distilled water to prepare a 1x working solution. e.g. 100 ml of Wash Solution Concentrate (10x) + 900 ml of distilled water for use of all 96-wells.

Stability and storage:

The diluted Wash Solution is stable 1 month when stored at 2–8°C. Opened Wash Solution Concentrate (10x) is stable 3 months when stored at 2–8°C.

## 10. PREPARATION OF SAMPLES

---

The kit measures CYPA in serum and plasma (EDTA, citrate, heparin).

Samples can be assayed immediately after collection, or after long-term storage at -70°C. Thoroughly mix thawed samples just prior to the assay and avoid repeated freeze/thaw cycles, which may cause erroneous results. Avoid using hemolyzed or lipemic samples.

Dilute samples **3x** with Dilution Buffer just prior to the assay, e.g. 50 µl of sample + 100 µl of Dilution Buffer for singlets, or preferably 100 µl of sample + 200 µl of Dilution Buffer for duplicates. **Mix well** (not to foam). Vortex is recommended.

Stability and storage:

Samples should be stored at -20°, or preferably at -70° for long-term storage. Avoid repeated freeze/ thaw cycles.

**Do not store the diluted samples.**

*Note: It is recommended to use a precision pipette and a careful technique to perform the dilution in order to get precise results.*

## 11. ASSAY PROCEDURE

1. Pipet **100 µl** of diluted Standards, Dilution Buffer (=Blank) and samples, preferably in duplicates, into the appropriate wells. See *Figure 1* for example of work sheet.
2. Incubate the plate at room temperature (ca. 25°C) for **1 hour**, shaking at ca. 300 rpm on an orbital microplate shaker.
3. Wash the wells 3-times with Wash Solution (0.35 ml per well). After final wash, invert and tap the plate strongly against paper towel.
4. Pipet **100 µl** of Biotin Labelled Antibody solution into each well.
5. Incubate the plate at room temperature (ca. 25°C) for **1 hour**, shaking at ca. 300 rpm on an orbital microplate shaker.
6. Wash the wells 3-times with Wash Solution (0.35 ml per well). After final wash, invert and tap the plate strongly against paper towel.
7. Pipet **100 µl** of Streptavidin-HRP Conjugate into each well.
8. Incubate the plate at room temperature (ca. 25°C) for **30 minutes**, shaking at ca. 300 rpm on an orbital microplate shaker.
9. Wash the wells 3-times with Wash Solution (0.35 ml per well). After final wash, invert and tap the plate strongly against paper towel.
10. Add **100 µl** of Substrate Solution into each well. Avoid exposing the microtiter plate to direct sunlight. Covering the plate with e.g. aluminium foil is recommended.
11. Incubate the plate for **10 minutes** at room temperature. The incubation time may be extended [up to 20 minutes] if the reaction temperature is below than 20°C. Do not shake the plate during the incubation.
12. Stop the colour development by adding **100 µl** of Stop Solution.
13. Determine the absorbance of each well using a microplate reader set to 450 nm, preferably with the reference wavelength set to 630 nm (acceptable range: 550-650 nm). Subtract readings at 630 nm (550-650 nm) from the readings at 450 nm. **The absorbance should be read within 5 minutes following step 12.**

*Note: If the microplate reader is not capable of reading absorbance greater than the absorbance of the highest standard, perform a second reading at 405 nm. A new standard curve, constructed using the values measured at 405 nm, is used to determine CYP A concentration of off-scale samples. The readings at 405 nm should not replace the reading for samples that were "in range" at 450 nm.*

*Note 2: Manual washing: Aspirate wells and pipet 0.35 ml Wash Solution into each well. Aspirate wells and repeat twice. After final wash, invert and tap the plate strongly against paper towel. Make certain that Wash Solution has been removed entirely.*

	<b>strip 1+2</b>	<b>strip 3+4</b>	<b>strip 5+6</b>	<b>strip 7+8</b>	<b>strip 9+10</b>	<b>strip 11+12</b>
<b>A</b>	<b>Standard 25</b>	Sample 1	Sample 9	Sample 17	Sample 25	Sample 33
<b>B</b>	<b>Standard 12.5</b>	Sample 2	Sample 10	Sample 18	Sample 26	Sample 34
<b>C</b>	<b>Standard 6.25</b>	Sample 3	Sample 11	Sample 19	Sample 27	Sample 35
<b>D</b>	<b>Standard 3.13</b>	Sample 4	Sample 12	Sample 20	Sample 28	Sample 36
<b>E</b>	<b>Standard 1.56</b>	Sample 5	Sample 13	Sample 21	Sample 29	Sample 37
<b>F</b>	<b>Standard 0.78</b>	Sample 6	Sample 14	Sample 22	Sample 30	Sample 38
<b>G</b>	<b>Standard 0.39</b>	Sample 7	Sample 15	Sample 23	Sample 31	Sample 39
<b>H</b>	<b>Blank</b>	Sample 8	Sample 16	Sample 24	Sample 32	Sample 40


*Figure 1: Example of a work sheet.*

## 12. CALCULATIONS

Most microplate readers perform automatic calculations of analyte concentration. The Standard curve is constructed by plotting the mean absorbance (Y) of Standards against the known concentration (X) of Standards in logarithmic scale, using the four-parameter algorithm. Results are reported as concentration of CYPA (ng/ml) in samples.

Alternatively, the *logit log* function can be used to linearize the standard curve (i.e. *logit* of absorbance (Y) is plotted against *log* of the known concentration (X) of standards).

**The measured concentration of samples calculated from the standard curve must be multiplied by their respective dilution factor, because samples have been diluted prior to the assay. e.g. 4.5 ng/ml (from standard curve) x 3 (dilution factor) = 13.5 ng/ml.**



*Figure 2: Typical Standard Curve for Human Cyclophilin A ELISA.*

## 13. PERFORMANCE CHARACTERISTICS

➤➤ Typical analytical data of BioVendor Human Cyclophilin A ELISA are presented in this chapter

- **Sensitivity**

Limit of Detection (LOD) (defined as concentration of analyte giving absorbance higher than mean absorbance of blank\* plus three standard deviations of the absorbance of blank:  $A_{\text{blank}} + 3 \times \text{SD}_{\text{blank}}$ ) is calculated from the real CYPA values in wells and is 0.28 ng/ml.

\*Dilution Buffer is pipetted into blank wells.

- **Limit of assay**

Samples with absorbances exceeding the absorbance of the highest standard should be measured again with higher dilution. The final concentration of samples calculated from the standard curve must be multiplied by the respective dilution factor.

- **Specificity**

The antibodies used in this ELISA are specific for human Cyclophilin A. No crossreactivity has been observed for other human recombinant Cyclophilins such as Cyclophilin F and Cyclophilin H at 25 ng/ml.

Sera of several mammalian species were measured in the assay. See results below.

For details please contact us at [info@biovendor.com](mailto:info@biovendor.com).

<i>Mammalian serum sample</i>	<i>Observed crossreactivity</i>
Bovine	no
Cat	no
Dog	no
Goat	no
Hamster	no
Horse	no
Monkey	no
Mouse	no
Pig	yes
Rabbit	yes
Rat	no
Sheep	no

## ➤ Presented results are multiplied by respective dilution factor

### • Precision

Intra-assay (Within-Run) (n=8)

<i>Sample</i>	<i>Mean</i> (ng/ml)	<i>SD</i> (ng/ml)	<i>CV</i> (%)
1	4.83	0.17	3.6
2	31.15	0.77	2.5

Inter-assay (Run-to-Run) (n=6)

<i>Sample</i>	<i>Mean</i> (ng/ml)	<i>SD</i> (ng/ml)	<i>CV</i> (%)
1	7.43	0.40	4.9
2	33.57	1.20	3.4

### • Spiking Recovery

Serum samples were spiked with different amounts of human CYPA and assayed.

<i>Sample</i>	<i>Observed</i> (ng/ml)	<i>Expected</i> (ng/ml)	<i>Recovery O/E</i> (%)
1	3.47	-	-
	8.07	8.15	99.0
	12.20	12.84	95.0
	20.17	22.22	90.8
2	2.30	-	-
	6.87	6.99	98.4
	10.79	11.68	92.4
	19.73	21.05	93.7

### • Linearity

Serum samples were serially diluted with Dilution Buffer and assayed.

<i>Sample</i>	<i>Dilution</i>	<i>Observed</i> (ng/ml)	<i>Expected</i> (ng/ml)	<i>Recovery O/E</i> (%)
1	-	31.31	-	-
	2x	15.13	15.66	96.7
	4x	7.75	7.83	99.0
	8x	3.66	3.91	93.5
2	-	23.05	-	-
	2x	11.88	11.52	103.1
	4x	5.48	5.76	95.1
	8x	2.47	2.88	85.8

- **Effect of sample matrix**

Citrate, heparin and EDTA plasmas were compared to respective serum samples from the same 10 individuals.

Results are shown below:

Volunteer No.	Serum (ng/ml)	Plasma (ng/ml)		
		EDTA	Citrate	Heparin
1	46.46	46.96	42.28	49.58
2	7.91	22.88	6.46	13.64
3	31.98	41.80	23.35	28.35
4	12.43	15.36	8.40	19.29
5	7.26	18.67	10.93	6.39
6	8.20	23.49	2.71	4.84
7	12.84	26.02	9.10	9.94
8	7.47	10.13	2.47	2.23
9	33.90	35.15	27.83	32.87
10	28.42	32.64	21.26	22.20
<b>Mean (ng/ml)</b>	<b>19.7</b>	<b>27.3</b>	<b>15.5</b>	<b>18.9</b>
<b>Mean Plasma/Serum (%)</b>	-	<b>138.7%</b>	<b>78.6%</b>	<b>96.2%</b>
<b>Coefficient of determination R<sup>2</sup></b>	-	<b>0.81</b>	<b>0.95</b>	<b>0.91</b>




Figure 3: CYPA levels measured using Human Cyclophilin A ELISA from 10 individuals using serum, EDTA, citrate and heparin plasma, respectively.

## 14. DEFINITION OF THE STANDARD

The Standard used in this kit is recombinant protein. The recombinant human CYPA, produced in *E. coli*, is 18 kDa protein consisting of 165 amino-acid residues.

## 15. PRELIMINARY POPULATION AND CLINICAL DATA

The following results were obtained when serum samples from 168 unselected donors (95 men + 73 women) 21 – 65 years old were assayed with the Biovendor Human Cyclophilin A ELISA in our laboratory.

Sex	Age (years)	n	Mean	Median	SD	Min	Max
			CYPA (ng/ml)				
Men	21-29	18	18.40	11.88	15.36	4.50	57.71
	30-39	28	13.11	11.68	7.29	3.96	34.14
	40-49	32	15.29	12.29	8.58	4.50	36.33
	50-65	10	22.00	16.91	15.22	7.86	54.02
Women	22-29	13	23.27	19.20	13.54	3.89	54.47
	30-39	28	14.01	13.86	6.04	3.94	25.67
	40-49	23	15.03	13.68	9.14	5.03	46.89
	50-61	8	17.75	13.01	14.49	6.51	54.02




Figure 4: Human CYPA concentration plotted against donor age and sex.

- **Reference range**

It is recommended that each laboratory include its own panel of control sample in the assay. Each laboratory should establish its own normal and pathological reference ranges for CYPA levels with the assay.

## 16. METHOD COMPARISON

The BioVendor Human Cyclophilin A ELISA has not been compared to any commercial immunoassay.

## 17. TROUBLESHOOTING AND FAQs

---

### »» Weak signal in all wells

Possible explanations:

- Omission of a reagent or a step
- Improper preparation or storage of a reagent
- Assay performed before reagents were allowed to come to room temperature
- Improper manual washing
- Improper wavelength when reading absorbance

### »» High signal and background in all wells

Possible explanations:

- Improper or inadequate washing
- Overdeveloping; incubation time with Substrate Solution should be decreased before addition of Stop Solution
- Incubation temperature over 30°C

### »» High coefficient of variation (CV)

Possible explanation:

- Improper or inadequate washing
- Improper mixing Standards or samples

## 18. REFERENCES









---

1. Liu MCH, Lee YW, Chang CHS, Tai YL, Yu JR, Su XT, Hsu LW, Lin SH, Wu CHH, Liu PY: cyclophilin A is Associated with Peripheral Artery Disease and Chronic Kidney Disease in Geriatrics: The Tianliao Old People (TOP) Study. *J of Scientific Reports*. April; 5 (9937): (2015)
2. Huang CHH, Chang CHCH, Kuo CHL, Huang CHS, Lin CHS, Liu CHS: Decrease in Plasma Cyclophilin A Concentration at 1 Month after Myocardial Infarction Predicts Better Left Ventricular Performance and Synchronicity at 6 Months: A Pilot Study in Patients with ST Elevation Myocardial Infarction. *J of Biological Sciences*. Jan; 11 (1): 38-47. (2015)
3. Qing J, Wang Y, Sun Y, Huang J, Yan W, Wang J, Su D, Ni CH, Li J, Rao Z, Liu L, Lou Z: Cyclophilin A Associates with Enterovirus-71 Virus Capsid and Plays an Essential Role in Viral Infection as an Uncoating Regulator. *J of PLOS Pathogens*. Oct; 10 (10): (2014)
4. Qu X, Wang CH, Zhang J, Qie G, Zhou J: The Roles of CD147 and/or Cyclophilin A in Kidney Diseases. *J of Mediators of Inflammation*. Dec; (2014)
5. Ramachandran S, Venugopal A, Kutty VR, Vinitha A, Divya G, Chitrasree V, Mullassari A, Pratapchandran NS, Santosh KR, Pillai MR, Kartha CC: Plasma level of cyclophilin A is increased in patients with type 2 diabetes mellitus and suggests presence of vascular disease. *J of Cardiovascular Diabetology*. 13(38): (2014)
6. Tian-tian Z, Jun-feng Z, Heng G: Functions of cyclophilin A in atherosclerosis. *J of Experimental and Clinical Cardiology*. 18(2): 118-124. (2013)
7. Nigro P, Pompilio G, Capogrossi MC: Cyclophilin A: a key player for human disease. *J of Cell Death and Disease*. Oct; 4: (2013)
8. Elvers M, Herrmann A, Seizer P, Münzer P, Beck Schönberger T, Borst O, Martin-Romero FJ, Lang F, May AE, Gawaz M: Intracellular cyclophilin A is an important  $\text{Ca}^{2+}$  regulator in platelets and critically involved in arterial thrombus formation. *J of The American Society of Hematology*. Aug; 120(6): 1317-1326. (2012)
9. Soe NN, Bradford MD, Berk C: Cyclophilin A: A Mediator of Cardiovascular Pathology. *J of Korean Society Hypertens*. Dec; 17(4): 133-147. (2011)
10. Obchoei S, Weakley SM, Wongkham S, Wongkham Ch, Sawanyawisuth K, Yoa O, Chen Ch: Cyclophilin A enhances cell proliferation and tumor growth of liver fluke-associated cholangiocarcinoma. *J of Molecular Cancer*. Aug; 10(102): (2011)
11. Nigro P, Satoh K, O'Dell MR, Soe NN, cui Z, Mohan A, Abe J, Alexis JD, Sparks JD, Berk BC: Cyclophilin A is an inflammatory mediator that promotes atherosclerosis in apolipoprotein E-deficient mice. *J of Experimental Medicine*. Dec; 208(1): 53-66. (2010)
12. Satoh K, Shimokawa H, Berk BC: Cyclophilin A. *J of Cardiovascular Therapy*. Nov; 74(11): 2249-2256. (2010)
13. Colgan J, Asmal M, Yu B, Luban J: Cyclophilin A-Deficient Mice Are resistant to Immunosuppression by Cyclosporine. *J of Immunology*. 174: 6030-6038. (2005)
14. Wang P, Heitman J: The cyclophilins. *J of Genome Biology*. Jun; 6(7): 226. (2005)


15. Jin ZG, Lungu AO, Xie L, Wang M, Wong Ch, Berk BC: Cyclophilin A Is a Proinflammatory cytokine that Activates Endothelial Cells. J of Arteriosclerosis, Thrombosis, and Vascular Biology. May; 24: 1186-1191. (2004)
16. Satoh K, Satoh T, Kikuchi N, Omura J, Kurosawa R, Suzuki K, Sugimura K, Aoki T, Nochioka K, Tatebe S, Miyamichi-Yamamoto S, Miura M, Shimizu T, Ikeda S, Yaoita N, Fukumoto Y, Minami T, Miyata S, Nakamura K, Ito H, Kadomatsu K, Shimokawa H: Basigin mediates pulmonary hypertension by promoting inflammation and vascular smooth muscle cell proliferation. Circ Res. Sep; 115(8): 738-750. (2014)

**»» For more references on this product see our WebPages at [www.biovendor.com](http://www.biovendor.com)**


## 19. EXPLANATION OF SYMBOLS

	Catalogue number
	Content
	Lot number
	Attention, see instructions for use
	Potential biological hazard
	Expiry date
	Storage conditions
	Name and registered office of the manufacturer

## Assay Procedure Summary



1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
	A	B	C	D	E	F	G	H



---

**BioVendor – Laboratorní medicína a.s.**

Karasek 1767/1, 621 00 Brno, Czech Republic

Phone: +420-549-124-185, Fax: +420-549-211-460

E-mail: [info@biovendor.com](mailto:info@biovendor.com), [sales@biovendor.com](mailto:sales@biovendor.com)

Web: [www.biovendor.com](http://www.biovendor.com)

There are BioVendor branches and distributors near you.

To find the office closest to you, visit [www.biovendor.com/contact](http://www.biovendor.com/contact)

