

## PRODUCT DATA SHEET

### L-Asparaginase

**Cat. No.:** RP1792872500

**Size:** 2500 IU

**Source:** E. coli

#### Other Names:

#### Description:

L-Asparaginase was purified from E. coli ASI.357

#### Purity:

Greater than 96.0% as determined by:  
(a) Analysis by RP-HPLC.

#### Formulation:

The enzyme was lyophilized with no additives.

#### Reconstitution:

It is recommended to reconstitute the lyophilized L-Asparaginase in 18MΩ water at 1 mg/ml.

#### Shipping:

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

#### Storage/Stability:

Two years when stored at -20°C, 2 weeks at 4°C.

#### Note:

BioVendor's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.