

## PRODUCT DATA SHEET

### **SARS-CoV-2 (2019-nCoV) Nucleocapsid Protein, Mosaic E.coli Recombinant**

**Cat. No.:** RP9720150050

**Size:** 0.05 mg

**Source:** E. coli

**Type:** Recombinant protein

**Species:** SARS

#### **Other Names:**

Severe acute respiratory syndrome coronavirus 2 nucleocapsid, 2019 novel coronavirus nucleoprotein, SARS-CoV-2 NP, SARS-CoV-2 N protein, COVID-19

#### **Description:**

The E.Coli derived recombinant protein contains the SARS-CoV-2 full length nucleocapsid Mosaic immunodominant regions, fused to His tag at the C-terminal.

#### **Purity:**

Greater than 90.0% as determined by Analysis by SDS-PAGE.

#### **Formulation:**

SARS-CoV-2 Nucleocapsid-Mosaic Protein 1mg/ml solution is supplied in 1x PBS.

#### **Storage/Stability:**

Shipped on ice packs. Store at -20°C upon arrival.

#### **Applications:**

COVID-19

#### **Note:**

BioVendor's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.