

## PRODUCT DATA SHEET

### SARS-CoV-2 (2019-nCoV) Spike-E-M Mosaic E.coli Recombinant

**Cat. No.:** RP9720140050

**Size:** 0.05 mg

**Source:** E. coli

**Type:** Recombinant protein

**Species:** SARS

**Other Names:**

Severe acute respiratory syndrome coronavirus 2 spike-envelope-membrane mosaic protein, 2019 novel coronavirus S-E-M mosaic protein, SARS-CoV-2 SEM mosaic protein, COVID-19

**Description:**

The E.Coli derived recombinant protein contains contains the SARS-CoV-2 spike (S), membrane (M), and envelope (E) immunodominant regions, fused to His tag at C-terminal.

**Purity:**

Greater than 90.0% as determined by Analysis by SDS-PAGE.

**Formulation:**

SARS-CoV-2 Spike Envelope-Mosaic Protein 1mg/ml solution is supplied in 1x PBS.

**Storage/Stability:**

Shipped on ice packs. Store at -20°C upon arrival.

**Applications:**

COVID-19

**Note:**

BioVendor's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.