

## MATERIAL SAFETY DATA SHEET

SARS-CoV-2 S1RBD IgG ELISA Kit      Version 2.0      Revision Date 04/15/2020

<b>SECTION 1. IDENTIFICATION</b>	
Product:	SARS-CoV-2 S1RBD IgG ELISA Kit
Company:	ImmunoDiagnostics Limited 5/F, Biotech Centre 2, No. 11 Science Park West Avenue Hong Kong Science Park, Sha Tin, Hong Kong
Tel:	(+852) 3502 2780/ Fax: (+852) 3502 2781
Email:	info@immunodiagnostics.com.hk
Emergency phone number:	Please contact your local Poison Control Centre Australia – 131126 New Zealand – (0)34747000 United States – 1800 222 1222
Recommended use:	This product is for laboratory research use only and is not intended for human/animal diagnostics, therapeutic or any other clinical uses.

### 2. INFORMATION ON PRODUCT/PREPARATION COMPOSITION

The product contains the following hazardous substances and those with the following highest permissible concentrations in the working environment:

#### 2.1. Microtiter strips (96 wells, coated with recombinant protein, sealed)

Identification numbers	Chemical name of substance	Concentrations	Risk and safety statements
CAS: None EC (EINECS): None	None		

#### 2.2. Negative control (mouse anti-S1RBD antibody, liquid)

Identification numbers	Chemical name of substance	Concentrations	Risk and safety statements
CAS: None EC (EINECS): None	None		

#### 2.3. HRP Labeled Antibody Solution (anti-human IgG antibody, HRP conjugate, in phosphate buffer solution)

Identification numbers	Chemical name of substance	Concentrations	Risk and safety statements
CAS: None EC (EINECS): None	None		

#### 2.4. Assay buffer concentrate

Identification numbers	Chemical name of substance	Concentrations	Risk and safety statements
CAS: None EC (EINECS): None	None		

#### 2.5. Wash buffer concentrate

Identification numbers	Chemical name of substance	Concentrations	Risk and safety statements
CAS: None EC (EINECS): None	None		

#### 2.6. Substrate solution

Identification numbers	Chemical name of substance	Concentrations	Risk and safety statements
CAS: 54827-17-7 EC (EINECS): 259-364-6 S26- 36/37	3,3',5,5'- Tetramethylbenzidine	0.4 g/L	R20/21/22-36/37/38- 40, S26-36/37

#### 2.7. Stop solution

Identification numbers	Chemical name of substance	Concentrations	Risk and safety statements
CAS: 7664-93-9 EC (EINECS):231- 639-5	Sulfuric acid	2 N	C, R35, S26-30-45

#### Risk symbols

C Corrosive

#### Risk phrases

R20/21/22 Harmful by inhalation, in contact with skin and if swallowed

R21/22 Harmful in contact with skin and if swallowed

R22 Harmful if swallowed

R35 Cause severe burns

R36/37/38 Irritating to eyes, respiratory system and skin

R40 Possible risk of irreversible effects

#### Safety phrases

S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice

S30 Never add water to this product

S36/37 Wear suitable protective clothing and gloves

S45 In case of accident or if you feel unwell, seek medical advice immediately (show label where possible)

---

### **3. SUBSTANCE/PREPARATION HAZARDS DATA**

The amount of dangerous substance is under the limit on which the kit is considered dangerous as a whole.

**When used, the most adverse impacts of the substance/preparation on human health include:**

The followed components can damage health if ingested and/or can irritate eyes and skin: Stop solution, Substrate solution.

---

### **4. FIRST AID MEASURES**

**Eye Contact:**

Check for and remove contact lenses. Flush eyes with water for at least 15 minutes. Get medical attention.

**Skin contact:**

Wash skin thoroughly with soap and water. Remove and wash contaminated clothing. Should irritation occur get medical attention.

**Inhalation:**

If inhaled remove to fresh air. If breathing is affected seek medical attention.

**Ingestion:**

Flush out mouth and drink large amounts of water. Seek medical attention if symptoms appear.

**Wounds:**

Allow to bleed freely. Wash thoroughly with soap and water. Get medical attention.

---

### **5. FIRE FIGHTING MEASURES**

**Suitable fire-extinguishing media:** Carbon dioxide, dry powder, foam, water.

**Thermal decomposition:** No thermal decomposition degradation products are expected.

**Special hazards:** None.

**Special protective means for firemen:** None.

---

### **6. ACCIDENTAL RELEASE MEASURES**

**Safety measures to protect humans:** Avoid contact with skin and eyes.

**Environmental safety measures:** Avoid penetration into sewerage systems, surface and ground water. Avoid soil pollution.

**Recommended cleaning and disposal methods:** Cover with suitable absorbing material. After removing the substance, rinse the spot of spilling thoroughly with water and soap.

---

### **7. HANDLING AND STORAGE**

**Handling instructions:** Avoid contact with skin, eyes and clothing. Use suitable protective means to work with the substance.

**Storage instructions:** Store at temperatures between + 2 and + 8°C in a dry and dark place.

---

## 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

**Technical measures:** Do not eat, drink and smoke when working with the kit. Use the kit only in rooms enabling good ventilation. Local exhaustion is necessary, general (forced) exhaustion is recommended.

**Personal protective means – protection of respiratory organs:** None

**Personal protective means – eye protection:** None

**Personal protective means – hand protection:** Protective gloves (wash your hands before and after work)

**Personal protective means – body protection:** Protective clothing

---

## 9. PHYSICAL AND CHEMICAL PROPERTIES

**State (at 20°C):**

Solid: Microtiter strips

Liquids: negative control, HRP labeled antibody Solution, Assay buffer concentrate, Wash buffer concentrate, Substrate solution, Stop solution

**Color:**

Yellow: HRP labeled antibody Solution, Assay buffer concentrate.

Colorless: negative control, Stop solution, Wash buffer concentrate, Substrate solution.

**pH value (at 25°C):** Stop solution: < 1

Others: 6.8 - 7.4

---

## 10. STABILITY AND REACTIVITY

**Conditions to be avoided:** Heat

**Substances and materials with which the product is not allowed to get in touch:**

Acids

**Hazardous decomposition products:** Not known

---

## 11. TOXICOLOGICAL INFORMATION

**Toxicological information**

The product contains sulphuric acid in total concentration of 2N.

**Acute toxicity-LD50 orally, rat (mg/kg):** 2140mg/kg - sulphuric acid

**Irritability:** No data available

**Mutagenicity:** No data available

**Reproduction toxicity:** No data available

**Tests on animals:** No data available

---

## 12. ECOLOGICAL INFORMATION

Water hazard class: 2 (water-damaging substances)

---

## 13. DISPOSAL CONSIDERATIONS

**The manner of disposing the substance/preparation:**

Mix or dissolve the material in a combustible solvent and burn up in a facility whose equipment matches all regulations in effect. Every waste disposal must be carried out in coincidence with national and local legislation or administrative regulations respectively.

---

#### **14. TRANSPORT INFORMATION**

**Overland transportation (ADR/RID):** As a mixture, the substance is subject to no limitations.

**Transatlantic transportation (IMDG):** As a mixture, the substance is subject to no limitations.

**Air transportation (ICAO/IATA):** As a mixture, the substance is subject to no limitations

---

#### **15. INFORMATION ON LEGAL REGULATIONS**

According to the Act No. 356/2003 Coll. on chemical substances and chemical preparations and on amendments in some other laws and acts, components contained in the SARS-CoV-2 S1RBD IgG ELISA Kit are non-hazardous substances in said concentrations.

SARS-CoV-2 S1RBD IgG ELISA Kit is subject to special marking regulations according to the EC regulations (28/10/1999).

SARS-CoV-2 S1RBD IgG ELISA Kit:

R20/21/22-21/22-22-35-36/37/38-40

S26-30-36/37-45

---

#### **16. OTHER INFORMATION**

*The material safety data sheet contains data necessary to ensure safety and health and environmental protection in working with chemical substances. The above-stated data match the contemporary state of knowledge and experience and are in coincidence with legal regulations currently in effect. This product is a chemical substance and can be solely used by persons with chemical education at their own risk.*

*SARS-CoV-2 S1RBD IgG ELISA Kit is designed for biomedical research. The manufacturer has no responsibility for damage caused by unsuitable use and by disrespecting the enclosed working instructions.*

*The above-stated information cannot be considered as complete and must be understood to be only a methodical instruction.*