

Prolactin Murine E. coli

Product Data Sheet

Type: Recombinant	Cat. No.:	
Source: E. coli	RBG20112010	(10 µg)
Species: Murine	RBG20112050	(50 µg)
	RBG20112100	(100 µg)

Description

Prolactin is a neuroendocrine hormone secreted by the pituitary gland. Its primary function is to promote and maintain lactation during pregnancy and suckling. In addition, prolactin plays an immunoregulatory role by stimulating the activities of ornithine decarboxylase and protein kinase C, which are important for the proliferation, differentiation, and function of lymphocytes. Recombinant Murine Prolactin is a 22.5 kDa globular protein containing 198 amino acid residues.

Introduction to the Molecule

Prolactin (about 24 kDa, 199 amino acids) is a peptide hormone primarily associated with lactation. It is synthesised and secreted by lactotrope cells in the adenohypophysis. It is also produced in other tissues including the breast and the decidua. Prolactin stimulates the mammary glands to produce milk (lactation).

Research topic

Animal studies

Amino Acid Sequence

MLPICSAGDC QTSRLRELFDR VVILSHYIHT LYTDMFIEFD KQYVQDREFM VKVINDCPTS SLATPEDKEQ ALKVPPEVLL
NLILSLVQSS SDPLFQLITG VGGIQEAPFY ILSRAKEIEE QNKQLLEGVE KIISQAYPEA KNGIYFVWS QLPSLQGVDE
ESKILSLRNT IRCLRRDSHK VDNFLKVLRC QIAHQNNC

Source

E. coli

Purity

98%

Biological Activity

Determined by its ability to induce the proliferation of rat Nb2-11 cells in the concentration range of 0.1-1.0 ng/ml.

Endotoxin

Endotoxin level is <0.1 ng/µg of protein (<1EU/µg).

Reconstitution

Centrifuge the vial prior to opening. Reconstitute in water to a concentration of 0.1-1.0 mg/ml. Do not vortex. For extended storage, it is recommended to further dilute in a buffer containing a carrier protein (example 0.1% BSA) and store in working aliquots at -20°C to -80°C.

Storage, Stability/Shelf Life

HEADQUARTERS:				
BioVendor Laboratorní medicína, a.s.	Karasek 1767/1	621 00 Brno CZECH REPUBLIC	Phone: +420-549-124-185 Fax: +420-549-211-460	E-mail: info@biovendor.com sales@biovendor.com Web: www.biovendor.com
AUSTRIA: BioVendor GesmbH	Nußdorfer Straße 20/10	1090 Vienna AUSTRIA	Phone: +43-1-89090-25 Fax: +43-1-89090-2515	E-mail: infoAustria@biovendor.com
GERMANY, SWITZERLAND: BioVendor GmbH	Otto-Hahn-Straße 16	34123 Kassel GERMANY	Phone: +49-6221-433-9100 Fax: +49-6221-433-9111	E-mail: infoEU@biovendor.com
USA, CANADA AND MEXICO: BioVendor LLC	128 Bingham Rd. Suite 1300	Asheville, NC 28806 USA	Phone: +1-828-575-9250 +1-800-404-7807 Fax: +1-828-575-9251	E-mail: infoUSA@biovendor.com