

Rabbit IgG Anti SARS-CoV-2 Spike (S1) Antibody, Clone: CR3022

Product Data Sheet

Host: Rabbit

Clone: CR3022

Isotype: IgG, kappa

Cat. No.:

RM812424100

(0.1 mg)

Other names: Severe acute respiratory syndrome coronavirus 2 spike glycoprotein S1, 2019 novel coronavirus S1 protein, SARS-CoV-2 S1 subunit, COVID-19

Research topic

Immune Response, Infection and Inflammation, COVID-19

Preparation

This is a recombinant monoclonal antibody. It is a reformatted antibody made using the variable domain sequences of the original Human IgG1 format. The original monoclonal antibody was generated by sequencing peripheral blood lymphocytes of a patient exposed to the SARS-CoV.

Amino Acid Sequence

Antibody binds to both SARS-CoV and SARS-CoV-2 (COVID-19) with high affinity at amino acids 318-510 in the S1 domain of the Spike protein.

Species Reactivity

SARS-CoV and SARS-CoV-2 (COVID-19)

Purification Method

Affinity chromatography on a column with immobilized protein A.

Antibody Content

0.1 mg; concentration: 1.0 mg/ml

Formulation

Supplied as a liquid in PBS, 0.02% Proclin 300.

Shipping

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

Storage/Stability

Short Term Storage: Up to 3 months at +4°C

Long Term Storage: -20°C

Expiration

See vial label.

Lot Number

See vial label.

Applications

ELISA, NTRL, SPR, Crystallography

Note

Products are for Research Use or for Further Manufacturing Use only. Not for Diagnostic or Therapeutic Use.

HEADQUARTERS:				
BioVendor Laboratorní medicína, a.s.	Karasek 1767/1	621 00 Brno CZECH REPUBLIC	Phone: +420-549-124-185 Fax: +420-549-211-460	E-mail: info@biovendor.com sales@biovendor.com Web: www.biovendor.com
AUSTRIA: BioVendor GesmbH	Nußdorfer Straße 20/10	1090 Vienna AUSTRIA	Phone: +43-1-89090-25 Fax: +43-1-89090-2515	E-mail: infoAustria@biovendor.com
GERMANY, SWITZERLAND: BioVendor GmbH	Otto-Hahn-Straße 16	34123 Kassel GERMANY	Phone: +49-6221-433-9100 Fax: +49-6221-433-9111	E-mail: infoEU@biovendor.com
USA, CANADA AND MEXICO: BioVendor LLC	128 Bingham Rd. Suite 1300	Asheville, NC 28806 USA	Phone: +1-828-575-9250 +1-800-404-7807 Fax: +1-828-575-9251	E-mail: infoUSA@biovendor.com