

Mouse Anti Influenza A Antibody, Clone: Y77

Product Data Sheet

Host: Mouse

Clone: Y77

Isotype: Mouse IgG₁

Cat. No.:

RM812163100

(0.1 mg)

Species Reactivity

This antibody is specific for *Influenza A* detection. The antibody does not cross react with *Adenovirus* or *Respiratory Syncytial virus*.

Purification Method

The antibody was purified by affinity chromatography on protein A.

Antibody Content

0.1 mg; concentration: 1.0 mg/ml

Formulation

Supplied as a liquid in Phosphate Buffered Saline pH 7.4, 0.09% Sodium Azide.

Shipping

At ambient temperature.

Storage/Stability

Store at +4°C for up to three months, or at -20°C for longer. Avoid repeated freeze/thaw cycles.

Applications

ELISA

Note

For research use only.

| | | | | |
|---|-------------------------------|-------------------------------|---|---|
| HEADQUARTERS: BioVendor Laboratorní medicína, a.s. | Karasek 1767/1 | 621 00 Brno CZECH REPUBLIC | Phone: +420-549-124-185 Fax: +420-549-211-460 | E-mail: info@biovendor.com sales@biovendor.com Web: www.biovendor.com |
| AUSTRIA: BioVendor GesmbH | Nußdorfer Straße 20/10 | 1090 Vienna AUSTRIA | Phone: +43-1-89090-25 Fax: +43-1-89090-2515 | E-mail: infoAustria@biovendor.com |
| GERMANY, SWITZERLAND: BioVendor GmbH | Otto-Hahn-Straße 16 | 34123 Kassel GERMANY | Phone: +49-6221-433-9100 Fax: +49-6221-433-9111 | E-mail: infoEU@biovendor.com |
| USA, CANADA AND MEXICO: BioVendor LLC | 128 Bingham Rd. Suite 1300 | Asheville, NC 28806 USA | Phone: +1-828-575-9250 +1-800-404-7807 Fax: +1-828-575-9251 | E-mail: infoUSA@biovendor.com |