

## Mouse Anti Respiratory Syncytial Virus Antibody, Clone: RV4

### Product Data Sheet

**Host:** Mouse

**Clone:** RV4

**Isotype:** Mouse IgG<sub>1</sub>

**Cat. No.:**

RM812148100

(0.1 mg)

### Species Reactivity

This antibody is specific for respiratory syncytial virus. The antibody does not cross react with Influenza Type A, Influenza Type B or Adenovirus.

### Purification Method

The antibody was purified by affinity chromatography on protein A.

### Antibody Content

0.1 mg; concentration: 1.0 mg/ml

### Formulation

Supplied as a liquid in Phosphate Buffered Saline pH 7.2, <0.1% Sodium Azide.

### Shipping

At ambient temperature.

### Storage

Short term storage at 2-8°C.

Longer term storage at -20°C.

Avoid repeated freeze/thaw cycles.

### Applications

ELISA

### Note

For research use only.

#### HEADQUARTERS:

BioVendor Laboratorní  
medicína, a.s.

Karasek 1767/1

621 00 Brno  
CZECH REPUBLIC

Phone: +420-549-124-185  
Fax: +420-549-211-460

E-mail: info@biovendor.com  
sales@biovendor.com

Web: www.biovendor.com

AUSTRIA:  
BioVendor GesmbH

Nußdorfer Straße 20/10

1090 Vienna  
AUSTRIA

Phone: +43-1-89090-25  
Fax: +43-1-89090-2515

E-mail: infoAustria@biovendor.com

GERMANY, SWITZERLAND:  
BioVendor GmbH

Otto-Hahn-Straße 16

34123 Kassel  
GERMANY

Phone: +49-6221-433-9100  
Fax: +49-6221-433-9111

E-mail: infoEU@biovendor.com

USA, CANADA AND MEXICO:  
BioVendor LLC

128 Bingham Rd.  
Suite 1300

Asheville, NC 28806  
USA

Phone: +1-828-575-9250  
+1-800-404-7807  
Fax: +1-828-575-9251

E-mail: infoUSA@biovendor.com