

## Mouse Anti Rhinovirus VP3 Antibody, Clone: 4967

### Product Data Sheet

<b>Host:</b> Mouse	<b>Cat. No.:</b>	
<b>Clone:</b> 4967	RM812384500	(0.5 mg)
<b>Isotype:</b> Mouse IgG <sub>1</sub>		

### Species Reactivity

Specific for the VP3 of numerous Rhino types (see table below). Not reactive with other Enteroviruses, Influenza A, Influenza B, RSV, Adenovirus, Para 1, Para 2, Para 3, Chlamydia pneumoniae, Mycoplasma pneumoniae & Hepatitis A virus.

### Purification Method

Purified preparations consist of >90% pure mouse monoclonal antibody which has been purified from ascites fluid or culture medium by protein A chromatography or sequential differential precipitations.

### Antibody Content

0.5 mg; concentration: 1.0 mg/ml

### Formulation

Supplied as a liquid in Phosphate Buffered Saline pH 7.2, 0.1% Sodium Azide.

### Shipping

At ambient temperature.

### Storage

This product contains no stabilizing proteins and should be stored at 2-8°C until ready for use.

### Applications

ELISA, IFA. Working dilution must be determined by the user. Suggested starting ranges are 1:10-1:50 for IFA and 1:20-1:200 for ELISA.

### Note

For research use only.

**Data:** Reactivity testing was carried out by indirect fluorescent-antibody (FA) at 10 µg/ml.

Product	Cat. No.:	Rhinovirus type									
		1a	2	14	16	17	27	42	70	80	
Mouse Anti Rhinovirus VP3 Antibody, Clone: 4965	RM812382100	+	+	+	+	+	-	0	+	0	
Mouse Anti Rhinovirus VP3 Antibody, Clone: 4965	RM812382500	+	+	+	+	+	-	0	+	0	
Mouse Anti Rhinovirus VP3 Antibody, Clone: 4967	RM812384100	+	+	+	+	+	0	0	+	+	
Mouse Anti Rhinovirus VP3 Antibody, Clone: 4967	RM812384500	+	+	+	+	+	0	0	+	+	

**HEADQUARTERS:**

BioVendor Laboratorní medicína, a.s.	Karasek 1767/1	621 00 Brno CZECH REPUBLIC	Phone: +420-549-124-185 Fax: +420-549-211-460	E-mail: info@biovendor.com sales@biovendor.com Web: www.biovendor.com
AUSTRIA: BioVendor GesmbH	Nußdorfer Straße 20/10	1090 Vienna AUSTRIA	Phone: +43-1-89090-25 Fax: +43-1-89090-2515	E-mail: infoAustria@biovendor.com
GERMANY, SWITZERLAND: BioVendor GmbH	Otto-Hahn-Straße 16	34123 Kassel GERMANY	Phone: +49-6221-433-9100 Fax: +49-6221-433-9111	E-mail: infoEU@biovendor.com
USA, CANADA AND MEXICO: BioVendor LLC	128 Bingham Rd. Suite 1300	Asheville, NC 28806 USA	Phone: +1-828-575-9250 +1-800-404-7807 Fax: +1-828-575-9251	E-mail: infoUSA@biovendor.com