



Endotoxin-free Water

Ultra-pure water for reconstitution and dilution of test samples and reagents (<0.001 EU/mL)

EndoGrade[®] Water is specially processed water for sensitive applications where low endotoxin levels are required. **EndoGrade[®] Water** is certified to contain less than 0.001 EU/mL and is ideal for preparing solutions and reagents for bacterial endotoxin testing and endotoxin removal.



EndoGrade[®] Water at a Glance

- Ultra-purified water ideal for endotoxin-free reconstitution and dilution of samples and reagents used in bacterial endotoxin testing and endotoxin removal
- Certified endotoxin content of less than 0.001 EU/mL
- Convenient plastics bottles with screw caps
- Shelf life of 5 years from date of certification

Distributed By:

BioVendor
Research
and Diagnostic Products





EndoGrade® Water - Specifications

Cat. No. 607100

- 16 bottles/pack
- Certified endotoxin-free water, <0.001 EU/ml



EndoLISA®

Cat. No. 609033

Endotoxin Detection Assay based on ELISA-technology and Recombinant Factor C (rFC)



EndoZyme® II

Cat. No. 890030

Recombinant Factor C (rFC) Endotoxin Detection Assay



Endo-RS®

Cat. No. 609065

Endotoxin Recovery Kit - A Toolbox for Demasking Endotoxin in Pharmaceutical Formulations exhibiting Low Endotoxin Recovery (LER)



EndoGrade® Glass Test Tubes

Cat. No. 800050

Endotoxin-free borosilicate glass test tubes with aluminium screw cap

Additionally Hyglos offers Endotoxin Detection Service including Hold Time Studies and Demasking Protocol Development applying state-of-the-art methods.

Contact BioVendor for more information!

Distributed By:

BioVendor

Research
and Diagnostic Products

BioVendor, LLC
128 Bingham Rd., Suite 1300
Asheville, NC 28806
Tel: 828-575-9250
Toll-free: 1800-404-7807
E-mail: infoUSA@BioVendor.com
www.BioVendor.com



A BIOMÉRIEUX COMPANY

