

- Quantitative measurement of both active and total human MMP-2
- 96-well plate format
- Colorimetric read-out at 405 nm
- Ability to detect high and low ranges in one assay
- Easy and sensitive
- Quantitative replacement for zymography
- Now based on monoclonal antibody

## Assay specifications

**Dynamic range:** 0.02 - 16 ng/ml

**Sensitivity:** 0.04 ng/ml (2 hr incubation)  
0.02 ng/ml (6 hr incubation)  
4 pg/ml (o/n incubation)


**Cross reactivity:** MMP-1, -3, -8, -9 < 0.05%,  
-13 < 1%, -7 < 7%

**Sample volume:** 10 - 100  $\mu$ l


**Format:** 96 - well plate

**Intended use:** Cell culture media, cell and tissue extracts, serum, urine

## Quantitative measurement



## Assay principle



Active MMP-2: direct measurement

Pro MMP-2: measurement after activation with APMA

## What's in the box?

**All standards, buffers and plates are included:** Pre-coated 96-well microtiter plate, assay buffer, APMA, pro-detection enzyme, peptide substrate, wash buffer, human MMP-2 standard.

## Related assays

- Human and mouse MMP-9 activity assay
- Human MMP-7, -8, -14 activity assays
- Soluble collagen assay
- Total collagen assay
- Hydroxyproline assay

## Distributor information

BioVendor, LLC

128 Bingham Road, Suite 1300

Asheville, NC 28806, United States

Tel: 828-575-9250 / Toll-free: 1800-404-7807

www.BioVendor.com

E-mail: infoUSA@BioVendor.com

**BioVendor**

Research  
and Diagnostic Products