

SAFETY INFORMATION

Product name: HUMAN INTERLEUKIN-8 ELISA
Catalogue number: RD194558200R
Characteristics: The RD194558200R Human Interleukin-8 ELISA is a sandwich enzyme immunoassay for the quantitative measurement of human IL-8.
Intended use: For professional use only. Users should have a thorough understanding of the IFU prior to their use of this kit.

Manufacturer: BioVendor – Laboratorní medicína a.s.
Registered office: Karásek 1767/1, 621 00 Brno, Czech republic
Company ID: 63471507

Kit Components	Contents hazardous substance
Antibody Coated Microtiter Strips	No
Biotin Labelled Antibody Conc. (100x)	No
Streptavidin-HRP Conjugate	StabilZyme™ HRP Conjugate Stabilizer
Master Standard	No
Biotin-Ab Diluent	No
Dilution Buffer	No
Wash Solution Conc. (10x)	No
Substrate Solution	No
Stop Solution	Sulphuric Acid < 5%

Safety Data Sheets in accordance with the current (EC/EU) Regulations as amended are attached.

SECTION 1	IDENTIFICATION OF THE SUBSTANCE / MIXTURE AND OF THE COMPANY / UNDERTAKING
------------------	---

1.1 Product identifier

Trade name: StabilZyme™ HRP Conjugate Stabilizer

Code: SZ02

Product form: mixture

1.2 Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses Use of the substance/mixture

Use of the substance/mixture: For further manufacture

Uses advised against

No additional information available

1.3 Details of the supplier of the safety data sheet

BioVendor - Laboratorní medicína a.s.

Karásek 1767/1

621 00 Brno

Czech Republic

Identification number: 63471507

Tel: +420 549 124 185

E-mail: info@biovendor.com

Website: www.biovendor.com

1.4 Emergency telephone number

European Chemicals Agency. National helpdesks contact details

<https://echa.europa.eu/support/helpdesks>

Links to Poison Centers and Clinical Toxicologists all over the World:

<https://www.eapcct.org/index.php?page=links>

SECTION 2	HAZARDS IDENTIFICATION
------------------	-------------------------------

2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin sensitisation, Category 1 H317

Full text of H statements: see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

2.2 Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms



GHS07

Signal word (CLP): Warning

Hazardous ingredients:

Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7]; and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1); 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6]

Hazard statements:

H317 May cause an allergic skin reaction

Precautionary statements:

P261 Avoid breathing mist, vapours

P272 Contaminated work clothing should not be allowed out of the workplace

P280 Wear eye protection, protective gloves, protective clothing

P302+P352 IF ON SKIN: Wash with plenty of soap and water

P321 Specific treatment (see supplemental first aid instruction on this label)

P333+P313 If skin irritation or rash occurs: Get medical advice/attention

P362+P364 Take off contaminated clothing and wash it before reuse

P501 Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation

EUH-statements:

EUH210 Safety data sheet available on request

EUH208 Contains CMIT/MIT. May produce an allergic reaction

2.3 Other hazards

Other hazards not contributing to the classification: None under normal conditions.

SECTION 3 COMPOSITION / INFORMATION ON INGREDIENTS
3.1 Substance

Not applicable

3.2 Mixtures

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
CMIT/MIT: Mixture, 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7]; and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1)	(CAS-No.) 55965-84-9 (EC Index-No.) 613-167-00-5	0.0024	Acute Tox. 3 (Oral), H301 Acute Tox. 2 (Dermal), H310 Acute Tox. 3 (Inhalation), H330 Eye Dam. 1, H318
2-methyl-2H-isothiazol-3-one [EC no. 220-239-6]	(CAS-No.) 2682-20-4 (EC Index-No.) 613-167-00-5	0.0126	Skin Corr. 1C, H314 Skin Sens. 1A, H317 Aquatic Acute 1, H400 Aquatic Chronic 1, H410

Specific concentration limits:

Name	Product identifier	Specific concentration limits: CLP
CMIT/MIT: Mixture, 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7]; and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1)	CAS-No.) 55965-84-9 (EC Index-No.) 613-167-00-5	(C ≥ 0.0015) Skin Sens. 1A, H317 (C ≥ 0.06%) Eye Dam. 1, H318

2-methyl-2H-isothiazol-3-one [EC no. 220-239-6]	(CAS-No.) 2682-20-4 (EC Index-No.) 613-167-00-5	(C >= 0.6) Skin Corr. 1C, H314
---	--	--------------------------------

Full text of R- and H-statements: see section 16

SECTION 4 FIRST AID MEASURES

4.1 Description of first aid measures

First-aid measures general:

If exposed or concerned, get medical attention/advice. Show this safety data sheet to the doctor in attendance. Wash contaminated clothing before re-use. Never give anything to an unconscious person

First-aid measures after inhalation: IF INHALED: Remove to fresh air and keep at rest in a comfortable position for breathing

First-aid measures after skin contact: IF ON SKIN (or clothing): Remove affected clothing and wash all exposed skin with water for at least 15 minutes. If irritation develops or persists, get medical attention.

First-aid measures after eye contact:

IF IN EYES: Immediately flush with plenty of water for at least 15 minutes. Remove contact lenses if present and easy to do so. Continue rinsing if pain, blinking, or irritation develops or persists, get medical attention. Continue rinsing.

First-aid measures after ingestion:

IF SWALLOWED: rinse mouth thoroughly. Do not induce vomiting without advice from poison control center. Get medical attention if you feel unwell.

4.2 Most important symptoms and effects, both acute and delayed

Symptoms/effects:

May cause an allergic skin reaction.

Symptoms/effects after inhalation:

May cause respiratory irritation.

Symptoms/effects after skin contact:

May cause an allergic skin reaction.

Symptoms/effects after eye contact:

Direct contact with eyes is likely to be irritating.

Symptoms/effects after ingestion:

May cause gastrointestinal irritation.

4.3 Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5 FIREFIGHTING MEASURES

5.1 Extinguishing media

Suitable extinguishing media : Foam. Dry powder. Carbon dioxide. Water spray. Sand.

5.2 Special hazards arising from the substance or mixture

Fire hazard:

No data available.

Explosion hazard:

No data available.

5.3 Advice for firefighters

Firefighting instructions:

Use water spray or fog for cooling exposed containers. Exercise caution when fighting any chemical
Do not dispose of fire-fighting water in the environment.

Protection during firefighting:

Do not enter fire area without proper protective equipment, including respiratory protection.

SECTION 6 ACCIDENTAL RELEASE MEASURES**6.1 Personal precautions, protective equipment, and emergency procedures**

General measures: Evacuate area. Keep upwind. Ventilate area. Spill should be handled by trained clean-up crews properly equipped with respiratory equipment and full chemical protective gear (see Section 8).

6.1.1 For non-emergency personnel

Protective equipment: Wear Protective equipment as described in Section 8.

Emergency procedures: Evacuate unnecessary personnel.

6.1.2 For emergency responders

Protective equipment: Wear suitable protective clothing, gloves and eye or face protection. Approved supplied-air respirator, in case of emergency.

6.2 Environmental precautions

Prevent entry to sewers and public waters. Notify authorities if liquid enters sewers or public waters. Avoid release to the environment.

6.3 Methods and material for containment and cleaning up**For containment:**

Contain any spills with dikes or absorbents to prevent migration and entry into sewers or streams.

Methods for cleaning up

Place in a suitable container for disposal in accordance with the waste regulations (see Section 13). Wash spill area thoroughly with plenty of soap and water.

6.4 Reference to other sections

No additional information available

SECTION 7 HANDLING AND STORAGE**7.1 Precautions for safe handling**

Precautions for safe handling: Do not handle until all safety precautions have been read and understood. Use appropriate personal protection equipment (PPE). Keep away from sources of ignition - No smoking. Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work.

7.2 Conditions for safe storage, including any incompatibilities

Storage conditions: Keep the container tightly closed. Store in a dry, cool and well-ventilated place. Avoid elevated temperatures.

Storage temperature: Refrigerated or room temperature

7.3 Specific end use(s)

No additional information available

SECTION 8 EXPOSURE CONTROLS / PERSONAL PROTECTION
8.1 Control parameters

Mixture , 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7]; and 2-methyl-2H -isothiazol-3-one [EC no. 220-239-6] (3:1) (CAS-No. 55965-84-9)		
Germany	TRGS 900 Occupational exposure limit value (mg/m ³)	0.2 mg/m ³ inhalable fraction
Germany	TRGS 900 Limitation of exposure peaks (mg/m ³)	0.4 mg/m ³ inhalable fraction

8.2 Exposure controls
8.2.1 Appropriate engineering controls:

Provide adequate general and local exhaust ventilation. Use process enclosures, local exhaust ventilation, or other engineering controls to control airborne levels below recommended exposure limits. Use explosion-proof equipment with flammable materials. Ensure adequate ventilation, especially in confined areas.

8.2.2 Individual protection measures, such as personal protective equipment

Gloves. Protective goggles. Protective clothing.

a) Eye protection:


Wear eye protection, including chemical splash goggles and a face shield when possibility exists for eye contact due to spraying liquid or airborne particles. END 166.

b) Skin protection

Body protection:

Wear long sleeves, and chemically impervious PPE/coveralls to minimize bodily exposure


Hand protection:

Use gloves chemically resistant to this material when prolonged or repeated contact could occur. Gloves should be classified under Standard EN 374 or ASTM

F1296. Suggested glove materials are: Neoprene, Nitrile/butadiene rubber, Polyethylene, Ethyl vinyl alcohol laminate, PVC or vinyl.

c) Respiratory protection:


Use NIOSH (or other equivalent national standard) -approved dust/particulate respirator. Where vapor, mist, or dust exceed PELs or other applicable OELs, use NIOSH-approved respiratory protective equipment. EN 250.

SECTION 9 PHYSICAL AND CHEMICAL PROPERTIES**9.1 Information on basic physical and chemical properties**

Physical state:	Liquid
Appearance:	Clear.
Colour:	Amber.
Odour:	None.
Odour threshold:	No data available
pH:	Neutral
Relative evaporation rate (butylacetate=1):	No data available
Melting point:	No data available
Freezing point:	No data available
Boiling point:	100 °C
Flash point:	No data available
Auto-ignition temperature:	No data available
Decomposition temperature:	No data available
Flammability (solid, gas):	No data available
Vapor pressure:	No data available
Relative vapor density at 20 °C:	No data available
Relative density:	No data available
Solubility:	Water: Complete
Log Pow:	No data available
Viscosity, kinematic:	No data available
Viscosity, dynamic:	No data available
Explosive properties:	Product is not explosive.
Oxidising properties:	Not an Oxidizer.
Explosive limits:	No data available

9.2 Other information

No additional information available

SECTION 10 STABILITY AND REACTIVITY**10.1 Reactivity**

No dangerous reactions known under normal conditions of use.

10.2 Chemical stability

Stable under recommended handling and storage conditions (see section 7).

10.3 Possibility of hazardous reactions

None known.

10.4 Conditions to avoid

None known.

10.5 Incompatible materials

None known.

10.6 Hazardous decomposition products

Thermal decomposition generates: Carbon oxides (CO, CO₂). Nitrogen oxides. Other toxic fumes.

SECTION 11 TOXICOLOGICAL INFORMATION
11.1 Information on toxicological effects

Acute toxicity:	Not classified
Skin corrosion/irritation:	Not classified pH: Neutral
Serious eye damage/irritation:	Not classified pH: Neutral
Respiratory or skin sensitisation:	May cause an allergic skin reaction.
Germ cell mutagenicity:	Not classified
Carcinogenicity:	Not classified
Reproductive toxicity:	Not classified
STOT-single exposure:	Not classified
STOT-repeated exposure:	Not classified
Aspiration hazard:	Not classified

SECTION 12 ECOLOGICAL INFORMATION
12.1 Toxicity

Ecology - general: No data available.

Mixture , 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7]; and 2-methyl-2H - isothiazol-3-one [EC no. 220-239-6] (3:1) (CAS-No. 55965-84-9)	
LC50 fish 1	0.19 mg/l Salmo Gairdneri
EC50 Daphnia 1	0.028 mg/l Crassostrea virg.
EC50 72h algae (1)	0.018 mg/l Selenastrum cap.

12.2 Persistence and degradability

Mixture , 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7]; and 2-methyl-2H - isothiazol-3-one [EC no. 220-239-6] (3:1) (CAS-No. 55965-84-9)	
Persistence and degradability	Not rapidly degradable

12.3 Bioaccumulative potential

Mixture , 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7]; and 2-methyl-2H - isothiazol-3-one [EC no. 220-239-6] (3:1) (CAS-No. 55965-84-9)	
Log Kow	> 5 Significant bioaccumulation

12.4 Mobility in soil

No additional information available

12.5 Results of PBT and vPvB assessment

No additional information available

12.6 Other adverse effects

No additional information available

SECTION 13 DISPOSAL CONSIDERATIONS
13.1 Disposal methods

Waste treatment methods: Obtain the consent of pollution control authorities before discharging to wastewater treatment plants.

Waste disposal recommendations:

Dispose in a safe manner in accordance with local/national regulations. Do not allow the product to be released into the environment.

SECTION 14 TRANSPORT INFORMATION**14.1 UN number**

Not hazardous for transport

14.2 UN proper shipping name

Proper Shipping Name (ADR): Not applicable

Proper Shipping Name (IMDG): Not applicable

Proper Shipping Name (IATA): Not applicable

Proper Shipping Name (ADN): Not applicable

Proper Shipping Name (RID): Not applicable

14.3 Transport hazard class(es)

ADR Transport hazard class(es) (ADR): Not applicable

IMDG Transport hazard class(es) (IMDG): Not applicable

IATA Transport hazard class(es) (IATA): Not applicable

ADN Transport hazard class(es) (ADN): Not applicable

RID Transport hazard class(es) (RID): Not applicable

14.4 Packaging group

Packing group (ADR) : Not applicable

Packing group (IMDG) : Not applicable

Packing group (IATA) : Not applicable

Packing group (ADN) : Not applicable

Packing group (RID) : Not applicable

14.5 Environmental hazards

Dangerous for the environment: No

Marine pollutant : No

Other information: No supplementary information available

14.6 Special precautions for user**Overland transport**

Not applicable

Transport by sea

Not applicable

Air transport

Not applicable

Inland waterway transport

Carriage prohibited (ADN): No

Not subject to AND: No

Rail transport

Carriage prohibited (RID): No

14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

SECTION 15 REGULATORY INFORMATION**15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture****15.1.1 EU-Regulations**

Contains no substances with Annex XVII restrictions StabilZyme® HRP Conjugate Stabilizer is not on the REACH Candidate List Contains no substance on the REACH Candidate List

Contains no REACH Annex XIV substances.

15.1.2 National regulations

All chemical substances in this product are listed in the EPA (Environment Protection Agency) TSCA (Toxic Substances Control Act) Inventory

Germany

VwVwS Annex reference: Water hazard class (WGK) 3, severe hazard to waters (Classification according to VwVwS, Annex 4)

12th Ordinance Implementing the Federal Immission : Is not subject of the 12. BImSchV (Hazardous Incident Ordinance)
Control Act - 12.BImSchV

15.2 Chemical Safety Assessment

No additional information available

SECTION 16 OTHER INFORMATION**Warning:**

The safety data sheet contains data needed for securing safety and health protection during work and environ-mental protection. The stated data correspond to the current state of knowledge and experience and is in accord-ance with valid legal regulations. It cannot be deemed as a guarantee of the properties, suitability, and usefulness of the product for specific application and therefore no contractual legal relationships are hereby created.

The safety data sheet is the property of the physical or legal entity stated in Section 1 and is protected by copy-right. All copying, distribution or sales without the consent of the owner is forbidden.

Relevant phrases:

H301 Toxic if swallowed
H310 Toxic in contact with skin
H314 Causes severe skin burns and eye damage
H317 May cause an allergic skin reaction
H318 Causes serious eye damage
H330 Fatal if inhaled
H400 Very toxic to aquatic life
H410 Very toxic to aquatic life with long lasting effects
EUH208 Contains . May produce an allergic reaction
EUH210 Safety data sheet available on request

Training hints:

Pursuant to article No 35 of the European Parliament and Council Regulation (ES) No 1907/2006, the employer must allow employees or their representatives access to information from the safety data sheet of the substance or preparation, which the employees use or to the effects of which they may be exposed during their work.

Physical entities performed individual activities within the scope of handling of this hazardous product are trained and regularly, at least once a year, retrained.

Product information sources: safety data sheet, product or technical information, safety instructions, and other ex-pert documents for the product, issued by the supplier.

Recommended restriction of use:

The product is designed only for professional purposes. It must not be used in households. The product can only be handled by a person older than 18 years, who is sufficiently informed about the work procedures, hazardous properties of the product, and also about the necessary safety measures.

The product is to be used only for the purpose, for which it is designed. It is up to the user's responsibility to adhere to the product usage conditions and to respect the safety instructions for health and environmental protection.

Further information: This product must be stored, sold, and used in accordance with valid hygienic regulations.

Classification according to Regulation (EC) No 1272/2008:	
Skin Sens. 1	Calculation method

First issue of SDS: 09.12.2022

Documents used to prepare SDS:

The original documents provided by the supplier or manufacturer related to the product (mixture), eventually to individual substances contained.

Sulphuric Acid < 5 %**SECTION 1 IDENTIFICATION OF THE SUBSTANCE / MIXTURE AND OF THE COMPANY / UNDERTAKING****1.1 Product identifier**

Trade name: Sulphuric Acid < 5 %

Other name: Sulfuric Acid < 5 %

1.2 Relevant identified uses of the substance or mixture and uses advised against

Stop solution for the ELISA kit.

1.3 Details of the supplier of the safety data sheet

BioVendor - Laboratorní medicína a.s.

Karásek 1767/1

621 00 Brno

Czech Republic

Identification number: 63471507

Tel: +420 549 124 185

E-mail: info@biovendor.com

Website: www.biovendor.com

1.4 Emergency telephone number

European Chemicals Agency. National helpdesks contact details

<https://echa.europa.eu/support/helpdesks>

Links to Poison Centers and Clinical Toxicologists all over the World:

<https://www.eapcct.org/index.php?page=links>

SECTION 2 HAZARDS IDENTIFICATION**2.1 Classification of the substance or mixture**

Classification according to Regulation 1272/2008/EC:

Corrosive to metals, Met. Corr. 1, H290

This material is considered hazardous by the OSHA Hazard Communication Standard (29 CFR 1910.1200). (USA)

2.2 Label elements

Hazard pictogram:



GHS05

Signal word: Warning

Hazard statement:

H290 May be corrosive to metals.

Precautionary statement:

not assigned

2.2.1 Labelling of packages where the contents does not exceed 125 ml:

Hazard pictogram:

none

Signal word:

none

Sulphuric Acid < 5 %**Hazard statement:**

none

Precautionary statement:

none

2.3 Other hazards

This mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1 % or higher.

SECTION 3 COMPOSITION / INFORMATION ON INGREDIENTS**3.1 Mixtures**

Ingredient	Conc. in w/w %	EINECS	CAS-Nr.	Index-Nr.
Sulphuric acid, 98 %	< 5	231-639-5	7664-93-9	016-020-00-8
REACH RN: 01-2119458838-20-xxxx				

Classification according to regulation 1272/2008/EC:

Met. Corr. 1, H290

Skin Corr. 1A, H314

Eye Dam. 1, H318

Specific concentration limits: Skin Corr. 1A, H314: C_{≥15} %; Skin Irrit. 2, H315: 5 % ≤ C < 15 %; Eye Irrit. 2, H319: 5 % ≤ C < 15 %.

For full text of and H-phrases see section 16.

SECTION 4 FIRST AID MEASURES**4.1 Description of first aid measures****General information:**

If feeling unwell or unsure, inform a doctor and show information from this Safety Data Sheet. If unconscious, place the affected person in the recovery position on their side with their head slightly backwards and take care there is no obstruction in their respiratory tract, NEVER induce vomiting. If the affected person is vomiting, ensure that they do not choke on their own vomit.

If inhaled: Immediately stop the exposure, take the affected person to fresh air. Keep the affected person warm. Get medical attention if irritation, shortness of breath or other symptoms persists.

In case of skin contact: Take off contaminated clothing. Wash the affected area with plenty of lukewarm water, if possible. If skin is not irritated, soap, soap solution or shampoo can be used. Get medical attention if irritation persists.

In case of contact with eyes: Rinse the eyes immediately with plenty of running water keeping eyelids open; if the affected person wears contact lenses, remove them immediately. Continue rinsing for minimum 10 minutes. Get medical attention.

If ingested: DO NOT INDUCE VOMITING – even the act of inducing vomiting may cause complications. Rinse mouth, if possible, administer a small dose of activated carbon (1-2 crushed tablets); do not attempt neutralization. Get medical attention.

4.2 Most important symptoms and effects, both acute and delayed

Inhalation: Possible respiratory tract irritation, cough, headache.

Skin contact: Painful reddening, irritation.

Eye contact: May cause irritation.

Ingestion: May cause irritation, sickness.

4.3 Indication of any immediate medical attention and special treatment needed

Immediate medical attention is not necessary.

Sulphuric Acid < 5 %**SECTION 5 FIREFIGHTING MEASURES****5.1 Extinguishing media**

Suitable extinguishing media: Alcohol resistant foam, carbon dioxide, powder, dispersed water, water spray.

Unsuitable extinguishing media: none known.

5.2 Special hazards arising from the substance or mixture

In case of fire sodium oxides can be released.

5.3 Advice for firefighters

Keep the contaminated fire-extinguishing media away from drains, surface and ground water. If necessary, wear a suitable self-contained breathing apparatus and full protective clothing.

SECTION 6 ACCIDENTAL RELEASE MEASURES**6.1 Personal precautions, protective equipment and emergency procedures**

Ensure sufficient ventilation. Wear proper protective equipment, refer to Section 8.

6.2 Environmental precautions

Avoid soil contamination and keep away from surface and ground water. Avoid spills and keep away from drains.

6.3 Methods and material for containment and cleaning up

Cover spills with suitable (non-flammable) absorbent (sand, diatomaceous earth, soil, and other suitable absorbent material); collect in a tightly closed container and dispose of in compliance with Section 13. Dispose of the contained material in accordance with local regulations. After containment, wash the affected area with plenty of water or other suitable detergent.

6.4 Reference to other sections

See section 7 for handling, section 8 for personal protective equipment, section 13 for disposal.

SECTION 7 HANDLING AND STORAGE**7.1 Precautions for safe handling**

Do not inhale gases and vapours. Avoid contact with the skin and the eyes. Wear suitable personal protective equipment according to Section 8. Comply with local and national regulations on occupational health and safety. Avoid long-term and repeated exposure.

7.2 Conditions for safe storage, including any incompatibilities

Keep in tightly closed containers and in cool, dry and well-ventilated designated places. Open packaging must be carefully closed and kept in a vertical position to avoid release or spill. Keep away from alkaline materials.

7.3 Specific end use(s)

See Section 1.2., no other uses are identified.

SECTION 8 EXPOSURE CONTROLS / PERSONAL PROTECTION**8.1 Control parameters**

Czech Republic (Government Decree no. 361/2007 Sb. as amended): Sulphuric acid, SO₃: long term 1 mg/m³, short term 2 mg/m³, note I (causes irritation of mucous membranes (eyes, respiratory tract) and skin irritation).

Sulphuric Acid < 5 %
8.2 Exposure controls
8.2.1 Appropriate engineering controls

Use in a well-ventilated area.

8.2.2 Individual protection measures, such as personal protective equipment

Take usual occupational health and safety measures. If it is impossible to meet exposition limits, suitable respiratory protective equipment must be used. When using, do not eat, drink and smoke. Wash hands carefully with water and soap after work and before breaks.

a) Eye/face protection:


Use protective glasses (according to EN 166).

b) Skin protection:

Hand protection:

Use suitable protective gloves resistant to the product. Follow the recommendations of the gloves manufacturer when selecting the appropriate thickness, material and permeability. In the event of long-term or repeated dermal exposure, apply suitable hand cream to the parts of the skin in direct contact with the mixture. Follow other recommendations of the manufacturer.


Body protection:

Protective cotton clothing. Carefully wash the skin exposed to the product.

c) Respiratory protection:


If necessary, wear a face respirator (according to EN 14387).

d) Thermal hazards: No information available.

8.2.3 Environmental exposure controls

Refer to Sections 6 and 12.

SECTION 9 PHYSICAL AND CHEMICAL PROPERTIES

9.1 Information on basic physical and chemical properties

Physical state:	Liquid
Colour:	Colourless
Odour:	Odourless
Odour threshold:	Not determined
pH:	0.8
Melting point/freezing point (°C):	No data available

Sulphuric Acid < 5 %

Initial boiling point and boiling range (°C):	No data available
Flash point (°C):	No data available
Evaporation rate:	No data available
Flammability:	No data available
Upper/lower flammability or explosive limits:	No data available
Vapour pressure:	No data available
Vapour density:	
Relative density:	No data available
Solubility:	Soluble in water
Partition coefficient n-octanol/water:	
Auto-ignition temperature (°C):	
Decomposition temperature:	
Viscosity:	No data available
Explosive properties:	No data available
Oxidizing properties:	No data available

9.2 Other information

None known.

SECTION 10 STABILITY AND REACTIVITY
10.1 Reactivity

No hazardous reactions occur under common conditions of use.

10.2 Chemical stability

Stable under normal conditions.

10.3 Possibility of hazardous reactions

Reacts with strong bases.

10.4 Conditions to avoid

No data available.

10.5 Incompatible materials

Not known under recommended conditions of use.

10.6 Hazardous decomposition products

They are not produced under normal use. In the event of fire: sulphur oxides.

SECTION 11 TOXICOLOGICAL INFORMATION
11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008

Acute toxicity:	The classification criteria are not met.
Sulphuric acid:	
LD50, oral, rat:	3140 mg/kg
LD50, inhalation, rat:	510 mg/m ³ /2 h
Skin corrosion/ irritation:	The classification criteria are not met. May cause mild irritation.
Serious eye damage/ irritation:	The classification criteria are not met. May cause mild irritation.
Respiratory or skin sensitisation:	The classification criteria are not met.
Germ cell mutagenicity:	The classification criteria are not met.
Carcinogenicity:	The classification criteria are not met.
Reproductive toxicity:	The classification criteria are not met.
STOT – single exposure:	The classification criteria are not met.
STOT – repeated exposure:	The classification criteria are not met.

Sulphuric Acid < 5 %

Aspiration hazard:
 Reproduction toxicity:

The classification criteria are not met.
 No data available.

11.2 Information on other hazards

The substance/mixture does not contain components considered to have endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at levels of 0.1 % or higher.

SECTION 12 ECOLOGICAL INFORMATION
--

12.1 Toxicity

No data available on the product itself.

Sulphuric acid:

LC50, fish, Gambusia affinis, 96 h:	42 mg/l
LC50, crustaceans, Pandalus montagui, 48 h:	42,5 mg/l

12.2 Persistence and degradability

No data available.

12.3 Bioaccumulative potential

No data available.

12.4 Mobility in soil

No data available.

12.5 Results of PBT and vPvB assessment

Does not meet the criteria for PBT and vPvB.

12.6 Endocrine disrupting properties

The substance/mixture does not contain components considered to have endocrine disrupting properties according to REACH Article 57(f) or Commission Delegated regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at levels of 0.1 % or higher.

12.7 Other adverse effects

No data available.

SECTION 13 DISPOSAL CONSIDERATIONS

13.1 Waste treatment methods

Dispose of in accordance with applicable national regulations, to regulate packaging and to amend some statutes (the Packaging Act), as amended and in accordance with implemented regulations on waste disposal.

Appropriate methods of waste treatment of both the substance or the mixture and any contaminated packaging:

The indicated waste, including the waste identification sheet, shall be handed over to a company authorized to treat and dispose of wastes according to the Waste Act and that the company producing waste has entered into a contract with.

Both completely empty and not completely empty packaging shall be placed in designated containers for waste collection and the indicated waste shall be handed over for disposal to a person authorized to dispose of the waste.

Waste type code	060101
Waste type	sulphuric acid and sulphurous acid *
Waste subgroup	waste from production, processing, distribution and use (VZDP) of acids
Waste group	WASTE CREATED BY ANORGANIC CHEMICAL PROCESSES

Waste type code for packaging 150110

Sulphuric Acid < 5 %

Waste type	packaging containing residues of dangerous substances or packaging contaminated with these substances*
Waste subgroup	Packaging (including separately collected household packaging material)
Waste group	WASTE PACKAGING; ABSORPTION AGENTS; CLEANING CLOTHS; FILTRATION MATERIALS AND PROTECTIVE CLOTHING NOT SPECIFIED IN ANY OTHER MANNER

(*) – hazardous waste according to Directive No. 91/689/EEC on hazardous waste

SECTION 14 TRANSPORT INFORMATION**14.1 UN number or ID number**

3264.

14.2 UN proper shipping name

CORROSIVE LIQUID, ACIDIC, INORGANIC, N.O.S. (sulphuric acid).

14.3 Transport hazard class

8.

14.4 Packing group

III.

14.5 Environmental hazards

No.

14.6 Special precautions for user

No data available.

14.7 Maritime transport in bulk according to IMO instruments

No data available.

SECTION 15 REGULATORY INFORMATION**15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**

Act No. 350/2011 Coll., to regulate chemical substances and chemical mixtures and to amend some statutes, as amended. Implemented regulations to Act No. 350/2011 Coll., as amended The Waste Act as amended. Government Decree No. 361/2007 Coll., to regulate the conditions of occupational health and safety, as amended. Regulation of the European Parliament and the Council (EC) No. 1907/2006 (REACH). Regulation of the European Parliament and the Council (EC) No. 1272/2008 (CLP). Commission Regulation (EU) No. 878/2020.

U.S. Federal regulations:

TSCA 8(a) CDR Exempt/Partial exemption: Not determined
United States inventory (TSCA 8b): All components are listed or exempted.
Clean Water Act (CWA) 311: sulphuric acid.

15.2 Chemical safety assessment

Chemical safety assessment has not been carried out.

SECTION 16 OTHER INFORMATION**Warning:**

The safety data sheet contains data needed for securing safety and health protection during work and environmental protection. The stated data correspond to the current state of knowledge and

Sulphuric Acid < 5 %

experience and is in accordance with valid legal regulations. It cannot be deemed as a guarantee of the properties, suitability, and usefulness of the product for specific application and therefore no contractual legal relationships are hereby created.

The safety data sheet is the property of the physical or legal entity stated in Section 1 and is protected by copy-right. All copying, distribution or sales without the consent of the owner is forbidden.

Relevant phrases:

- H290 May be corrosive to metals.
- H314 Causes severe skin burns and eye damage.
- H315 Causes skin irritation.
- H318 Causes serious eye damage.
- H319 Causes serious eye irritation.

Hazardous Material Information System (U.S.A.)

Health	1
Chronic Health Hazard	
Flammability	0
Physical Hazards	0

National Fire Protection Association (U.S.A.)

Health	1
Flammability	0
Instability/Reactivity	0
Special	

Training hints:

Pursuant to article No 35 of the European Parliament and Council Regulation (ES) No 1907/2006, the employer must allow employees or their representatives access to information from the safety data sheet of the substance or preparation, which the employees use or to the effects of which they may be exposed during their work.

Physical entities performed individual activities within the scope of handling of this hazardous product are trained and regularly, at least once a year, retrained.

Product information sources: safety data sheet, product or technical information, safety instructions, and other expert documents for the product, issued by the supplier.

Recommended restriction of use:

The product is designed only for professional purposes. It must not be used in households. The product can only be handled by a person older than 18 years, who is sufficiently informed about the work procedures, hazardous properties of the product, and also about the necessary safety measures.

The product is to be used only for the purpose, for which it is designed. It is up to the user's responsibility to adhere to the product usage conditions and to respect the safety instructions for health and environmental protection.

Further information: This product must be stored, sold, and used in accordance with valid hygienic regulations.

Date of issue: 30.5.2022

Revision notes: Commission Regulation (EU) No. 830/2015 was replaced by Commission Regulation (EU) No. 878/2022. Change of classification of the mixture. Exception for labelling of packaging not exceeding 125 ml, addition of REACH RN of the component, change of section 11 structure, addition of transport information in section 14.