

PRODUCT DATA SHEET

FAM3B Human, Rabbit Polyclonal Antibody

Cat. No.: RD181259100**Size:** 0.1 mg**Source of Antigen:** E. coli**Type:** Polyclonal Antibody**Host:** Rabbit**Other Names:**

PANDER, Cytokine-like protein 2-21, Pancreatic-derived factor, C21orf11, C21orf76, PRED44, UNQ320/PRO365

Preparation:

The antibody was raised in rabbits by immunization with the recombinant Human FAM3B.

Amino Acid Sequence of Immunogen:

Total 216 AA. MW: 24,10 kDa (calculated). UniProtKB acc.no. P58499. N-Terminal His-tag (10 extra AA)

MKHHHHHHAS ELIPDAPLSS AAYSIRSIGE RPVLKAPVPK RQKCDHWTPC PSDTYAYRLL
SGGGRSKYAK ICFEDNLLMG EQLGNVARGI NIAIVNYVTG NVTATRCFDM YEGDNSGPMT
KFIQSAAPKS LLFMVTYDDG STRLNNDAKN AIEALGSKEI RNMKFRSSWV FIAAKGLELP
SEIQREKINH SDAKNNRYSG WPAEQIEGC IPKERS**Purification Method:**

Immunoaffinity chromatography on a column with immobilized recombinant Human FAM3B.

Species Reactivity:

Human. Not yet tested in other species.

Antibody Content:

0.1 mg (determined by BCA method, BSA was used as a standard)

Formulation:

The antibody is lyophilized in 0.05 M phosphate buffer, 0.1 M NaCl, pH 7.2.

Reconstitution:

Add 0.2 ml of deionized water and let the lyophilized pellet dissolve completely. Slight turbidity may occur after reconstitution, which does not affect activity of the antibody. In this case clarify the solution by centrifugation.

Shipping:

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

Storage/Stability:

The lyophilized antibody remains stable and fully active until the expiry date when stored at -20°C. Aliquot the product after reconstitution to avoid repeated freezing/thawing cycles and store frozen at -80°C. Reconstituted antibody can be stored at 4°C for a limited period of time; it does not show decline in activity after one week at 4°C.

Quality Control:

Indirect ELISA – to determine titer of the antibody

SDS PAGE – to determine purity of the antibody

BCA - to determine quantity of the antibody

Note:

This product is for research use only.