

## PRODUCT DATA SHEET

### S100A6 Human, Rabbit Polyclonal Antibody

**Cat. No.:** RD181215100

**Size:** 0.1 mg

**Source of Antigen:** E. coli

**Type:** Polyclonal Antibody

**Host:** Rabbit

**Other Names:**

S100 calcium-binding protein A6, Calcyclin, Prolactin receptor-associated protein, PRA, Growth factor-inducible protein 2A9, MLN 4, CACY

**Preparation:**

The antibody was raised in rabbits by immunization with the recombinant Human S100-A6.

**Amino Acid Sequence of Immunogen:**

The immunization antigen (11.42 kDa) is a protein containing 100 AA of recombinant Human S100-A6 and N-Terminal His-tag 10 extra AA (highlighted).

MKHHHHHHAS MACPLDQAIG LLVAIFHKYS GREGDKHTLS KKELKELIQK ELTIGSKLQD  
AEIARLMEDL DRNKDQEVNF QEYVTFLGAL ALIYNEALKG

**Purification Method:**

Immunoaffinity chromatography on a column with immobilized recombinant Human S100-A6.

**Species Reactivity:**

Human. Not yet tested in other species.

**Antibody Content:**

0.1 mg (determined by BCA method, BSA was used as a standard)

**Formulation:**

The antibody is lyophilized in 0.05 M phosphate buffer, 0.1 M NaCl, pH 7.2. **\*\*AZIDE FREE\*\***.

**Reconstitution:**

Add 0.2 ml of deionized water and let the lyophilized pellet dissolve completely. Slight turbidity may occur after reconstitution, which does not affect activity of the antibody. In this case clarify the solution by centrifugation.

**Shipping:**

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

**Storage/Stability:**

The lyophilized antibody remains stable and fully active until the expiry date when stored at  $-20^{\circ}\text{C}$ . Aliquot the product after reconstitution to avoid repeated freezing/thawing cycles and store frozen at  $-80^{\circ}\text{C}$ . Reconstituted antibody can be stored at  $4^{\circ}\text{C}$  for a limited period of time; it does not show decline in activity after one week at  $4^{\circ}\text{C}$ .

**Quality Control:**

Indirect ELISA – to determine titer of the antibody

SDS PAGE – to determine purity of the antibody

BCA - to determine quantity of the antibody

**Applications:**

ELISA, Western blotting

**Note:**

This product is for research use only.