

PRODUCT DATA SHEET

Visinin-Like Protein 1 Human, Rabbit Polyclonal Antibody

Cat. No.: RD181119100

Size: 0.1 mg

Source of Antigen: E. coli

Type: Polyclonal Antibody

Host: Rabbit

Other Names:

Visinin-Like Protein 1, Hippocalcin-like protein 3, HLP3, VILIP, VLP-1, VSNL1, VISL1

Preparation:

The antibody was raised in rabbits by immunization with the recombinant Human Visinin-like protein 1.

Amino Acid Sequence of Immunogen:

Recombinant Human Visinin-like protein 1, total 201 AA. MW: 23.4 kDa (calculated).
N-terminal His-tag, 10 extra AA (highlighted). The AA sequence is identical to
UniProtKB/Swiss-Prot entry P62760.

MKHHHHHHAS MGKQNSKLAP EVMEDLVKST EFNEHELKQW YKGFLKDCPS
GRLNLEEFQQ LYVKFFPYGD ASKFAQHAFR TFDKNGDGTI DFREFICALS ITSRGSFEQK
LNWAFNMYDL DGDGKITRVE MLEIIEAIYK MVGTVIMMKM NEDGLTPEQR VDKIFSKMDK
NKDDQITLDE FKEAAKSDPS IVLLLQCDIQ K

Purification Method:

Immunoaffinity chromatography on a column with immobilized recombinant Human Visinin-like protein 1.

Species Reactivity:

Human. Not yet tested in other species.

Antibody Content:

0.1 mg (determined by BCA method, BSA was used as a standard)

Formulation:

The antibody is lyophilized in 0.05 M phosphate buffer, 0.1 M NaCl, pH 7.2. ****AZIDE FREE****.

Reconstitution:

Add 0.2 ml of deionized water and let the lyophilized pellet dissolve completely. Slight turbidity may occur after reconstitution, which does not affect activity of the antibody. In this case clarify the solution by centrifugation.

Shipping:

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

Storage/Stability:

The lyophilized antibody remains stable and fully active until the expiry date when stored at -20°C . Aliquot the product after reconstitution to avoid repeated freezing/thawing cycles and store frozen at -80°C . Reconstituted antibody can be stored at 4°C for a limited period of time; it does not show decline in activity after one week at 4°C .

Quality Control:

Indirect ELISA – to determine titer of the antibody

SDS PAGE – to determine purity of the antibody

BCA - to determine quantity of the antibody

Applications:

ELISA, Western blotting

Note:

This product is for research use only.