

PRODUCT DATA SHEET

Uromodulin Human, Rabbit Polyclonal Antibody

Cat. No.: RD181163100

Source of Antigen: Human urine

Host: Rabbit

Size: 0.1 mg

Type: Polyclonal Antibody

Other Names:

Tamm-Horsfall urinary glycoprotein, THP, UMOD

Preparation:

The antibody was raised in rabbits by immunization with the native Human Uromodulin (urine).

Purification Method:

Immunoaffinity chromatography on a column with immobilized native Human Uromodulin (urine).

Species Reactivity:

Human. Not yet tested in other species.

Antibody Content:

0.1 mg (determined by BCA method, BSA was used as a standard)

Formulation:

The antibody is lyophilized in 0.05 M phosphate buffer, 0.1 M NaCl, pH 7.2.

Reconstitution:

Add 0.2 ml of deionized water and let the lyophilized pellet dissolve completely. Slight turbidity may occur after reconstitution, which does not affect activity of the antibody. In this case clarify the solution by centrifugation.

Shipping:

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

Storage/Stability:

The lyophilized antibody remains stable and fully active until the expiry date when stored at -20°C. Aliquot the product after reconstitution to avoid repeated freezing/thawing cycles and store frozen at -80°C. Reconstituted antibody can be stored at 4°C for a limited period of time; it does not show decline in activity after one week at 4°C.

Quality Control:

Indirect ELISA – to determine titer of the antibody

SDS PAGE – to determine purity of the antibody

BCA - to determine quantity of the antibody

Applications:

ELISA, Immunohistochemistry, Immunoprecipitation, Western blotting

Note:

This product is intended for research use only.