

PRODUCT DATA SHEET

Fibroblast Growth Factor 21 Human, Sheep Polyclonal Antibody

Cat. No.: RD184108100

Source of Antigen: E. coli

Host: Sheep

Size: 0.1 mg

Type: Polyclonal Antibody

Other Names:

FGF-21, UNQ3115/PRO10196

Preparation:

The antibody was raised in sheep by immunization with the recombinant Human FGF-21.

Amino Acid Sequence of Immunogen:

The immunization antigen (20.1 kDa) is a protein containing 195 AA of recombinant Human FGF-21. N-terminal His-tag, 14 extra AA.

MRGSHHHHHH GMASHPIPDS SPLIQFGGQV RQRYLYTDDA QQTEAHLEIR EDGTVGGAAD
QSPESLLQLK ALKPGVIQIL GVKTSRFLCQ RPDGALYGSL HFDPEACSFH ELLLEDGYNV
YQSEAHGLPL HLPGNKSPHR DPAPRGPARF LPLPGLPPAP PEPPGILAPQ PPDVGSSDPL
SMVGPSQGRS PSYAS

Purification Method:

Immunoaffinity chromatography on a column with immobilized recombinant Human FGF-21.

Species Reactivity:

Human. Not yet tested in other species.

Antibody Content:

0.1 mg (determined by BCA method, BSA was used as a standard)

Formulation:

The antibody is lyophilized in 0.05 M phosphate buffer, 0.1 M NaCl, pH 7.2.

Reconstitution:

Add 0.2 ml of deionized water and let the lyophilized pellet dissolve completely. Slight turbidity may occur after reconstitution, which does not affect activity of the antibody. In this case clarify the solution by centrifugation.

Shipping:

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

Storage/Stability:

The lyophilized antibody remains stable and fully active until the expiry date when stored at -20°C. Aliquot the product after reconstitution to avoid repeated freezing/thawing cycles and store frozen at -80°C. Reconstituted antibody can be stored at 4°C for a limited period of time; it does not show decline in activity after one week at 4°C.

Quality Control:

Indirect ELISA – to determine titer of the antibody

SDS PAGE – to determine purity of the antibody

BCA - to determine quantity of the antibody

Applications:

ELISA, Immunohistochemistry, Western blotting

Note:

This product is for research use only.