

PRODUCT DATA SHEET

L-Asparaginase

Cat. No.: RP1792872500

Size: 2500 IU

Source: E. coli

Other Names:

Description:

L-Asparaginase was purified from E. coli ASI.357

Purity:

Greater than 96.0% as determined by:
(a) Analysis by RP-HPLC.

Formulation:

The enzyme was lyophilized with no additives.

Reconstitution:

It is recommended to reconstitute the lyophilized L-Asparaginase in 18MΩ water at 1 mg/ml.

Shipping:

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

Storage/Stability:

Two years when stored at -20°C, 2 weeks at 4°C.

Note:

BioVendor's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.