

PRODUCT DATA SHEET

Fibronectin Human, E.coli Recombinant

Cat. No.: RP1722621100

Size: 100 µg

Source: E. coli

Type: Recombinant protein

Species: Human

Other Names:

Cold insoluble globulin, FINC, LETS, MSF, FN

Description:

Fibronectin Human Recombinant produced in E. coli is a single, non-glycosylated polypeptide chain containing 574 amino acids and having a molecular mass of 62.6kDa. The Fibronectin is purified by proprietary chromatographic techniques.

Fibronectin is useful for the induction of cell attachment to a variety of surfaces including plastic and glass tissue culture labware, petri dishes, coverslips, microcarrier beads, etc. Fibronectin is useful for growth and maintenance of cells in low serum conditions. In general, the reconstituted fibronectin should be diluted with sterile physiological saline or serum-free medium to a concentration of 10-50µg/ml.

Purity:

Greater than 95.0%

Formulation:

Lyophilized from a 0.2µm filtered concentrated solution in 20 mM Tris-HCl, pH 8.0, 150 mM NaCl, with 5 % Trehalose and 0.02 % Tween-20.

Reconstitution:

It is recommended to reconstitute the lyophilized Fibronectin in sterile 18MΩ-cm H₂O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.

Shipping:

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

Storage/Stability:

Lyophilized Fibronectin although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Fibronectin should be stored at 4°C between 2-7 days and for future use below -18°C.

Please prevent freeze-thaw cycles.

Applications:

Cell culture and/or animal studies, In vitro

Note:

This BioVendor product is furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.