



**NGS**



## **fastGEN**

**A breakthrough system  
for a one-step NGS  
library preparation**




**GATGCTAC  
AGATGGTA  
GTTCATA  
CCGATGGT  
GATGCTAC**

**010000110  
000100110  
100010010  
1110010110  
000100110**

# It's faster. Way faster!

The fastGEN technology is based on a rapid and easy-handling workflow that significantly lowers the probability of mistakes.



# 4 reasons why fastGEN worths attention

**BioVendor Group developed user-friendly diagnostic NGS kits with excellent analytical parameters and high-speed processing.**

## Unique

Simple, ultra-sensitive,  
specific and effective technology

## Reliable

Significantly lower probability of mistakes,  
uniform coverage for different amplicons  
and different samples

## Fast

Excellent speed (<30 min hands-on time)

## Robust


Robustness achieved by using short amplicons  
obtained by a single polymerase chain reaction  
with special tagged hybrid primers

# GENOVESA

## Cloud-based database bioinformatic and interpretation system

### Key features

- User-friendliness
- Automatized bioinformatic analysis
- Advanced quality control of sequencing data
- Easy variants filtration
- Local clinical variants database
- Visualization of NGS data
- Internal interpretation of variants and sharing between clinics
- Clinical report generation
- Custom baseline creation



Find out more about  
the software on our [website](#).



# Quick and easy detection of mutations in genes

*NRAS | KRAS | BRAF*

*PIK3CA*

*TERT*

*TP53*

*BCR::ABL1*

*POLE | CTNNB1*

*H3F3A | IDH1 | IDH2*

*MSI*

*EGFR | HER2*

*JAK2 | CALR | MPL*

*ALDOB | AOC1 | MCM6*



# Contact our NGS team



**Iveta Tóthová, Ph.D.**

Product & Scientific Manager  
tothova@biovendor-mdx.com



**Anna Kolouchová, Ph.D.**

Business Development Manager  
kolouchova@biovendor-mdx.com



**MSc. Adam Novotný**

Application Specialist  
novotny@biovendor-mdx.com



**MSc. Tereza Hubinková**

Product Specialist  
hubinkova@biovendor-mdx.com



**Petr Fajkus, Ph.D.**

Application & Product Specialist GENOVESA  
fajkus@biovendor-mdx.com

PRODUCER:



DISTRIBUTOR:



**BioVendor MDx a.s.**

Karásek 1767/1, 621 00 Brno

Czech Republic

info@biovendor-mdx.com

www.biovendor.com