

## PRODUCT DATA SHEET

### Mouse Anti Respiratory Syncytial Virus Antibody, Clone: RV3

**Cat. No.:** RM812147100

**Size:** 0.1 mg

**Type:** Monoclonal Antibody

**Host:** Mouse

**Isotype:** IgG1

#### **Purification Method:**

The antibody was purified by affinity chromatography on protein A.

#### **Species Reactivity:**

This antibody is specific for respiratory syncytial virus. The antibody does not cross react with Influenza Type A, Influenza Type B or Adenovirus.

#### **Antibody Content:**

0.1 mg; concentration: 1.0 mg/ml

For actual concentration see the Certificate of Analysis.

#### **Formulation:**

Supplied as a liquid in Phosphate Buffered Saline pH 7.2, < 0.1% Sodium Azide.

#### **Shipping:**

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

#### **Storage/Stability:**

Short term storage at 2–8°C.

Longer term storage at -20°C.

Avoid repeated freeze/thaw cycles.

#### **Applications:**

ELISA

#### **Note:**

For research use only.