

PRODUCT DATA SHEET

Mouse Anti Respiratory Syncytial Virus Antibody, Clone: RV11

Cat. No.: RM812398100

Size: 0.1 mg

Type: Monoclonal Antibody

Host: Mouse

Isotype: IgG1

Purification Method:

Purified in a multi-step procedure to $\geq 95\%$ purity by SDS-PAGE.

Species Reactivity:

This antibody is specific for respiratory syncytial virus. The antibody does not cross react with Influenza Type A, Influenza Type B or Adenovirus.

Antibody Content:

0.1 mg; concentration: 1.0 mg/ml

For actual concentration see the Certificate of Analysis.

Formulation:

Supplied as a liquid in Phosphate Buffered Saline pH 7.4, 0.09% Sodium Azide.

Shipping:

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

Storage/Stability:

Short term storage at 2–8°C.

Longer term storage at -20°C.

Avoid repeated freeze/thaw cycles.

Applications:

ELISA, LFA

Note:

For research use only.