

## PRODUCT DATA SHEET

### Mouse Anti Influenza B Antibody, Clone: YB91

**Cat. No.:** RM812164100

**Size:** 0.1 mg

**Type:** Monoclonal Antibody

**Host:** Mouse

**Isotype:** IgG2

#### **Purification Method:**

The antibody was purified by affinity chromatography on protein A.

#### **Species Reactivity:**

This antibody is specific for Influenza B detection. The antibody does not cross react with Adenovirus or Respiratory Syncytial virus.

#### **Antibody Content:**

0.1 mg; concentration: 1.0 mg/ml

For actual concentration see the Certificate of Analysis.

#### **Formulation:**

Supplied as a liquid in Phosphate Buffered Saline pH 7.4, 0.09% Sodium Azide.

#### **Shipping:**

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

#### **Storage/Stability:**

Store at +4°C for up to three months, or at -20°C for longer. Avoid repeated freeze/thaw cycles.

#### **Applications:**

ELISA

#### **Note:**

For research use only.