

PRODUCT DATA SHEET

Mouse Anti Rhinovirus VP3 Antibody, Clone: 4967

Cat. No.: RM812384500

Size: 0.5 mg

Type: Monoclonal Antibody

Host: Mouse

Isotype: IgG1

Purification Method:

Purified preparations consist of >90% pure mouse monoclonal antibody which has been purified from ascites fluid or culture medium by protein A chromatography or sequential differential precipitations.

Species Reactivity:

Specific for the VP3 of numerous Rhino types (see table below). Not reactive with other Enteroviruses, Influenza A, Influenza B, RSV, Adenovirus, Para 1, Para 2, Para 3, Chlamydia pneumoniae, Mycoplasma pneumoniae & Hepatitis A virus.

Antibody Content:

0.5 mg; concentration: 1.0 mg/ml

Formulation:

Supplied as a liquid in Phosphate Buffered Saline pH 7.2, 0.1% Sodium Azide.

Shipping:

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

Storage/Stability:

This product contains no stabilizing proteins and should be stored at 2–8°C until ready for use.

Applications:

ELISA, IFA. Working dilution must be determined by the user. Suggested starting ranges are 1:10-1:50 for IFA and 1:20-1:200.

Note:

For research use only.