

Activin A Human Rat Murine E. coli (E.coli-derived)

Product Data Sheet

Type: Recombinant	Cat. No.:	
Source: E. coli	RBG10052E002	(2 µg)
Species: Murine, Rat, Human	RBG10052E010	(10 µg)
Other names: Inhibin beta-1, FRP, FSH (Follicle-stimulating hormone)-releasing protein	RBG10052E100	(100 µg)

Description

Activin A is a TGF-beta family member that exhibits a wide range of biological activities, including regulation of cellular proliferation and differentiation, and promotion of neuronal survival. Elevated levels of Activin A in human colorectal tumors and in postmenopausal women have been implicated in colorectal and breast cancers, respectively. The biological activities of Activin A can be neutralized by inhibins and by the diffusible TGF-beta antagonist, follistatin. Activin A binds to the two forms of activin receptor type I (Act RI-A and Act RI-B) and two forms of activin receptor type II (Act RII-A and Act RII-B). Activins are homodimers or heterodimers of different beta subunits. They are produced as precursor proteins with an amino terminal propeptide that is cleaved to release the C-terminal bioactive ligand. Human/Murine/Rat Activin A is a 26.0 kDa disulfide-linked homodimer of two betaA chains, each containing 116 amino acid residues.

Research topic

Apoptosis, Cardiovascular disease, Neural tissue markers, Oncology

Amino Acid Sequence

GLECDGKVINI CCKKQFFVSF KDIGWNDWII APSGYHANYC EGECPSHIAG TSGSLSLFHS TVINHYMRG HSPFANLKSC
CVPTKLRPMS MLYYDDGQNI IKKDIQNMIV EECGCS

Source

E. coli

Purity

97

Biological Activity

The ED50 as determined by its ability to inhibit the proliferation of murine MPC-11 cells is 2.0 ng/ml, corresponding to a specific activity of 5×10^5 units/mg.

Storage, Stability/Shelf Life

- 20 °C

HEADQUARTERS: BioVendor Laboratorní medicína, a.s.	Karasek 1767/1	621 00 Brno CZECH REPUBLIC	Phone: +420-549-124-185 Fax: +420-549-211-460	E-mail: info@biovendor.com sales@biovendor.com Web: www.biovendor.com
AUSTRIA: BioVendor GesmbH	Nußdorfer Straße 20/10	1090 Vienna AUSTRIA	Phone: +43-1-89090-25 Fax: +43-1-89090-2515	E-mail: infoAustria@biovendor.com
GERMANY, SWITZERLAND: BioVendor GmbH	Otto-Hahn-Straße 16	34123 Kassel GERMANY	Phone: +49-6221-433-9100 Fax: +49-6221-433-9111	E-mail: infoEU@biovendor.com
USA, CANADA AND MEXICO: BioVendor LLC	128 Bingham Rd. Suite 1300	Asheville, NC 28806 USA	Phone: +1-828-575-9250 +1-800-404-7807 Fax: +1-828-575-9251	E-mail: infoUSA@biovendor.com