

Prolactin Receptor Human E. coli

Product Data Sheet

Type: Recombinant
Source: E. coli
Species: Human

Cat. No.:
RP1764690002 (2 µg)
RP1764690005 (5 µg)
RP1764690010 (10 µg)

Description

Cytoplasmatic Prolactin Receptor Human Recombinant produced in E.coli is a non-glycosylated, Polypeptide chain amino acids 432-623 and having a molecular mass of 45 kDa, the PRLR is fused with a GST tag. The Prolactin Receptor is purified by proprietary chromatographic techniques.

Source

E. coli

Biological Activity

Not Known.

Formulation

PRLR Human at 0.1 mg/ml in 50 mM Tris-Acetate, pH-7.5, 1 mM EDTA and 20% Glycerol.

Shipping

At ambient temperature. Upon receipt, store the product at the temperature recommended below.

Storage, Stability/Shelf Life

Store vial at -20°C to -80°C. When stored at the recommended temperature, this protein is stable for 12 months **Please prevent freeze-thaw cycles.**

Applications

ELISA, Inhibition Assays, Western blotting

Physical Appearance

Sterile filtered liquid.

HEADQUARTERS: BioVendor Laboratorní medicina, a.s.	Karasek 1767/1	621 00 Brno CZECH REPUBLIC	Phone: +420-549-124-185 Fax: +420-549-211-460	E-mail: info@biovendor.com sales@biovendor.com Web: www.biovendor.com
AUSTRIA: BioVendor GesmbH	Nußdorfer Straße 20/10	1090 Vienna AUSTRIA	Phone: +43-1-89090-25 Fax: +43-1-89090-2515	E-mail: infoAustria@biovendor.com
GERMANY, SWITZERLAND: BioVendor GmbH	Otto-Hahn-Straße 16	34123 Kassel GERMANY	Phone: +49-6221-433-9100 Fax: +49-6221-433-9111	E-mail: infoEU@biovendor.com
USA, CANADA AND MEXICO: BioVendor LLC	128 Bingham Rd. Suite 1300	Asheville, NC 28806 USA	Phone: +1-828-575-9250 +1-800-404-7807 Fax: +1-828-575-9251	E-mail: infoUSA@biovendor.com