

BioVendor new products

August 31, 2015

Dear customer, we would like to introduce our new products and hope you will find them interesting. Below is a list of events, in which we plan to participate in 2015:

4th European Congress of Immunology (06 \ 09 \ 2015 - 09 \ 09 \ 2015 - Vienna)

14th International Conference on Lyme Borreliosis and other Tick-Borne Diseases (27 \ 09 \ 2015 - 30 \ 09 \ 2015 - Vienna)

MEDICA 2015 (16 \ 11 \ 2013 - 19 \ 11 \ 2015 - Düsseldorf)

XIX. Lipid Meeting Leipzig (10 \ 12 \ 2015 - 12 \ 12 \ 2015 - Leipzig)

ENDOTOXIN DETECTION PRODUCTS

EndoTrap[®], affinity chromatography for highly efficient endotoxin removal, manufactured by Hyglos GmbH.

EndoTrap[®] reliably removes endotoxins (lipopolysaccharides / LPS) from aqueous solutions of e.g. proteins, antibodies, vaccines, nucleic acids, buffers and various other substances with excellent sample recovery rates.

» FEATURES

- Highly specific and efficient LPS removal from complex mixtures
- Quick and easy to use in column (gravity flow) or batch mode
- Upscaling on fully automated liquid-chromatographic systems possible
- Average sample/protein recovery > 95%
- Broad pH tolerance and broad pI range of target proteins
- Compatible with commonly used buffers and high salt conditions
- Reusable at least 10 times (EndoTrap[®] HD)
- Not Polymyxin B based, non toxic

The EndoTrap[®] product family consists of two members:

EndoTrap[®] red

especially for the use with PBS based samples.

EndoTrap[®] HD

is a further development of EndoTrap EndoTrap[®] blue for challenging samples and large scale endotoxin removal in e.g. biopharmaceutical production processes. For EndoTrap[®] HD additionally is available the EndoTrap[®] leakage ELISA.

CAT. NO.	NAME	SIZE
311053	EndoTrap [®] blue 1/1	1x1 ml column
311063	EndoTrap [®] blue 5/1	5x1 ml column
311064	EndoTrap [®] blue 10/1	10 ml resin
311075	EndoTrap [®] blue 50/1	50 ml resin
321053	EndoTrap [®] red 1/1	1x1 ml column
321063	EndoTrap [®] red 5/1	5x1 ml column
321064	EndoTrap [®] red 10	10 ml resin
321075	EndoTrap [®] red 50	50 ml resin
800001	EndoTrap [®] HD 5x Regeneration Buffer 500 ml (equals 2.5 L buffer)	500 ml
800033	EndoTrap [®] Leakage ELISA	96-wells
800034	EndoTrap [®] HD 10	10 ml
800035	EndoTrap [®] HD 50	50 ml
800036	EndoTrap [®] HD 250	250 ml
800037	EndoTrap [®] HD Bulk	> 1 l

CAT. NO.	NAME	SIZE
800044	EndoTrap® HD 5x Equilibration Buffer 500 ml (equals 2.5 L buffer)	500 ml
800046	EndoTrap® HD 5x Storage 500 ml (equals 2.5 L buffer)	500 ml
800053	EndoTrap HD 1/1	1x1 ml column
800063	EndoTrap HD 5/1	5x1 ml column
311067	EndoTrap® blue Regeneration Buffer	125 ml
311108	EndoTrap® blue Equilibration Buffer	125 ml
321080	EndoTrap® red Regeneration Buffer	125 ml
321081	EndoTrap® red Equilibration Buffer	125 ml
321082	EndoTrap® red Low Salt Equilibration Buffer	125 ml
800003	EndoTrap® red Regeneration Buffer, 10X	1000 ml
800004	EndoTrap® blue Regeneration Buffer, 10X	1000 ml
800006	Disposable Polypropylene Funnel	5 pack
800007	Empty Columns, Large	5 pack
800008	Empty Columns, Small	5 pack

ENDOTOXIN DETECTION PRODUCTS

Highly specific and robust detection of endotoxin (LPS) helps you make confident in process control (IPC) and product release decisions:

- Selecting the ideal endotoxin (LPS) test method for your samples
- Removing inhibitory effects due to buffer and sample components
- Improving reproducibility of results
- Optimizing endotoxin recovery

CAT. NO.	NAME	SIZE
609050	EndoZyme® Recombinant Factor C Assay	192 tests
609033	EndoLISA® Endotoxin Detection Assay based on ELISA-technology	192 tests
609065	Endo-RS® (Endotoxin Recovery Kit)	1 kit
137101	Plate for EndoZyme assay	1 plate

BioVendor – Laboratorni medicina a.s.
 Phone: +420 549 124 185
 Fax: +420 549 211 460
 E-mail: info@biovendor.com
www.biovendor.com

Distributed by:



BioVendor
 Research
 and Diagnostic Products